

"PHASE-5-COMBO" from The TENS Experts!

TENS Program	P01	P02	P03	P04	P05	P06	P07	P08	P09
Type >	Conventional	Burst	Modulated	"HAN"	Lo-R + Hi-W	Hi-Rate	Mid-Range	Lo-Rate	Phases = 2
Rate	80Hz	2Hz	70-100Hz	100-2Hz	10Hz	110Hz	50Hz	2Hz	100/2 Hz
Width	180uS	150uS	200-100uS	200-350uS	250uS	175uS	100uS	175uS	100/200uS
Tx Time	30 Mins	30 Mins	30 Mins	30 Mins	30 Mins	30 Mins	30 Mins	30 Mins	30/15 Min

** Tx time adjustable in Single Phase programs -- See operating manual for more program details.*

** Program slots P17, 18, 19 available for custom regimens!*

NMES Program	P10	P11	P12	P13	P14	P15	P16
Type >	Warm Up	Endurance	Strength	Rehab. Stim	Resistance	Toning	Relaxation
Rate	5Hz	35Hz	50Hz	35Hz	5-45Hz	5-75Hz	2-15Hz
Width	300u/sec	300u/sec	300u/sec	450u/sec			
Ramp	0 / 0	2 / 1.5	1 / 1	1 / 1	Phases = 5	Phases = 5	Phases = 3
On/Off	Cont.	10/20	10/10	3/15			
Tx Time	5 Mins	40 Mins	20 Mins	20 Mins	33 Mins	16 Mins	55 Mins



*** Contact JOHN REEVES : 613.800.2525**

Regular Kit comes complete with both small & large size Electrodes

or TENSexperts@gmail.com