



Analysis of Preservatives (1): Paraoxybenzoic Acid Esters

Preservatives are added to maintain the freshness of food products by suppressing bacterial growth in food products and prevent rotting¹⁾ In this study, paraoxybenzoic acid esters (known as parabens) were separated by using a phenyl-hexyl column and detected by an UV detector and Chromaster 5610 MS Detector, and the results are presented here. Isomers have the same m/z, but the retention times are different due to the differences in their structures. The confirmation using both the chromatogram and mass spectrum increases the reliability for the component identification.



5610 MS Detector

¹⁾ Bunseki, Masakazyu Horie p.124 (2009)

Analysis of Paraoxybenzoic Acid Ester by LC-MS

Analytical Conditions

Table 1 MS Detector Settings

Ionization method	ESI
Ionization mode	Negative
Ionization voltage	2000 V
Measurement mode	Scan (m/z 100-200)

Table 2 Conditions for HPLC Analysis

Column	Xselect® CSH Phenyl-Hexyl (5 μm) 2.1 mm I.D. x 150 mm (Waters)
Mobile phase	CH ₃ CN / 10 mmol/L Acetic acid buffer (pH4.0) = 1.8 / 1
Flow rate	0.2 mL/min (split ratio = 1:50)
Injection vol.	5 μL (100 mg/L each)
Detection wavelength	230 nm

LC-MS Analysis

Highly reliable analysis is possible by confirming the retention time and mass spectrum of the target component.




Figure 1 UV Chromatogram (Top) and Mass Spectra (Bottom) of 6 Paraoxybenzoic Acid Esters

* Methylparaben is not approved as a food additive in Japan.

The data introduced here were provided by Kita-ku Public Health Center, Tokyo.

<Main system configuration> Chromaster 5110 Pump, 5210 Autosampler, 5310 Column Oven, 5410 UV Detector, 5610 MS Detector

NOTE: These data are an example of measurement; the individual values cannot be guaranteed.