

RECOMMENDATIONS

Indications for Use:

For the adjunctive treatment of musculoskeletal pain

Note: Treatment effects of device use were clinically assessed for up to 4 weeks. Pain relief results may vary for each user. Always read the directions. Use only as directed. If symptoms persist see your healthcare professional

Any serious incident that has occurred in relation to the device should be reported to the manufacturer.


Recommended Treatment Duration (Use Time):

Use the device for a minimum of 12 hours per day, up to 24 hours per day.

HOW TO TURN DEVICE ON & OFF

How to turn the Device On:

Step 1: To activate the device, remove the white tab from the back of the device and push the silver on/off button for 1-2 seconds. Release the button.




Step 2: Once the device is activated, the green LED light on the front of the device will turn on. If the green LED light does not turn on, please repeat Step 1.



How to turn the Device Off:

To deactivate the device, press the silver button and hold it down for 1-2 seconds. Once the device is deactivated, the green LED light will turn off.



For Best Results:

The device loop area should be placed directly over the source of the pain. For maximum pain relief, wear continuously in one area until pain diminishes. The device should be placed as close to the skin as possible. The therapy will also be effective through light clothing.

You can use any type of wraps, adhesives, bandages or clothing to help hold the device in place.

MAINTENANCE AND STORAGE

- Use a damp cloth and mild soap to gently wipe clean after each use, when the device is soiled, or to remove any buildup of residue from medical adhesives.
- Device Operation: a temperature range of +5°C to +40°C; a relative humidity range of 15% to 90%, non-condensing, but not requiring a water vapor partial pressure greater than 50hPa
- Device Transport/Storage: a temperature range of -25°C to +5°C, and +5°C to +35°C at a relative humidity up to 90%, non-condensing; >35°C to 70°C at a water vapor pressure up to 50hPa
- The device should be operated, stored and transported at an atmospheric pressure between 50 kPa and 106 kPa (0.5 atm and 1.04 atm), up to an altitude of 5,575 m above sea level

Note: Consult your local electronics store or waste management company for guidance on proper disposal of the device.


APPLICATION

One or more of the following ACCESSORY may be provided in your Therapy Kit: Adhesives

Note: Accessories are only intended to secure the ActiPatch device on the body. Accessory use is optional.

Adhesives

For best results and maximum hold, clean and thoroughly dry the skin around the area of application to remove dirt, oils, and lotions.



Bend the tape in the center with the paper backer facing you. (Figure 1).

Peel back the raised backer tabs exposing the middle section of tape without removing the backer from the ends. Do not touch tape adhesive. (Figure 2).

Place the Adhesive tape on a flat surface with the exposed adhesive facing you, and place the activated ActiPatch device so the wire is positioned in the center of the tape (Figure 3).

Position the loop of the ActiPatch directly over the area of pain. Using the ends with the backer, pick the tape up and apply it to the skin securing the loop to the skin. (Figure 4).

Once the center of the tape is holding on the ActiPatch, remove the backer from the ends of the tape and finish securing the tape in place using gentle pressure. (Figure 5).

Note: Adhesives should be changed every 3 days or when the stickiness wears off, whichever may occur first.



/ActiPatch



@ActiPatch



/ActiPatch



@tryActiPatch

Find ActiPatch on Social Media

CHECK OUT

www.ActiPatch.com

FOR ACTIPATCH PRODUCTS, ATTACHMENTS, TUTORIALS, & MORE!

Adverse Reactions:

- If pain persists within 3-4 days of use or worsens with use, discontinue use and seek medical attention.

Precautions¹:

- There are no user-serviceable parts inside the unit. Do not attempt to modify or break open the device.
- Do not wear the device in the shower or bath: the device is not waterproof
- Keep this unit out of the reach of children
- If the LED light does not come on, it indicates that the device is no longer operational and can be disposed of according to local regulations.
- The device should not be used by/on children under the age of 17.
- The device is not intended for use on multiple patients.
- The IP (Ingress Protection) rating for the device is IP22 and therefore offers protection from touch by fingers and objects greater than 12 millimeters. Additionally, the device is protected from water spray less than 15 degrees from vertical.
- The device is non-sterile. Avoid exposing the device to lint, dust and light (including sunlight) to prevent discoloration and to prevent build up of residue.
- The time required for the ME EQUIPMENT to warm from the minimum storage temperature, or cool from the maximum storage temperature, is 1 hour.

Contraindications²:

- Do not use this device directly over a cardiac pacemaker, implanted defibrillator, deep brain stimulator and nerve stimulators or other active implantable device.
- Do not use this device if you are experiencing sudden, unexplained pain. Sudden, unexplained pain can be an indicator of a serious medical condition and may require immediate medical attention.
- Do not use the device if you do not know the cause of your musculoskeletal pain. Contact your doctor to know more about the source of your pain.
- ActiPatch is a therapy for the adjunct treatment of musculoskeletal pain. Do not use for pain which is located deeper in the body, for example in the chest or stomach. This device is not intended to treat pain deep in the body.
- Do not use this device if you are pregnant or think you are pregnant.
- Do not use this device to treat cancer related pain. This device is not intended to treat cancer related pain.

Warnings³:

- If you are in the care of a doctor, consult your doctor before using this device.
- If your pain does not improve after using the device for 7 days, stop using the device and consult your doctor.
- ActiPatch is not a sterile device, so it should not come in direct contact with open wounds or irritable spots.
- Choking hazard: do not swallow the unit.
- Keep out of the reach of children.
- Before using, check for damage to the module and cable insulation and use if there is no problem.

¹ A precaution is used to identify a hazard that may result in minor or moderate injury to the user or patient or damage to the equipment or other property.

² Contraindications are known and reasonably foreseeable conditions under which the device should not be used because the risk of use clearly outweighs any possible benefit.

³ A warning is used to identify a hazard that may lead to death or serious injury.

EXPLANATION OF SYMBOLS

Symbol	Description	Location	Other Information Location
	Manufacturer: This Symbol is accompanied by the name and address of the manufacturer.	Box	On Box: Barcode, Part#, Rev#, Warnings, Contents, Patent#
	Symbol for Authorized Representative In The European Community	Box	On Label: Quantity, Description, Model number
	Upper and Lower limit of temperature	Box	
	Attention, see warning statement	Box	
	Symbol for Not Sterile Product	Box	
	Upper and Lower limit of humidity	Box	
	Symbol for Follow instructions for use	Box	
	Type BF Applied Part	Box	
	Conformity marking for medical devices sold in the EU	Box	
	Non-ionizing radiation	Box	
	Upper and Lower limit for operation, storage and transport for atmospheric pressure	Box	
	Not made with natural rubber latex	Box	