



The VISION™ DSX open gantry extra large FOV Rectangular gamma-camera design combines mechanics, detection and control software optimally; thus providing easy to learn and swift operation, patient comfort and safety, clinical value, reliability, maintainability.

Its open gantry gives clinical scan versatility, patient contouring in SPECT and Whole Body scans, AUTO-LOOP™ automatic acquisition of anterior and posterior whole body scan and better patient access. It also enables the detector to be calibrated or checked in proper physical conditions.

VISION™ DSX uses a high PM-tube per cm^2 density detector, which, combined with SMV's Digital Deconvolution Detection System, possesses of the best intrinsic characteristics on the market.

VISION DSX acquisition console is used for acquisitions and system quality check. It may also be used for image hardcopy, transfer over the network or archiving onto optical disk. With The Processing workstation option,

Gantry Motions

The gantry has 5 motorized independent movements that can be performed simultaneously. All rotational movements have an accuracy control of 0.1° , whereas rectilinear movements of 0.1mm .

Orbit motion:

Clockwise and counter clockwise, used for SPECT acquisitions.

Range: -35° to $+565^\circ$

Speed: 0.4 turn/min to 1.5 turn/min

Radial motion

This motion is performed perpendicularly to the axis of rotation, allowing Patient contour following in SPECT and Whole body scan.

Range: -115 to 343 mm.

Speed: 16 to 40 mm/second

Vision™ DSX TECHNICAL DATA

it offers a cost effective, yet comprehensive fully integrated nuclear imaging system.

The DSX offers the operational simplicity and ease of use given by the sophia touchscreen graphic console and dialogue software developed specially for nuclear medicine use. Acquisition protocols are created in an instant entering parameters with the help of the touch-screen and the exclusive Sophia medical "point the box" system. Protocols of different acquisition types may be chained, thus meeting the needs of the most elaborate acquisition protocols. Patient care is improved through shorter preparation and scanning times with optimized gantry control software.

Highly automated acquisition procedures insure consistently reliable results, yet non-routine, experimental use of the camera is possible.

The VISION™ DSX open gantry has been designed to ensure maximum flexibility, precision and stability in all its movements. The microprocessor controlled gantry motions are all independent, quiet and safe.



VISION™ DSX

Lateral motion

Motor driven, used for whole body scan , Organ centered SPECT and imaging on hospital bed.

Range 230 cm

speed: 2 cm/minute to 200 cm/minute

Detector Angulation

Angle -100° to $+200^\circ$, $2^\circ/\text{sec}$ or $10^\circ/\text{s}$. Used for cranial/caudal view.

Patient bed

The scintibed I is a patient bed designed for use in SPECT (patient lies perpendicular to the gantry) as well as whole body imaging (the patient lies parallel to the gantry). Self-locking floor mounts in both positions guarantee total stability and accuracy. The panel supporting the patient is made of low-absorption carbon fiber.

Bed motions

Vertical range


72 cm to 97 cm




speed 50 and 100 cm/min
Manual Longitudinal range 110 cm

Patient Support Physical characteristics


Length 200 cm
With Head Extension 232 cm
Width (SPECT section) 350 mm
Width (whole body section) 580 mm
Thickness: 30 mm
Absorption: 5.5% (140 keV)
Panel flexion (160 kg load) 2 cm (SPECT)
Maximum load¹ 160 kg



Scintibed I



VISION DSX SPET body contouring



Whole Body Contouring & Autoloop™

VISION DSX Detector

| | |
|---|--|
| Overall dimensions: | 580x710x290 mm |
| Shielding | Up to 400 KeV |
| Sodium iodine crystal doped with thallium: NaI (Tl) | |
| Cristal size | 465 x 600 mm |
| thickness | 9.5 mm (3/8") 12.5 mm (1/2") ² |
| UFOV | 540 mm x 400 mm |
| number of PMTs | 94 |
| Hexagonal | 86 (2,3" across flats) |
| circular | 8 (1,5") |

Collimators

The collimator storage unit can contain five collimators. The associated handling cart is parked within it, providing a very compact storage area.

Collimator changing position is programmable on the gantry.

| | |
|-----------|---|
| LEAP140 | Low energy All Purpose |
| LEUHR140 | Pediatric, Bone, Brain SPECT |
| LEHS140 | First Pass |
| MEAP 250 | Indium General purpose imaging |
| MEAP 300 | Gallium General purpose imaging |
| HEAP 360 | Thyroid, Whole Body with I131 |
| LE-PH 360 | Low energy pin-hole for THYROID scans or Pediatrics scans |
| HE-PH 360 | Highenergy pin-hole for THYROID scans |
| HR-FB 140 | High resolution Fan Beam |
| UHR-FB | Ultra High resolution fan Beam |

VISION DSX: Available collimator types³

¹ with a safety coefficient of 6, (i.e. 960 kg distributed over total surface area > 2.5)

²Option

³see characteristics in appendix



DSX NEMA Detection Performance

| Intrinsic Spatial (mm) | | CFOV | UFOV |
|---|----------|--------|--------|
| Intrinsic spatial resolution | FWHM | 3.4 | 3.5 |
| | FWTM | 6.5 | 6.6 |
| Intrinsic spatial resolution @75Kcps | FWHM | n/a | n/a |
| | FWTM | | |
| Intrinsic spatial linearity @ 20Kc/s | Diff. | 0.06 | 0.07 |
| | Integral | 0.16 | 0.16 |
| Intrinsic flood field uniformity | | CFOV % | UFOV % |
| Without correction @ 140Kev, 20Kcps | Diff. | 1.4 | 1.5 |
| | Integral | 1.7 | 1.8 |
| With correction @ 140Kev, 20Kcps | Diff. | 1.4 | 1.5 |
| | Integral | 1.5 | 1.6 |
| With correction @ 140Kev, 75Kcps | Diff. | 1.4 | 2.0 |
| | Integral | 1.9 | 2.2 |
| Energy resolution @ 140keV , 20Kcps | FWHM | 10.3% | |
| Energy resolution @ 140keV, 75Kcps | FWHM | 10.5% | |
| Maximum Count Rate | | | |
| (w/o scatter) | | 250 | Kcps |
| 20% loss | Incident | 270 | Kcps |
| | Observ. | 220 | Kcps |
| Energy Range | 40 | 400 | Kev |
| Multiple window spatial registration ⁴ | < | 1.0 | mm |
| Reconstructed SPET LEUHR @ 150mm without scatter | Tang | 8.5 | mm |
| | radial | 9.3 | mm |
| | Central | 9.3 | mm |
| Reconstructed SPET LEUHR @ 150mm with scatter | Tang | 9.9 | mm |
| | radial | 12.3 | mm |
| | central | 11.3 | mm |

Vision DSX performances⁵(Typical)

⁴Nema maximum displacement of x+, x-, y+ ,y-, measurements Ga67, Ba133

⁵Preliminary data

Room Requirements


For more details refer to the installation manual.

Room size: recommended 490 x 460cm .

Minimum Height: 2,3 m

Gantry weight: 1700Kg with a high energy collimators. (1400 Kg on the front rail)

Required power: 2,5KVA nominal.



VISION DSX Recommended Room Requirements