


[Medit Help Center](#) > [Medit Link](#) > [Medit Link App](#)

 Search

Articles in this section ▾

Export scan data

9 months ago · Updated

Follow

The **Export** feature allows you to export scan files in various formats: meditmesh, OBJ, PLY, STL.

1. Go to **Work Box**.
2. Select the case that you wish to export files from.
3. Click the **Export** button in Work File List.



Export files in case

OptionsFile List

Export Files To

C:/Users/MEDIT/Desktop/

Folder Name

2021-07-05-Complete

File Name

2021-07-05-Complete

CAD Type

exocad

Mesh Format

Convert mesh files to the below formats

☐ meditMesh
☐ OBJ
☒ PLY
☐ STL

☐ Do not convert the attached files

Image Format

☒ Convert images to the below format

☒ JPG
☐ JPEG
☐ PNG
☐ BMP

☐ Compress Files
☐ Combine Individual Mesh

Close

Export

(1) CAD program for which you want to use your exported data.

- CAD program options: exocad, 3Shape, 3Shape 2020, Maestro, Kuzler, DWOS & Cares Visual, DTX Studio, CEREC, TiZian Creative RT
- Exported files will be compatible with the 3D axis of a program you chose.(It does not mean that the program has been integrated.)

💡 If the program is set to exocad, a dentalProject file is automatically created when exporting a case.

(2) Export File to

: Choose a storage location for exported files.

(3) Preferred Mesh Format

: Choose your preferred mesh format: meditmesh, OBJ(color data), STL (grayscale), PLY

: You can choose multiple formats. If you do so, the program will export scan files by multiple formats to the designated storage location.

💡 A "meditmesh" format — the data post-processed from its raw data — has drastically smaller file size than the size of its original data, which is helpful when transmitting data over internet. Furthermore, meditmesh data includes extra information which can be utilized for various Medit Apps.

(4) Combine Individual Mesh

: You can select whether to combine scanbodies with jaw data or export them individually.

: If the option is disabled, scanbody and jaw will be exported separately.

: If the option is enabled, the following data will be combined.

- Scanbody and its jaw data
- Abutment and its jaw data
- Occlusion data



Was this article helpful?

Yes

No

53 out of 117 found this helpful

Have more questions? [Submit a request](#)

Return to top ^

Recently viewed articles

[Intraoral Scanner \(PDF\)](#)

Comments

0 comments

Please [sign in](#) to leave a comment.

Medit Help Center

Terms of Use | Privacy Policy

English (US) ▾