



Weekly Check List

Month/Year:

OWNER NAME:

MODEL:

Next Month:

MAKE:

SERIAL NUMBER:

LOCATION:

Weekly Checklist for N2O/O2 Equipment	
TESTS:	Weekly Comments
Ensure blender is turned on	
Ensure dial for O2/flow control is turned off	
Ensure dial for N2O flow control is turned off	
Ensure Sure Seal washer is in place	
Turn on Oxygen and Nitrous oxide	
New O2: 1800-2200 PSI, N2O 600-700 PSI	
Ensure cylinders are full or have spare on cart	
Ensure no audible leaks at gauges or blender	
Attach your circuit, without nasal hood attached	
Taking the circuit block delivery tubes use thumbs	
Use the flush button on your blender, inflate rebreathing bag,	
Inflate rebreathing bag till firm, release flush button.	
Watch and ensure the bag stays inflated to ensure no leaks	
Remove thumbs from delivery tube ends, allow bag to deflate	
Once this is done the system is ready for patient use.	
Additional Comments:	

**At end of procedure ensure all oxygen and nitrous oxide cylinders are turned off.
 Ensure the blender is turned off, and returned to its storage location
 Ensure all equipment and supplies are cleaned between each patient, and at end of day.**

Requirement	Company	Contact Information:
Gas Supplier:		
Technical Support:	Vanguard Medical	(403) 828.0846
Parts/Equipment:		