
# Microstream<sup>®</sup>


ETCO<sub>2</sub> SAMPLING

**Products for  
Intubated Patients  
Infant and Neonatal Patients  
Directions for Use**

## FilterLine<sup>®</sup> H Sets

## VitaLine<sup>™</sup> H Sets

US Patents: [www.covidien.com/patents](http://www.covidien.com/patents)

© 2015 Oridion Medical 1987 Ltd.  
All rights reserved.

PN 006147N

**Oridion<sup>®</sup>**

# Microstream<sup>®</sup> EtCO<sub>2</sub> Sampling for Intubated Infant and Neonatal Patients

**FilterLine<sup>®</sup> H Sets** for humidified ventilation  
**VitaLine<sup>™</sup> H Sets** for use in humidity-controlled  
incubators and very high humidity environments

**Airway Adapters for Infant and Neonatal patients:**

- Added dead space <0.5 cc.
- To be used with an endotracheal tube ≤ 4.5mm bore.

Single patient use, not to be re-processed.

Non-sterile / non-pyrogenic.

Does not contain natural rubber latex

Intended use: to measure CO<sub>2</sub> concentration in exhaled breath for intubated infant/neonatal patients

**Warning:** When used with closed suction system, do not place the airway adapter between the suction catheter and endotracheal tube. This is to ensure that the airway adapter does not interfere with the functioning of the suction catheter.

**Warning:** Loose or damaged connections may compromise ventilation or cause an inaccurate measurement of respiratory gases. Securely connect all components and check connections for leaks according to standard clinical procedures.

**Warning:** Do not use the FilterLine<sup>®</sup> H Set/VitaLine<sup>™</sup> H Set Infant/Neonatal during magnetic resonance imaging (MRI) scanning. Using the products during MRI scanning could create an artefact on the MRI image.

**Warning:** Check CO<sub>2</sub> tubing regularly during use to ensure that no kinks are present. Kinked tubing may cause inaccurate CO<sub>2</sub> sampling.

**Warning:** Do not use the FilterLine<sup>®</sup> H Set or the VitaLine<sup>™</sup> H Set with products with protruding inner connector pieces (i.e., the Hamilton flow sensor), as this may cause breakage of the airway adapter.

**Caution:** Non-standard dimensioning of some fittings on ventilator or anesthesia equipment may make secure mating with the airway adapter difficult. Use only with equipment having standard 15 mm fittings. The sampling lines are to be used with an inlet ET tube connector with an inner diameter of greater than 10 mm.

**Caution:** Ensure that tubing is not stretched during use.

**Caution:** Do not attempt to clean, disinfect, sterilize or flush any part of the sampling line as this can cause damage to the monitor.

**Caution:** Dispose of sampling lines according to standard operating procedures or local regulations for the disposal of contaminated medical waste.

**Caution:** Incorrect placement of the Nafion<sup>®</sup>\*\* humidity filters may affect the duration of the VitaLine<sup>™</sup> (placement is described below).

**Note:** During nebulization or suction, in order to avoid moisture buildup and sampling line occlusion, remove the sampling line luer connector from the monitor.

**Note:** Replace the sampling line according to hospital protocol or when a blockage is indicated by the device. Excessive patient secretions or a build-up of liquids in the airway tubing may occlude the sampling line, requiring more frequent replacement.

**Note:** When connecting a sampling line to the monitor, insert the sampling line connector clockwise into the monitor CO<sub>2</sub> port and turn the connector until it can no longer be turned, to ensure that it is connected securely to the monitor. This will assure that there is no leak of gases during measurement at the connection point and that measurement accuracy is not compromised.

**Note:** Following connection of the CO<sub>2</sub> sampling line, check that CO<sub>2</sub> values appear on the monitor display.

**Note:** To remove the airway adapter from the breathing circuit, twist the airway adapter outwards by a quarter turn while firmly gripping the closest section of the ventilation source or patient's endotracheal tubing, and disconnect.

**Note:** During setup, ensure that airway adapter can be easily attached and detached from the breathing circuit before proceeding.

**Note:** Expect delay times for the following products to be longer than the standard by the given times: Long (4m) FilterLines - 3 seconds.

**Note:** Sampling lines with H in their names include a moisture reduction component (Nafion<sup>®</sup> or its equivalent) for use in higher humidity environments where long duration use of CO<sub>2</sub> sampling is required.




Federal law (U.S.) restricts this device to sale by or on the order of a physician.

\*\* or alternate moisture reduction technology

# Microstream® EtCO<sub>2</sub> Sampling for Intubated Patients

## Instructions for use:

1. Firmly connect the small-end (male) of the Microstream® Airway Adapter ① to the female-end of the Wye piece (ventilation source).
2. Firmly connect the patient's endotracheal tube connector into the large-end (female) ② of the Microstream® Airway Adapter.
3. When using a VitaLine™ with an incubator, string the VitaLine™ tubing through the tubing access door of the incubator such that one Nafion®\* humidity filter ③ is positioned inside the incubator, while the second Nafion®\* humidity filter ④ is positioned outside the incubator, close to the incubator tubing access door. Ensure that the tubing is not twisted or crimped.
4. Twist the female luer connector ⑤ securely but gently into the Microstream® input port ⑥.
5. Following connection of the CO<sub>2</sub> sampling line, check that CO<sub>2</sub> values appear on the monitor display.
6. To remove the airway adapter from the breathing circuit, twist the airway adapter outwards by a quarter turn while firmly gripping the closest section of the ventilation source or patient's endotracheal tubing, and disconnect.
7. For use with a closed suction system, see sketches on the following pages.

## Legend


1. Male-end of airway adapter
2. Female end of airway adapter
3. Humidity filter in incubator
4. Humidity filter outside incubator
5. Female luer connector
6. Microstream® input port

# Microstream® EtCO<sub>2</sub> Sampling with a Closed Suction System

## **Instructions for use with a closed suction system:**

1. Set up the parts of the system as seen below. See legend.
2. Connect the system as seen below, as follows:
3. Firmly connect the small-end (male) of the Microstream® Airway Adapter 8 to the female-end of the Wye piece (ventilation source) 9.
4. Firmly connect the large-end (female) of the Microstream® airway adapter 8 to the ventilator connector on the closed suction system 10, so that the airway adapter does not interfere with the functioning of the suction tubing 11.
5. The setup seen at the bottom right is incorrect; if set up in this way, the airway adapter will interfere with the functioning of the suction tubing.


## **Setup of Microstream EtCO<sub>2</sub> Sampling with a closed suction system**


## **Legend**

1. Male-end of airway adapter
2. Female end of airway adapter
3. Humidity filter in incubator
4. Humidity filter outside incubator
5. Female luer connector
6. Microstream® input port
7. Endotracheal tube
8. FilterLine Set with airway adapter
9. Wye piece (ventilation source)
10. Ventilator connector on closed suction system, with lavage
11. Suction tubing

**Correctly Connected  
Microstream EtCO<sub>2</sub> Sampling  
with a closed suction system**

**Incorrectly Connected  
Microstream EtCO<sub>2</sub> Sampling  
with a closed suction system**

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 Single patient use	For use with a single patient
<b>Rx ONLY</b>	Rx Only
 Consult instructions for use	Consult instructions for use
 Caution, consult accompanying documents	Caution, consult accompanying documents
 -4°F -20°C 158°F 70°C Temperature limitation	Storage temperature
 This product does not contain natural rubber latex	Does not contain natural rubber latex
 MR Unsafe	Not safe for use in MR environments
	CE Mark

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Oridion®, Microstream®, FilterLine®, Smart CapnoLine®, CapnoLine®, Smart CapnoLine Guardian™ and VitaLine™ are trademarks or registered trademarks of Oridion Medical 1987 Ltd. ™\* and ®\* are trademarks of their respective owners.

Oridion Medical 1987 Ltd.  
7 Hamarpe St.  
P.O. Box 45025  
Jerusalem 9777407, Israel  
Tel: 972-2-589-9111  
Fax: 972-2-582-5873  
US Toll-free: 1-888-ORIDION  
(1-888-674-3466)

Obelis S.A.  
Bd Général Wahis 53  
1030 Brussels, Belgium  
Tel: +(32) 2 732-59-54

[www.oridion.com](http://www.oridion.com)