







# Microstream<sup>®</sup>

## ETCO<sub>2</sub> SAMPLING

### Products for Non-Intubated Patients

### Directions for Use

	 <p>CO<sub>2</sub> Nasal Sampling with O<sub>2</sub></p>  <p>CO<sub>2</sub> Nasal Sampling</p>
	 <p>CO<sub>2</sub> Oral/ Nasal Sampling with O<sub>2</sub></p>  <p>CO<sub>2</sub> Oral/ Nasal Sampling</p>
	 <p>CO<sub>2</sub> Oral/Nasal Sampling with O<sub>2</sub></p>
 <p>O<sub>2</sub></p>	 <p>≤ 5 L/min</p> <p>O<sub>2</sub></p> <p>CO<sub>2</sub></p>

US Patents: [www.covidien.com/patents](http://www.covidien.com/patents).

© 2015 Oridion Medical 1987 Ltd.

All rights reserved.

PN 009809M

Oridion<sup>®</sup>

## Microstream® EtCO<sub>2</sub> Sampling for Non-Intubated Patients

For use with all Microstream-enabled monitors

Single patient use, not to be re-processed

Non-sterile / non-pyrogenic

Does not contain natural rubber latex

**Warning:** Loose or damaged connections may compromise ventilation or cause an inaccurate measurement of respiratory gases. Securely connect all components and check connections for leaks according to standard clinical procedures.

**Warning:** Do not cut or remove any part of the sample line. Cutting the sample line could lead to erroneous readings.

**Warning:** Check CO<sub>2</sub> and O<sub>2</sub> tubing regularly during use to ensure that no kinks are present. Kinked tubing may cause inaccurate CO<sub>2</sub> sampling or affect O<sub>2</sub> delivery to patient.

**Caution:** Ensure that tubing is not stretched during use.

**Caution:** Do not attempt to clean, disinfect, sterilize or flush any part of the sampling line as this can cause damage to the monitor.

**Caution:** Dispose of sampling lines according to standard operating procedures or local regulations for the disposal of contaminated medical waste.

**Caution:** Use of a CO<sub>2</sub> sampling line with H in its name (indicating that it is for use in humidified environments) during MRI scanning may cause interference. These sampling lines include CapnoLine H/Long, CapnoLine H O<sub>2</sub>, Smart CapnoLine H/Long, Smart CapnoLine H O<sub>2</sub>, and Smart CapnoLine H Plus O<sub>2</sub>/Long. The use of non H sampling lines is advised.

**Note:** Replace the sampling line according to hospital protocol or when a blockage is indicated by the device. Excessive patient secretions or a build-up of liquids in the airway tubing may occlude the sampling line, requiring more frequent replacement.

**Note:** When connecting a sampling line to the monitor, insert the sampling line connector clockwise into the monitor CO<sub>2</sub> port and turn the connector until it can no longer be turned, to ensure that it is connected securely to the monitor. This will assure that there is no leak of gases during measurement at the connection point and that measurement accuracy is not compromised.

**Note:** Following connection of the CO<sub>2</sub> sampling line, check that CO<sub>2</sub> values appear on the monitor display.

**Note:** Expect delay times for the following products to be longer than the standard by the given times: Long (4m) FilterLines/CapnoLines - 3 seconds; Nasal FilterLine Infant/Neonatal - 1 second.

**Note:** Sampling lines with H in their names include a moisture reduction component (Nafion®\* or its equivalent) for use in higher humidity environments where long duration use of CO<sub>2</sub> sampling is required.

**Note:** Sampling lines are recommended for use with oxygen provided at up to 5l/min. At higher levels of oxygen provision, dilution of CO<sub>2</sub> readings may occur, leading to lower CO<sub>2</sub> values.

Federal law (U.S.) restricts this device to sale by or on the order of a physician.

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 Single patient use	For use with a single patient
	Rx Only
 Consult instructions for use	Consult Instructions for Use
 Caution, consult accompanying documents	Caution, consult accompanying documents
 -4°F -20°C 131°F 55°C Temperature limitation	Storage temperature
	DEHP
 This product does not contain natural rubber latex	Does not contain natural rubber latex
 MR Unsafe	Not safe for use in MR environments
 Fire hazard	Fire Hazard
 0482	CE Mark

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Oridion®, Microstream®, FilterLine®, Smart CapnoLine®, CapnoLine®, and Smart CapnoLine Guardian™ are trademarks or registered trademarks of Oridion Medical 1987 Ltd.™\* and ®\* are trademarks of their respective owners.

Oridion Medical 1987 Ltd.  
7 Hamarpe St.  
P.O. Box 45025  
Jerusalem 9777407, Israel  
US Toll Free: 1-888-674-3466  
1-888-ORIDION

Obelis S.A.  
Bd Général Wahis 53  
1030 Brussels  
Belgium

[www.oridion.com](http://www.oridion.com)