

Importance of Properly Fit TF Prosthetics & Microprocessor Knees

Scott Sabolich, CP, LP 12-November-2025
Senior Director Ottobock Patient Care



Scott Sabolich, CP, LP

Senior Director, Ottobock Patient Care


- **Certified & Licensed Prosthetist since 1995, in the field since 1989 and before**
- **Specializes in difficult to fit TF & challenging cases**
- **Instructor of the patented Sabolich Socket through Ottobock**
- **In 1999, opened Scott Sabolich Prosthetics and Research in OKC as a destination facility with a state-of-the-art Prosthetic Patient care center**
- **In 2012, he opened a second facility to extend exceptional care and expertise in Dallas TX**
- **In 2019, Scott Sabolich Prosthetics and Research joined Ottobock, a global leader in prosthetic and orthotic technology as Ottobock entered into Patient Care**



50+


Clinics Across the U.S.

ottobock.



The vicious cycle of falling in TF prosthetics


ottobock.



MacKay et al. Fear of Falling in Older Adults: A Scoping Review of Recent Literature. Can Geriatr J. 2021;24(4):379-394.

The vicious cycle of falling in TF prosthetics

ottobock.



Socket Stability (Handcast)

ottobock.



Socket Stability (Pre Fab)

ottobock.

Most Popular Adjustable Prefab Sockets

- Ottobock Varos
- Martin Bionics
- LIM Innovations
- Ifit Immediate Fit
- Click Medical Revo Design



Insert Content on Gait Deviations

Insert Content on Socket Fittings


Socket Stability

ottobock.



Knee Stability

ottobock.



Impact of MPKs (and Sockets) on Falls Research

The Prevalence and Risk Factors of Falling and Fear of Falling Among Lower Extremity Amputees

William C. Miller, PhD, MSc, OT, Mark Speechley, PhD, Barry Deathe, MD

Pre-MPK era data (2001)!

- 52% fell in previous 12 months (19% sought medical attention)
- TF **2.78x more likely** to fall vs. TT

Impact of Vascular Disease, Amputation Level, and the Mismatch Between Balance Ability and Balance Confidence in a Cross-Sectional Study of the Likelihood of Falls Among People With Limb Loss

Post-MPK era data (2018/2019)!

Christopher Kevin Wong, PT, PhD, OCS, and Stanford T. Chihuri, MPH

- 66% fell in previous 12 months
- After adjusting for vascular comorbidity and balance, TF **less likely to fall** vs. TT
- TT's **2.32x more likely** to have a fall-related injury vs. TF (Chihuri and Wong 2018)

With the advent of MPK adoption for AK amputees:

- TF level alone no longer associated with higher fall risk vs TT
- TF amputation associated with a lower **fall injury** risk vs TT

MPKs and Cognitive Load

ottobock.

Where are the eyes looking?

Stair descent & multi-tasking



Courtesy of Dale Berry, CPO

Lower Mobility



Higher Mobility



Kenevo

K2 level solution
Lightest and shortest MPK
Progresses with the user



C-Leg 4

K3 level solution
Most studied knee in the world
Fits all daily activities



Rheo Knee

Low to Moderate Impact
Dependable swing initiation
Weatherproof



Genium X4

K4 level solution
Advanced gait algorithms
Rugged for impact activity

Default Stance


Default Swing

Default Stance




- ✓ Activity Modes
- ✓ Supported sitting/standing
- ✓ Wheelchair function
- ✓ Donning function
- ✓ Bicycle function
- ✓ Battery life: 1 day


Activity Modes



Functions as a
locked knee



Functions as a weight-
activated knee



Functions as a
yielding knee


Activity modes are set by the prosthetist according to the end-user's functional abilities. The modes can be changed when the end-user progresses through the recovery process.




Basic Functions



Donning Function
(A/B/B+)



Supported sitting
down & standing
up



Sitting Function



Wheelchair
Function



Stance Function
(C)



Bicycle Function

These underlying functions are designed with the lower mobility individual in mind. Their purpose is safety and stability during functional activities.

C-Leg 4

ottobock.



- ✓ Basic Mode
- ✓ Stumble Recovery Plus
- ✓ Stance function
- ✓ Training feedback
- ✓ MyModes (5)
- ✓ Battery life: 2 days


C-Leg 4

Basic Mode

ottobock.



Knee flexion angle
65°




Two basic adjustments are required: stance flexion resistance and the swing flexion angle. In addition, two MyModes can be set for specific activities.



C-Leg 4

Training function



An acoustic feedback signal is emitted when stance release is successfully triggered

(User enabled)

C-Leg 4

ottobock.

Stance Function (intuitive or deliberate)

- ✓ Stability on uneven surfaces
- ✓ Rest for the sound limb
- ✓ Positioning for tasks



**Can hamper training efforts if not trained*

Rheo / Rheo XC

ottobock.



- ✓ Mimicks natural gait
- ✓ Effortless Swing Initiation
- ✓ Weight activated stance
- ✓ Weather proof
- ✓ Manual extension lock
- ✓ Step over and step backwards features
- ✓ XC Stair over stair and walk to run

Genium X4

ottobock.



- ✓ Start to Walk
- ✓ PreFlex
- ✓ Added support on ramps
- ✓ Optimized slope ascent
- ✓ Stair and obstacles
- ✓ Intuitive cycling
- ✓ Dynamic backward movement
- ✓ Battery life: 5 days

Genium X4 Features

ottobock.



Start-to-walk



On Set value

This function makes it easier to go from standing to walking and to walk in tight spaces.

The user can turn this function on and off in the Cockpit app.

[? More information](#)



PreFlex



On Set value

This function facilitates the initiation of stance phase flexion and allows for more harmonious forward movement.

The user can turn this function on and off in the Cockpit app.

[? More information](#)



Dynamic backward movement

20

Minimum

Maximum



The user can adjust this value between 15 and 30 in the Cockpit app.

This function enables dynamic backward movement due to increased flexion resistance.

[? More information](#)



Intuitive cycling



On Set value


This function enables automatic switching into and out of cycling mode.

The user can turn this function on and off in the Cockpit app.

[? More information](#)

Genium X4 Features

ottobock.



Early stance phase support

Minimum Maximum

0 3

☐ The user can adjust this value between 0 and 1 in the Cockpit app.

Set value

This function defines the level of support in early stance phase when descending flat ramps. If the user feels too strongly restricted, they can start at setting 0.

[More information](#)



Late stance phase support

Dynamic Supported


1 7

☐ The user can adjust this value between 1 and 7 in the Cockpit app.

Set value

This function defines the level of support in late stance phase when descending ramps and stairs.

[More information](#)




Optimised ascent ☒

☐ On

This function makes it easier to walk up hills and ramps.

The user can turn this function on and off in the Cockpit app.

[More information](#)



Stairs and obstacles ☒

☐ On

This function makes it possible to climb stairs step-over-step and to overcome obstacles.

The user can turn this function on and off in the Cockpit app.

[More information](#)

Water & MPKs

Moisture protection for the electronics

ottobock.



IP 22

Dripping water



IP 68

Continuous water immersion
3m water up to 1 hour



IP 68

Continuous water immersion
3m water up to 1 hour
Corrosion-resistant

ASCENT K2 Study Findings

Randomized study including over 100 subjects

MPK

Fear of falling:	Reduced
Falls:	1.3
Near-Falls:	1.7
Walking speed:	Faster
Fall risk (TUG):	Lower
Fall impact:	-2.21 ± 1.30
Quality of Life:	No change

Green: significant improvement or difference

*Fear of falling activity avoidance; **Consequences of falling; ***Quality of life



NMPK

No change
2.7
5.7
No change
No change
No change
Decline

Red: significant decline

Conclusion

~~“Why would this person benefit from an MPK?”~~

“Why would an MPK not be appropriate for this person?”



- MacKay S, Ebert P, Harbidge C, Hogan DB. Fear of Falling in Older Adults: A Scoping Review of Recent Literature. *Can Geriatr J*. 2021;24(4):379-394. Published 2021 Dec 1. doi:10.5770/cgj.24.521
- Miller WC, Speechley M, Deathe B. The prevalence and risk factors of falling and fear of falling among lower extremity amputees. *Arch Phys Med Rehabil*. 2001;82(8):1031-1037. doi:10.1053/apmr.2001.24295
- Blumentritt S: Biomechanische Aspekte zur Indikation von Prothesenkniegelenken. [Biomechanical aspects of the indications of prosthetic knee joints.] *Orthopädie-Technik*. 2004;55(6):508-524. Article in German.
- Köhler TM, Bellmann M, Blumentritt S. Polycentric exoprosthetic knee joints – extent of shortening during swing phase. *Can Prosthet Orthot J* 2020;Volume, Issue1, No.5. <https://doi.org/10.33137/cpoj.v3i1.33768>
- Wong CK, Chihuri ST. Impact of Vascular Disease, Amputation Level, and the Mismatch Between Balance Ability and Balance Confidence in a Cross-Sectional Study of the Likelihood of Falls Among People With Limb Loss: Perception Versus Reality. *Am J Phys Med Rehabil*. 2019;98(2):130-135. doi:10.1097/PHM.0000000000001034
- Hahn A, Bueschges S, Prager M, Kannenberg A. The effect of microprocessor controlled exo-prosthetic knees on limited commClinical Researchy ambulators: systematic review and meta-analysis. *Disabil Rehabil* 2021 Oct 25:1-19. Download
- Highsmith MJ, Kahle JT, Bongiorno DR, Sutton BS, Groer S, Kaufman KR. Safety, energy efficiency, and cost efficacy of the C-leg for transfemoral amputees. *Prosth Orthot Int* 2010;34(4):362-377. Download (open access)
- Hafner BJ, Smith DG. Differences in function and safety between Medicare Functional Classification Level-2 and -3 transfemoral amputees and influence of prosthetic knee joint control. *J Rehabil Res Dev* 2009;46(3):417-434.
- Blumentritt S, Schmalz T, Jarasch R. The safety of C-leg: Biomechanical tests. *J Prosthet Orthot* 2009;21(1):2-17. Download (open access for AAOP members only)
- Kahle JT, Highsmith MJ, Hubbard SL. Comparison of Non-microprocessor Knee Mechanism versus C-Leg on Prosthesis Evaluation Questionnaire, Stumbles, Falls, Walking Tests, Stair Descent, and Knee Preference; *J Rehabil Res Dev* 2008;45(1):1-14.
- Hafner BJ, Willingham LL, Buell NC, Allyn KJ, Smith DG: Evaluation of Function, Performance, and Preference as Transfemoral Amputees Transition from Mechanical to Microprocessor Control of the Prosthetic Knee. *Arch Phys Med Rehabil* 2007;88(2):207-217. Download (available for purchase)
- Campbell JH, Stevens PM, Wurdeman SR. OASIS I. Retrospective analysis of four different microprocessor knee types. *Journal Rehabil Assist Technol Eng* 2020;7: 1-10. Download (open access)
- Highsmith MJ, Kahle JT, Miro RM, Mengelkoch, MJ. Ramp descent performance with the C-leg and interrater reliability of the Hill Assessment Index. *Prosthet Orthot Int* 2013;37(5):362-368. Download (open access)
- Seymour R, Engbretson B, Kott K, Ordway N, Brooks G, Crannell J, Hickernell E, Wheller K. Comparison between the C-leg(R) microprocessor-controlled prosthetic knee and non-microprocessor control prosthetic knees: A preliminary study of energy expenditure, obstacle course performance, and quality of life survey. *Prosthet Orthot Int* 2007;31(1):51 – 61. Download (open access)
- Wurdeman SR, Hafner BJ, Sawers A, England DL, Lundstrom R, Kannenberg A. ASsessing Clinical outcomes with microprocEssor kNee uTilization in a K2 population (ASCENT K2): randomized controlled trial results for above-knee prosthesis users over age 65., *Disabil Rehabil* 2025 Jul1-18

Q&A

Thank you