

## Typical results

Before



Photo courtesy of Diana Demidion, M.D.

After




Photo courtesy of Hiroki Kanda, M.D.



Photo courtesy of Melanie D. Palm, M.D.



Photo courtesy of Mariela Nazar, M.D.



To start your own  
before and after journey,  
ask about **Stellar M22™**  
skin treatment, today.

Romero Advanced Dermatology  
24953 Paseo de Valencia Suite 16b  
Laguna Hills, CA 92653  
(949) 447-2099

[AESTHETIPEDIA.COM](http://AESTHETIPEDIA.COM)



**Stellar M22**

Shine on

## Photorejuvenation

Photo-ready skin every day





# Why Photorejuvenation?

A fast, yet gentle treatment that visibly improves:

- ✓ Pigmentation
- ✓ Age/sun spots
- ✓ Broken capillaries and
- ✓ Overall skin appearance by stimulating collagen and elastic fiber production.

## What is it?

**The Stellar M22™ Photorejuvenation treatment** enables your physician to offer you remarkable results on a range of skin concerns such as age and sun spots, freckles, birthmarks, rosacea and broken capillaries.

Our unique IPL (Intense Pulsed Light) technology gently and effectively eliminates spots, discolorations, and unwanted blood vessels.

### Is It right for me?

Stellar M22™ treatments are not suitable for everyone and may carry some risks. Photorejuvenation is not suitable when you have active infections, viral, fungal or bacterial diseases, inflammatory skin conditions or skin cancer. Risks may include: redness, swelling, scarring or a change of pigmentation. Be sure to consult with your treatment provider before choosing this treatment.

## What should I expect?

### When will I see results?

The Stellar M22™ Photorejuvenation treatment works with your skin to remove your imperfections and stimulate collagen and elastic fiber production over a sequence of treatments. Some skin concerns such as shallow pigmentation and age/sun spots can be cleared with as few as 1-2 treatments. Typically 2-6 sessions with 2-4 weeks in between are needed to achieve optimal results for most skin concerns.

### What can I expect after treatment?

Immediately following treatment, you may experience some redness, depending on your customized treatment settings. The redness will usually disappear within a few hours. In most cases, makeup may be immediately applied, and daily activities can be resumed the very same day. Typically your physician should advise to stay out of direct sun for a few weeks following the treatment and to apply sun screen.

## How does it work?

### Photorejuvenation technology


The IPL (Intense Pulsed Light) technology allows treatment of pigmentation and vascular imperfections, while improving overall skin appearance. Stellar M22™ IPL with Optimal Pulse Technology (OPT™) enables your clinician to customize the treatment according to your skin type and desired results, with maximum safety. Light pulses penetrate the tissue and create heat, targeting the skin imperfection in a controlled manner. The body's natural process then removes the treated tissue, while collagen and elastic fiber stimulation occurs, resulting in a visible improvement in skin tone and overall appearance.

### The treatment

To ensure utmost comfort, your physician should typically apply a topical numbing cream before treatment and/or use air cooling during treatment. Depending on your specific concern, you can expect the session to last 15-30 minutes. You may experience a warm or prickling sensation as the laser is applied to the skin, but the treatment is gentle and should not feel uncomfortable.



Filtered light specifically targets blood vessels below skin surface.



Intense pulses of light are delivered through the lightguide with integrated cooling to ensure maximum comfort.

