

»» **For multi-unit restorations -**
conventional or digital ««




MPS
Multi Platform Systems

**MULTI-UNIT
ABUTMENTS**


» Multi-unit abutment «

The Multi-unit abutment supports a variety of prosthetic restorations. Thus it is ideal for creating patient oriented individual hybrid restorations or being the base for a customized full arch restoration according to The Full Arch principle.

Implant with conventional shoulder and straight abutment



Implant with conventional shoulder and 30° angled abutment




»» The Full Arch principle ««

The Full Arch principle for atrophic ridges allows for a comprehensive treatment plan for edentulous patients, of full-arch immediate restoration, using just two straight and two 30° angulated abutments.

The temporary restoration is immediately placed over the implants with the Multi-unit abutments.

The straight and 30° angled Multi-unit abutments allow for a common insertion direction for the prosthetic restoration.



All Multi-unit abutments are supplied sterile including the insertion aid.

»» There are various possibilities for prosthetic restorations ««

- Cement or bond the titanium base into the denture to create a temporary restoration.
- The framework production of superstructures by means of CoCr caps or gold caps.
- CAD/CAM milled bars or bridges.
- The Full Arch treatment concept.



Full denture

A provisional full denture has been sufficiently blocked out/relieved in the area of the titanium bases and then inserted passively into the denture on the model and / or in the mouth.



Final bridge restoration on two implants in the posterior region

A metal framework is blocked out/relieved sufficiently in the area of the titanium bases and then orally glued without tension. Next exists the possibility of making the framework by means of directly cast on the gold caps or CoCr caps. After final processing, the bridge will be screwed occlusally.

This procedure is also suitable for the preparation of a final full arch restoration on at least four implants in the lower jaw or six implants in the upper jaw.



Bar with removable full denture

Onto the Multi-unit abutments a individually milled or CAD/CAM made bar is screwed on. Afterwards, screw the cast bar or milled bar onto the multi-unit abutments. Depending on the indication and patient request, the final removable denture can be manufactured.

The Multi-unit product range is being continually extended to include other abutment series.

Available for a large number of implant systems:

EV-Series	compatible with	Dentsply Implants	ASTRA TECH OsseoSpeed EV
F-Series	compatible with	Nobel Biocare	NobelActive ^{®*}
			NobelReplace [®] Conical*
H-Series	compatible with	BIOMET 3i	Certain ^{®*}
N-Series	compatible with	Straumann	Tissue Level
R-Series	compatible with	Zimmer Dental	Tapered Screw-Vent ^{®*}

All Multi-unit abutments including inserting tool are supplied in sterile packaging.


* is a registered trademark of an independent third party


»» Multi-unit abutments ««







»» Determination of the gingiva height ««

The gingiva height stated in the catalogue is marked here as **GH1**.
With the two angled abutments **GH2** is calculated as shown:





Angulation	Variant		Gingiva height (GH)*
straight			low
			medium
			high
17°	Type 1: over flat 	Type 2: over corner 	low
			medium
			high
30°	Type 1: over flat 	Type 2: over corner 	low
			medium
			high

*Exact lengths depend on the implant system

Multi-unit prosthetic components



Castable gold cap
to cast
gold alloys



CoCr cap
to cast
non-precious metals



Titanium base
to glue into the bridge/denture
(passive fit)



Titanium cap
to shorten individually and glue into the
bridge or denture as well as to create
temporary restorations



**Modelling aid for modelling on the
titanium cap or titanium base**
only usable in combination with the
titanium base or titanium cap



Bridge screw
included with every cap
and the titanium base



Multi-unit prosthetic components for digital restoration

Scanbody Multi-unit straight
for digital impression of the
straight Multi-unit abutments



Scanbody Multi-unit angled
for digital impression of the
angled Multi-unit
abutments 17° and 30°



3shape
exocad
dental wings

Titanium base
to glue into the bridge/denture
(passive fit)



Titanium cap
to shorten individually and glue into the
bridge or denture and to create
temporary restorations



Laboratory implant CAD/CAM
for use with digital model
production (print model)



Bridge screw
included with every cap
and the titanium base



»» Laboratory implants Multi-unit abutments ««



Laboratory implant
for use with conventional
model production.



Laboratory implant angled
for use with conventional model
production. Angled 17° & 30°.



Laboratory implant CADCAM
for use with digital model
production (print model).



Laboratory implant CADCAM angled
for use with digital model production
(print model). Angled 17° & 30°.

» Multi-unit accessories «

Implant pick up
for impression of the
Multi-unit abutments.



Cover cap
to protect the abutments
during the healing phase.




Placement instrument
used with the straight Multi-unit abutments.



» Multi-unit divergency compensation «

The Multi-unit abutments have a uniform platform to the prosthetic components with a conicity of 21° (one-sided).


When using the 30° angled Multi-unit abutments divergencies of up to 102° between implants can be covered. This results out of the 30° angulation of the multi-unit abutments and the conicity of the platform of 21°.



»» How to use ««

Multi-unit abutment straight

After the implant placement, the Multi-unit abutment straight is connected to the implant with its special insertion aid, and tightened with the aid according to the torque value indicated in the catalog.




All Multi-unit abutments are supplied sterile including the insertion aid.

»» How to use ««

Multi-unit abutment angled

After the implant placement, the angulated Multi-unit abutment is connected to the implant with its special insertion aid, and tightened with the aid according to the torque value indicated in the catalog.



All Multi-unit abutments are supplied sterile including the insertion aid.



>> Portfolio overview
Multi-unit abutments <<

Multi-unit scanbody

- Titanium specially coated
- incl. bridge screw



Version	straight	angled
Article No.	EV 8500	EV 8513

Please note: The Scanbody is sterilisable and for intra-oral scanning.

Multi-unit Laboratory implant

- Titanium Grade 5 CF



Version	straight	angled	angled
Type		17°	30°
Article No.	8504	8514-17	8514-30

Multi-unit Laboratory implant CAD/CAM

- Titanium Grade 5 CF



Version	straight	angled	angled
Type		17°	30°
Article No.	8509	8515-17	8515-30

Abutment screw

- Titanium Grade 5 CF
- Hex 1,26



Article No.	EV 64	EV 65	EV 66
-------------	-------	-------	-------

Please note: The correct choice of the abutment screw is stated with the respective abutment in the catalogue.

Multi-unit abutment straight

- Titanium Grade 5 CF
- Sterile packaged
- Recommended torque: 25 Ncm



Implant connection	D 3,6	D 3,6	D 3,6
Gingiva height	1,5 mm	3,5 mm	5,5 mm
Article No.	EV 8010	EV 8011	EV 8012

Multi-unit abutment straight

- Titanium Grade 5 CF
- Sterile packaged
- Recommended torque: 25 Ncm



Implant connection	D 4,2	D 4,2	D 4,2
Gingiva height	1,5 mm	3,5 mm	5,5 mm
Article No.	EV 8020	EV 8021	EV 8022

Multi-unit abutment straight

- Titanium Grade 5 CF
- Sterile packaged
- Recommended torque: 25 Ncm



Implant connection	D 4,8	D 4,8	D 4,8
Gingiva height	1,5 mm	3,5 mm	5,5 mm
Article No.	EV 8030	EV 8031	EV 8032

Multi-unit abutment angled 17°

- Titanium Grade 5 CF
- Sterile packaged
- incl. abutment screw
- Type 1 = angled over flat
- Recommended torque: 25 Ncm



Implant connection	D 3,6	D 3,6	D 3,6
Gingiva height (mm)	1,1/2,5	2,1/3,5	4,1/5,5
Article No.	EV 8110	EV 8111	EV 8112
Abutment screw	EV 64	EV 64	EV 64

Multi-unit abutment angled 30°

- Titanium Grade 5 CF
- Sterile packaged
- incl. abutment screw
- Type 1 = angled over flat
- Recommended torque: 25 Ncm



Implant connection	D 3,6	D 3,6	D 3,6
Gingiva height (mm)	0,6/3,0	1,6/4,0	3,1/5,5
Article No.	EV 8310	EV 8311	EV 8312
Abutment screw	EV 64	EV 64	EV 64

Multi-unit abutment angled 17°

- Titanium Grade 5 CF
- Sterile packaged
- incl. abutment screw
- Type 1 = angled over flat
- Recommended torque: 25 Ncm



Implant connection	D 4,2	D 4,2	D 4,2
Gingiva height (mm)	1,1/2,5	2,1/3,5	4,1/5,5
Article No.	EV 8120	EV 8121	EV 8122
Abutment screw	EV 65	EV 65	EV 65



Multi-unit abutment angled 30°

- Titanium Grade 5 CF
- Sterile packaged
- incl. abutment screw
- Type 1 = angled over flat
- Recommended torque: 25 Ncm



Implant connection	D 4,2	D 4,2	D 4,2
Gingiva height (mm)	0,6/3,0	1,6/4,0	3,1/5,5
Article No.	EV 8320	EV 8321	EV 8322
Abutment screw	EV 65	EV 65	EV 65



Multi-unit abutment angled 17°

- Titanium Grade 5 CF
- Sterile packaged
- incl. abutment screw
- Type 1 = angled over flat
- Recommended torque: 25 Ncm



Implant connection	D 4,8	D 4,8	D 4,8
Gingiva height (mm)	1,1/2,5	2,1/3,5	4,1/5,5
Article No.	EV 8130	EV 8131	EV 8132
Abutment screw	EV 66	EV 66	EV 66



Multi-unit abutment angled 30°

- Titanium Grade 5 CF
- Sterile packaged
- incl. abutment screw
- Type 1 = angled over flat
- Recommended torque: 25 Ncm








Implant connection	D 4,8	D 4,8	D 4,8
Gingiva height (mm)	0,6/3,0	1,6/4,0	3,1/5,5
Article No.	EV 8330	EV 8331	EV 8332
Abutment screw	EV 66	EV 66	EV 66



Multi-unit prosthetic components

- EV 8501 Multi-unit bridge screw Hex 1,26
- EV 8508 Multi-unit titanium base incl. bridge screw
- EV 8505 Multi-unit titanium cap incl. bridge screw
- EV 8506 Multi-unit gold cap, castable incl. bridge screw
- EV 8507 Multi-unit CoCr cap incl. bridge screw
- Recommended torque: 15 Ncm

					
Material	Titanium Grade 5 CF	Titanium Grade 5 CF	Titanium Grade 5 CF	"[AU 60%; Pd 20%; Pt 19%; Ir 1%]"	CrCo alloy / CTE 14.1
Article No.	EV 8501	EV 8508	EV 8505	EV 8506	EV 8507





Multi-unit modelling sleeve

- Tecanat (PC)

	
Article No.	8510
Please note:	The Multi-unit modelling sleeve can be used with the Multi-unit titanium base and Multi-unit titanium cap.

Multi-unit accessories

- EV 8502 Multi-unit cover cap Hex 1,26
- EV 8503 Multi-unit implant pick-up incl. retention screw
- M 11-11 Placement instrument Multi-unit abutment Manual and ratchet
- 8515 T Placement instrument laboratory implant CAD/CAM Multi-unit

				
Article No.	EV 8502	EV 8503	M 11-11	8515 T

Multi-unit scanbody

- Titanium specially coated
- incl. bridge screw



Version	straight	angled
Article No.	F 8500	F 8513

Please note: The Scanbody is sterilisable and for intra-oral scanning.

Multi-unit Laboratory implant

- Titanium Grade 5 CF



Version	straight	angled	angled
Type		17°	30°
Article No.	8504	8514-17	8514-30

Multi-unit Laboratory implant CAD/CAM

- Titanium Grade 5 CF



Version	straight	angled	angled
Type		17°	30°
Article No.	8509	8515-17	8515-30

Abutment screw

- Titanium Grade 5 CF
- UG




Article No.	F 600	F 601
-------------	-------	-------

Please note: The correct choice of the abutment screw is stated with the respective abutment in the catalogue.

Multi-unit abutment straight

- Titanium Grade 5 CF
- Sterile packaged
- Recommended torque: 35 Ncm



Implant connection	NP 3,5	NP 3,5	NP 3,5
Gingiva height	1,5 mm	3,5 mm	5,5 mm
Article No.	F 8000	F 8001	F 8002

Multi-unit abutment straight

- Titanium Grade 5 CF
- Sterile packaged
- Recommended torque: 35 Ncm



Implant connection	RP 4,3/5,0	RP 4,3/5,0	RP 4,3/5,0
Gingiva height	1,5 mm	3,5 mm	5,5 mm
Article No.	F 8010	F 8011	F 8012

Multi-unit abutment angled 17°

- Titanium Grade 5 CF
- Sterile packaged
- incl. abutment screw
- Type 1 = angled over flat
- Recommended torque: 35 Ncm



Implant connection	NP 3,5	NP 3,5	NP 3,5
Gingiva height (mm)	1,1/2,5	2,1/3,5	4,1/5,5
Article No.	F 8100	F 8101	F 8102
Abutment screw	F 600	F 600	F 600

Multi-unit abutment angled 30°

- Titanium Grade 5 CF
- Sterile packaged
- incl. abutment screw
- Type 1 = angled over flat
- Recommended torque: 35 Ncm



Implant connection	NP 3,5	NP 3,5	NP 3,5
Gingiva height (mm)	0,6/3,0	1,6/4,0	3,1/5,5
Article No.	F 8300	F 8301	F 8302
Abutment screw	F 600	F 600	F 600

Multi-unit abutment angled 17°

- Titanium Grade 5 CF
- Sterile packaged
- incl. abutment screw
- Type 1 = angled over flat
- Recommended torque: 35 Ncm



Implant connection	RP 4,3/5,0	RP 4,3/5,0	RP 4,3/5,0
Gingiva height (mm)	1,1/2,5	2,1/3,5	4,1/5,5
Article No.	F 8110	F 8111	F 8112
Abutment screw	F 601	F 601	F 601

Multi-unit abutment angled 30°

- Titanium Grade 5 CF
- Sterile packaged
- incl. abutment screw
- Type 1 = angled over flat
- Recommended torque: 35 Ncm



Implant connection	RP 4,3/5,0	RP 4,3/5,0	RP 4,3/5,0
Gingiva height (mm)	0,6/3,0	1,6/4,0	3,1/5,5
Article No.	F 8310	F 8311	F 8312
Abutment screw	F 601	F 601	F 601



Multi-unit prosthetic components

- F 8501 Multi-unit bridge screw UG
- F 8508 Multi-unit titanium base incl. bridge screw
- F 8505 Multi-unit titanium cap incl. bridge screw
- F 8506 Multi-unit gold cap, castable incl. bridge screw
- F 8507 Multi-unit CoCr cap incl. bridge screw
- Recommended torque: 15 Ncm

Material	Titanium Grade 5 CF	Titanium Grade 5 CF	Titanium Grade 5 CF	"(AU 60%; Pd 20%; Pt 19%; Ir 1%)"	CrCo alloy / CTE 14.1
Article No.	F 8501	F 8508	F 8505	F 8506	F 8507



Multi-unit modelling sleeve

- Tecanat (PC)

Article No.	8510
-------------	------

Please note: The Multi-unit modelling sleeve can be used with the Multi-unit titanium base and Multi-unit titanium cap.



Multi-unit accessories

- F 8502 Multi-unit cover cap UG
- F 8503 Multi-unit implant pick-up incl. retention screw
- M 11-11 Placement instrument Multi-unit abutment Manual and ratchet
- 8515 T Placement instrument laboratory implant CAD/CAM Multi-unit

Article No.	F 8502	F 8503	M 11-11	8515 T
-------------	--------	--------	---------	--------



Multi-unit scanbody

- Titanium specially coated
- incl. bridge screw



Version	straight	angled
Article No.	H 8500	H 8513

Please note: The Scanbody is sterilisable and for intra-oral scanning.

Multi-unit Laboratory implant

- Titanium Grade 5 CF



Version	straight	angled	angled
Type		17°	30°
Article No.	8504	8514-17	8514-30

Multi-unit Laboratory implant CAD/CAM

- Titanium Grade 5 CF



Version	straight	angled	angled
Type		17°	30°
Article No.	8509	8515-17	8515-30

Abutment screw

- Titanium Grade 5 CF
- Hex 1,2



Article No.	H 62
-------------	------

Please note: The correct choice of the abutment screw is stated with the respective abutment in the catalogue.

Multi-unit abutment straight

- Titanium Grade 5 CF
- Sterile packaged
- Recommended torque: 20 Ncm



Implant connection	D 3,4	D 3,4	D 3,4
Gingiva height	1,0 mm	3,0 mm	5,0 mm
Article No.	H 8000	H 8001	H 8002

Multi-unit abutment straight

- Titanium Grade 5 CF
- Sterile packaged
- Recommended torque: 20 Ncm



Implant connection	D 4,1	D 4,1	D 4,1
Gingiva height	1,0 mm	3,0 mm	5,0 mm
Article No.	H 8010	H 8011	H 8012

Multi-unit abutment angled 17°

- Titanium Grade 5 CF
- Sterile packaged
- incl. abutment screw
- Type 1 = angled over flat
- Recommended torque: 20 Ncm



Implant connection	D 3,4	D 3,4	D 3,4
Gingiva height (mm)	0,6/2,0	1,6/3,0	3,6/5,0
Article No.	H 8100	H 8101	H 8102
Abutment screw	H 62	H 62	H 62

Multi-unit abutment angled 30°

- Titanium Grade 5 CF
- Sterile packaged
- incl. abutment screw
- Type 1 = angled over flat
- Recommended torque: 20 Ncm



Implant connection	D 3,4	D 3,4	D 3,4
Gingiva height (mm)	0,6/3,0	1,6/4,0	2,4/5,0
Article No.	H 8300	H 8301	H 8302
Abutment screw	H 62	H 62	H 62

Multi-unit abutment angled 17°

- Titanium Grade 5 CF
- Sterile packaged
- incl. abutment screw
- Type 1 = angled over flat
- Recommended torque: 20 Ncm



Implant connection	D 4,1	D 4,1	D 4,1
Gingiva height (mm)	0,6/2,0	1,6/3,0	3,6/5,0
Article No.	H 8110	H 8111	H 8112
Abutment screw	H 62	H 62	H 62

Multi-unit abutment angled 30°

- Titanium Grade 5 CF
- Sterile packaged
- incl. abutment screw
- Type 1 = angled over flat
- Recommended torque: 20 Ncm



Implant connection	D 4,1	D 4,1	D 4,1
Gingiva height (mm)	0,6/3,0	1,6/4,0	2,4/5,0
Article No.	H 8310	H 8311	H 8312
Abutment screw	H 62	H 62	H 62



Multi-unit prosthetic components

- H 8501 Multi-unit bridge screw Hex 1,2
- H 8508 Multi-unit titanium base incl. bridge screw
- H 8505 Multi-unit titanium cap incl. bridge screw
- H 8506 Multi-unit gold cap, castable incl. bridge screw
- H 8507 Multi-unit CoCr cap incl. bridge screw
- Recommended torque: 15 Ncm

Material	Titanium Grade 5 CF	Titanium Grade 5 CF	Titanium Grade 5 CF	"(AU 60%; Pd 20%; Pt 19%; Ir 1%)"	CrCo alloy / CTE 14.1
Article No.	H 8501	H 8508	H 8505	H 8506	H 8507



Multi-unit modelling sleeve

- Tecanat (PC)

Article No.	8510
--------------------	------



Please note: The Multi-unit modelling sleeve can be used with the Multi-unit titanium base and Multi-unit titanium cap.

Multi-unit accessories

- H 8502 Multi-unit cover cap Hex 1,2
- H 8503 Multi-unit implant pick-up incl. retention screw
- M 11-11 Placement instrument Multi-unit abutment Manual and ratchet
- 8515 T Placement instrument laboratory implant CAD/CAM Multi-unit

Article No.	H 8502	H 8503	M 11-11	8515 T
--------------------	--------	--------	---------	--------



Multi-unit scanbody

- Titanium specially coated
- incl. bridge screw



Version	straight
Article No.	N 8500

Please note: The Scanbody is sterilisable and for intra-oral scanning.

Multi-unit Laboratory implant

- Titanium Grade 5 CF



Version	straight
Type	
Article No.	8504

Multi-unit Laboratory implant CAD/CAM

- Titanium Grade 5 CF



Version	straight
Type	
Article No.	8509

Multi-unit abutment straight

- Titanium Grade 5 CF
- Sterile packaged
- Recommended torque: 35 Ncm



Implant connection	RN 4,8
Gingiva height	0,9 mm
Article No.	N 8010

Multi-unit abutment straight




- Titanium Grade 5 CF
- Sterile packaged
- Recommended torque: 35 Ncm



Implant connection	WN 6,5
Gingiva height	1,2 mm
Article No.	N 8020


Multi-unit prosthetic components

- N 8501 Multi-unit bridge screw Torx T6
- N 8508 Multi-unit titanium base incl. bridge screw
- N 8505 Multi-unit titanium cap incl. bridge screw
- N 8506 Multi-unit gold cap, castable incl. bridge screw
- N 8507 Multi-unit CoCr cap incl. bridge screw
- Recommended torque: 15 Ncm

					
Material	Titanium Grade 5 CF	Titanium Grade 5 CF	Titanium Grade 5 CF	"(AU 60%; Pd 20%; Pt 19%; Ir 1%)"	CrCo alloy / CTE 14.1
Article No.	N 8501	N 8508	N 8505	N 8506	N 8507





Multi-unit modelling sleeve

- Tecanat (PC)

	
Article No.	8510
Please note:	The Multi-unit modelling sleeve can be used with the Multi-unit titanium base and Multi-unit titanium cap.

Multi-unit accessories

- M 11-11 Placement instrument Multi-unit abutment Manual and ratchet
- N 8502 Multi-unit cover cap Torx T6
- N 8503 Multi-unit implant pick-up incl. retention screw
- 8515 T Placement instrument laboratory implant CAD/CAM Multi-unit

				
Article No.	M 11-11	N 8502	N 8503	8515 T

Multi-unit scanbody

- Titanium specially coated
- incl. bridge screw



Version	straight	angled
Article No.	R 8500	R 8513
Please note: The Scanbody is sterilisable and for intra-oral scanning.		

Multi-unit Laboratory implant

- Titanium Grade 5 CF



Version	straight	angled	angled
Type		17°	30°
Article No.	8504	8514-17	8514-30

Multi-unit Laboratory implant CAD/CAM

- Titanium Grade 5 CF



Version	straight	angled	angled
Type		17°	30°
Article No.	8509	8515-17	8515-30

Abutment screw

- Titanium Grade 5 CF
- Hex 1,26



Article No.	R 64
Please note: The correct choice of the abutment screw is stated with the respective abutment in the catalogue.	

Multi-unit abutment straight







- Titanium Grade 5 CF
- Sterile packaged
- Recommended torque: 30 Ncm



Implant connection	D 3,5	D 3,5	D 3,5
Gingiva height	1,0 mm	3,0 mm	5,0 mm
Article No.	R 8000	R 8001	R 8002

Multi-unit abutment straight

- Titanium Grade 5 CF
- Sterile packaged
- Recommended torque: 30 Ncm







			
			
Implant connection	D 4,5	D 4,5	D 4,5
Gingiva height	1,0 mm	3,0 mm	5,0 mm
Article No.	R 8010	R 8011	R 8012

Multi-unit abutment angled 17°

- Titanium Grade 5 CF
- Sterile packaged
- incl. abutment screw
- Type 1 = angled over flat
- Recommended torque: 30 Ncm

Type 1









			
			
Implant connection	D 3,5	D 3,5	D 3,5
Gingiva height (mm)	1,1/2,5	2,1/3,5	4,1/5,5
Article No.	R 8100	R 8101	R 8102
Abutment screw	R 64	R 64	R 64

Multi-unit abutment angled 30°

- Titanium Grade 5 CF
- Sterile packaged
- incl. abutment screw
- Type 1 = angled over flat
- Recommended torque: 30 Ncm

Type 1









			
			
Implant connection	D 3,5	D 3,5	D 3,5
Gingiva height (mm)	0,6/3,0	1,6/4,0	3,1/5,5
Article No.	R 8300	R 8301	R 8302
Abutment screw	R 64	R 64	R 64

Multi-unit abutment angled 17°

- Titanium Grade 5 CF
- Sterile packaged
- incl. abutment screw
- Type 1 = angled over flat
- Recommended torque: 30 Ncm

Type 1



			
			
Implant connection	D 4,5	D 4,5	D 4,5
Gingiva height (mm)	1,1/2,5	2,1/3,5	4,1/5,5
Article No.	R 8110	R 8111	R 8112
Abutment screw	R 64	R 64	R 64

Multi-unit abutment angled 30°

- Titanium Grade 5 CF
- Sterile packaged
- incl. abutment screw
- Type 1 = angled over flat
- Recommended torque: 30 Ncm



Implant connection			
Gingiva height (mm)	0,6/3,0	1,6/4,0	3,1/5,5
Article No.	R 8310	R 8311	R 8312
Abutment screw	R 64	R 64	R 64

Multi-unit prosthetic components

- R 8501 Multi-unit bridge screw Hex 1,26
- R 8508 Multi-unit titanium base incl. bridge screw
- R 8505 Multi-unit titanium cap incl. bridge screw
- R 8506 Multi-unit gold cap, castable incl. bridge screw
- R 8507 Multi-unit CoCr cap incl. bridge screw
- Recommended torque: 15 Ncm

Material					
Article No.	R 8501	R 8508	R 8505	R 8506	R 8507

Multi-unit modelling sleeve

- Tecanat (PC)

Article No.		8510
Please note:	The Multi-unit modelling sleeve can be used with the Multi-unit titanium base and Multi-unit titanium cap.	

Multi-unit accessories

- M 11-11 Placement instrument Multi-unit abutment Manual and ratchet
- R 8502 Multi-unit cover cap Hex 1,26
- R 8503 Multi-unit implant pick-up incl. retention screw
- 8515 T Placement instrument laboratory implant CAD/CAM Multi-unit

Article No.				
	M 11-11	R 8502	R 8503	8515 T

**Placement instrument
Multi-unit abutment straight**

- Stainless steel



Version

Manual and ratchet

Article No.

M 11-11

**Placement instrument
laboratory implant CAD/CAM
Multi-unit**

- Stainless steel



Article No.

8515 T

Placement instrument Hex 1.26

• Hardened stainless steel



Version	Contra-angle	Contra-angle	Manual and ratchet	Manual and ratchet
Type	long	short	long	short
Article No.	M 03-1	M 09-1	M 10-1	M 11-1
Suitable for:	I-Series, R-Series, S-Series, Medentiloc abutment straight			

Placement instrument Hex 1.2

• Hardened stainless steel



Version	Contra-angle	Contra-angle	Manual and ratchet	Manual and ratchet
Type	long	short	long	short
Article No.	M 03-2	M 09-2	M 10-2	M 11-2
Suitable for:	T-Series, H-Series			

Placement instrument Hex 1.0

• Hardened stainless steel



Version	Contra-angle	Contra-angle	Manual and ratchet	Manual and ratchet
Type	long	short	long	short
Article No.	M 03-3	M 09-3	M 10-3	M 11-3
Suitable for:	Y-Series			

Placement instrument Torx T6

• Hardened stainless steel



Version	Contra-angle	Contra-angle	Manual and ratchet	Manual and ratchet
Type	long	short	long	short
Article No.	M 03-4	M 09-4	M 10-4	M 11-4
Suitable for:	N-Series			

Placement instrument UG

• Hardened stainless steel



Version	Contra-angle	Contra-angle	Manual and ratchet	Manual and ratchet
Type	long	short	long	short
Article No.	M 03-5	M 09-5	M 10-5	M 11-5
Suitable for:	E-Series, F-Series, K-Series			

Placement instrument MedentiLOC abutment, angled

- Hardened stainless steel
- Ball Hex



Version	Contra-angle	Manual and ratchet
Article No.	M 03-7	M 10-7

Placement instrument MedentiBASE abutment

- Hardened stainless steel



Version	Manual and ratchet
Article No.	M 11-6

Torque ratchet

- with infinitely variable torque setting
- Hardened stainless steel



Article No.	M 12
-------------	------

ISO shank adapter

- Ratchet M12
- Hardened stainless steel



Article No.	M 14
-------------	------

Altus gingiva gauge



Article No.	2013.901
-------------	----------

MedentiBASE finisher Connection geometry

- Tungsten carbide



Version	Connection geometry	Screw head	Connection geometry
Article No.	M 17	M 18	M 19

Please note: This is used to finish the fitting surfaces of the MedentiBASE plastic cap 4710 after casting. For the MedentiBASE plastic cap N 4710 the Finisher M 19 is required. The Finishers are used in the Processing holder prosthetic M 04.

**Screw channel block-out pin
Titanium base**

•Aluminium



Version	blue coated
Article No.	M 05-1
Suitable for:	H 1000 - H 1020, H 1100 - H 1120, H 1500 - H 1520, H 800 - H 820, K 1100, K 1500, K 800, K 1000, N 1100, N 800, N 1000, S 1030, S 1130, S 1530, T 1000 - T 1020, T 1100 - T 1120, T 1500 - T 1520, T 800 - T 820

**Screw channel block-out pin
Titanium base**

•Aluminium



Version	black coated
Article No.	M 05-2
Suitable for:	F 1100, F 1200, F 1300, F 1500, F 800, F 1000, S 1100, S 1200, S 1300, S 1500, S 800, S 1000

**Screw channel block-out pin
Titanium base**

•Aluminium



Version	red coated
Article No.	M 05-3
Suitable for:	S 1020, S 1220, S 1120, S 1320, S 1520, S 820

**Screw channel block-out pin
Titanium base**

•Aluminium



Version	yellow coated
Article No.	M 05-4
Suitable for:	E 1000 - E 1030, E 1100 - E 1130, E 1500 - E 1520, E 800 - E 830, F 1110, F 1210, F 1310, F 1510, F 810, F 1010, I 1000 - I 1020, I 1100 - I 1120, I 1500 - I 1520, I 800 - I 820, K 1110, K 1510, K 810, K 1010, N 1010, N 1020, N 1110, N 1120, N 1510, N 1520, N 810, N 820, R 1000 - R 1020, R 1100 - R 1120, R 1500 - R 1520, R 800 - R 820

**Screw channel block-out pin
Titanium base**

•Aluminium



Version	green coated
Article No.	M 05-5
Suitable for:	EV 1020 - EV 1040, EV 1120 - EV 1140, EV 1220 - EV 1240, EV 1320 - EV 1340

Compartment box Abutments and scanbodies

- 30 compartments
- Polystyrene



Version

without contents

Article No.

MB 04

Tray Eco

- sterilisable
- incl. silicone pad



Article No.

M 76

Tool Box

*without contents



Article No.

M 89



Publisher: MEDENTiKA® GmbH
Hammweg 8-10
76549 Hügelshheim
Tel: +49 (0)7229 69912-10
info@medentika.de
www.medentika.com

Design: The WeberFink GbR
graphic design studio
www.weberfink.de

Date: February 2019

We are certified:
DIN EN ISO 13485
Medizinprodukterichtlinie 93/42/EEC, Annex II

CE0483

Technical changes and errors reserved.

You can find the Instructions for use and warranty conditions on the website www.medentika.com

More information on the warranty can also be requested directly from the manufacturer

Distributed by Straumann USA

MEDENTiKA®

A Straumann Group Brand

60 Minuteman Road

Andover, MA 01810

Phone 800/448 8168 (US)

Fax 978/747 2490

www.medentika.us

USLIT.2066 4/20 V1 PMR
medentika.us/ifu

» Passion for
Precision «

© MEDENTiKA® 2020. All rights reserved.

MEDENTiKA® and/or other trademarks and logos from MEDENTiKA® that are mentioned herein are the trademarks or registered trademarks of MEDENTiKA® and/or its affiliates. All rights reserved.

Products indicated with ® are registered brand names of respective manufacturers.