


Mackworth Easy Fit Gantry System






User Manual

Contents

1. Safety Instructions and Warnings	Page 3
• Symbols Used	Page 3
• Contraindications	Page 4
• Intended Use/ Overview of the Easy Fit System	Page 4
• Caution	Page 5
2. Components of the Easy Fit System	Page 6
3. Specifications	Page 7
4. Placement of the Easy Fit System in a room	Page 8
5. Assembly Instructions	Page 9-12
6. Completed Easy Fit System track	Page 13
7. Disassembly Instructions	Page 14-16
8. Troubleshooting	Page 17-18
9. Pre-use inspection and maintenance	Page 19
10. Cleaning Instructions	Page 20
• Dilution Chart	Page 21
• Safety Precautions	Page 22
11. Standards Applied	Page 22
18. Service record history	Page 23-26
Warranty	Page 27

1. Safety Instructions and Warnings

Symbols used

	Warning – Consult instructions for use
	Manufacturer
	Please observe local laws on recycling

CAUTION: DO NOT ATTEMPT TO USE THIS EQUIPMENT WITHOUT FIRST UNDERSTANDING THE CONTENTS OF THIS MANUAL

Before using this equipment, and to ensure the safe operation of your Easy Fit System, carefully read this entire manual, especially the section on “Cautions”.

Failure to comply with warnings in this manual may result in injury to the operator, or the individual being lifted/transferred. Damage to the lift and/or related components may occur. Be sure that the contents of this manual are completely understood prior to using this piece of equipment.

Store this manual with the documents included with the lift system and sling(s). Contents of this manual are subject to change without prior written notice.

Contraindications

There are no known “contraindications” associated with the use of the Easy Fit System and its accessories, provided they are used as per manufacturer’s recommendations and guidelines. However, it is recommended that a client specific assessment is completed by a trained and knowledgeable health care professional to determine the method of transfer.

Mackworth does not recommend a required number of caregivers for the use of our products. This information and recommendation can only be provided after a thorough personalized, case specific assessment, as there are many factors that can influence these decisions.

It is however, “obligatory” that a client that is assessed as being an independent user of our ceiling hoist technology have the ability to receive assistance, during the transfer, in the event of a hoist malfunction or personal concern. This assistance can be provided in the form of; a nearby qualified caregiver, a phone, a communication device etc.

Intended use/ Overview of the Easy Fit Gantry System

The Easy Fit System is a lifting aid used by health care professionals and those providing care in the home to lift, position and transfer clients or a disabled family member.

The Easy Fit System is part of what is termed, overhead lift technology, which takes advantage of lifting from above and not from below or the side. The Easy Fit System makes it possible to move mobility impaired individuals with minimal strain or risk to the caregiver, while providing complete safety, dignity and comfort for the client or family member.

The Easy Fit System is designed to work with a portable hoist system such as the P-200. Use of other portable lift units with this system may be possible. However, please contact your local authorized dealers in order to obtain approval before use.

The Easy Fit System can be used to lift and transfer individuals from a bed, chair, wheelchair or similar. The length is adjustable and thus can be assembled to suit a wide range of applications. The Pressure Fit System is easy to assemble and can be completed by just one person in a short period of time. Additionally, no tools are required. It is also light in weight and once dismantled, can be moved to another location such as a hotel, apartment, or home.

Please review the following pages that outline the parts included with your package. Should you have any questions about this product or its use contact your local authorized dealer.

Caution




The Easy Fit System must be assembled prior to use. Should you have any questions during assembly contact your local authorized dealer.

- Under no circumstance should the track, lift and sling(s) or entire system be put in control of a person who has not been properly trained in the use and care of this equipment. Failure to adhere to this warning may result in serious injury to the operator, and/or the individual being lifted/transferred.
- Never install the Bottom foot plate on the top portion of post section or vice versa. See section ‘**Assembly**’ instructions.
- The Easy Fit System and associated lift and sling(s) are not toys. Do not use them for unsafe practices. Do not allow children to play with this equipment or any of its components.
- The manufacturer’s warranty is void if persons unauthorized to perform work undertake this on the Easy Fit System.
- Ensure rubber feet are cleaned and dry prior to each use. See section “**Cleaning Rubber Base Instructions.**”
- Ensure that the pressure gauge is in the safe zone during set up.
- Check the posts during installation are level both visually and with master spirit level.
- In places where more than one operator will be responsible for using the Easy Fit System, associated lift and sling(s), it is imperative that all operators be trained in its proper use.
- Never expose the Easy Fit directly to water. Warranty does not cover any misuse or abuse of the Easy Fit system.
- To maintain optimum function, the Easy Fit should be inspected and maintained on a regular basis. See the section titled “**Pre-Use Inspection and Maintenance**”.
- Any accessories used with the Easy Fit System including lift and sling(s), should be checked to ensure that they are in good working order. Check for signs of wear or fraying prior to use. Report any unusual wear or damage immediately to your local authorized dealer.
- The Easy Fit System and associated lift, and sling(s) are intended only for lifting and transferring of a person. The manufacturer will not be responsible for any damage caused by the misuse, neglect or purposeful destruction of the lift and/or its associated components.
- Do not in any circumstance exceed the maximum load for this piece of equipment. Refer to the “**Specifications of the Easy Fit System**” section of this manual.
- The installation of the Easy Fit System, lift, accessories, and sling are certified to a maximum load. Do not exceed the maximum rated load of any of the components.
- Ensure that a clear space is maintained around the Easy Fit System. Remove all furniture and other obstacles out of the way before performing a transfer.
- The Easy Fit System has been designed to lift vertically at its maximum load and at maximum height. Do not attempt to lift an individual at an angle to the track.

2. Components

Easy Fit Components




Component List

Description	Quantity
Bottom plate	2
Top plate	2
Supporting post	2
Adjustable track	1
Post pins	2
Pin 'R' Clips (for post pins)	2
Owner's manual	1

3. Specifications

- High-strength lightweight Aluminium
- Portable or semi-permanent installation
- Eliminates the mounting of a connection system in the ceiling
- Pressure gauge indicator to help ensure proper installation
- Quick attachment with secure fixings
- Addresses a number of different bedroom & bathroom situations
- Weight capacity 440 lbs
- Adjustable in length
- Built-in easy-slide trolley
- Exerting pressure between Floor and Ceiling is 80 to 100 lbs.
- Specially designed and engineered rubber moulded feet for ultimate grip
- Levelling indicators on posts
- One (1) year warranty


Maximum & Minimum dimensions (fit)



Model table for Easy Fit

Code	Description
540100	Easy Fit – 2 post version

4. Placement of the Easy Fit System in a room



Drawing showing front view of a common Easy Fit System set up in a room. In this drawing arrangement, a general placement of a bed is shown in relationship to the system. Notice that there is still room for a wheelchair to be placed beside the bed.



DO NOT GO OUTSIDE OF THIS RANGE OF PLACEMENT OTHERWISE INJURY MAY OCCUR TO THE OPERATOR OF THE LIFT AND/OR THE INDIVIDUAL BEING TRANSFERRED, AND/OR THE EASY FIT SYSTEM

5. Assembly



Before using the Easy Fit System, the top plate, bottom plate, supporting posts and adjustable track must be visually checked to ensure that there are no missing parts or unusual wear and tear. Top plate and Bottom plate should be cleaned and dry with wet cloth before installation.

Should anything look unusual contact your local dealer prior to use. Failure to comply with this caution could result in serious injury to the operator, the individual being lifted and/or damage to the Easy Fit System and/or port-able lift unit.


1. Place the circular Bottom Plate on the floor in required position for use




2. Insert the bottom of the Support Post to the Bottom Plate



3. Insert the Top Plate onto the top of the Support Post




4. Depress the red button and release the Safety Lock



Assembly continued

5. Depress the Ratchet Button to release the Lever, then push the Lever up until it meets the Magnet Symbol. Note: Pushing down on the top part of the Support Post may help to release the Lever



6. Push the top part of the Support Post upwards until the Top Plate compresses against the ceiling. Push until you hear two clicks on the ratchet



7. Open the Post Level (opposite the Safety Lock) and check that the bubble is within the central circle by adjusting the position of the Support Post




8. Pull the lever down until it flat to the support leg and clicks into place



Assembly (continued).....


9. Check to ensure that the blue marker on Pressure Gauge (located at the top of the Support Post) is comfortably into the green 'safe' zone on both sides. If it is not, it will be necessary to release the Lever again and apply more than two clicks before closing it.



10. Push the Safety Lock back into place until it clicks and rotate the post to ensure that the lever is facing 90 degrees to the position of the second post.



11. Repeat the procedure to built the second Support Post and set it perpendicular to the first.



12. Lift the Track (N.B. You may need assistance from another person and steps to reach) and hook one end on to the Support Post Pins.



(Assembling (continued ...))


13. Extend or contract the Track to fit between the two Support Posts and hook it on to the Support Post pins at the other side.



14. Push the Post Pins



15. Push the Safety Pin



7. Disassembly Instructions



Failure to read and follow these instructions carefully may result in serious injury.

Track disassembly

1. Unhook the track from one post end by releasing the locking pin and clip then proceeding to lift the track from the hook as indicated by the arrow in (fig 3a.)



fig 3a

2. Squeeze the track and unhook the track from other post end (fig 3b.) below. Once the track has been removed, store the track in a safe place.




fig 3b.

3. Press the red button to release the safety lock as shown in (fig 4a.) with the red directional arrow, indicated, to show the direction of the opening mechanism for the lock




4. Press the lever in all the way to release the pressure from the ratchet as shown below



Disassembly Instructions (continued...)

5. Push in and hold the ratchet




7. Hold ratchet and pull out the lever as shown in the image below.



Lever may swing out – always keep one hand on the lever during this step

8. Push the lever up to engage with magnet




9. Once this has been completed, push in and hold the ratchet while pushing the top post down



Disassembly Instructions (continued...)

10. Make sure bottom pin touches the plastic cover.



11. Bring back the level and safety lock to original position and lock in place shown in the image opposite.



12. Take off top foot.
Store the components in a safe place for further use.



13. Take off the bottom foot.







NOTE - Repeat steps 1 to 7 for the other post disassembly.

8. Troubleshooting

Should problems arise with the use of the Easy Fit System review the following chart. Find the fault and complete the recommended solution. If the fault is not found and/or the solution does not correct the problem, contact your local authorized dealer for service immediately.

Fault	Reason	Recommended Solution
Visually the system is not aligned straight	<ol style="list-style-type: none"> 1. Level is not working. 2. Floor is not levelled. 3. Level is damaged. 4. Post is knocked with impact 	<ol style="list-style-type: none"> 1. Set up the post following "Post Assembly Instructions". Check visually and with master level to ensure the post is aligned straight. 2. Check the floor with master level for flatness. 3. Visually check the level is fine. 4. Dismantle the post following "Post Disassembly Instructions" and set up again following "Post assembly instructions". If problem persists, contact your local Authorized dealer immediately
Feet are slipping after installation	<ol style="list-style-type: none"> 1. Feet are not cleaned. 2. Floor is not cleaned. 	<ol style="list-style-type: none"> 1. Clean the feet as per "Cleaning the Rubber Base" Instructions. 2. Ensure floor area is dry and free of obstacles. If problem persists, contact your local authorized dealer immediately .
Lever is not engaged with magnet.	<ol style="list-style-type: none"> 1. Post Assembly instructions not followed 2. Clutch bearing function fails. 3. Lever not cleaned. 	<ol style="list-style-type: none"> 1. Press in the ratchet and hold it. Grab the lever and push up to the magnet label. 2. Contact local authorized dealer immediately. 3. Ensure lever is cleaned and free of grease or dirt. <p>If problem persists, contact your local authorized dealer immediately</p>
Pressure gauge is not in safe zone.	<ol style="list-style-type: none"> 1. Post assembly instructions not followed. 2. Barrel nut is loose or not secured tightly. 3. Top plate Spring broken or bent. 	<ol style="list-style-type: none"> 1. Set up the post following post assembly instructions. 2. Tight the barrel nut with flat screw driver. <p>If problem persists, contact your local authorized dealer immediately</p>
The adjustable track does not extend smoothly or doesn't extend up to required limit.	<ol style="list-style-type: none"> 1. Physical damage to track. 2. Track slots are not cleaned. 	<ol style="list-style-type: none"> 1. Check for physical damage to track. If so, contact local authorized dealer immediately 2. Clean the track slots with cloth and wipe out the dirt. Again try extending the tracks. <p>If problem persists, contact your local authorized dealer immediately .</p>

Fault	Reason	Recommended Solution
Trolley does not move smoothly along the adjustable track.	<ol style="list-style-type: none"> 1. Physical damage to track or trolley wheels. 2. Track slots or trolley wheels are not cleaned. 	<ol style="list-style-type: none"> 1. visually make sure tracks and wheels are fine and are free of obstacles. 2. Ensure that tracks are cleaned with cloth and free of dirt <p>If problem persists, contact local authorized dealer immediately</p>
The adjustable track is not locking on the Post Pins.	<ol style="list-style-type: none"> 1. One or more track connector bracket(s) are damaged or broken. 2. Track lock pin(s) is(are) broken or damaged. 	 <p>Do not use the system. Contact your local authorized dealer immediately.</p>
Misalignment of adjustable track on the post pins. (for example; track is not sitting straight on the post assembly)	<ol style="list-style-type: none"> 1. Track connector hooks are not Sitting right on the post pins. 	<ol style="list-style-type: none"> 1. Check Mounting Adjustable Track Assembly Instructions in the user manual 2. Check visually to ensure track connector hooks are sitting at same height <p>If problem persists, contact local authorized dealer immediately.</p>
Post disassembly problem for example; cannot retract post.	<ol style="list-style-type: none"> 1. Ratchet is not held during disassembly. 2. Lever arm is not engaged with magnet. 	<ol style="list-style-type: none"> 1. Follow Post Disassembly instructions to dismantle the post 2. Follow Post Disassembly instructions to dismantle the post
Safety lock not working.	<ol style="list-style-type: none"> 1. Red button is not pressed in enough. 2. Physical damage to safety lock cover. 3. Lock latch Spring breaks 	<ol style="list-style-type: none"> 1. Press in the Red button all the way so that safety lock can be engaged 2. Contact local authorized dealer immediately 3. Contact local authorized dealer immediately
Difficult to insert top and bottom foot into post assembly.	<ol style="list-style-type: none"> 1. Top post spacer damaged 2. Bottom post plug damaged. 3. Bottom plate foot post damaged or bend. 4. Top plate foot post damaged or bend. 	 <p>Do not use the system. Contact your local authorized dealer immediately.</p>
Track connector bracket are loose.	<ol style="list-style-type: none"> 1. Rivets are broken or loose. 	 <p>Do not use the system. Contact your local authorized dealer immediately.</p>
Parts missing	<ol style="list-style-type: none"> 1. Locking pin missing. 2. Locking clip missing 	 <p>Do not use the system. Contact your local authorized dealer immediately.</p>



Please ensure any accessory used with this system is stored safely when not in use to prevent any damage or any unintended use which would introduce an inadvertent risk to minors who may be present around the system

9. Pre-Use Inspection and Maintenance



Prior to each lift, this will be completed by the user, the Easy Fit system and associated hoist, accessories and sling(s), must be visually inspected. Refer to the hoist, accessory and sling user guides for specific details regarding their inspection.

Should any of these items below fail the inspection do not use the Easy Fit system - Contact your local authorized dealer for service

Visually check for the following:

- The Easy Fit system parts should not be loose, broken or bent
- The bases are on a stable, level surface and the side support assemblies are secure and in the green safe zone (see (fig 5a.) below)
- Ensure that both side support assemblies are completely vertical from all directions (see (fig 5b.) below)
- The trolley moves easily along the track
- The top plate and bottom plate must be cleaned with wet cloth before installation in order to remove dirt on rubber moulds. See “cleaning the rubber base” section on page 26
- The adjustment track is placed in the correct location for transferring
- The portable hoist has been inspected as outlined in its owner’s manual
- The sling (s) has been inspected as outlined in its owner’s manual.




fig 5a



fig 5b

1. See (fig 5a.) - Ensure Pressure guide is in the Safe Zone. **Note: Green section must be visible on both sides of the guide.**

2. See (fig 5b.) - Ensure both sides are vertical from the front and the side

10. Cleaning



Prior to each installation, top and bottom rubber feet need to be cleaned with a wet cloth soaked with water.

1. Take bottom rubber foot and wipe out the dirt attached to its bottom side of foot.
2. Take top rubber foot and wipe out the dirt attached to its bottom side of foot. As shown in the picture

Note: After cleaning, make sure rubber feet are cleaned and dry with no residue of dirt or dust left on them. Also, ensure floor area is cleaned and free of obstacles.

Recommended Cleaning Instruction

The exterior of the Easy Fit System can be cleaned using a damp soapy cloth for general cleaning duties. Please ensure the cloth is damp and not wet. Ensure the exterior of the device is dry after cleaning, dry using a clean dry cloth.



Care should always be taken when cleaning around electrical components where these are installed.

Disinfecting (if necessary)

Should the Easy Fit System require a more thorough clean, the use of the Actichlor disinfectant product, which is widely available in tablet form and used throughout the healthcare industry, is recommended for use to ensure a thorough clean.

- Please follow the manufacturer's safety instruction for the use of the cleaning product before use to ensure safe use for the operator and the patient.



- Ensure the cloth is **damp** before the cleaning process. **DO NOT** use a wet cloth over electrical systems.
- Be careful not to let water ingress into the device as although the device is IP rated, it is not water tight.

Application is through a clean soaked (but damp) cloth applied to wipe the device down

Used in the following dilutions to ensure an effective clean:

- Actichlor dissolvable chlorine tablets provide a concentration of 1000 ppm of available chlorine (0.1%) per 1 tablet
- 1 tablet (1.7g formed tablet (x1)) will create a virucidal solution, diluted in 1 litre of water to provide effective means to clean a "Dirty" device. This is also ideal for use after an outbreak of the Norovirus / winter vomiting and can be used as a precaution against C.Diff. It is effective against viruses, bacteria, spores, yeasts and moulds.
- The contact time against the outer components of the device should be for 5 minutes to prevent any virucidal infections without a degradation to the functionality of the device. 5 minutes is a recommended contact time. The device can withstand a longer contact period but the 5 minute recommendation as a minimum must be followed to provide an effective cleaning regime.
- Blood spills should be dealt with by an increased concentration of the solution – please refer to the instructions on the manufacturer's product labelling.

Concentration limits for differing cleans will be shown on the manufacturer's Actichlor tablet container, however, this is reflected in the table below:

Dilution chart

Product used as:	Device Condition	Concentration (ppm)	Dilution qty*	Tablets per 0.26gal (1l)	Contact time
Bactericidal	Clean	200	1.32gal (5l)	1	1 minute
	Dirty	1000	0.26gal (1l)	1	5 minutes
Yeasticidal	Clean	200	1.32gal (5l)	1	1 minute
	Dirty	1000	0.26gal (1l)	1	5 minutes
Fungicidal	Clean	2000	0.26gal (1l)	2	15 minutes
	Dirty	5000	0.26gal (1l)	5	15 minutes
Mycrobactericidal	Clean	1000	0.26gal (1l)	1	15 minutes
	Dirty	5000	0.26gal (1l)	5	15 minutes
Virucidal	Clean	500	0.53gal (2l)	1	5 minutes
	Dirty	1000	0.26gal (1l)	1	5 minutes
Sporcidal (C.Diff)	Clean	1000	0.26gal (1l)	1	10 minutes
	-	-	-		
Sporcidal	Clean	5000	0.26gal (1l)	5	10 minutes
	-	-	-		

* Dilution is made within water

- When diluted in water, one tablet gives 1000 ppm of available chlorine DO NOT dilute within any other medium
- The concentration of the solution depends upon whether the device being cleaned is noticeably dirty or note (indicated in the table by “Device Condition”)

Safety precautions when using this cleaning agent

Handling and Storage:

Advice on Safe Handling



Avoid contact with skin and eyes. Do not breathe dust / fumes / gas / mist / vapors / spray. Use only with adequate ventilation. Wash hands thoroughly after handling. Mixing this product with acid or ammonia releases chlorine gas

Hygiene Measures

Handle in accordance with good industrial hygiene and safety practice. Remove and wash contaminated clothing before re-use. Wash face, hands and any exposed skin thoroughly after handling.

Conditions for safe storage, including and incompatibilities

- Keep out of reach of children
- Keep container tightly closed
- Store in suitable labelled containers
- Storage temperature: 32°F to 77°F (0 to 25°C)

Individual protective measures:

- Hand protection: Gloves

Dissolve

Dissolve in cold water – With no agitation, 1 tablet will take approx. 10 minutes to fully dissolve in the water used.

The information above has been extracted from the Actichlor MSDS (Manufacturers Safety Data Sheet). For a full review of the data please follow the link below:

<http://www.nhsggc.org.uk/media/236215/msds-actichlor-plus.pdf>

11. Standards Applied

The standards that have been applied to the device are as follows:

- BS EN ISO 10535

Hoists for the transfer of disabled persons. Requirements and test methods

- BS EN ISO 12182

Assistive products for persons with disability. General requirements and test methods

12. Service Record History

18.1 Initial Information

- Complete the following section on **Purchase and Service Information** as soon as this equipment is supplied.
- Use the service record history to record to any completed service and repairs.
- Ensure that the service record is signed and dated each time it is used.
- Be sure to have this piece of equipment serviced on a regular basis.

Purchase Information

Product name: Easy Fit System

Date of purchase: **Serial No.:**

Purchased from:

Address:

City: **State:** **ZIP code:**

Telephone number:

Comments:

Contact the following company for service

Company:

Address:

City: **State:** **ZIP code:**

Telephone number:

Comments:

Service Record History

Complete this section after each service, repair inspection and/or maintenance.

* Photocopy additional pages as required*

Date:		Time:	
Service Type: <input type="checkbox"/> Periodic inspection <input type="checkbox"/> Monthly inspection <input type="checkbox"/> 6 inspection <input type="checkbox"/> Repair <input type="checkbox"/> Yearly inspection <input type="checkbox"/> Other			
Completed by:		(signature)	
Company:			
Remarks & Action Taken:			
Device left in a safe usable condition: YES <input type="checkbox"/> NO <input type="checkbox"/> if 'NO' explain in remarks the action) taken			

Date:		Time:	
Service Type: <input type="checkbox"/> Periodic inspection <input type="checkbox"/> Monthly inspection <input type="checkbox"/> 6 inspection <input type="checkbox"/> Repair <input type="checkbox"/> Yearly inspection <input type="checkbox"/> Other			
Completed by:		(signature)	
Company:			
Remarks & Action Taken:			
Device left in a safe usable condition: YES <input type="checkbox"/> NO <input type="checkbox"/> if 'NO' explain in remarks the action) taken			

Date:		Time:	
Service Type: <input type="checkbox"/> Periodic inspection <input type="checkbox"/> Monthly inspection <input type="checkbox"/> 6 inspection <input type="checkbox"/> Repair <input type="checkbox"/> Yearly inspection <input type="checkbox"/> Other			
Completed by:		(signature)	
Company:			
Remarks & Action Taken:			
Device left in a safe usable condition: YES <input type="checkbox"/> NO <input type="checkbox"/> if 'NO' explain in remarks the action) taken			

Date:		Time:	
Service Type: <input type="checkbox"/> Periodic inspection <input type="checkbox"/> Monthly inspection <input type="checkbox"/> 6 inspection <input type="checkbox"/> Repair <input type="checkbox"/> Yearly inspection <input type="checkbox"/> Other			
Competed by:		(printed name) (signature)	
Company:			
Remarks & Action Taken:			
Device left in a safe usable condition: YES <input type="checkbox"/> NO <input type="checkbox"/> if 'NO' explain in remarks the action) taken			

Date:		Time:	
Service Type: <input type="checkbox"/> Periodic inspection <input type="checkbox"/> Monthly inspection <input type="checkbox"/> 6 inspection <input type="checkbox"/> Repair <input type="checkbox"/> Yearly inspection <input type="checkbox"/> Other			
Competed by:		(printed name) (signature)	
Company:			
Remarks & Action Taken:			
Device left in a safe usable condition: YES <input type="checkbox"/> NO <input type="checkbox"/> if 'NO' explain in remarks the action) taken			

Date:		Time:	
Service Type: <input type="checkbox"/> Periodic inspection <input type="checkbox"/> Monthly inspection <input type="checkbox"/> 6 inspection <input type="checkbox"/> Repair <input type="checkbox"/> Yearly inspection <input type="checkbox"/> Other			
Competed by:		(printed name) (signature)	
Company:			
Remarks & Action Taken:			
Device left in a safe usable condition: YES <input type="checkbox"/> NO <input type="checkbox"/> if 'NO' explain in remarks the action) taken			

Date:		Time:	
Service Type: <input type="checkbox"/> Periodic inspection <input type="checkbox"/> Monthly inspection <input type="checkbox"/> 6 inspection <input type="checkbox"/> Repair <input type="checkbox"/> Yearly inspection <input type="checkbox"/> Other			
Competed by:		(printed name)	
Company:		(signature)	
Remarks & Action Taken:			
Device left in a safe usable condition: YES <input type="checkbox"/> NO <input type="checkbox"/> if 'NO' explain in remarks the action) taken			

Date:		Time:	
Service Type: <input type="checkbox"/> Periodic inspection <input type="checkbox"/> Monthly inspection <input type="checkbox"/> 6 inspection <input type="checkbox"/> Repair <input type="checkbox"/> Yearly inspection <input type="checkbox"/> Other			
Competed by:		(printed name)	
Company:		(signature)	
Remarks & Action Taken:			
Device left in a safe usable condition: YES <input type="checkbox"/> NO <input type="checkbox"/> if 'NO' explain in remarks the action) taken			

Date:		Time:	
Service Type: <input type="checkbox"/> Periodic inspection <input type="checkbox"/> Monthly inspection <input type="checkbox"/> 6 inspection <input type="checkbox"/> Repair <input type="checkbox"/> Yearly inspection <input type="checkbox"/> Other			
Competed by:		(printed name)	
Company:		(signature)	
Remarks & Action Taken:			
Device left in a safe usable condition: YES <input type="checkbox"/> NO <input type="checkbox"/> if 'NO' explain in remarks the action) taken			

WARRANTY

This Warranty does not affect or in any way limit your Statutory Rights

1. All equipment supplied as new against failure within the period of 1 year from date of purchase by virtue of defects in material or workmanship.
2. This guarantee does not apply to failure attributable to normal wear and tear, damage by natural forces, user neglect or misuse or to deliberate destruction, or to batteries more than 90 days after original purchase.
3. This guarantee shall be void if the equipment is not serviced by Mackworth or its authorized service agents in accordance with the manufacturers recommendations or if any unauthorized person carries out works on the equipment.
4. The liability of Mackworth under the terms of this guarantee shall be limited to the replacement of defective parts and in no event shall Mackworth incur liability for any consequential or unforeseeable losses.



T: 437-317-4453

HighStar Healthcare Inc.
1060 Stacey Court,
Unit 1,
Mississauga
ON
L4W 2X8

Disclaimer

While every effort has been made to ensure the accuracy of information contained in this user manual, no liability can be accepted by Mackworth for any errors or omissions. Mackworth operates a policy of continuous improvement. Specifications and other data are subject to change without notice.



Mackworth