

Endocrinology



Comprehensive Assessment of Neuroendocrine Status

- Salivary Sex Steroid Hormones and Adrenal/HPA Axis Function
- Hormone and Urinary Metabolites Assessment Profile (HuMap™)
- Neurotransmitter and Metabolites Profiles
- Thyroid Profile



SCIENCE + INSIGHT

Endocrinology

Doctor's Data offers analysis of hormones, neurotransmitters, and metabolites that may impact symptoms, metabolic function, mood or behavior.

Salivary Sex Steroid Hormones and HPA Axis/Adrenal Function



Saliva Testing

Saliva contains only the active, bioavailable portion of the body's hormones, which are secreted in a pulsatile manner. Doctor's Data salivary testing accounts for the normal pulsatile secretion of hormones by averaging 4 collections to give a better idea of how best to treat. This method is noninvasive, and easy to collect, which is particularly beneficial for diurnal cortisol, CAR (Cortisol Awakening Response), and luteal surge testing. Many providers consider saliva testing for determining clinical diagnosis and to monitor any route of hormone replacement therapy. Saliva testing is the matrix of choice to monitor transdermal and sublingual hormone therapies.

Comprehensive Plus Hormone Profile; saliva

Order: SAMPLE REPORT

Client #: 12345

Doctor: Sample Doctor
Doctor's Data, Inc.
3755 Illinois Ave.
St. Charles, IL 60174

Patient: Sample Patient

Age: 39

Sex: Female

Menopausal Status: Pre-menopausal

Sample Collection Date/Time

Date Collected 04/05/2022

AM30 04/05/2022 06:00

Noon 04/05/2022 12:00

Evening 04/05/2022 17:00

Night 04/05/2022 19:00

Date Received 04/06/2022

Date Reported 04/07/2022

Analyte	Result	Unit	L	WRI	H	Reference Interval	Supplementation Range**
Estrone (E1)*	29.0	pg/mL				0 – 35	
Estradiol (E2)	0.60	pg/mL				0.6 – 4.5	1.0 – 6.0
Estriol (E3)*	<5.0	pg/mL				7.5 – 66	45 – 680
EQ (E3 / (E1 + E2)) Ratio	0.17					≥ 1.0	
Progesterone (Pg)	26	pg/mL				127 – 446	400 – 4000
Pg/E2 Ratio*	43.3					>200	≥ 200
Testosterone	7	pg/mL				6 – 49	25 – 60
DHEA*	15	pg/mL				106 – 300	

Hormone Comments

- Low estriol levels are often associated with vaginal dryness.
- Henry Lemon MD developed the Estrogen Quotient (EQ), a simple ratio of the cancer protective E3 relative to the proliferative estrogens E1 and E2, to assess breast cancer risk. A lower number (<1.0) indicates increased risk, and a higher number (>1.0) signifies lower risk. Dr. Lemon stated that for maximum protection, an optimal EQ is >1.5.
- The Estrogen Quotient (EQ) is low. Estriol supplementation is a consideration to balance this quotient and reduce associated risks.
- Progesterone to estradiol (Pg/E2) ratio is consistent with progesterone insufficiency (estrogen dominance). Supplementation with progesterone to correct this relative deficiency is a consideration depending on the clinical picture. Note: The progesterone level is suggestive of an anovulatory cycle or luteal phase defect. Query BCP usage.
- DHEA levels typically decline with age and the level measured here is below the reference range. Note: Supplementation with DHEA may increase testosterone and/or estradiol levels.
- Supplementation reference ranges are based on adherence to proper dosage interval(s). Please visit <https://www.DoctorsData.com/Resources/BestPractices.pdf> for more information.

Notes:
 The current samples are routinely held three weeks from receipt for additional testing.
R = Reference Interval, *L* (blue) = Low (below *R*), *WRI* (green) = Within *R* (optimal), *WRI* (yellow) = Within *R* (not optimal), *H* (red) = High (above *R*)
 **This test was developed and its performance characteristics determined by Doctor's Data Laboratories in a manner consistent with CLIA requirements. The U. S. Food and Drug Administration (FDA) has not approved or cleared this test; however, FDA clearance is not currently required for clinical use. The results are not intended to be used as a sole means for clinical diagnosis or patient management decisions.
 †The Pg/E2 ratio is an optimal range established based on clinical observation. Reference intervals for Pg/E2 ratio have not been established in males and post-menopausal women who are not supplementing with progesterone and/or estrogens.
 **If supplementation is reported then the supplementation ranges will be graphed. The supplementation ranges depicted are for informational purposes only and were derived from a cohort of adult men and women utilizing physiologic transdermal bioidentical hormone therapy.
 Methodology: Enzyme Immunoassay

Analyzed by DOCTOR'S DATA, INC. • 3755 Illinois Avenue, St. Charles, IL 60174-2420 USA • LAB DIR: Erto Roth, MD • CLIA ID: 14D0846470

HPA Axis/Adrenal Function Testing

Cortisol profiles from Doctor's Data require four saliva samples taken at 30 minutes post-awakening, noon, dinner and bed time. These points capture the diurnal rhythm of cortisol secretion and can identify dysregulation in the natural cortisol circadian rhythm.


The cortisol awakening response (CAR) is the natural rise in cortisol that is seen minutes after awakening, followed by a noticeable drop within 60 minutes. CAR can be utilized as a biomarker for assessment of the HPA axis function in routine clinical practice—it can be influenced by overall HPA reactivity, as well as a person's anticipation of stress. Combined with the 30-minute post-awakening collection, Doctor's Data measures CAR with two additional saliva collections—one immediately upon waking and the other at 60 minutes post-awakening. The CAR Profile can be ordered alone or added to any cortisol profile with an AM30 collection.

Secretory IgA (sIgA)


sIgA is most often measured in stool or saliva. Measuring sIgA in stool can reveal information about gut immunity, inflammation, recent or current infections, and potential acute or chronic stress generally associated with GALT (gut-associated lymphoid tissue). sIgA measured in the saliva primarily provides insight into the body's stress response. However, there is some evidence that activated B cells can migrate from GALT to salivary glands, which could potentially demonstrate systemic inflammation and possibly link GI pathology via salivary sampling.

Melatonin

Melatonin is a hormone that also acts as a neurotransmitter. It is the major indole compound synthesized by the pineal gland and is converted from serotonin. Melatonin levels follow a diurnal rhythm



Comprehensive Adrenal Function Profile + CAR; saliva



Order: SAMPLE REPORT

Client #: 12345

Doctor: Sample Doctor
Doctor's Data, Inc.
3755 Illinois Ave.
St. Charles, IL 60174

Patient: Sample Patient

Age: 39

Sex: Female

Menopausal Status: Pre-menopausal

Sample Collection **Date/Time**

Date Collected 04/05/2022

AM30 04/05/2022 06:00

Noon 04/05/2022 12:00

Evening 04/05/2022 17:00


Night 04/05/2022 19:00

Date Received 04/06/2022

Date Reported 04/07/2022

Analyte	Result	Unit	L	WRI	H	Optimal Range	Reference Interval
Cortisol Waking	7.0	nmol/L	▲	▲	▲	12–17	6.0–20
Cortisol AM30	15	nmol/L	▲	▲	▲	14.0–25.0	7.0–30.0
Cortisol AM60	12	nmol/L	▲	▲	▲	10.0–18.0	4.9–22.0
CAR Rise (Calculated)	114.3	%	▲	▲	▲		35–60
CAR Decline (Calculated)	71.43	%	▲	▲	▲		-33–0
Cortisol Noon	1.8	nmol/L	▲	▲	▲	5.0–10.0	2.1–14.0
Cortisol Evening	1.2	nmol/L	▲	▲	▲	2.0–5.0	1.5–8.0
Cortisol Night	0.44	nmol/L	▲	▲	▲	1.0–4.0	0.33–7.0
DHEA*	15	pg/mL	▲	▲	▲		106–300
Secretory IgA*	125	pg/mL	▲	▲	▲		75–330


Cortisol Graph




Adrenal Phase: 1

Hormone Comments

- AM cortisol level appears adequate, although the suboptimal diurnal cortisol pattern is suggestive of early (Phase 1) HPA axis (adrenal gland) dysfunction.
- Cortisol Awakening Response (CAR) is the expected maximum rise in cortisol levels observed at 30 minutes (AM30) post awakening (approximately 35-60% above the waking value). This is followed by an expected decline sixty minutes after waking. The behavior of cortisol is a critical marker to understand the HPA axis' physiologic responsiveness and is a key indicator of HPA axis adaptability and reactivity.
- This patient's CAR Rise (the percent change from waking to AM30) exceeds the expected degree of increase.
- This patient's CAR Decline (the percent change from waking to AM60) is observed as an unexpected increase.
- DHEA levels typically decline with age and the level measured here is below the reference range. Note: Supplementation with DHEA may increase testosterone and/or estradiol levels.
- Every mucosal membrane surface such as the eyes, nose, throat, and gastrointestinal system represent a large portal of entry for pathogenic bacteria, viruses, and yeasts. Secretory IgA (sIgA) is the predominant antibody found on these mucosal membranes, and represents the body's first line of defense. sIgA levels change in response to stress.



Melatonin Report



Order: SAMPLE REPORT

Client #: 12345

Doctor: Sample Doctor
Doctor's Data, Inc.
3755 Illinois Ave.
St. Charles, IL 60174

Patient: Sample Patient

Age: 39

Sex: Female

Menopausal Status: Pre-menopausal

Sample Collection **Date/Time**

Date Collected 04/05/2022

Date Received 04/06/2022

Date Reported 04/07/2022

Analyte	Result	Unit	L	WRI	H	Reference Interval
Melatonin Morning*	5.3	pg/mL	▲	▲	▲	3.0–25
Melatonin Evening*	2.0	pg/mL	▲	▲	▲	1.4–19
Melatonin Night*	<1.4	pg/mL	▲	▲	▲	4.3–25

in response to the light/dark cycle, with melatonin and light occurring at opposite times. Endogenous melatonin production begins rising approximately two hours before bedtime, provided light

is dim. Melatonin and cortisol share an inverse relationship—when melatonin levels are low, cortisol levels should be high, and vice versa.

Clinical Relevance and Patient Compliance

Free, unbound sex steroid hormones can fluctuate throughout the day and thus have many peaks and troughs. In order to obtain the most clinically accurate and truly representative results of an individual's hormone status, Doctor's Data averages four saliva specimens,

collected throughout the day, to minimize the risk of reporting a peak or trough. We go an extra step to deliver the most clinically relevant results—a fifth tube of saliva is created by pooling an equal amount of saliva from each of the four submitted samples. This pooled tube is mixed thoroughly to provide homogenization and becomes the saliva

source from which estriol, estradiol, estrone, progesterone, testosterone and DHEA are measured. This additional effort provides a far superior reflection of each patient's hormonal status. Convenient, easy and low-stress sample collection, with no needles or 24-hour urine collection, greatly increases patient compliance.

Doctor's Data pooled tube methodology equalizes hormone peaks and troughs for
a more precise illustration of your patients' sex hormone status

Salivary Hormone and HPA Axis/Adrenal Function Profiles



	Comprehensive Plus Profile	Comprehensive Hormone Profile	Basic Hormone Profile	Comprehensive Adrenal Function Profile	Adrenal Function Profile	Diurnal Cortisol Profile	CAR Profile	Melatonin Profile
Estrone	✓							
Estradiol	✓	✓	✓					
Estriol	✓							
Progesterone	✓	✓	✓					
Testosterone	✓	✓	✓					
DHEA	✓	✓	✓	✓	✓			
Cortisol	x4	x4	x2	x4	x4	x4	x3	
Melatonin								x3
slgA				✓				
PG/E2 Ratio	✓	✓	✓					
Estrogen Quotient	✓							


CAR, slgA and melatonin can be added to profiles. Smaller profiles and single-analyte tests are also available. Call Doctor's Data for assistance in selecting the tests that will maximize value for your patients.

Hormone and Urinary Metabolites Assessment Profile (HuMap™)



The Hormone and Urinary Metabolite Assessment Profile (HuMap™) provides a complete overview of steroid hormones, their metabolites, and the efficiency of the enzymes that metabolize these hormones. This non-invasive test requires only 4 or 5 separate urine collections. Because the breakdown of hormones relies so heavily on processes within the liver, this test can also elucidate areas of interest as it pertains to conjugation of each metabolite. Additionally, testing urinary hormone metabolites can contribute to further understanding of endogenous hormone secretion, supplemental hormone utilization, enzyme activity, oxidative stress, and insight into whether your body is safely metabolizing hormones.

- Color-coded results and expanded metabolic pathways specific to each patient
- Summary page to highlight important findings
- Key enzyme activity
- Full commentary



Comprehensive Hormone and Metabolite Markers Available in the HuMap™

- Extensive estrogen metabolites available: E1, E2, E3, 2-OH-E1/E2, 4-OH-E1/E2, 16-OH-E1, 2-Methoxy-E1/E2, 4-Methoxy-E1/E2
- Comprehensive androgens and metabolites
- Progesterone and metabolites
- Daily cortisol (5), Daily cortisone (5)
- Metabolized cortisol, cortisone, and corticosterone
- Biomarker for Oxidative DNA Damage: 8-hydroxy-2'-deoxyguanosine (8-OHdG)
- Key enzyme activity: COMT, 11βHSD2, aromatase, and 5α-reductase



The HuMap™ is the most comprehensive urinary metabolite profile offered by Doctor's Data and is designed to provide deep insight into the secretion and metabolism of unconjugated sex hormones and their metabolites, and the enzymes that fuel

these conversions. Doctor's Data also offers smaller profiles that assess a more select set of hormones and metabolites, allowing for more targeted assessment for initial assessment as well as follow-up testing.

	 HuMap™	Sex Hormones Profile	Androgens and Progesterones Profile	Adrenal Corticoids Profile	Estrogen Metabolites Profile
Estrogens Fractionated: E1, E2, E3, 2-OH-E1, 4-OH-E1, 16-OH-E1, 2-M-E1, 4-M-E1, 2-OH-E2, 4-OH-E2, 2-M-E2, and 4-M-E2	✓	✓			✓
Androstenedione (A4)	✓	✓	✓		
EPI-Testosterone (EPI-T)	✓	✓	✓		
Testosterone (T)	✓	✓	✓		
Androsterone (AN)	✓	✓	✓		
11-hydroxy-Androsterone (OHAN)	✓	✓	✓		
5a-Androstenediol (5A-AD)	✓	✓	✓		
5-Dihydrotestosterone (5A-DHT)	✓	✓	✓		
Etiocholanolone (ET)	✓	✓	✓		
11-hydroxy-Etiocholanolone (OHET)	✓	✓	✓		
5-Androstenediol (5B-AD)	✓	✓	✓		
Dehydroepiandrosterone (DHEA)	✓	✓	✓	✓	
Dehydroepiandrosterone Sulfate (DHEAS)	✓	✓	✓	✓	
Progesterone	✓	✓	✓		
5-Pregnanediol (5A-PD)	✓	✓	✓		
5-Pregnanediol (5B-PD)	✓	✓	✓		
Allopregnanolone (ALLOP)	✓	✓	✓		
21-Hydroxyprogesterone (21-OHP)	✓	✓	✓		
17-Hydroxyprogesterone (17-OHP)	✓	✓	✓		
5-pregnanetriol (metabolite of 17-OHP)	✓	✓	✓		
Cortisol and Cortisone	x5			x5	
Corticosterone (B)	✓			✓	
Tetrahydrodehydrocorticosterone (THA)	✓			✓	
5-Tetrahydrocorticosterone (5B-THB)	✓			✓	
5a-Tetrahydrocorticosterone (5A-THB)	✓			✓	
11-Deoxycortisol (11-DOC)	✓			✓	
5a-Tetrahydrocortisol (5A-THF)	✓			✓	
5-Tetrahydrocortisol (5B-THF)	✓			✓	
Tetrahydrocortisone (THE)	✓			✓	
8-hydroxy-2'-deoxyguanosine (8-OHdG)	✓				
Key Enzyme Activity	✓				

Expanded Metabolic Pathways and Full Commentary included in the HuMap™

Estrogen Metabolites; urine

Order: SAMPLE REPORT
Client #: 12345
Doctor: Sample Doctor
 Doctor's Data, Inc.
 3755 Illinois Ave.
 St. Charles, IL 60174

Patient: Sample Patient
Age: 49
Sex: Female
Menopausal Status: Pre-menopausal

Sample Collection **Date/Time**
Dinnertime 05/01/2022 17:00
Bedtime 05/01/2022 21:00
Waking 05/02/2022 06:00
2 Hr. Post Waking 05/02/2022 08:00
Collection Period Multipoint daily
Date Received 05/03/2022
Date Reported 05/04/2022

Estrogens	Result	Unit	L	WRI	H	Reference Interval
Estrone (E1)	21.7	ng/mg Creat/Day				9.41 – 34.73
2-Hydroxyestrone	0.446	ng/mg Creat/Day				13.4 – 33.7
4-Hydroxyestrone	0.348	ng/mg Creat/Day				0 – 2.85
16α-Hydroxyestrone	44.6	ng/mg Creat/Day				4.2 – 15
2-Methoxyestrone	0.196	ng/mg Creat/Day				2 – 5.9
4-Methoxyestrone	0.777	ng/mg Creat/Day				0.051 – 0.28
Estradiol (E2)	12.6	ng/mg Creat/Day				2.58 – 13.9
2-Hydroxyestradiol	0.446	ng/mg Creat/Day				0.8 – 3.9
4-Hydroxyestradiol	0.098	ng/mg Creat/Day				0.144 – 0.93
2-Methoxyestradiol	0.134	ng/mg Creat/Day				0.174 – 0.502
4-Methoxyestradiol	0.250	ng/mg Creat/Day				0.049 – 0.106
Estrone (E1)	1.61	ng/mg Creat/Day				8.57 – 17.8
Ratios and Calculations						
2-OH-E1 %	0.983	%				50 – 81
4-OH-E1 %	0.786	%				2 – 5

Analyzed by DOCTOR'S DATA, INC. • 3755 Illinois Avenue, St. Charles, IL 60174-2420 USA • LAB DIR: Erlo Roth, MD • CLIA ID: 14D0646470

Free Cortisol and Cortisone	Result
Cortisol Waking	16.8
Cortisol Waking+2hrs	55.6
Cortisol Dinnertime	25.9
Cortisol Bedtime	16.7
Cortisol/day	31.3 (F)
Cortisone Waking	32.7
Cortisone Waking+2hrs	106
Cortisone Dinnertime	104
Cortisone Bedtime	63.7
Cortisone/day	49.1 (E)
Creatinine Waking	113 mg/dL
Creatinine Waking+2hrs	99.0 mg/dL
Creatinine Dinnertime	85.0 mg/dL
Creatinine Bedtime	102 mg/dL
Creatinine/day	112 mg/dL/Day

Analyzed by DOCTOR'S DATA, INC. • 3755 Illinois Avenue, St. Charles, IL 60174-2420 USA • LAB DIR: Erlo Roth, MD • CLIA ID: 14D0646470

Much like a human fingerprint, metabolism is inherently unique to everyone. Understanding how the major hormones of the body are metabolized, as well as the activity of the enzymes responsible for metabolism, is a rare glimpse into a patient's health.

Convenience, Clinical Utility and Researched Methodology of the HuMap™

Enhanced Specificity and Specificity

Liquid chromatography (LC) tandem mass spectrometry (MS) is an extremely sensitive and specific way to measure many substances. In the case of urinary steroid hormone testing, LC-MS is used to determine both free hormone and their metabolite amounts in urine. Put simply, both parts work together to figure out the exact amount of each hormone component in a patient sample. LC separates hormones in a liquid urine sample, which are then injected at various times for MS analysis. MS technology monitors the injected sample for specific hormones based on their molecular weights and expected injection times (commonly called retention times). Tandem MS verifies each hormone identity based on fragmentation and determines its amount. The combination of LC and tandem MS allows for extremely sensitive and specific hormone measurements, even in samples containing similar substances that would interfere with other methods.

Ability to Test Neurotransmitters from HuMap Urine Specimens

Many of the symptoms that would drive one to test urinary metabolite imbalance in a patient (fatigue, sleep difficulties, stress, mood concerns, cognitive concerns, vasomotor symptoms) can also be influenced by neurotransmitter imbalance. Additionally, the COMT (Catechol-O-Methyltransferase) enzyme plays an essential role in estrogen metabolism as well as catecholamine metabolism. Issues with COMT activity can result in both estrogen metabolite and neurotransmitter imbalance. Adding neurotransmitter testing to HuMap provides a deeper dive into the biochemistry that may be contributing to a patient's symptom picture. Providers can choose to add on the NeuroBasic or Comprehensive Neurotransmitter Profile to the same urine samples for additional convenience to the patient.

Advantages of Liquid Urine Collection

The main advantage of liquid urine collection is enhanced sensitivity, especially for low concentration metabolites. Dried urine must be reconstituted from the filter paper once the sample arrives. This reconstitution can lead to loss of polar steroid metabolites or creatinine for some patient samples. With liquid urine, samples can be shipped after being frozen for 4-6 hours, can be processed faster, and concentrated further to enhance the detection of low-level analytes. Steroids are also quite stable in properly preserved liquid urine.



When to Consider the HuMap™ for your Patients

Women

- Breast Health
- Endometriosis
- PCOS
- PMS/PMDD
- Menopausal Symptoms

Men

- Erectile Dysfunction
- Prostate Health
- Loss of Muscle Mass
- Breast Health

General

- Metabolic Syndrome
- Thyroid Pathologies
- Inflammation
- Oxidative Stress
- Fatigue/Insomnia
- Libido
- Mood or Cognitive Concerns
- Family History of Hormone-Driven Cancers
- HRT/BHRT Utilization
- Weight Gain

Urine and Saliva Testing for BHRT

Urinary hormone and metabolite testing is uniquely suited to provide insight into how hormones and their metabolites are moving through the body as well as risk assessment from the generation of certain metabolites. When using urinary metabolite testing, the level of detectable unconjugated sex hormone levels is not representative of circulating or bioavailable tissue hormone levels, because urine is not reflective of tissue uptake. When the goal of testing is to understand bioavailable levels of hormones, salivary testing may be a better option, especially if utilizing topical hormones.

Neurotransmitter Profiles



Urinary neurotransmitter profiles provide insight into a patient's overall ability to synthesize and metabolize neurotransmitters. Neurotransmitters are secreted from pre-synaptic neurons into the synapse between nerve cells to stimulate receptors on post-synaptic neurons. These chemical messengers regulate many physical and emotional processes including movement, stress response, cognition, emotions, energy, cravings, pain and more. Doctor's Data Neurotransmitter profiles help identify alterations in urinary neurotransmitter status which may be associated with a variety of conditions from metabolic to behavioral disorders.

Value in Clinical Practice

Analysis of neurotransmitters and their metabolites provides a non-invasive assessment of neurotransmitter status. A review of current scientific literature indicates that neurotransmitter testing may be useful in a variety of areas:

Identification of Imbalances—Research indicates that urinary neurotransmitter levels may correlate with conditions such as depression and PTSD. Further clinical associations have also been made with imbalances in neurotransmitters and:

- Other mood disorders
- HPA axis (adrenal) dysfunction: fatigue, insomnia
- Loss of mental focus: ADD, ADHD, cognitive concerns
- Cravings, addiction and dependency
- Hormonal imbalances: estrogen dominance (progesterone insufficiency), estrogen deficiency, androgen imbalance
- Loss of appetite control: obesity and insulin resistance

Functional Testing - Neurotransmitter metabolism is mediated by catechol-O-methyltransferase (COMT), monoamine

Comprehensive Neurotransmitter, urine

DOCTOR'S DATA

Order: 999999-9999

Client #: 12345
Doctor: Sample Doctor
3755 Illinois Ave.
St. Charles, IL 60174 U.S.A.

Patient: Sample Patient
Id: 999999
Age: 30
Sex: Female
Body Mass Index (BMI): 29

Sample Collection
Date Collected: 01/08/2022
Wake Up Time: 08:50
Collection Period: 2nd morning void
Date Received: 01/10/2022
Date Reported: 01/12/2022

Analyte	Result	Unit per Creatinine	L	WRI	H	Reference Interval
Phenethylamine (PEA)	32	nmol/g		▲		32 – 84
Tyrosine	94	µmol/g		▲	▲	32 – 80
Tyramine	5.3	µmol/g		▲	▲	2.0 – 4.0
Dopamine	40300	µg/g		▲	▲	125 – 250
3,4-Dihydroxyphenylacetic acid (DOPAC)	392600	µg/g		▲	▲	390 – 1500
3-Methoxytyramine (3-MT)	7320	nmol/g		▲	▲	90 – 210
Norepinephrine	22.7	µg/g		▲		22 – 50
Normetanephrine	144	µg/g		▲		85 – 300
Epinephrine	5.8	µg/g		▲		1.6 – 8.3
Metanephrine	72	µg/g		▲		45 – 119
Norepinephrine / Epinephrine ratio	3.9			▲		< 13
Tryptamine	0.63	µmol/g		▲		0.20 – 0.90
Serotonin	429	µg/g		▲	▲	60 – 125
5-Hydroxyindoleacetic acid (5-HIAA)	12420	µg/g		▲	▲	2000 – 8000
Glutamate	28	µmol/g		▲		12.0 – 45.0
Gamma-aminobutyrate (GABA)	4.8	µmol/g		▲		2.0 – 5.6
Glycine	875	µmol/g		▲		450 – 2200
Histamine	22	µg/g		▲		14 – 44
Taurine	2499	µmol/g		▲		320 – 1000
Creatinine	12.3	mg/dL	▲			30 – 225

NT Neurotransmitter Comments:

- Urinary neurotransmitter levels provide an overall assessment of the body's ability to make and break down neurotransmitters and are representative of whole body levels. Neurotransmitters are secreted all through the body, in neurons of both the central and peripheral nervous systems. The enzymes, cofactors and precursors in neurotransmitter metabolism in general are the same in the periphery and in the central nervous system. Therefore, alterations in urinary neurotransmitter levels assessed in urine provide important clinical information, and may be associated with many symptoms including cognitive and mood concerns, diminished drive, fatigue and sleep difficulties, cravings, addictions and pain.
- Tyrosine is the non-essential amino acid precursor for dopamine, norepinephrine and epinephrine. Increased tyrosine may exacerbate migraine headaches and hyperthyroid conditions. Elevated tyrosine levels may occur due to supplementation (phenylalanine or tyrosine), heritable enzyme defects, or liver disease. Tyrosine hydroxylase converts tyrosine into the dopamine precursor L-DOPA; BH4, Vitamin D and iron are cofactors for that enzymatic activity.

Notes:
Results are creatinine corrected to account for urine dilution variations. Creatinine is not meant to be used as an indicator of renal function.
RI= Reference Interval, L (blue)= Low (below RI), WRI (green)= Within RI (optimal), WRI (yellow)= Within RI (not optimal), H (red)= High (above RI)
Methodology: LCMS-GCQ-Creatinine by Jaffe Reaction Analyzed by DOCTOR'S DATA, INC. • 3755 Illinois Avenue, St. Charles, IL 60174-2420 USA • LAB DIR: Eric Roth, MD • CLIA ID: 14D0646470

oxidase (MAO), and other enzymes. Test results may provide functional information about these important enzymes.

Response to Therapy—Neurotransmitters (serotonin, for example) may be altered by the addition of precursor amino acids such as 5-hydroxytryptophan (5-HTP). These changes may be apparent in the urine.

Toxicology Risk Assessment—Changes in urinary neurotransmitters like serotonin, dopamine, and glutamate may be clinically relevant outcomes for neurobehavioral toxicology resulting from chemical or environmental exposures.

Which Profile to Consider

The NeuroBasic Profile offers a broad-brush assessment that covers most clinical needs and which many providers find highly clinically relevant

and useful. For providers interested in a more comprehensive look at neurotransmitter secretion and the metabolism of these markers, the

Comprehensive Neurotransmitter Profile includes catecholamine and serotonin metabolites, trace amines, and additional markers that expand the clinical view.

	Comprehensive Neurotransmitter Profile 	NeuroBasic Profile 	NeuroAdrenal Profile 
Serotonin	✓	✓	✓
GABA	✓	✓	✓
Dopamine	✓	✓	✓
Norepinephrine	✓	✓	✓
Epinephrine	✓	✓	✓
Glutamate	✓	✓	✓
Histamine	✓	✓	✓
PEA (Phenethylamine)	✓	✓	✓
5-HIAA (5-Hydroxyindolacetic acid)	✓		
Glycine	✓	✓	✓
Taurine	✓		
Tryptamine	✓		
Tyroamine	✓		
Tyrosine	✓		
DOPAC (3,4 Dihydroxyphenylacetic acid)	✓		
3-MT (3-Methoxytyramine)	✓		
Metanephrine Fractionation (Metanephrine, Normetanephrine)	✓		
Creatinine	✓	✓	✓
Cortisol			x4
DHEA			✓

Thyroid Profile



The analysis of thyroid hormones and auto-antibodies together may improve the accuracy diagnosis and clinical success. The American Thyroid Association estimates that approximately 20 million Americans have thyroid disease, and approximately 60% of those with thyroid disease are unaware of their condition. Many patients with thyroid disorders may remain undiagnosed in many patients with asymptomatic or non-specific clinical presentations. The recognition of auto-immunity as a leading cause of thyroid dysfunction has led to the evaluation of auto-antibodies in thyroid testing.

Measuring only thyroid stimulating hormone (TSH) may be misleading in a variety of circumstances, including the recent treatment of thyrotoxicosis, pituitary disease, non-thyroid illness, thyroid hormone resistance or rare, TSH-secreting tumors. Under these circumstances, and in many other cases, the evaluation of thyroid hormones and thyroid antibodies may clarify the diagnosis of thyroid conditions and improve clinical success.

This test is useful for

- Hypothyroid conditions
- Hyperthyroid conditions
- Autoimmune conditions
- Arrhythmia
- Infertility
- Cholesterol disorders
- Fatigue
- Pituitary disorders




LAB #: Sample Report
 PATIENT: Sample Patient
 ID:
 SEX: Female
 DOB: 01/01/1951 AGE: 67

CLIENT #: 12345
 DOCTOR: Sample Doctor
 Doctor's Data, Inc.
 3755 Illinois Ave.
 St. Charles, IL 60174 U.S.A.

Thyroid Profile; serum

	RESULT / UNIT	REFERENCE INTERVAL	PERCENTILE				
			2.5 th	16 th	50 th	84 th	97.5 th
Free T3	3.9 pg/mL	2.2 - 4.0					
Free T4	0.8 ng/dL	0.6 - 1.3					
Thyroid Stimulating Hormone (TSH)	0.01 µIU/mL	0.30 - 4.5					
Thyroglobulin Antibody (Anti-TG)	1.1 IU/mL	< 4.0				95 th	99 th
Thyroid Peroxidase Antibody (Anti-TPO)	250 IU/mL	< 9.0					



This diagnostic algorithm is intended for baseline assessments only and may not be accurate for patients on thyroid medications.

SPECIMEN DATA	
Comments:	
Date Collected: 01/18/2022	Time Collected:
Date Received: 01/20/2022	Fasting:
Date Completed: 01/21/2022	
Methodology: Chemiluminescent Immunoassay	

©DOCTOR'S DATA, INC. • ADDRESS: 3755 Illinois Avenue, St. Charles, IL 60174-2420 • CLIA ID NO: 14D0646470 • MEDICARE PROVIDER NO: 148453
 0001912

OUR MISSION:

To research, develop and offer innovative specialty tests that help doctors identify health risks and improve outcomes for patients with chronic conditions.

To educate and support healthcare professionals.

To improve lives through science.



SCIENCE + INSIGHT



3755 Illinois Avenue
St. Charles, IL 60174-2420

800.323.2784 (US AND CANADA)
0871.218.0052 (UK)
+1.630.377.8139 (GLOBAL)

doctorsdata.com

©2022 Doctor's Data, Inc. All rights reserved.

About Doctor's Data

Doctor's Data, Inc. has provided innovative specialty testing to healthcare practitioners around the world from our advanced, CLIA-licensed clinical laboratory since 1972.

A specialist and pioneer in essential and toxic elemental testing, the laboratory provides a wide array of functional testing to aid in decision making and better patient outcomes. Choose Doctor's Data to help you assess and treat heavy metal burden, nutritional deficiencies, gastrointestinal function, hormone status, cardiovascular risk, liver and metabolic abnormalities, and more.