

Eubacterium limosum* Competent Cells Prep*1) Procedure Objective**

Prepare electrocompetent cells from *E. limosum* cultures to perform genetic transformation by electroporation.

2) Health and Safety

Lab coat, closed toed shoes, gloves, safety glasses should be worn for this procedure.

3) Materials List

- 50mL tubes
- Eppendorf tubes
- Large refrigerated centrifuge
- Sucrose
- Glycerol
- DSMZ 135 medium (see media recipes)
- *E. limosum*
- Ice/ice bucket
- Autoclave and/or 0.22um filter cups
- Anaerobic chamber

4) Pre-lab Steps

- Ensure that all plastics, reagents, and media have been in the anaerobic chamber for at least 48 hours to become totally anoxic. The ice is an exception to this rule.
- Prepare 270mM sucrose solution and pass into chamber
- Prepare 270mM sucrose solution + 20% v/v glycerol and pass into chamber

5) Procedure Steps

- 1) Inside the anaerobic chamber, transfer 50mL sterile DSMZ 135 medium without sulfide to a sterile 50mL centrifuge tube (with o-ring cap) and inoculate with a stab from an *E. limosum* cryostock.
 - a. Grow up culture overnight at 37C in the anaerobic chamber incubator.
 - b. Should grow up for 18-19 hours so that the culture is just starting to become opaque. Ideal procedure is to inoculate at 1 or 2pm and the perform competent cells prep at 9am the following day.
- 2) Pass a large ice bucket filled with ice into the chamber.
- 3) Place the culture in the tube on the ice and fully submerge.
 - a. Leave in the ice for ~20 minutes or until it becomes ice cold.

- 4) Pass the tube out of the chamber through the glove hole
- 5) Spin down culture at 15,000xg for 10 minutes at 4C in a pre-chilled centrifuge and pass back into the chamber.
- 6) Decant supernatant and resuspend fully in 10mL ice cold 270mM sucrose solution
- 7) Repeat centrifugation step and resuspend fully in 10mL ice cold 270mM sucrose solution+20% v/v glycerol
- 8) Centrifuge again and decant supernatant.
- 9) Resuspend cells in the leftover liquid – this should be about ~1mL total volume
- 10) Aliquot 50uLs into pre-chilled Eppendorf tubes and store at -80C for future use.

6) References

[1] "Genome Engineering of *Eubacterium limosum* Using Expanded Genetic Tools and the CRISPR-Cas9 System" Shin et.al. 2019

[a] "A genetic System for *Clostridium ljungdahlii*: a chassis for Autotrophic Production of Biocommodities and a Model Homoacetogen" Leang et.al. 2013

7) Revision Table

Revision No.	Date of Change	Change Description	Person Responsible
0	1/21/21	Initial Development	PAS
1	2/1/22	Modified to grow culture in 50mL tube from start and increased speed to 15xkg	PAS
