




# LEAK TESTING



PRESENTED BY:  
RGN.Ibtisam almasruri  
RGN.Noora alzadjali

# Contents:

---

- 01** Definition and Purpose of leak test.
- 02** Know when to perform a leak test, how to interpret pass/fail results, and how to escalate appropriately.
- 03** Scope cleaning steps.



- **Definition & Purpose of leak test**

---

**Leak testing is performed under pressure to detect any internal or external damage to the endoscope immersion in liquids.**






- **When to perform a leak test**

---

- Before every manual cleaning of endoscope.
- After every procedure to check for hidden damage.
- If the scope was dropped and mishandled.
- If you see any functional issues "suction, air/water not working".
- Any time moisture is suspected inside the channels.



# • Equipment needed for leak test




**Leakage tester (MB-155)**  
(Sold separately. Refer to its instruction manual.)



**Maintenance unit (MU-1)**  
(Sold separately. Refer to its instruction manual.)


- Clean, large basins  
(size: 40 (W) × 40 (H) × 25 (D) cm or more)
- Water (for reprocessing)  
(Refer to Section 3.5, “Water”)

- Clean lint-free cloths
- Detergent solution  
(Refer to Section 3.3, “Detergent solution for manual cleaning”)





# Leakage tester




- **How to interpret pass/fail results.**

---

## PASS

- 1.No bubbles appear during wet leak testing.
- 2.Pressure remains stable on the leak tester.
- 3.No visible external defects.



## FAIL

- 1.Continueus bubbles appear from any part of the scope.
- 2.Pressure drops on the leak tester.
- 3.Visible cracks, dents, torn bending rubber or fluid invasion signs.



- **Correct escalation**

---

- 1.Immediately stop reprocessing.**
- 2.Label the scope “leak test fail”.**
- 3.Remove from clinical use.**
- 4.Document the failure in tracking system.**
- 5.Send to biomedical engineering or service center for evaluation /repair.**






# • **Mechanical “wet” leak testing steps**

---

- 1. Remove suction, air/water, and biopsy valves.**
- 2. Discard discard disposable parts and keep reusable ones the same scope.**
- 3. Attach the leak tester and pressure the scope before submerging it in clean water.**
- 4. Don't add detergent during leak testing.**
- 5. Submerge the scope completely , flex the distal tip, and observe for bubbles from any part of the scope.**



- **Scope Cleaning steps**

---

**1. Cleaning & Brushing : to remove all debris and bioburden from the endoscope.**

- **Wash the exterior of the scope while submerged in detergent to prevent splashing.**
- **Use a small, soft brush to clean removable parts and openings(suction, air /water, biopsy ports).**
- **Brush all channels and external surfaces thoroughly using properly sized brushes.**



- **Scope Cleaning steps**

---


## **1.Cleaning &Brushing :**

- **Special scopes( duodenoscope, EUS) need extra steps due to complex channels.**
- **Rinse brushes frequently in detergent continue brushing until no debris is visible.**
- **Clean and high level disinfect reusable brushes between cases replace damage ones.**





# Brushing scope




- **Scope Cleaning steps**


---

## **2. Flushing and Soaking.**

- **Attach manufacturer's cleaning adapters for suction, air/water and biopsy channels.**
- **Use special adapters for elevator or auxiliary channels if needed.**
- **Flush all channels with detergent solution to remove debris.**
- **Soak the endoscope and channels for the time specified on the detergent lable.**



# Flushing & Soaking





- **Scope Cleaning steps**

---

## **NOTE**

- 1. All steps must be performed sequentially and immediately after the procedure.**
- 2. Always follow the manufacturer's recommendations for delayed re-cleaning and reprocessing.**



# REFERENCE


---

- Rutala, W. A., & Weber, D. J. (2018). Guideline for disinfection and sterilization in healthcare facilities. Centers for Disease Control and Prevention (CDC).
- World Gastroenterology Organisation. (2024). Global Guidelines: Endoscope Disinfection. World Gastroenterology Organisation.



Thank You for your attention  
Visit our online platform and learn more

---



[www.muscatendoscopyacademy.com](http://www.muscatendoscopyacademy.com)

