ARTERIAL EMBOLCTOMY CATHETER

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NOTE: Not all symbols are used in the labeling of this product.

CE 0459

EC REP

Manufactured By
LUCAS MEDICAL, INC.
1554 North Case Street
Orange, CA 92867 USA

P/N: 10024 Rev. J

CMC MEDICAL DEVICES
C/HORACIO LENGU 18
CP 29006
MALAGA, SPAIN
(TEL) +34 951 214 054

U.S. PATENT NOS. 5,762,996 & 5,795,332

English.....	1
French.....	1
Deutsch.....	2
Espanol.....	3
Italiano.....	3
Nederlands.....	4
Dansk.....	5
Svenska.....	5
Turkish.....	6

ARTERIAL EMBOLCTOMY CATHETER

ENGLISH

INDICATIONS FOR USE: The device is intended for intra-operative removal of arterial emboli from peripheral vessels.

DESCRIPTION: The catheter consists of single lumen, stainless wire reinforced silicone tubing with a silicone balloon attached at the distal end. Syringe attachment is to a poly-carbonate female luer lock connector. A recommended size syringe is provided.

PRINCIPLE OF OPERATION: The affected blood vessel is accessed surgically proximal to the region containing the emboli. An appropriately sized Embolectomy catheter for the blood vessel is selected. The catheter is inserted into the artery and advanced distally through the clot/embolus with the balloon deflated. The catheter balloon is then inflated using a sterile physiological fluid (i.e. sterile saline), and pulled back toward the arterial incision. The balloon pushes the loosened embolus out to the incision sight where the clot is removed.

PRECAUTIONS: As with all catheter procedures, complications may occur. The medical judgment of the physician must be exercised at all times.

WARNINGS:

- Use of highly viscous or liquid suspension media is not recommended due to the possibility of catheter lumen occlusion.
- Sterile only if package is not opened, damaged or broken. Do not re-sterilize. Catheters are intended for single patient use only.
- Exceeding the maximum recommended volume greatly increases the possibility of balloon rupture and/or arterial damage.
- This catheter is not designed or recommended for use in the venous system, endarterectomy, dilation or occlusion procedures.
- Federal Law (United States of America) restricts this device to sale by or on the order of a physician.
- The possibility of balloon rupture must be considered, particularly if air is used for the balloon inflation media. Air should not be used for balloon inflation in instances where balloon rupture could produce a dangerous air embolus.
- The physician must use proper judgment in limiting balloon inflation diameter and pull forces to avoid vessel and intimal lining injury or dissection.

SPECIFICATIONS: 3F (1.00mm), 5F (1.67mm), 6F (2.0mm) & 7F (2.3mm) supplied in 40 & 80cm lengths; 2F (0.67mm) supplied in 40 & 60 cm lengths; 4F (1.33mm) in 40, 80 & 100 cm lengths.

Catheter French Size	Catheter Color Ident	Recom- mended Syringe Size	Maximum Liquid Inflation Volume (MLV)	Approx. Balloon Diameter @ MLV
2F (0.67mm)	Purple	1.0 cc	0.05 cc	4 mm
3F (1.00mm)	Green	3.0 cc	0.10 cc	5 mm
4F (1.33mm)	Red	3.0 cc	0.50 cc	9 mm
5F (1.67mm)	White	3.0 cc	0.75 cc	11 mm
6F (2.0mm)	Blue	3.0 cc	1.25 cc	13 mm
7F (2.3mm)	Yellow	3.0 cc	2.00 cc	14 mm

PRIME MEDIA: The catheter should be primed with sterile, blood-compatible liquid prime (e.g. isotonic saline) before use. Gas or air is not recommended.

TESTING: A visual inspection for leaks should be performed prior to use with the catheter fully primed and balloon inflated. **NOTE: Do not exceed the Maximum Liquid Inflation Volume (MLV) for balloon inflation.**

LIQUID PRIME PROCEDURE: Fill the recommended-size syringe with the Maximum Liquid Inflation Volume (MLV) of a sterile, blood-compatible liquid (e.g. isotonic saline). Expel air bubbles and attach the syringe to the catheter Luer connector. Draw back on the syringe plunger to aspirate air from the catheter body and balloon. Inflate and deflate the balloon with a small amount of fluid. Gently tap the catheter to facilitate air removal. Repeat this procedure until all air is purged from the catheter. Before fully inflating the balloon, follow the recommended balloon exercise procedure.

BALLOON EXERCISE: For optimum balloon performance, the balloon should be exercised in the following sequence after the catheter has been primed. Adjust the liquid volume in the syringe to the Maximum Liquid Inflation Volume. Inflate and deflate the balloon in three successive steps using an increasingly larger volume of inflation media with each inflation, (i.e. 50, 75% and 100% of the Maximum Liquid Inflation Volume). Aspirate air bubbles, if present, during the exercise procedure.

PATIENT POPULATION: Adult and Pediatric

NOTE: Drawing and artwork under review. Awaiting translation for the changes of the indications for use and patient population.

REFER TO THE SYMBOL LEGEND AT THE END OF THIS DOCUMENT

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1. Billig, D.M., et al. "Arterial Embolism," Arch. Surg., 95:1-6, 1967.
2. Cranley, J.J. et al. "A Complication with the Use of the Fogarty Balloon Catheter for Arterial Embolectomy," J. Cardiovasc. Surg., 10:407-409, 1969.
3. Cranley, J.J. et al. "Catheter Technique for Arterial Embolectomy: A Seven Year Experience," J. Cardiovasc. Surg., 11:44-51, 1970.
4. Dainko, E.A. "Complications of the Use of the Fogarty Balloon Catheter," Arch. Surg., 105:79-82, 1972.
5. Fogarty, T.J., et al. "Experience with Balloon Catheter Technic for Arterial Embolectomy," Am. J. Surg., 122:231-237, 1971.
6. Foster, J.H., et al. "Arterial Injuries Secondary to the Use of the Fogarty Catheter," Ann. Surgery, 171:971-978, 1970.
7. Krause, R.J. et al. "Further Experience with a New Embolectomy Catheter," Surgery, 59:81-87, 1966.
8. Oehlert, W.H. "A Complication of the Fogarty Arterial Embolectomy Catheter," Am. Heart J., 84:484-486, 1972.
9. Prolo, D.J., et al. "Intraarterial Occlusion of a Carotid-Cavernous Sinus Fistula with a Balloon Catheter," J. Neurosurg., 35:237-242, 1971.
10. Provan, J.L., A.O. Ransford. "The Role of the Fogarty Embolectomy Catheter in the Treatment of Arterial Embolism of the Limbs," Br. J. Surg., 57:59-62, 1970.
11. McCaughan, J.J., J.M.Young. "Intra-Arterial Occlusion in Vascular Surgery," Ann. of Surg., 695-703, May 1970.

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