DIRECTORY

1. ENGLISH Carotid shunt
2. ITALIANO Shunt carotideo
3. DEUTSCH Carotis shunt
4. ESPAÑOL Desviación Carótida
5. NEDERLANDS Hoofdslagader shunt
6. PORTUGUÊS Shunt carotídeo
7. FRANÇAIS Shunt carotidien
8. TÜRKÇE Karotis şant


SILICONE CAROTID ENDARTERECTOMY SHUNT CATHETER (DOUBLE BALLOON)

SILICONE CAROTID ENDARTERECTOMY SHUNT CATHETER INSTRUCTION FOR USE (DOUBLE BALLOON)

INDICATIONS FOR USE: The LMI silicone carotid endarterectomy shunt catheter is intended for use during a carotid endarterectomy procedure to provide temporary carotid artery (blood) flow for cerebral circulation.

A. TYPICAL SHUNT POSITION

DESCRIPTION: The LMI carotid shunt is a single use, sterile device that is used on by the order of a physician. The shunt consists of a multi-lumen straight or tempered (on various models) silicone tube with or without stainless steel wire reinforcement, that is formed to create a rounded leading edge. Silicone balloons are formed on the distal ends of the shunt that can be inflated and deflated independently. The proximal ends of the inflation tubes consist of a female luer connector that is attached to the silicone shaft. A stopcock is also affixed to the luer connector. The central communicating tube is also fitted with a female luer connector and a stopcock. Syringe attachment is to the female luer lock connector. Markings are provided to identify the common carotid end of the shunt.

INSTRUCTIONS FOR USE:

1. The affected blood vessel is accessed surgically and an appropriately sized carotid shunt for the carotid artery is selected.
2. The carotid shunt is first inserted into the internal carotid side of the carotid artery with the balloon in the deflated position.
3. The balloon is then inflated, using saline, to occlude arterial flow and to position the shunt to prevent movement with the artery. Volumes of less than $\frac{3}{4}$ cc of saline are usually sufficient to fill the balloon properly in the internal carotid artery.
4. With the balloon inflated and positioned in the internal carotid artery blood is allowed to flow through the central lumen of the shunt until all of the air has been expelled.
5. The remaining end of the shunt is then inserted into the common carotid artery and the balloon is then inflated to a size that will occlude the vessel and prevent any movement of the shunt. Usually a volume of saline of approximately $\frac{3}{4}$ to $1\frac{1}{4}$ cc is adequate to achieve this desired result. The decision is left to the discretion of the surgeon.

A safety balloon has been developed, and is utilized on several models to reduce the possibility of damaging the internal carotid artery resulting from excessive balloon pressure due to over-inflation of the balloon. As artery characteristics vary from patient to patient, care should be taken not to force a shunt into a smaller vessel, or stretching the vessel by over-inflating the balloon. It has been reported that human internal carotid arteries can rupture at pressures of 20 psi. The safety balloon is intended to prevent the pressure from reaching this level by inflating at a level that is less than 20 psi. In the event that the internal balloon is inadvertently over-inflated, the safety balloon will inflate thereby reducing the pressure in the internal balloon. In the event that this occurs, deflate both the safety balloon and the internal balloon and repeat the inflation of the internal balloon without inflating the safety balloon. All that is required of the internal balloon is that it will stay in its position and prevent back bleeding during the procedure. After satisfactory inflation of the internal balloon, slide the sheath that is provided over the safety balloon to assure that any back pressure from the internal balloon that occurs during the procedure will not inflate the safety balloon and cause the internal balloon to deflate and become dislodged.

CAUTION: Care must be taken at all times to **prevent over-inflation** of the balloons, as it may result in obstruction of blood flow and injury to the vessel.

PRECAUTIONS: As with all catheter procedures, complications may occur. The medical judgment of the physician must be exercised at all times.

B. POTENTIAL COMPLICATIONS

1. INTIMAL TRAUMATISM
2. RUPTURE OF THE ARTERIAL WALL
3. EMBOLISATION OF THE ATHEROMATOUS PLAQUE
4. THROMBOSIS
5. PROLIFERATION OF THE SUB-ENDOTHELIAL LAYER
6. OCCLUSION OF THE LUMEN OF A VESSEL

WARNINGS:

1. Sterile only if package is not opened, damaged or broken. **Do not re-sterilize.** The LMI carotid shunt catheter is intended for single patient use only. **DO NOT REUSE**
2. Exceeding the recommended inflation volume greatly increase the possibility of balloon rupture and/or arterial damage.
3. The recommended inflation media is saline. Air or any gas should not be used for balloon inflation as it could produce dangerous air embolus if balloon rupture occurs.
4. The LMI carotid shunt is a temporary device and cannot and should not be implanted.
5. Federal Law (United States of America) restricts this device to sale by or on the order of a physician.
6. Aspirate the proximal and distal balloons of all air prior insertion.

7. Ensure that connectors between syringe and luer connectors are secure to prevent introduction of air.
8. Do not grasp the balloons with traumatic instruments during insertion. Care must also be taken to prevent the use of sharp instruments in the carotid shunt shaft.
9. In the internal carotid artery, do not attempt to insert a shunt that is larger than is needed as this can result in vessel damage. The internal carotid artery varies in length from approximately 7 to 10 cm; therefore care must be taken to prevent insertion of the carotid shunt more than 4 to 5 cm.
10. Deflate both of the balloons prior to withdrawing the shunt.

INSPECTION AND TESTING: Prior to use of the LMI carotid shunt, the package should be inspected for any evidence of damage. The shunt should also be inspected for any evidence of damage. Both balloons should be inflated with air and immersed in sterile water. Do not use the shunt if there is any evidence of damage or air bubbles from around the balloons.

WARRANTY: LMI, Inc. warrants the reasonable care has been used in the manufacture of this device. This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether expressed or implied by operation of law or otherwise, including but not limited to any implied warranties of MERCHANTABILITY OR FITNESS, since handling, storage, cleaning and sterilization of this device as well as factors relating to the patient, his diagnosis, treatment, surgical procedures and other matters beyond the control of LMI, Inc. directly affect the device and the CONSEQUENTIAL LOSS, DAMAGE OR EXPENSE directly or indirectly arising from the use of this device other than replacement of all part of it. LMI, Inc. neither assumes, nor authorizes any other persons to assume for it, any other or ADDITIONAL LIABILITY OR RESPONSIBILITY in connection with this device

Sterilization Method: 100% Ethylene Oxide Gas

	Internal Balloon		Common Balloon	
Catheter French Size Inflation Volume (MLV)	Maximum Liquid Inflation Volume (MLV)	Approximate Balloon Diameter @ MLV	Maximum Liquid	Approximate Liquid
8 Fr	0,75 cc	11,0 mm	1,25 cc	12,0 mm
9 Fr	0,75 cc	12,0 mm	1,25 cc	13,0 mm
10 Fr	0,75 cc	12,0 mm	1,25 cc	13,5 mm
11 Fr	0,75 cc	12,5 mm	1,25 cc	13,5 mm
12 Fr	0,75 cc	12,5 mm	1,25 cc	14,0 mm

Models: CES, CESB, CBS

PRIME MEDIA: The catheter should be printed with sterile, blood-compatible liquid prime (e.g. isotonic saline) before use. Gas or air is not recommended

LIQUID PRIME PROCEDURE: Fill the recommended-size syringe with the Maximum Liquid Inflation Volume (MLV) of a sterile, blood-compatible liquid (e.g. isotonic saline). Expel air bubbles and attach the syringe to the catheter luer connector. Draw back on the syringe plunger to aspirate air from the catheter body and balloon. Inflate and deflate the balloon with a small amount of fluid. Gently tap the catheter to facilitate air removal. Repeat this process until all air is purged from the catheter.

PROCEDURAL RISKS: Complications can occur with any surgical procedure. Some possible complications that may occur with carotid endarterectomy, but are not limited to, include: (1). Stroke or transient ischemic attack. (TIA), (2). Heart attack, (3). Wound hematoma, (4). Cranial nerve deficits, (5). Intracerebral hemorrhage, (6) Seizures, (7). Bleeding, (8) Infection, (9). Hypertension and (10). Arrhythmias.

REFERENCES:

1. Billing, D.M., et al. "Arterial Embolism" Arch. Surg., 95:1-6, 1967.
2. Cranley, J.J. et al. "A complication With The Use of the Fogarty Balloon Catheter For Arterial Embolectomy" J. Cardiovasc. Surg., 10:407-409, 1969.
3. Dainko, E.A., "Complications Of The Use Of The Fogarty Balloon Catheter", Arch. Surg. 105:79-82, 1972.
4. Foster, J.H. et al. "Arterial Injuries Secondary To The Use Of The Fogarty Catheter", Ann. Surgery, 171:971-978, 1970.
5. Ochler, W.H. "A Complication of the Fogarty Arterial Embolectomy Catheter", Am Heart J. 84:484-486, 1972.
6. McCaughan, J.J., J.M. Young, "Intra-Arterial Occlusion In Vascular Surgery", Ann. Of Surg., 695-703, May 1970.
7. Pruitt, J. Crayton, "1099 Consecutive Carotid Endarterectomies Using Local Anesthesia, EEG and Selective Shunting with Pruitt-Inahara Carotid Shunt, "Contemporary Surgery, Vol. 23, Sept. 1993".

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BIBLIOGRAFIA:

- Billing, D.M., et al. "Arterial Embolism" Arch. Surg., 95:1-6, 1967.
- Cranley, J.J. et al. "A complication With The Use of the Fogarty Balloon Catheter For Arterial Embolectomy" J. Cardiovasc. Surg., 10:407-409, 1969.
- Dainko, E.A., "Complications Of The Use Of The Fogarty Balloon Catheter", Arch. Surg. 105:79-82, 1972.
- Foster, J.H. et al. "Arterial Injuries Secondary To The Use Of The Fogarty Catheter", Ann. Surgery, 171:971-978, 1970.
- Ochert, W.H. "A Complication of the Fogarty Arterial Embolectomy Catheter", Am Heart J. 84:484-486, 1972.
- McCaughan, J.J., J.M. Young, "Intra-Arterial Occlusion In Vascular Surgery", Ann. Of Surg., 695-703, May 1970.
- Pruitt, J. Crayton, "1099 Consecutive Carotid Endarterectomies Using Local Anesthesia, EEG and Selective Shunting with Pruitt-Inahara Carotid Shunt, "Contemporary Surgery, Vol. 23, Sept. 1993".

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REFERENZEN:

- Billing, D.M., et al. "Arterial Embolism" Arch. Surg., 95:1-6, 1967.
- Cranley, J.J. et al. "A complication With The Use of the Fogarty Balloon Catheter For Arterial Embolectomy" J. Cardiovasc. Surg., 10:407-409, 1969.
- Dainko, E.A., "Complications Of The Use Of The Fogarty Balloon Catheter", Arch. Surg. 105:79-82, 1972.
- Foster, J.H. et al. "Arterial Injuries Secondary To The Use Of The Fogarty Catheter", Ann. Surgery, 171:971-978, 1970.
- Ochlerl, W.H. "A Complication of the Fogarty Arterial Embolectomy Catheter", Am Heart J. 84:484-486, 1972.
- McCaughan, J.J., J.M. Young, "Intra-Arterial Occlusion In Vascular Surgery", Ann. Of Surg., 695-703, May 1970.
- Pruitt, J. Crayton, "1099 Consecutive Carotid Endarterectomies Using Local Anesthesia, EEG and Selective Shunting with Pruitt-Inahara Carotid Shunt, "Contemporary Surgery, Vol. 23, Sept. 1993".

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REFERENCIAS:

1. Billing, D.M., et al. "Arterial Embolism"Arch. Surg., 95:1-6, 1967.
2. Cranley, J.J. eat al. "A complication With The Use of the Fogarty Balloon Catheter For Arterial Embolectomy"J. Cardiovasc. Surg., 10:407-409, 1969.
3. Dainko, E.A., "Complications Of The Use Of The Fogarty Balloon Catheter", Arch. Surg. 105:79-82, 1972.
4. Foster, J.H. et al. "Arterial Injuries Secondary To The Use Of The Fogarty Catheter", Ann. Surgery, 171:971-978, 1970.
5. Ochert, W.H. "A Complication of the Fogarty Arterial Embolectomy Catheter", Am Heart J. 84:484-486, 1972.
6. McCaughan, J.J., J.M. Young, "Intra-Arterial Occlusion In Vascular Surgery", Ann. Of Surg., 695-703, May 1970.
7. Pruitt, J. Crayton, "1099 Consecutive Carotid Endarterectomies Using Local Anesthesia, EEG and Selective Shunting with Pruitt-Inahara Carotid Shunt, "Contemporary Surgery, Vol. 23, Sept. 1993".

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REFERENTIES:

- Billing, D.M., et al. "Arterial Embolism" Arch. Surg., 95:1-6, 1967.
- Cranley, J.J. et al. "A complication With The Use of the Fogarty Balloon Catheter For Arterial Embolectomy" J. Cardiovasc. Surg., 10:407-409, 1969.
- Dainko, E.A., "Complications Of The Use Of The Fogarty Balloon Catheter", Arch. Surg. 105:79-82, 1972.
- Foster, J.H. et al. "Arterial Injuries Secondary To The Use Of The Fogarty Catheter", Ann. Surgery, 171:971-978, 1970.
- Ochert, W.H. "A Complication of the Fogarty Arterial Embolectomy Catheter", Am Heart J. 84:484-486, 1972.
- McCaughan, J.J., J.M. Young. "Intra-Arterial Occlusion In Vascular Surgery", Ann. Of Surg., 695-703, May 1970.
- Pruitt, J. Crayton, "1099 Consecutive Carotid Endarterectomies Using Local Anesthesia, EEG and Selective Shunting with Pruitt-Inahara Carotid Shunt, "Contemporary Surgery, Vol. 23, Sept. 1993".

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REFERÊNCIAS:

1. Billing, D.M., et al. "Arterial Embolism" Arch. Surg., 95:1-6, 1967.
2. Cranley, J.J. et al. "A complication With The Use of the Fogarty Balloon Catheter For Arterial Embolectomy" J Cardiovasc. Surg., 10:407-409, 1969.
3. Dainko, E.A., "Complications Of The Use Of The Fogarty Balloon Catheter", Arch. Surg. 105:79-82, 1972.
4. Foster, J.H. et al. "Arterial Injuries Secondary To The Use Of The Fogarty Catheter", Ann. Surgery, 171:971-978, 1970.
5. Ochler, W.H. "A Complication of the Fogarty Arterial Embolectomy Catheter", Am Heart J. 84:484-486, 1972.
6. McCaughan, J.J., J.M. Young, "Intra-Arterial Occlusion In Vascular Surgery", Ann. Of Surg., 695-703, May 1970.
7. Pruitt, J. Crayton, "1099 Consecutive Carotid Endarterectomies Using Local Anesthesia, EEG and Selective Shunting with Pruitt-Inahara Carotid Shunt, "Contemporary Surgery, Vol. 23, Sept. 1993".

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RÉFÉRENCES :

- Billing, D.M., et al. "Arterial Embolism" Arch. Surg., 95:1-6, 1967.
- Cranley, J.J. et al. "A complication With The Use of the Fogarty Balloon Catheter For Arterial Embolectomy" J. Cardiovasc. Surg., 10:407-409, 1969.
- Dainko, E.A., "Complications Of The Use Of The Fogarty Balloon Catheter", Arch. Surg. 105:79-82, 1972.
- Foster, J.H. et al. "Arterial Injuries Secondary To The Use Of The Fogarty Catheter", Ann. Surgery, 171:971-978, 1970.
- Ochler, W.H. "A Complication of the Fogarty Arterial Embolectomy Catheter", Am Heart J. 84:484-486, 1972.
- McCaughan, J.J., J.M. Young, "Intra-Arterial Occlusion In Vascular Surgery", Ann. Of Surg., 695-703, May 1970.
- Pruitt, J. Crayton, "1099 Consecutive Carotid Endarterectomies Using Local Anesthesia, EEG and Selective Shunting with Pruitt-Inahara Carotid Shunt, "Contemporary Surgery, Vol. 23, Sept. 1993".

SİLİKON KAROTİS ENDARTEREKTOMİ ŞANT KATETERİ KULLANMA TALİMATLARI (ÇİFT BALONLU)

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REFERANSLAR:

1. Billing, D.M., et al. "Arterial Embolism"Arch. Surg., 95:1-6, 1967.
2. Cranley, J.J. et al. "A complication With The Use of the Fogarty Balloon Catheter For Arterial Embolectomy"J. Cardiovasc. Surg., 10:407-409, 1969.
3. Dainko, E.A., "Complications Of The Use Of The Fogarty Balloon Catheter", Arch. Surg. 105:79-82, 1972.
4. Foster, J.H. et al. "Arterial Injuries Secondary To The Use Of The Fogarty Catheter", Ann. Surgery, 171:971-978, 1970.
5. Ochler, W.H. "A Complication of the Fogarty Arterial Embolectomy Catheter", Am Heart J. 84:484-486, 1972.
6. McCaughan, J.J., J.M. Young, "Intra-Arterial Occlusion In Vascular Surgery", Ann. Of Surg., 695-703, May 1970.
7. Pruitt, J. Crayton, "1099 Consecutive Carotid Endarterectomies Using Local Anesthesia, EEG and Selective Shunting with Pruitt-Inahara Carotid Shunt, "Contemporary Surgery, Vol. 23, Sept. 1993".

**Symbol Legend • Legenda dei simboli • Symbollegende • Significado de los símbolos •
 Legende van symbolen • Legenda de símbolos • Légende des symboles • Sembol Açıklaması**

	English	Italiano	Deutsch	Español	Nederlands	Português	Français	Türkçe
	Usable Length	Lunghezza utilizzabile	Nutzbare Länge	Longitud utilizable	Bruikbare lengte	Comprimento Útil	Longueur utilisable	Kullanılabilir Uzunluk
	Catalogue Number	Numero di catalogo	Katalognummer	Número de catálogo	Catalogusnummer	Número de catálogo	Numéro de catalogue	Katalog Numarası
	Batch code	Codice lotto	Chargencode	Código de lote	Batchcode	Código de lote	Code de lot	Toplu kod
	Use-by date	Data di scadenza	Verfallsdatum	Fecha de caducidad	Verbruiksdatum	Uso por data	Date limite de consommation	Son kullanma tarihi
	Balloon Capacity	Capacità del palloncino	Ballonkapazität	Capacidad del balón	Balloncapaciteit	Capacidade do Balão	Capacité du ballon	Balon Kapasitesi
	Size	Misurare	Größe	Tamaño	Maat	Diâmetro	Taille	Boyut
	Manufacturer	Produttore	Hersteller	Fabricante	Fabrikant	Fabricante	Fabricant	Üretici
	Date of Manufacture	Data di produzione	Herstellungsdatum	Fecha de fabricación	Datum van vervaardiging	Data de Fabrico	Date de fabrication	Üretim Tarihi
	Sterilized Using Ethylene Oxide	Sterilizzato utilizzando ossido di etilene	Mit Ethylenoxid sterilisiert	Esterilizado con óxido de etileno	Gesteriliseerd met behulp van ethyleenoxide	Esterilizado con óxido de etileno	Stérilisé à l'oxyde d'éthylène	Etilen Oksit Kullanılarak Sterilize Edilmiştir
	Do not use if package is damaged.	Non utilizzare se la confezione è danneggiata.	Nicht verwenden, wenn die Verpackung beschädigt ist.	No usar si el paquete está dañado.	Niet gebruiken als de verpakking beschadigd is.	Não utilizar se a embalagem estiver danificada.	Ne pas utiliser si l'emballage est endommagé.	Paket hasarlıysa kullanmayın.
	Not made with natural rubber latex.	Non realizzato con lattice di gomma naturale.	Nicht aus Naturkautschuklatex hergestellt.	No contiene látex de caucho natural.	Niet gemaakt van natuurlijk rubberlatex.	Não fabricado com látex de borracha natural.	Non fabriqué avec du latex de caoutchouc naturel.	Doğal kauçuk lateksten yapılmamıştır.
	Do not re-use	Non riutilizzare	Nicht wiederverwenden	No reutilizar	Niet hergebruiken	Não reutilizar	Ne pas réutiliser	Tekrar kullanmayın
	Do not resterilize	Non risterilizzare	Nicht erneut sterilisieren	No reesterilizar	Niet opnieuw steriliseren	Não voltar a esterilizar	Ne pas restériliser	Tekrar sterilize etmeyin
	Store in a cool, dry place.	Conservare in un luogo fresco e asciutto.	An einem kühlen, trockenen Ort aufbewahren.	Guárdese en un lugar fresco y seco.	Op een koele, droge plaats bewaren.	Guarde num local fresco e seco.	Conserver dans un endroit frais et sec.	Serin ve kuru bir yerde saklayın.
	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Attenzione: La legge federale (USA) limita la vendita di questo dispositivo ai medici o su prescrizione di un medico.	Vorsicht: Laut Bundesgesetz (USA) darf dieses Gerät nur von einem Arzt oder auf Anweisung eines Arztes verkauft werden.	Advertencia: Las leyes federales de Estados Unidos establecen que solo se puede vender este dispositivo a personal médico o a petición de este.	Voorzichtigheid: De federale wet (VS) beperkt de verkoop van dit apparaat door of op voorschrift van een arts.	Aviso: A lei federal (dos EUA) limita a venda deste dispositivo aos médicos ou por ordem dos mesmos.	Prudence: La loi fédérale (États-Unis) restreint la vente de cet appareil par ou sur ordonnance d'un médecin.	Dikkat: Federal (ABD) kanunları, bu cihazın bir doktor tarafından veya bir doktor emriyle satışına izin vermemektedir.
	Caution	Attenzione	Vorsicht	Precaución	Voorzichtigheid	Aviso	Prudence	Dikkat
	Consult instructions for use	Consultare le istruzioni per l'uso	Gebrauchsanweisung beachten	Consulte las instrucciones de uso	Raadpleeg de gebruiksaanwijzing	Consultar as instruções de utilização	Consulter les instructions d'utilisation	Kullanım talimatlarına bakın
	Non-pyrogen	Apirogeno	Nicht pyrogen	No pirogénico	Niet-pyrogen	Não-pirógeno	Approgène	Pirojenik olmayan
	Balloon Diameter	Diámetro del palloncino	Ballondurchmesser	Diámetro del balón	Ballondiameter	Diâmetro do balão	Diamètre du ballon	Balon Çapı
	Distributor	Distributore	Verteiler	Distribuidor	Distributeur	Distribuidor	Distributeur	Distribütör
	Medical Device	Dispositivo medico	Medizinisches Gerät	Producto sanitario	Medisch apparaat	Dispositivo Médico	Dispositif médical	Tıbbi Cihaz
	Single sterile barrier system with protective packaging inside	Sistema a barriera sterile singola con imballo protettivo interno	Einzelnes Sterilbarriersystem mit Schutzverpackung im Inneren	Sistema de barrera estéril con embalaje protector interno	Enkelvoudig steriel barrièresysteem met beschermende verpakking aan de binnenkant	Sistema de barreira única com embalagem protetora no interior	Système de barrière stérile unique avec emballage de protection à l'intérieur	İçinde koruyucu ambalaj bulunan tek steril bariyer sistemi
	Keep away from sunlight	Tenere lontano dalla luce solare	Von Sonnenlicht fernhalten	Proteger de la luz solar	Uit de buurt van zonlicht houden	Manter longe da luz solar	Tenir à l'écart du soleil	Güneş ışığından uzak tutun
	Unique Device Identifier	Identificatore univoco del dispositivo	Eindeutige Gerätekennung	Identificador inequívoco de un producto sanitario	Unieke apparaat-ID	Identificador exclusivo do dispositivo	Identifiant unique de l'appareil	Benzersiz Cihaz Tanımlayıcı

Note: Not all symbols may be included in the labeling of this product. • **Nota:** Non tutti i simboli potrebbero essere inclusi nell'etichetta di questo prodotto. • **Hinweis:** Möglicherweise sind nicht alle Symbole in der Kennzeichnung dieses Produkts enthalten. • **Nota:** Es posible que no todos los símbolos aparezcan en el etiquetado de este producto. • **Let op:** Het is mogelijk dat niet alle symbolen op de etikettering van dit product voorkomen. • **Nota:** Poderão não estar incluídos todos os símbolos na rotulagem deste produto. • **Remarque:** Tous les symboles peuvent ne pas être inclus dans l'étiquetage de ce produit. • **Not:** Bu ürünün etiketinde tüm semboller bulunmayabilir.

Manufactured by
LUCAS MEDICAL, INC.
1554 North Case Street
Orange, CA 92867
Made in the USA
TEL: (714) 938-0233

U.S. PATENT NUMBERS 5,762,996 & 5,796,332
 PN 10020 – REV – J