

Case Study Example for a Mobile Application

ETHICON
A part of the Johnson & Johnson family of companies

Courses Podcasts 9:25 AM Videos Case Studies 57%

Umbilical Hernia Case
This case study app explains the anatomy and causes of umbilical and small ventral hernias and how to repair them using the PROCEED™ Ventral Patch.

CLEAR
PLAY

CLINICAL EXPERTISE.COM™ Expand your skills with e-modules & tutorials of procedures, webinars from industry leaders, simulations of surgical scenarios, and connect with peers around the world. REGISTER NOW AT CLINICALEXPERTISE.COM™

ETHICON
A part of the Johnson & Johnson family of companies

9:26 AM Umbilical Hernia Case: Treatment 86% EP-434-12

Table of Contents

- Introduction
 - 1 - Objectives
- Patient
 - 2 - General Patient History
 - 3 - Prior Medical History
 - 4 - Physical Examination
 - 5 - Abdominal Examination
 - 6 - Lower Body Examination
- Diagnosis
 - 7 - Diagnostic Tests


Page 10 of 20
Tap for Index

VIRTUAL EXPERT

Case Study Example for a Mobile Application

9:25 AM **Umbilical Hernia Case: Introduction** 87% EP-434-12

ETHICON
Part of the Johnson & Johnson family of companies



Objectives

After finishing this interactive case study you will be able to:

- Describe the risk factors for and presentation of umbilical hernia
- Define umbilical hernia and severity criteria
- Compare the different repair techniques for umbilical hernia
- Describe the advantages of mesh-based hernia repair
- List the steps of a mesh-based hernia repair using the PROCEED™ Ventral Patch

Disclaimer:

This application includes a demonstration of the use of a surgical device; it is not intended to be used as a surgical training guide. Other surgeons may employ different techniques. The steps demonstrated may not be the complete steps of the procedure. Individual surgeon preference and experience, as well as patient needs, may dictate variation in procedure steps. Before using any medical device, including those demonstrated or referenced in this application, review all relevant package inserts, with particular attention to the indications, contraindications, warnings and precautions, and steps for use of the device.

Page 1 of 20
Tap for Index

9:25 AM **Umbilical Hernia Case: Patient** 86% EP-434-12

ETHICON
Part of the Johnson & Johnson family of companies



General Patient History

The patient is a 52 year old white male, who presents in mild acute distress secondary to mild-abdominal pain. He is a plumber and earlier in the day, he had a short period of severe periumbilical pain while straining to remove a pipe. After a few minutes of nausea and one episode of vomiting, most of his pain resolved. He states that he has had a golf-ball-sized bulge at his umbilicus for as long as he can remember, but never previously experienced periumbilical or severe abdominal pain.

Page 2 of 20
Tap for Index



Case Study Example for a Mobile Application

IPad ETHICON 9:26 AM Outcome and Follow-Up: Treatment EP-434-12



Outcome and Follow-Up

The most appropriate procedure for this patient was the mesh-based retrodefect device repair. Following this procedure, the patient did very well post-operatively.


He was discharged home from the recovery room on oral pain medication with codeine. He discontinued the codeine on Postoperative Day 2 and started on ibuprofen, which he subsequently discontinued on Postoperative Day 3. The patient returned to most of his routine activities on Postoperative Day 5 and returned to work on Postoperative Day 6.

The patient was seen again by the surgeon on Postoperative Day 10. The patient had no complaints and examination was normal. Then at 3 months after surgery, the patient was healing nicely and had experienced no complications.

Page 9 of 20
Tap for Index

VIRTUAL EXPERT

IPad ETHICON 9:26 AM Umbilical Hernia Case: Treatment EP-434-12



Treatment Options

Now that you've had a chance to review the results of diagnostic tests, what treatment option might best repair this umbilical hernia?

- Laparoscopic hernia repair - Intraoperative onlay mesh (POM)
- Mesh-based retrodefect device repair
- Open primary hernia repair
- Open repair with onlay mesh
- Open repair with sublay mesh

The patient required general anesthesia and tolerated the 80-minute procedure well. In the recovery room, the patient complained of a very sore throat which resolved on post-operative day 1. The patient's abdomen was too sore to return to work for 10 days.

Page 8 of 20
Tap for Index

VIRTUAL EXPERT