




## Form and Fit

---

The Gaia Photonics Surgical Light Pod has been designed as a drop-in replacement for the Light Pods of the Skytron Aurora and Aurora II surgical Lighting systems (Skytron part numbers B1-610-38, B1-710-0, B1-710-1).


The exact mechanical fit and compatibility with the existing Skytron Aurora electronics make the replacement of the Light Pods effortless and time efficient. No changes to the system are required.




## Performance Enhancements:

Optical performance starts with choosing the correct light source for the application. State-of-the-art ceramic LED packages from an industry leading Tier 1 LED supplier have been selected to power our Light Pods. This ensures not only the highest LED performance and reliability, but also optimum color rendering and minimal chromaticity drift over the life of the product.

A solid optomechanical design for mounting all critical components in exactly the right position and a consistent and repeatable assembly process result in a significantly improved beam uniformity, consistency from unit to unit and stability over time.



Combining the improvements in optical alignment with the state-of-the-art large format ceramic LEDs used in the Surgical Light Pods enables Aurora LED5 Lightheads to deliver more than 300,000 lux. This more than twice the flux specified in the original product.



Selection of three different types of LEDs with very high color rendering index allows for a chromaticity adjustment between 4100K and 4500K while maintaining exceptional color rendering Indices of CRI  $\geq 97$ , R9 (red)  $\geq 93$  and R13 (skin tone)  $\geq 99$ .

Metric	20-00001		OEM	
	4100 K	4500 K	4100 K	4500 K
CRI	97	96	95	96
R9 - deep red	93	97	87	97
R13 – skin	99	99	NR	NR
TM-30 Rf	94	94	NR	NR
TM-30 Rg	103	103	NR	NR

NR = not reported



## Reliability:

---





A thermal design that ensures the heat generated by the LEDs is efficiently spread and conducted to the carrier minimizes the temperature rise of the LED junction over ambient.



The efficient thermal design combined with the use of large format LED die capable of handling operating currents multiple times of the actual application warrants a L90 lifetime prediction of the Gaia Surgical Light pod of more than 50,000hrs.


## Process Control:

Selection of the best components, a solid design as well as repeatable manufacturing and assembly processes establish a very consistent product performance including center illuminance level, chromaticity coordinates, color rendering and beam distribution. The following charts display actual product variation with a 3 sigma tolerance band (4 Step MacAdam ellipse for chromaticity).



## Individual Performance Reports

Every Gaia Photonics Surgical Light Pod undergoes a final test and alignment process to warrant that the product meets and exceeds our specifications.



## Technical Specifications:

---

The Gaia Photonics Surgical Light Pod 20-00001 is a UL recognized component under UL file E548175.



The Gaia Photonics Surgical Light Pod 20-00001 is exclusively designed for installation in Skytron Aurora LED and Skytron Aurora II series Lightheads and is only to be operated under the following conditions:

Product Ratings	
Input current CC	410 mA
Input Power	10W
Input Type	Isolated Circuit
Environmental conditions	Damp or Dry

The typical operating characteristics of a single standalone unit at 410mA constant current operation and in the Skytron LED5 Lighthead with Skytron driver at dimming setting 5 providing the equivalent of 410mA constant current drive current are captured in the following table.

Configuration	20-00001 Single Light Pod <sub>1</sub>		20-00001 Surgical Light pod in LED5 Light Head <sub>2</sub>	
	Amber off	Amber on	Amber off	Amber on
<b>CCT setting</b>				
<b>Color temperature</b>	4400K	4100	4400K	4100K
<b>CRI</b>	96	97	97	96
<b>R9 (red)</b>	97	93	96	91
<b>R13 (skin tone)</b>	99	99	99	99
<b>Chromaticity x</b>	0.3631	0.3744	0.3635	0.3749
<b>Chromaticity y</b>	0.3644	0.3697	0.3622	0.3680
<b>Center Illuminance Ec [lx]</b>	62000	71000	260000	310000
<b>Light field diameter d50 [mm]</b>	87	87	90	90
<b>Light field diameter d10 [mm]</b>	138	138	140	140
<b>Center Irradiance [W/m<sup>2</sup>]</b>	230	270	980	1200
<b>UV center Irradiance [W/m<sup>2</sup>]</b>	0.03	0.035	0.24	0.28

Illuminance and Irradiance measurements performed at 1m distance

<sub>1</sub>Single Light Pod operated at 410mA DC current through both channels

<sub>2</sub> 20-00001 in LED5 Light Head operated with Aurora II wall control in dim5 setting

