

**FEATURES**

**DP40T**

- ◆ Measuring Specific Rotation Optical Rotation, Concentration and Sugar Degree.
- ◆ Polarimeter adopts the principle of photoelectric automatic measurement.
- ◆ 4 measurement modes to choose. It automatically retested 6 times and calculate the average value
- ◆ Touch screen Colour Display.
- ◆ USB communication Interface.



**SPECIFICATIONS**

**HSN : 90275030**

<b>PRODUCT NAME</b>	<b>AUTOMATIC POLARIMETER</b>
<b>MODEL</b>	<b>DP40T</b>
Display	Touch screen
Measuring Mode	Optical Rotation, Specific Rotation, Sugar Degree, Concentration
Light Source	LED + Interference Filter
Wavelength	589.44nm
Measuring Range	± 45° (Optical Rotation)
Measuring Range - Sugar	± 120°Z
Minimum Indication Value	0.001°
Repeatability- Optical Rotation	≤ 0.003°
Test Tube – Standard	100mm & 200mm
Transmission	10%
Interface	USB
Power Supply	AC 220V / 50Hz
Dimension & Gross Weight	730 x 480 x 470 mm & 36kg

\* Specifications are subject to change without Notice

**Functions**

**DP40T**




Mode: Optical Rotation  
Tube Length: 200  
Measure Times: 1  
 $\sigma_{n-1}$ :  
**α** 0.000

Mode  
Tube Length  
Measure Times

**Mode Selection**




Optical Rotation  
Specific Rotation  
Concentration  
Sugar level  
Mode: Specific Rotation  
TUBE SELECTION  
Optical Rotation Specific Rotation  
Concentration Sugar Scale  
 $\sigma_{n-1}$ :  
**α** 0.000



Mode: Optical Rotation  
Tube Length: 200  
TUBE SELECTION  
50 mm 100 mm 200 mm  
 $\sigma_{n-1}$ :  
**α** 0.000

Tube length.  
Press 'Tube length'  
to choose

**Measure times.**  
Press 'measure times'  
to choose 1~6.




Mode: Specific Rotation  
Tube Length: 200  
Measure Times: 1  
Concentration: 1.0000  
 $\sigma_{n-1}$ :  
**α** 0.000




Date & Time  
Meas.  
System Set  
**α** 0.000

System settings.  
Press 'System set'

**Measure specific rotation:**  
(Need input tube length,  
concentration and the  
measure times. a  
—value of specific rotation.  
a<sub>av</sub>—average value)




Mode: Specific Rotation  
Tube Length: 200  
Measure Times: 1  
Concentration: 1.0000  
 $\sigma_{n-1}$ :  
**[α]** 0.000



Mode: Concentration  
Tube Length: 200  
Measure Times: 1  
Concentration: 1.0000  
 $\sigma_{n-1}$ :  
**C** 0.000

Measure concentration:  
(Need input tube length,  
specific rotation and the  
measure times. )

**Measure sugar degree**  
(Need input tube length  
and the measure times.)



Mode: Sugar Scale  
Tube Length: 200  
Measure Times: 1  
Specific Rotation: 1.0000  
 $\sigma_{n-1}$ :  
**Z** 0.000