

Step 1. Obtain a prescription for the test

The process starts with the Healthcare Professional (HCP) ordering the test for you and scheduling an appointment for your specimen collection.

The HCP ordering the test will initiate one of two options:

- The HCP will provide a test kit and shipping supplies directly to you; or
- P23 Lab will ship a test kit overnight directly to your home after you authenticate your order (Step 2 below). A valid email address will be required.



Step 2. Authenticate your order

- Go to P23 Lab's website at www.P23Labs/covid19
- Verify your identity by uploading your driver license or a current photo of yourself.
- You will be prompted for a response to the following question:

Did you elect to receive a test kit via the HCP's office?

If the answer is **yes**, do you have a scheduled appointment for your supervised specimen collection? *The specimen collection and packaging must be supervised by a HCP.*

If the answer is **no**, contact your HCP to schedule your specimen collection appointment.

- If payment is required, you will be asked to pay at this time.




Step 3. Receive your kit

Once you have received your kit, initiate the following steps:

- Go to P23 Lab's website at www.P23Labs/register
- Inspect your test kit to ensure all components are included.



- Locate your unique 14-digit barcode on the collection tube and register your kit by entering the barcode on the secure patient portal.



IMPORTANT: Follow all directions carefully to ensure an adequate specimen is collected for testing.

Do not collect your specimen until your scheduled appointment with your HCP. Do NOT eat, drink, smoke or chew gum for 30 minutes before giving your saliva specimen. Do not remove the plastic film on the tube's funnel.

For questions at any time, please call the laboratory help line at **901-341-6464**. Our professionals are available to guide you.

This home collection kit has not been FDA cleared or approved.

This home collection kit has been authorized by FDA under an EUA.

This home collection kit has been authorized only for the home collection and maintenance of saliva specimens as an aid in detection of nucleic acid from SARS-CoV-2, not for any other viruses or pathogens.



IMPORTANT: Do NOT eat, drink, smoke or chew gum for 30 minutes before giving a saliva specimen.

Healthcare Professionals (HCP) trained according to 49CFR 172.700 can find additional guidance at:

Collection steps: <https://www.dnagenotek.com/US/support/collection-instructions/omnigene-oral/om-505.html>

Collection video: <https://www.youtube.com/watch?v=kwb4r2moy3g>

Step 1. Collection appointment


Under the **supervision** of your HCP during your appointment, you will collect your saliva specimen.

Note: This may be done in person at the HCP office or via online Telemedicine.

Step 2. Specimen collection

Complete the specimen collection with your HCP's guidance:


- IMPORTANT:** Wash your hands **thoroughly for 20 seconds** then dry your hands **before** starting the collection.
- Open your kit and place all of the contents on a clean, dry surface.
- Remove the saliva collection device from its packaging and collect your specimen.



- IMPORTANT:** Wash your hands **thoroughly for 20 seconds** again, then dry your hands **after** completing the collection.

Step 4. Specimen shipping preparation with your HCP

- Make sure the cap on your specimen tube is tightly closed.
- Place specimen into the bio hazard bag and press firmly across the closure to securely zip.
- Follow HCP instructions for placing the bagged specimen onto the protective bubble wrap in the shipping box. Place the absorbent sheet on top and close the shipping box lid.
- IMPORTANT:** Wash your hands **thoroughly for 20 seconds** again, then dry your hands **after** completing the shipping preparation.
- The shipping box goes into the labeled return overpack, remove the adhesive liner and seal overpack closed.
- Drop off the sealed return overpack containing your specimen at the closest Fedex or UPS depot as determined by the shipping label with P23 Lab's address.



IMPORTANT: Drop off the return overpack at the shipping depot within 24 hours of specimen collection.

Step 5. Results are ready within 72 hours

Your HCP will contact you when your results are available.

For questions at any time, please call the laboratory help line at **901-341-6464**. Our professionals are available to guide you.