

Extraction Free Portable COVID-19 RT-PCR Detection Solutions



Phone: 800-848-7240



Less than 90 minutes from sample to result
up to 16 sample a time

A New COVID-19 detection solution

RNA Extraction Free! ! !

01

Sample Collection

Collect swab sample in VTM

02

Release viral RNA

Vortex or shake for 2 minutes

03

Transfer to PCR Tube

Pipette 17 uL VTM into 8 uL PCRmastermix
Make 25 uL PCR reaction mix

04

PCR Detection

Load the reaction mix onto Fast-16 real time
PCR instrument, start experiment run.
Complete the PCR thermal cycling, read the
results.



Sample Collection

1-2 min



**RNA Extraction
Extra
Instrument
Extra labor**

35-45 min



PCR Detection

55 min Fast Protocol
80 min Standard Protocol



Results

The new generation PCR technology

Accuracy And Safety

sensitivity

97.9%



Manufacturing Capacity

100,000 tests/day

sample

Detection limit
1000 copies/mL (Standard protocol)
1500 copies/mL (Fast protocol)
Sample volume
up to 17uL

3 Gene


specificity

95.7%




3 Gene Detected
ORF1ab(FAM)
N gene (HEX)
Reference gene(Cy5)


FAM channel
ORF1ab gene



Cy5 channel
Internal control gene



HEX channel
N gene



Faster, Lower Cost
Same Performance

01

Same Ct value

The Extraction free 1 step RT-PCR detection kit generates same Ct value as regular RT-PCR detection kit with RNA extraction step.

02

Viral RNA released completely

Without additional RNA extraction step and instrumentation, releasing agent is able to release viral RNA completely and generate results comparable with regular RNA extraction required RT-PCR reaction.

Portable Instrument Size

Adapt to Mobile or On-board Testing Environment

- ✓ Small footprint PCR instruments for high throughput or portable market needs.
- ✓ Portable "shoe box size" PCR instrument process 16 samples per run.
- ✓ Minimal hands on time, only simple vortex and pipetting required.
- ✓ Gold Standard RT-PCR Technology

16 Sample per Run Portable RT PCR instrument



Total < 90 minutes for Processing 16 samples

Instrument Size:
(L x W x H) 320 x250 x 177mm
4.5 KG

Read Positive/Negative directly from the instrument touch screen. Actual unit physical design may vary slightly.

4.3.3.4 Data


Well	Sample Name	Channel	Quantitative	CT	TM	Result	Barcode
A4	n						
B4	Sample1			42.34	80.91	Negative	
A5	p			27.18	81.42	Positive	
B5	Sample1						

Criteria Setting Report

v1.0.3 User: Public 08:25 2017/02/14

Criteria setting: choose qualitative in detection type. Enter the qualitative threshold and click ok to view the analysis results in the data interface after clicking the analysis button on the amplification plot interface.


Extractionless COVID-19 RT-PCR Detection Kit Workflow



1. Drop swab to VTM (virus transport medium) containing COVID-19 viral release reagent
2. Vortex or shake VTM tubes for 2 mins
3. Spin down solution, pipette 17 ul VTM solution to 8 ul PCR mix, prepare 25 ul reaction volume
4. Load reaction mix onto the Mini RT PCR instrument, start experiment run
5. Complete PCR thermal cycling and read results in 55 (Fast protocol) minutes or 81 (standard protocol) minutes for 16 samples
6. Generate reports of SARS-CoV-2 positive or SARS-CoV-2 negative
7. Reports can be generated via Excel or PFD files
8. Export reports using USB disk. Exported reports can be printed or store on labtop.

448 Tests of Standard Process per day (2 shifts) just
 4 PCR instrument + 2 Technicians

Standard Protocol



4 Mini 16 Well Plate RT PCR + 2 Technicians



224 Samples in 6 hrs


X 2 Shifts (12 hours)

448 Samples/ Day

			6 Hours		
	Time/Run	Tests/Run	No of Runs	No of Tests	Note
Fast-16 Portable PCR instrument	81 min	14	4	56	16 well plate (14 specimens + 1 Negative control + 1 positive control)

672 Tests of Standard Process per day (2 shifts)
just 4 PCR instrument + 2 Technicians

Fast Protocol



4 Mini 16 Well Plate RT PCR instruments + 2 Technicians




336 Samples in 6 hrs

X 2 Shifts (12 hours)

672 Samples/Day

				6 Hours	
	Time/Run	Tests/Run	No of Runs	No of Tests	Note
Fast-16 Portable PCR instrument	55 min	14	6	84	16 well plate (14 specimens + 1 Negative control + 1 positive control)

Order Information



Product Name

Fast-16 (Portable qPCR)



**SARS-CoV-2 Detection
Kit (RT-PCR, RNA extraction free)
PCR assay+ Sample collection VTM**



Sample Vortexer (50 samples throughput) (Optional)



Portable Power Station (Optional)

Manufacturing Compliance



ISO 13485 Certified