

*Verdel Instruments: Total Correlation Mass Spectrometry™:
Fast, effective, true DIA*




VERDEL
INSTRUMENTS LTD



Innovate UK


Mass Spectrometry specialises in analysis of mixtures

- Mass Spectrometry gives the total number of species and overall mass.
- Further tandem mass spectrometry is needed to increase confidence of identification.



Mass Spectrometry specialises in analysis of mixtures


- Quadrupoles have a non-zero isolation width.
- Possibility of overlap producing mixed/chimeric spectra decreasing confidence of identification.
- Practical isolation wide $> 1 m/z$



Liquid Chromatography is often added

- DDA –LCMS systems operate on TopX structure.
- High resolution MS scan detects precursors and followed by MS/MS scans of precursors detected.
- Top3, Top5, and Top10 strategies are common but precursors may be ignored in the interest of speed.


Key targets often not high intensity




Is there another option?

- Isolation results in the loss of precursor information.

Tandem mass spectrometry




Total Correlation Mass Spectrometry™



Is there another option?

- How can we correlate precursors and fragments without separation: Total Correlation Mass Spectrometry TM.




- Total Correlation Mass Spectrometry is the name for the FTICR 2DMS technique applied using a QTOF

Is there another option?

- TC-MS™ – groups fragments and precursors together **without** prior isolation – no LC or quadrupole isolation needed.
- True DIA – the number of precursors does not affect total acquisition time.
- Developed on an FT-ICR– Verdel is moving the technology to a number of QTOF instruments from different manufacturers.

How does TC-MS™ outside of FT-ICR work?


- Modified RF power supply for quadrupole modulates ions.



*Repeated multiple times
with different modulation*


How does TC-MS™ outside of FT-ICR work?

- The radii of precursors are iterated multiple times producing an inverse relationship between precursors and their corresponding fragments:




How does TC-MS™ outside of FT-ICR work?

- Small modification of standard TOF instrument (Bruker MaXis).




How does TC-MS™ outside of FT-ICR work?

- Addition of Verdel RF power supply for the Quadrupole.
- Addition of external UV laser and mirror setup within TOF section.




How does TC-MS™ outside of FT-ICR work?

- Quadrupole modification can be carried out without quadrupole replacement.



How does TC-MS™ outside of FT-ICR work?


- Use of UVPD gives method with a wide array of applications, highly robust to sample type.



*Verdel's TC-MS™ technology
allows parallel analysis of
all molecules in a mixture
without prior separation,
removing chromatography.*


What results does TC-MS™ generate?

- TC-MS produces large TC spectra where fragments are correlated with their precursor ions within m/z space.




What results can TC-MS™ generate?

- The Autocorrelation line contains all the precursors.



TC-MS™ data

- The fragments corresponding to each precursor are present on horizontal fragment lines.



Facile Determination of Phosphorylation Sites in Peptides Using Two-Dimensional Mass Spectrometry

Johanna Paris, Tomos E. Morgan, Christopher A. Wootton, Mark P. Barrow, John O'Hara, and Peter B. O'Connor*

- TC-MSTM of a phosphopeptide mixture.

- Autocorrelation and fragment lines show high species characterisation.







Figure 1. (A) Full 2DMS spectrum of phosphopeptide mixture; color code in Figure 1A: purple, 2+ species; pink, 3+ species. (B) Zoom of [Ph01 + 2H]²⁺. Resolution in the vertical and the horizontal axes. (C) Autocorrelation line reveals all the fragmented precursors (similar to a 1D MS spectrum). (D) Extracted 2+ phosphate neutral (H₃PO₄) loss line. (E) Extracted 3+ phosphate neutral (H₃PO₄) loss line.

Differentiating Fragmentation Pathways of Cholesterol by Two-Dimensional Fourier Transform Ion Cyclotron Resonance Mass Spectrometry

Maria A. van Agthoven,¹ Mark P. Barrow,¹ Lionel Chiron,² Marie-Aude Coutouly,² David Kilgour,³ Christopher A. Wootton,¹ Juan Wei,¹ Andrew Soulbey,¹ Marc-André Delsuc,^{2,4} Christian Rolando,⁵ Peter B. O'Connor¹



Publications describing TC-MS™

analytical chemistry Article
pubs.acs.org/ac

Two-Dimensional Mass Spectrometry for Proteomics, a Comparative Study with Cytochrome c



Maria A. van Agthoven,^{||} Christopher A. Wootton,^{||} Lionel Chiron,[⊥] Marie-Aude Coutouly,^{§,¶} Andrew Soulby,^{||} Juan Wei,^{||} Mark P. Barrow,^{||} Marc-André Delsuc,^{‡,⊥} Christian Rolando,^{†,‡} and Peter B. O'Connor^{*,||}


European Biophysics Journal (2019) 48:213–229
<https://doi.org/10.1007/s00249-019-01348-5>

TECHNICAL PRIMER

Two-dimensional mass spectrometry: new perspectives for tandem mass spectrometry

Maria A. van Agthoven¹ · Yuko P. Y. Lam¹ · Peter B. O'Connor¹ · Christian Rolando² · Marc-André Delsuc^{3,4}

American Society for Mass Spectrometry, 2016  J. Am. Soc. Mass Spectrom. (2016) 27:1531–1538
DOI: 10.1007/s13361-016-1431-z

AS MS RESEARCH ARTICLE

2D FT-ICR MS of Calmodulin: A Top-Down and Bottom-Up Approach

Federico Floris,¹ Maria van Agthoven,¹ Lionel Chiron,² Andrew J. Soulby,¹ Christopher A. Wootton,¹ Yuko P. Y. Lam,¹ Mark P. Barrow,¹ Marc-André Delsuc,^{2,3} Peter B. O'Connor¹

analytical chemistry Article
pubs.acs.org/ac

Polymer Analysis in the Second Dimension: Preliminary Studies for the Characterization of Polymers with 2D MS

Federico Floris,[†] Claudio Vallotto,[†] Lionel Chiron,[‡] Alice M. Lynch,[†] Mark P. Barrow,[†] Marc-André Delsuc,^{‡,§} and Peter B. O'Connor^{*,†}


Analyst


ARTICLE

Uncoiling Collagen: A Multidimensional Mass Spectrometry Study

H. J. Simon,^a M. A. van Agthoven,^a P. Y. Lam,^a F. Floris,^a L. Chiron,^b M.-A. Delsuc,^{b,c} C. Rolando,^d M. P. Barrow,^a and P. B. O'Connor^a

Received 00th January 20xx, Accepted 00th January 20xx



American Society for Mass Spectrometry, 2015  J. Am. Soc. Mass Spectrom. (2015) 26:2105–2114
DOI: 10.1007/s13361-015-1226-7

AS MS RESEARCH ARTICLE

Differentiating Fragmentation Pathways of Cholesterol by Two-Dimensional Fourier Transform Ion Cyclotron Resonance Mass Spectrometry

Maria A. van Agthoven,¹ Mark P. Barrow,¹ Lionel Chiron,² Marie-Aude Coutouly,² David Kilgour,³ Christopher A. Wootton,¹ Juan Wei,¹ Andrew Soulby,¹ Marc-André Delsuc,^{2,4} Christian Rolando,⁵ Peter B. O'Connor¹



Huge array of applications

- Any analyte present within a complex mixture and fragmented by UVPD
- Proteins/peptides
 - Large m/z range accessible making top-down and bottom-up possible
- Polymers
- Lipids
- Agrochemicals
 - Present within complex matrices
- Small molecule pharmaceuticals

- Plus many more

One fast, parallel, and detailed experiment

- Removes the need for separation techniques.
 - Eliminates column bias or variation
- Analyses all the ions in a sample, every time; offers true DIA.
- No prior sample knowledge needed prior to analysis.
- Offers both high quality data and rapid analysis.
- Can be supplied as an upgrade to a range of QTOF instruments from different suppliers.
- No specialist training needed.



Please contact us with any further questions.

www.verdelinstruments.co.uk