

Determination of Water Content in Peanut Butter

1. Apparatus

- a. AQUACOUNTER Coulometric Karl Fischer Titrator Model AQ-2100S
- b. AQUACOUNTER Oil Evaporator EV-2000L
- c. Grease Sampler (P/N E730017-A)
- d. Evaporation Chamber Type A (P/N E323310-A)

2. Reagents

- a. Anode solution: HYDRANAL Coulomat AG
- b. Cathode solution: HYDRANAL Coulomat CG
- c. Azeotropic Solvent: Dehydrated Toluene

3. Test Procedure

- 1) Add 10mL of dehydrated toluene into the evaporation chamber of EV-2000L.
- 2) Lift the heating block up to eliminate the moisture contained in the toluene.
- 3) Take approx. 0.3 g of sample using the grease sampler.
- 4) After the background level becomes stable, press the SAMPLE key.
- 5) Remove the glass stopper of evaporation chamber and insert the grease sampler into the evaporation chamber.
- 6) Press the TITRATION key to start the measurement.
- 7) When the titration is completed, discard the toluene and the sample in the evaporation chamber.

Titration parameter settings:

Parameters	Values
Cal Mode	0
T-Timer	10 min
Current	Medium
Heating Temperature	120°C
Carrier Gas Flow	50mL/min

4. Result and Discussion

Table 1. Analytical Result of Peanut Butter

Sample		Sample Size (g)	Water (ugH ₂ O)	Water (%)	Statistics	
A03-068-01	1	0.3764	3944.2	1.0479	AVG(3)	1.04 %
	2	0.3213	3308.4	1.0291	SD	0.011%
	3	0.3135	3229.8	1.0302	CV	1.0 %
A03-068-02	1	0.3140	3554.8	1.1321	AVG(3)	1.12 %
	2	0.3800	4303.1	1.1324	SD	0.018%
	3	0.3104	3419.7	1.1017	CV	1.6 %
A03-068-03	1	0.3340	3727.9	1.1161	AVG(3)	1.10 %
	2	0.3296	3649.9	1.1074	SD	0.014%
	3	0.3124	3400.5	1.0885	CV	1.3 %

Table 1. shows the analytical result of peanut butter samples. Using the oil evaporator EV-2000L, good reproducibility and water content very close to the customer report was obtained.


On the other hand, the measurement using solid evaporator EV-2000 was examined and only 50% of the water was measured at heating temperature of 125 degree C, extraordinary long measurement time was taken. Figure 1 through 3 show the titration curve of each sample, Figure 4 shows the titration curve of the measurement using solid evaporator.

Sample should be stirred and mixed well prior to sampling.

```

***** AQ-2100 Result *****
Analysis Name  Peanst Butter
Date          2005/04/11 17:11
Sample No.    2
ID Code       A03-068-01
Total H2O     3229.8 ug
Tit Time      0:10:01
Back Ground   3 ug/min
H2O           3229.8 ug
Size          0.31350 g
*Conc         1.0302 %
Temp.         20.0 °C
Humid.        65.0 %

```




T-Data Time	ug	Conc
0:01:00	1002.0	0.3196 %
0:02:00	2027.6	0.6468 %
0:03:00	2656.4	0.8473 %
0:04:00	2842.6	0.9067 %
0:05:00	2961.8	0.9448 %
0:06:00	3046.1	0.9716 %
0:07:00	3109.5	0.9919 %
0:08:00	3164.0	1.0093 %
0:09:00	3203.4	1.0218 %
0:10:00	3229.8	1.0302 %

Fig 1. Titration curve for the sample A03-068-01

```

***** AQ-2100 Result *****
Analysis Name  Peanst Butter
Date          2005/04/12 11:00
Sample No.    2
ID Code       A03-068-02
Total H2O     4303.1 ug
Tit Time      0:10:01
Back Ground   3 ug/min
H2O           4303.1 ug
Size          0.38000 g
*Conc         1.1324 %
Temp.         -----
Humid.        -----

```




T-Data Time	ug	Conc
0:01:00	1048.3	0.2759 %
0:02:00	2249.7	0.5920 %
0:03:00	3451.0	0.9082 %
0:04:00	3875.0	1.0197 %
0:05:00	4000.2	1.0527 %
0:06:00	4100.5	1.0791 %
0:07:00	4167.9	1.0968 %
0:08:00	4222.4	1.1112 %
0:09:00	4272.8	1.1244 %
0:10:00	4303.1	1.1324 %

Fig 2. Titration curve for the sample A03-068-02

```

***** AQ-2100 Result *****
Analysis Name  Peanst Butter
Date          2005/04/12 12:38
Sample No.    3
ID Code       A03-068-03
Total H2O     3649.9 ug
Tit Time      0:10:01
Back Ground   3 ug/min
H2O           3649.9 ug
Size          0.32960 g
*Conc         1.1074 %
Temp.         -----
Humid.        -----


```



T-Data Time	ug	Conc
0:01:00	1044.6	0.3169 %
0:02:00	2246.0	0.6814 %
0:03:00	3088.4	0.9370 %
0:04:00	3302.7	1.0020 %
0:05:00	3402.0	1.0322 %
0:06:00	3486.4	1.0578 %
0:07:00	3540.8	1.0743 %
0:08:00	3590.3	1.0893 %
0:09:00	3620.5	1.0985 %
0:10:00	3649.9	1.1074 %

Fig.3 Titration curve for the sample A03-068-03

***** AQ-2100 Result *****
 Analysis Name Peanst Butter
 Date 2005/04/11 15:33
 Sample No. 1
 ID Code A03-068-01
 Total H2O 767.8 ug
 Tit Time 0:30:02
 Back Ground 4 ug/min
 H2O 767.8 ug
 Size 0.15350 g
 *Conc 0.5002 %
 Temp. 20.0 °C
 Humid. 65.0 %



T-Data Time	ug	Conc	°C
0:01:00	33.9	220.8 ppm	129
0:02:00	88.1	573.9 ppm	125
0:03:00	164.3	0.1070 %	124
0:04:00	231.4	0.1507 %	125
0:05:00	285.7	0.1861 %	126
0:06:00	333.9	0.2175 %	126
0:07:00	369.0	0.2404 %	126
0:08:00	399.0	0.2599 %	125
0:09:00	429.0	0.2795 %	125
0:10:00	454.0	0.2958 %	125
0:11:00	479.0	0.3121 %	125
0:12:00	502.9	0.3276 %	125
0:13:00	523.9	0.3413 %	125
0:14:00	543.9	0.3543 %	125
0:15:00	559.8	0.3647 %	125
0:16:00	578.7	0.3770 %	125
0:17:00	596.6	0.3887 %	125
0:18:00	613.6	0.3997 %	125
0:19:00	628.5	0.4094 %	126
0:20:00	643.4	0.4192 %	125
0:21:00	658.3	0.4289 %	125
0:22:00	668.1	0.4352 %	126
0:23:00	683.0	0.4450 %	126
0:24:00	697.9	0.4547 %	125
0:25:00	707.8	0.4611 %	126
0:26:00	722.7	0.4708 %	125
0:27:00	732.6	0.4773 %	126
0:28:00	747.5	0.4870 %	126
0:29:00	757.3	0.4934 %	126
0:30:00	767.8	0.5002 %	126

Fig.4 Titration curve for the sample A03-068-0 1 using Solid evaporator EV-2000.