



AQUACOUNTER Application Sheet	COM series	DATA No. D5	1st edition
Environmental	Measurement of alkalinity of mineral water		

1. Measurement outline

Natural water contains alkaline components such as hydroxides and bicarbonates. Such water is alkaline and alkalinity is used as an index. Alkalinity is expressed as mg/L of calcium carbonate (CaCO₃) equivalent for these alkaline components. Alkalinity is divided into phenolphthalein alkalinity (P alkalinity) and total alkalinity (T alkalinity or M alkalinity) by the pH value of neutralization point.



In the measurement method, 1/2 of the hydroxides and carbonates are measured when it is titrated to about pH8.3 with sulfuric acid titrant.



All of the bicarbonates are neutralized when it is titrated successively to about pH4.8.



This section introduces a measurement example using potentiometric titration for end point detection method in conformance to the test method for drinking water.


2. Reagents and Electrodes

(1) Reagents	Titrant	0.02mol/L sulfuric acid titrant
(2) Electrodes	Indicator electrode	*Glass electrode GE-101B
	Reference electrode	*Reference electrode RE-201 (4M potassium chloride inner solution used)
*standard accessories		

3. Measurement conditions example (for COM-1600S + 1unit of Buret B-2000-20)

Master File No.1					
Condition file : 1 + 2					
Parameters for Condition file 1 (For 1 st EP)		Parameter for Condition file 2 (For 2 nd EP)		Mode No. 21	
Method	Set	Method	Set	Pre Int	0 sec
Amp No.	1	Amp No.	1	Del K	5
Buret No.	1	Buret No.	1	Del Sens	0 mV
Meas Unit	pH	Meas Unit	pH	Int Time	3 sec
S Timer	10 sec	S Timer	0 sec	Int Sens	3 mV
CP pH	14.00 pH	CP pH	14.00 pH	Brst Speed	2
Direction	Down	Direction	Down	Pulse	8
DP pH	10.00 pH	DP pH	10.00 pH		
End Point	8.30 pH	End Point	4.80 pH		
Over mL	0.00 mL	Over mL	0.00 mL		
Max. Vol.	30 mL	Max. Vol.	20 mL		
Unit	mgEq/L	Unit	mgEq/L		
Size	100 g	Size	100g		
Blank	0	Blank	0		
Factor	Titer of the titrant	Factor	Titer of the titrant		
Molarity	0.02	Molarity	0.02		
K	0	K	0		
Formula	D×F×M×1000/S	Formula	D×F×M×1000/S		
Mode No.	21	Mode No.	21		

4. Measurement example



Initial	8.99 pH
Temp	25.0 °C
EP	8.30 pH
	0.378 ml
dE/dV	1
SIZE	100.0000 g
*Conc	0.0770 mgEq/L
EP	4.80 pH
	2.842 ml
dE/dV	1
SIZE	100.0000 g
*Conc	0.5789 mgEq/L

Measurement results on alkalinity of mineral water

Sample No.	Sample volume (mL)	P alkalinity		T (M) alkalinity	
		Titration value (mL)	P alkalinity	Titration value (mL)	T (M) alkalinity
1	100	0.380	0.0774	2.851	0.581
2	100	0.364	0.0741	2.852	0.581
3	100	0.378	0.0770	2.842	0.579
Avg.		0.0762 CaCO₃mg/L		0.580 CaCO₃mg/L	
Std. Dev.		0.0018 CaCO₃mg/L		0.0012 CaCO₃mg/L	
C.V.		2.40 %		0.20 %	

5. Note

Potentiometric titration was conducted for alkalinity measurement instead of titration with indicator, and measurement was possible with standard deviation 0.002 CaCO₃mg/L and coefficient of variation within 3%.

Key words

Alkalinity, P alkalinity, T (M) alkalinity

Hitachi High-Technologies Corporation

Head Office 1-24-14, Nishishinbashi, Minato-Ku, Tokyo 105-8717, Japan

Tel : 81-3-3504-7239 Fax : 81-3-3835-7302

<http://www.hitachi-hitech.com>

Hiranuma Sangyo Co., Ltd.

1739, Motoyoshidacho, Mito-City, Ibaraki 310-0836, Japan

Tel : 81-29-247-6411 Fax : 81-29-247-6942

<http://www.hiranuma.com>