

# Determination of Scrophularia root

Scrophularia root is the dried root of Scrophularia ningpoensis Hemsl. In this report, the analysis of harpagide and harpagoside in Scrophulariaceae is performed using a Hitachi Primaide HPLC system refers to the first part of the 2020 edition of the Chinese Pharmacopoeia, with a GL Sciences InertSustain AQ-C18 column as the analytical column, acetonitrile-0.03 % phosphoric acid solution as the mobile phase, and the detection is carried out with a UV detector at 210 nm.

The analysis of standard samples and the determination of commercially available samples is measured.

Liquid Chromatograph Primaide

High performance



System suitability

Concentration(mg/L)

Component	theoretical plate number in Pharmacopoeia	Measured value of theoretical plate number
harpagide	≥5000	33400
harpagoside	≥5000	386000

SD

RSD

0.007662

0.07%

413

0.23%

0.004676

0.02%

The good repeatability of the retention time and area was obtained. Harpagoside and harpagoside has a good linear relationship. The measurement results show that the method meets the requirement of the theoretical plate number specified in the Chinese Pharmacopoeia.

551

0.18%





## **Determination of Scrophulariacease sample**

Analysis of Scrophulariaceae sample



#### Result of Quantitative Analysis

Sample name	Component	Measurement result (total content %)	Labeled value (total content %)	Judgement
Scrophulariaceae	harpagide	0.71	≥0.45	ОК
	harpagoside			

The contents of harpagide and harpagoside in the commercially available Scrophulariaceae samples were determined, and the determination results met the requirements of the pharmacopoeia.

#### Sample preparation method

Scrophulariaceae sample	Weigh 0.5g medicine powder
$\bigvee$	add 50 mL of 50% methanol
Pass through No. 3 sieve after crushing	Weigh
$\bigvee$	$\bigvee$
Get powder	Soak for 1 hour
	$\bigvee$
	Ultrasonic treatment for 45 min
	$\downarrow$
	Cool down, weigh it again, and use 50% methanol to make up the lost weight
	$\downarrow$
	After shaking, filter with a 0.22 μm hydrophilic PTFE filter and analyze by HPLC

<System> Primaide 1110 Pump, 1210 Autosampler, 1310 Column oven, 1410 UV detector

NOTE: These data are an example of measurement; the individual vales can not be guaranteed.

### Hitachi High-Tech Science Corporation