

## Determination of Water Content in Aldehydes and Ketons

This application method is applicable for the determination of water in Aldehydes and Ketons. Aldehydes and Ketons reacts with the Karl Fischer reagent and subsequently generates water. To avoid this reaction, the reagent of special formulation is available.

<b>Titration</b>	AQV-2000/AQV-200 (Volumetric titrator) AQ-2000/AQ-200 (Coulometric titrator)	
<b>KF Oven</b>	N/A	
<b>Reagents for Volumetric Titrator</b>	Titrant	HYDRANAL <sup>®</sup> Composite 5K or equivalent
	Titration medium	HYDRANAL <sup>®</sup> KetoSolver or equivalent
<b>Reagents for Coulometric Titrator</b>	Anode solution	HYDRANAL <sup>®</sup> Coulomat AK or equivalent
	Cathode solution	HYDRANAL <sup>®</sup> Coulomat CG-K or equivalent

### PARAMETER SETTINGS for the Titrator

Volumetric Titrator		Coulometric Titrator	
End Mode	4	CAL Mode	0
CAL Mode	0	Interval	20 sec
Interval	20 sec	Current	Fast
MIN Feed	0.01 ml	S-timer	0 min
S-timer	0 min	T-timer	0 min
T-timer	0 min		

### PARAMETER SETTINGS for Oven

<b>Temperature</b>	-
<b>Carrier flow rate</b>	-

### PROCEDURE:

1. Fill a 5 ml glass syringe with the sample. Before filling, rinse the syringe with the sample 2 – 3 times.
2. Titrant the titration medium to zero. (Blanking)
3. Press SAMPLE key.
4. Inject approx. 1 ml of the sample into titration cell by piercing the rubber septum.
5. Press TITRATION key.
6. Press S.SIZE key and enter the sample size.

**OPERATING NOTES:**

1. The background level may increase after several injection of the sample. In this case, exchange the anode solution with the fresh one.
2. Stop the end of the syringe needle with a small silicon block. This prevents the loss of the sample caused by evaporation.
3. After injection of the sample, pull slightly the syringe plunger to suck the sample drops sticking on the end of syringe needle.
4. Injected sample size is calculated by subtracting AFTER injection syringe weight from BEFORE injection syringe weight.
5. Pierce different part of the rubber septum for each injection. If the same part is pierced for many times, atmospheric moisture will get into the titration cell, which will result in a high background or taking a long time for blanking.
6. To get the most accurate result, the titrant should be consumed around 5 ml. To make titrant consumption into this range, calculate the optimum sample size by this equation ;

$$\text{SIZE (g)} = 5 \times F / (C \times 10)$$

Where ;

F : Factor of titrant. ex. Factor of Composite 5 is 5 (mgH<sub>2</sub>O/ml).

C : Sample moisture in %

For example, sample moisture is 1%, and Factor is 5, then optimum sample size is,

$$5 \times 5 / (1 \times 10) = 2.5\text{g}$$

ⓘ When using this homepage, we recommend Chrome browser. Chrome can translate web

pages.

ABOUT US

PRODUCT

APPLICATIONS


FAQ

CONTACT

<b>HIRANUMA APPLICATION DATA</b>		Karl Fischer Titrator	Data No.	KF2	Jun.6. 2017
<b>Water contents</b>	<b>Alcohols</b>				

## 1. Abstract

Water content of Alcohols is determined by Karl Fischer coulometric titrator. In coulometric titration, iodine of Karl Fischer reagent is generated by electrolysis and generated iodine quantitatively reacts with water. Reaction formula is described below.



Alcohols do not interfere the Karl Fischer reaction and direct injection method could apply. Anode solution is selected in accordance with sample solubility. General-use anode solution contains methanol as solvent. Alcohols with long carbon chain have low solubility in methanol. In that case, use of anode solution for oil is appropriate. When fritless cell is used, cathode solution is not necessary.

## 2. Apparatus and Reagents

### (1) Apparatus

Titration	:	HIRANUMA Karl Fischer Coulometric titrator AQ-series or MOICO-A19
Electrolytic cell	:	Standard Cell Fritless Cell

### (2) Reagents

Anode solution	:	Hydranal coulomat AG (for general use, nonhalogenated)
Cathode solution	:	Hydranal coulomat CG

## 3. Procedure

- (1) Fill 100 mL of anode solution and one ampoule of cathode solution into the electrolytic cell as shown in Fig.3.1.
- (2) Start blanking to attain stable background.
- (3) Wash the syringe with sample.
- (4) Draw the sample into syringe and then weigh the syringe.
- (5) Inject sample from rubber septum of electrolytic cell as shown in Fig.3.2.
- (6) Start titration. Measurement parameter is shown in Table 4.1.
- (7) Weigh the syringe again and then set the difference of weight to sample size.

ⓘ When using this homepage, we recommend Chrome browser. Chrome can translate web pages.

ABOUT US

PRODUCT

APPLICATIONS

FAQ

CONTACT



Fig.3.1. Preparation of the reagents.




Fig.3.2. Injection of sample.

#### 4. Parameters and results

Table 4.1. Parameters

Condition File	
Cal Mode	o:Sample weight(net) X-(H <sub>2</sub> O-BLANK)/SIZE
Interval Time	20 sec
Current	SLOW
S-Timer	0 min
Blank Value	0 ug
Unit Mode	AUTO
Auto Interval	0 g
Minimum Count	5 ug
Back Ground	ON
Sample Size Input	Every Time
Cell Type	Standard/Fritless

Table 4.2. Results of water content measurement in alcohols

ⓘ When using this homepage, we recommend Chrome browser. Chrome can translate web pages.

				Sample (g)	Water (μg)	Water Content
ABOUT US	PRODUCT	APPLICATIONS		88	672.9	1124 ppm
			CG	0.8890	991.6	1115
				0.8710	973.6	1118
		Fritless	AG	0.9415	1047.4	1112 ppm
				0.8964	996.7	1112
				0.8739	972.3	1113
2-propanol	AQ	Standard	AG	1.5078	138.9	92.1 ppm
			CG	1.6015	145.7	91.0
				1.6574	152.6	92.1
		Fritless	AG	1.5623	142.7	91.3 ppm
				1.5496	139.4	90.0
				1.5832	142.9	90.3
2-methoxyethanol	AQ	Standard	AG	0.86583	156.4	180.6 ppm
			CG	1.36865	243.5	177.9
				0.83448	150.1	179.9
		Fritless	AG	1.0574	188.2	178.0 ppm
				0.95761	172.1	179.7
				1.02781	185.3	180.3
n-hexylalcohol	AQ	Standard	AG	0.87698	132.1	150.6 ppm
			CG	0.92031	138.0	149.9
				0.81105	124.0	152.9
		Fritless	AG	0.80510	122.9	152.7 ppm
				0.85949	131.1	152.5
				0.97448	148.7	152.6

## 5. Note

- (1) Use dried syringe and syringe vial in Fig.5.1, for preventive of contamination by atmospheric water.
- (2) Put appropriate anode solution in use according to the solubility of the sample. For example, Hydranal Coulomat AG-H and Oil are suitable for oils.

ⓘ When using this homepage, we recommend Chrome browser. Chrome can translate web pages.

[ABOUT US](#)

[PRODUCT](#)

[APPLICATIONS](#)

[FAQ](#)

[CONTACT](#)



Fig.5.1. Draw the sample from syringe vial.



[ABOUT US](#)

[PRODUCT](#)

[APPLICATIONS](#)

[FAQ](#)

[CONTACT](#)

[TERMS OF USE](#)

[PRIVACY POLICY](#)

