

ⓘ When using this homepage, we recommend Chrome browser. Chrome can translate web

pages.

ABOUT US

PRODUCT

APPLICATIONS

FAQ


CONTACT

<b>HIRANUMA APPLICATION DATA</b>		Automatic Titrator	Data No.	N2	Feb. 10,2021
<b>Cement Concrete</b>	<b>Determination of chloride ion in concrete</b>				

## 1. Abstract

This report introduces an example of the determination of chloride ion in concrete.

This measurement method is described in “Potentiometric titration method using chloride ion-selective electrode” of “JIS A 1154 Method of test for chloride ion content in hardened concrete”. Nitric acid is added to adjust the pH to 3 or less, and extract chloride ion with boiling. Cool to room temperature, filtrate to separate undissolved component and collect filtrated solution as sample. Chloride ion is determined by precipitation titration using a silver nitrate standard solution.



## 2. Configuration of instruments and Reagents

### (1) Configuration

Main unit : Automatic Titrator COM series

Electrode : Chloride ion-selective electrode CLi-o81 (Connect to IE-2)

Reference electrode MS-231Z (Connect to RE-2)

\*Remark

The general reference electrode (RE-201Z) cannot be used for this titration because KCl inner solution might come out to sample solution and it causes measurement error.

The inner electrodes of MS-231Z use mercury (I) sulfate. When these electrodes are disposed, please ask the specialized industrial waste disposal operator.

Filter paper : For quantitative analysis

### (2) Reagents

Titrant : 0.005 mol/L Silver nitrate standard solution

Additive : Diluted nitric acid

Prepared by mixing nitric acid (60 %) in a volume ratio of 1 to 6 in DI water

## 3. Measurement procedure

- (1) Take 10 g of sample into a 200 mL beaker and accurately weigh it.
- (2) Add 70 mL of diluted nitric acid slowly and stirrer it. Then measure the pH to confirm the pH is 3 or less.
- (3) Cover the beaker with watch glass and heat to boil for 5 minute. Then cool it to room temperature.
- (4) Filtrate the solution using filter paper with aspirating and wash the filter paper with water.
- (5) Collect filtrated solution into a volumetric flask to make it 200 mL with DI water.
- (6) Take the solution prepared at (5) with a 50 mL volumetric pipet into a 100 mL beaker.
- (7) Immerse electrodes and titrate with 0.005 mol/L silver nitrate standard solution.

ⓘ When using this homepage, we recommend Chrome browser. Chrome can translate web

pages.

[ABOUT US](#)


[PRODUCT](#)

[APPLICATIONS](#)

[FAQ](#)

[CONTACT](#)

Cndt No.	1	Constant No.	1	Mode No.	8
Method	Auto	Size	10.0314 g	Pre Int	0 sec
Buret No.	1	Blank	0 mL	Del K	5
Amp No.	2	Molarity	0.005 mol/L	Del Sens	0 mV
D. Unit	mV	Factor	1	Int Time	5 sec
S-Timer	5 sec	K	35.45	Int Sens	3 mV
C.P. mL	0 mL	L	50	Brst Speed	2
T Timer	0 sec	Unit	%	Pulse	40
D.P. mL	0.1 mL	Formula	$(D-B)*K*F*M/(S*10)*200/L$		
End Sens	30	Decimal Places	4		
Over mL	0.5 mL				
Max Vol.	20 mL				



Example of titration curve

#### Measurement results

Number of measurement	Sample size(g)	Titration volume(mL)	Chloride ion (%)
1		2.326	0.0164
2	10.0314	2.340	0.0165
3		2.296	0.0162
		Avg.	0.0164 %
		SD	0.0002 %
		RSD	0.93 %

#### 5.Note

· Indicator electrode

Chloride ion-selective electrode was used as an indicator electrode for this measurement.

In addition to the above electrode, a silver electrode coated with silver chloride (model: AG-311A) can also be used for this measurement. However, since the silver electrode coated with silver chloride deteriorates after long-term use, the electrode potential change near the end point becomes small and unclear. On the other hand, the chloride ion-selective electrode has the advantage that when the sensitivity deteriorates, the sensitivity can be easily restored by lightly polishing the sensitive membrane with a sandpaper.

Keywords: Chloride ion, Precipitation titration, Concrete, JIS A1154

 When using this homepage, we recommend Chrome browser. Chrome can translate web

pages.

[ABOUT US](#)

[PRODUCT](#)

[APPLICATIONS](#)

[FAQ](#)

[CONTACT](#)



[ABOUT US](#)

[PRODUCT](#)

[APPLICATIONS](#)

[FAQ](#)

[CONTACT](#)

[TERMS OF USE](#)

[PRIVACY POLICY](#)

