

Silver & Gold Pellets

Gold and Silver colored metal, non-invasive pellets, or disks, are applied to corresponding and hand micro-meridian points to provide a unique stimulation. Pellets can be used for both chronic and acute conditions. They are an effective, non invasive alternative to needles.

How Pellets Work

Pellets effect the body and acu-points in three ways:

1. **Bio-Electric Effect:**

The electrical properties of the metal interacts with the energy field of the acu-points. This causes a dispersing action which opens the blockage in the area. Impulses sent to the brain initiate a healing response in the body.

2. **Color Effect:**

The colors of gold and silver have unique electrical properties that affect the body energy. Either color can be used for pain syndromes.

Gold is used to tonify and silver to sedate the Five element points.

Energy moves from the Gold and Silver colored metal, non-invasive pellets, or disks, are applied to corresponding and hand micro-meridian points to provide a unique stimulation. Pellets can be used for both chronic and acute conditions. They are an effective, non invasive alternative to needles. meridian energy.

3. **Pressing Effect:**

Massaging or pressing the pellet into hand or fingers sends an electrical signal to the brain. This initiates the release of healing effects in the corresponding sites in the body.

Choosing Pellet Size

Small Silver or Gold pellets are used on:

- tender reflex points,
- acu-points and
- smaller areas such as the fingers.

Note: When treating corresponding points, Use all silver or all gold.

DO NOT MIX COLORS



Small Silver Color Pellets
1/8th in. Dia
100 / box



Small Gold Color Pellets
1/8th in. Dia
100 / box

Sedating and Tonifying Micro-Meridians

- **To Sedate:** Place the [silver](#) pellet on higher number and the [gold](#) pellet on the lower number.
- **To Tonify:** Place the [silver](#) on the lower number and the [gold](#) on the higher.
- **Location** Place the two pellets between the distal and the proximal finger joints

Note: Center the pellet on the point for best results.



Medium Gold Color Pellets
1/4th in. Dia
30 / box
Available in Silver & Gold

Medium Pellets

are used on large pain or reflex areas on the hands,
i.e: neck, stomach,
knee, spine, etc.

Large Pellets

are used on larger pain or reflex areas on the hands,
i.e: sacrum, hips, lower abdomen, chest, etc.



Large Gold Color Pellets
5/8th in. Dia
10 / box
Available in Silver & Gold

Procedure for Pellet Application

1. Wash your hands (see rationale below in notes)
2. Locate and mark the exact correspondence point(s)
3. Center the pellet precisely on the point
4. Press the pellet on the into the point
5. Re-check point location if # 4 does not produce a strong reflex response
6. Find the point with the strong reflex response
7. Secure the pellet with a secondary tape



Notes:

1. Oily skin prevents good pellet adhesion and electrical conductivity.
2. Greater accuracy is guaranteed with skin marking pen
3. Carefully center the pellet point on the acu-point .
4. Press the pellet into the point to elicit a reflex reaction to verify accuracy.

Reflex reactions may include:


- Giving a verbal response
- Sharp breath intake
- Jerk their hand away
- Grimacing

5. If there is no reflex reaction, stop & re-locate the precise point location.


Tape Pellets for Long Term Retention

Long Term Can Mean:

- Overnight retention
- 2 - 5 days
- To prevent the return of pain.



1. Leave the pellet on its original tape, then tape it on the acupoint.



2. Use a secondary tape to completely encircle the finger. Overlap to insure good retention.




3. A full circle wrap will keep the pellets in place during normal use of their hands.

Note: In our clinic and seminars, we routinely applied a secondary tape to insure the pellets remained in the correct location.

Taping Ideas for Low Adhesion Pellets

If your pellets have been stored long term or in a hot area, the adhesive may not adhere to the skin. The pellets will still produce the KHT effect. Therefore use a secondary tape over the original pellet tape to

1. Use a piece of tape long enough to encircle the finger with a small overlap.
2. Do not remove the pellet from the original tape. Place the pellet + tape in the center of the secondary tape
4. Center the protrusion of the pellet on the acu- point
5. Wrap the strip of tape around the finger with 1/4th inch overlap.



Notes:

1. If clients have moist or oily skin, clean the skin with an alcohol pad and dry.
2. Change pellets after 2-3 days to avoid skin issues.
3. Ask clients to avoid hand lotion before hand treatments

Tips for Storing Hy-Tape

Store the roll upright (Fig.1).

Hy-Tape adhesive is temperature sensitive, so if the roll is stored leaning or on an uneven surface, it can slide off the roll, (Fig.2).




Fig. 1: Correct Storage



Fig 2: Disaster Storage