

SARS-CoV-2 Ag Rapid Test kit

SARS-CoV-2 Ag Rapid Test Kit

Kat: 6973315390026

[PRODUKTNAME]

SARS-CoV-2 Ag Rapid Test Kit

[VERPACKEN]

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Hangzhou Jucheng Medical Products Co., Ltd
3rd floor, Building 1, No.8 Mingyuan Road, Gaohong Town,
Lin'an District, Hangzhou City, Zhejiang Province, China

CMC Medical Devices & Drugs S.L
C/Horacio Lengua N°18
CP 29006, Málaga-Spain

SARS-CoV-2 Ag Rapid Test Kit

Cat: 6973315390026

[PRODUCT NAME]

SARS-CoV-2 Ag Rapid Test Kit

[PACKING]

1 piece/bag, 5 pieces/box, 25 pieces/box.

[INTENDED USE]

This product is intended for the qualitative detection of novel coronavirus, or COVID-19, in Nasal swab samples. It aids in the diagnosis of infection with novel coronavirus. This product can be used for self-test.

[SUMMARY]

The novel coronaviruses (SARS-CoV-2) belong to the β genus. COVID-19 is an acute respiratory infectious disease. People are generally susceptible to infection. Currently, the patients infected by the novel coronavirus are the main source of infection; asymptomatic infected people can also be an infectious source. Based on the current epidemiological investigation, the incubation period is 1 to 14 days, particularly 3 to 7 days. The main symptoms include fever, fatigue, and dry cough. Nasal congestion, runny nose, sore throat, myalgia, and diarrhea are also found in some cases.

[PRINCIPLE]

The **SARS-CoV-2 Ag Rapid Test Kit** is an immunochromatographic membrane assay that uses highly sensitive monoclonal antibodies to detect nucleocapsid protein from SARS-CoV-2 in Nasal swab samples. The test strip is composed of the following parts: namely **sample pad**, **reagent pad**, **reaction membrane**, and **absorbing pad**. The **reagent pad** contains the colloidal gold conjugated with the monoclonal antibodies against the nucleocapsid protein of SARS-CoV-2; the **reaction membrane** contains the secondary antibodies for nucleocapsid protein of SARS-CoV-2. The whole strip is fixed inside a plastic device. When the sample is added into the sample well, conjugates dried in the **reagent pad** are dissolved and migrate along with the sample. If the SARS-CoV-2 antigen presents in the sample, a complex formed between the anti-SARS-2 conjugate and the virus will be captured by the specific anti-SARS-2 monoclonal antibodies coated on the test line region (T). The absence of the T line suggests a negative result. To serve as a procedural control, a red line will always appear in the control line region (C), indicating that proper volume of sample has been added, and membrane wicking has occurred.

[COMPOSITION]

- Disposable test device
- Disposable sample extraction tube
- Specimen collection swab (nasal swab)

[STORAGE AND STABILITY]

- Store as packaged in the hermetically-sealed bag at the temperature (2-30°C or 38-86°F) and avoid direct sunshine. The kit is stable within the expiration date printed on the labeling.
- Once the sealed bag is opened, the test should be used within one hour. Prolonged exposure to hot and humid environments will cause product deterioration.
- The lot number and the expiration date are printed on each sealed bag.

[TEST PROCEDURE]

Allow the test device and specimens to equilibrate to room temperature (15-30°C or 59-86°F) prior to testing.

NASAL SWAB SAMPLE:


- Use the nasal swab provided in the kit to wipe the anterior part of the nasal cavity.
- Open the cap of the disposable sample extraction tube, cut off the tip of the cotton swab inserted into the test tube and close the extraction tube.
- Shake the disposable sample extraction tube to dissolve the sample on the swab.
- Take the test cassette from the packaging bag and place it on a table.
- Cut off the protrusion of the collection tube.
- Add 3 drops of the sample into the sample hole vertically and wait 15 minutes for the appearance of the T line.

[INTERPRETATION OF RESULTS (AFTER 15 MINUTES)]

Positive (+): Both of T and C lines appear in 3-15 minutes.

Negative (-): C line appears while no T line appears in 15 minutes after the sample is added.

Invalid: If the C line does not appear, indicating that the test result is invalid, and you should retest specimen with another test device.

[NOTES]

- SARS-CoV-2 Ag Rapid Test Kit** is only applicable to Nasal swabs samples. Blood, serum, plasma, urine, and other samples may cause abnormal results if placed in device sample hole. If any sample is tested positive, please see your local healthcare authority for further clinical diagnosis and reporting of results.
- Please make sure that 3 drops of the sample are added to the sample hole. Too much or too little sample may cause deviations in results.
- For positive judgment, it can be confirmed as soon as both T and C lines have appeared, typically 3-15 minutes after the sample has been added. For negative judgment, please wait for 15 minutes after the sample has been added. The C line will appear alone while no T line will be present.
- The test device is a disposable product and will contain biohazards after use. Please properly dispose of the test devices, specimens, and all collection materials after use.
- The product must be used prior to the expiration date on product labeling. If part of the test membrane containing the reagents is out of the test window, or more than 2 mm of filter paper or latex pad is exposed in the test window, do not use it because the test results will be invalid. Use a new test kit instead.

PERFORMANCE CHARACTERISTICS:

1. Clinical Performance

A clinical evaluation was carried out to confirm that the sensitivity and specificity of the **SARS-CoV-2 Ag Rapid Test Kit** for SARS-CoV-2, and the results of SARS-CoV-2 Ag Rapid Test Kit were compared with those of RT-PCR. The results are as follows summarized:

Nasal Swab Sample		RT-PCR		Total
		Positive	Negative	
Test reagent	Positive	116	0	116
	Negative	5	105	110
Total		121	105	226

Nasal swab samples: The **SARS-CoV-2 Ag Rapid Test Kit** showed 95.9% sensitivity and 100% specificity in nasal swab samples.

Clinical sensitivity (%) = $[116 / (116 + 5)] \times 100\% = 95.9\%$

Clinical specificity (%) = $[105 / (0 + 105)] \times 100\% = 100\%$

Total agreement rate (%) = $[(116 + 105) / (116 + 5 + 0 + 105)] \times 100\% = 97.8\%$

3. Limit of Detection (LoD)

SARS-COV-2 nucleocapsid protein expressed in vitro and National Standard Reference sample of SARS-CoV-2 were used for Limit of Detection (LoD) tests. The LoD of the **SARS-CoV-2 Ag Rapid Test Kit** is 0.5 pg/mL SARS-COV-2 nucleocapsid protein. The LoD of the **SARS-CoV-2 Ag Rapid Test Kit** is 1×10^4 TCID50/mL SARS-COV-2.

Nasal Swab Sample	
10 pg/mL	30/30 (100%)
5 pg/mL	30/30 (100%)
1 pg/mL	30/30 (100%)
0.5 pg/mL	29/30 (96.7%)
0.1 pg/mL	7/30 (23.3%)
0.05 pg/mL	0/30 (0%)
0 pg/mL	0/30 (0%)

Nasal Swab Sample	
1×10^4 TCID50/mL	30/30 (100%)
1×10^3 TCID50/mL	30/30 (100%)
1×10^2 TCID50/mL	28/30 (93.3%)
1×10^1 TCID50/mL	6/30 (20%)
1×10^0 TCID50/mL	0/30 (0%)
0 TCID50/mL	0/30 (0%)

3. Cross-reactivity:

The cross-reactivity with the following organisms and viruses was examined. The following substances will not produce false positive or false negative reactions when tested with the **SARS-CoV-2 Ag Rapid Test Kit** for the SARS-CoV-2.

Organism	Concentration (TCID50/mL)	Organism	Concentration (TCID50/mL)
HKU1	1.5×10^5	EnterovirusD	4×10^5
Oc43	1.5×10^5	Epstein-Barr virus	2.5×10^5
NI63	1.5×10^5	Measles virus	3×10^5
229E	1.5×10^5	Human cytomegalovirus	3×10^5
MERS	1.5×10^5	Rotavirus	5×10^5
Influenza A H1N1	3×10^5	Norovirus	5×10^5
Seasonal Influenza H1N1	2×10^5	Mumps virus	5×10^5
Influenza A H3N2	3×10^5	Varicella-zoster virus	5×10^5
Influenza A H5N1	3×10^5	Human Metapneumovirus (hMPV)	4×10^5
Influenza A H7N9	3×10^5	Parainfluenza virus 1	4×10^5
Influenza B	5×10^5	Parainfluenza virus 2	2.5×10^5
Syncytial virus	4×10^5	Parainfluenza virus 3	3×10^5
Rhinovirus A	2.5×10^5	Parainfluenza virus 4	3×10^5
Rhinovirus B	2.5×10^5	Respiratory syncytial virus	3.5×10^5
Rhinovirus C	2.5×10^5	Haemophilus influenzae	5×10^5
Adenovirus type 1	5×10^5	Mycoplasma pneumoniae	6×10^4 cells/mL
Adenovirus type 2	5×10^5	Chlamydia pneumoniae	4.5×10^4 cells/mL
Adenovirus type 3	5×10^5	Legionella pneumophila	6×10^4 cells/mL
Adenovirus type 4	3.5×10^5	Staphylococcus aureus	6×10^4 cells/mL
Adenovirus 5	5×10^5	Streptococcus pneumoniae	5×10^4 cells/mL
Adenovirus type 7	3.5×10^5	Streptococcus pyogenes	5×10^4 cells/mL
Adenovirus 55	4×10^5	Candida albicans	5×10^4 cells/mL
Enterovirus A	4×10^5	Pooled human sampling site wash - representative of normal I respiratory microbial flora	4.5×10^4 cells/mL
Enterovirus B	4×10^5	Bordetella pertussis	4.5×10^4 cells/mL
Enterovirus C	4×10^5		

4. Endogenous/exogenous material interference test

The following substances, which occur naturally in breath samples or which can be artificially introduced into the airways, were evaluated as listed below. The **SARS-CoV-2 Ag Rapid Test Kit** does not report false positive or false negative.

Substance	Substance	Substance
Purified Mucin	alpha-interferon	Tobramycin
Bilirubin	Zanamivir	Histamine hydrochloride
Blood lipids	Ribavirin	Benfurin
Hemoglobin	Oseltamivir	Oxymetazoline
Rheumatoid factor	Paramivir	Sodium chloride
Antinuclear antibody	Lopinavir	Beclomethasone
Antimitochondrial antibody	Ritonavir	Dexamethasone
HAMA	Abidol	Flunisolone
Total IgG	Levofloxacin	Triamcinolone
Total IgM	Azithromycin	Budesonide
Hematocrit	Ceftriaxone	Momisson
Meropenem	Fluticasone	

5. Hook effect

The hook effect refers to the false-negative phenomenon caused by the incorrect ratio of antigen to antibody. Antibody excess is called proband effect.

For **SARS-CoV-2 Ag Rapid Test Kit**, even if the concentration of SARS-COV-2 nucleocapsid protein reaches $200 \mu\text{g/mL}$, the **SARS-CoV-2 Ag Rapid Test Kit** still has no hook effect.

[INDEX OF SYMBOLS]

	Consult instructions for use		Tests per kit		Authorized Representative
	For in vitro diagnostic use only		Use by		Do not reuse
	Store between 2-30°C		Lot Number		Catalogue number

Hangzhou Jucheng Medical Products Co., Ltd
3rd floor, Building 1, No.8 Mingyuan Road, Gaohong Town,
Lin'an District, Hangzhou City, Zhejiang Province, China

CMC Medical Devices & Drugs S.L
C/Horacio Lengua N°18
CP 29006, Málaga-Spain

EC REP CERTIFICATE

CMC MEDICAL DEVICES & DRUGS SL NO. CMC/CE/2020/07012021.20

CONFIRMED THAT CMC MEDICAL DEVICES & DRUGS S.L. Is the European Authorized Representative of

Hangzhou Jucheng Medical Products Co., Ltd
3rd floor, building 1, No.8 Mingyuan Road, Gaohong Town, Lin'an District,
Hangzhou City, Zhejiang Province

The certificate remains valid until the expiration agreement of EC REP, manufacturing conditions, the quality system or relevant legislation are changed. The validity is conditioned by positive results of periodic surveillance audits.

The product liability rests with the manufacturer in accordance with applicable directive and standard, after fulfilling of the relevant EU legislation requirements, the manufacturer shall affix relevant CE marking to all above mentioned models of the medical device.

Complies with the applicable essential requirements of the council directive 98/79/EEC in vitro diagnostics as amended.

The products in Annex I was registered in Spanish MOH with number **RPS/28/2021**


Issued on: 07/01/2021

Valid until: 06/01/2022

Authorized Signatory
CMC Medical Devices & Drugs SL

EC REP CERTIFICATE

ANNEX I Medical Device Products

SARS-CoV-2 Ag Rapid Test Kit

COVID-19 Neutralizing Antibody Detection Kit

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**THE SPANISH AGENCY OF MEDICINES AND
MEDICAL DEVICES (AEMPS)**

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CERTIFIES: On the basis of the declarations and/or
certificates provided:

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The company CMC MEDICAL DEVICES & DRUGS S.L., based in Calle HORACIO LENGU, nº 18, 29006 - MÁLAGA SPAIN, is the Authorised Representative in the European Union of the manufacturer Hangzhou Jucheng Medical Products Co., Ltd, based in 3rd floor, building 1, No.8 Mingyuan Road, Gaohong Town, Lin'an District, Hangzhou City, Zhejiang Province, PEOPLES REPUBLIC OF CHINA, of the products listed in **Annex I**.

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The products listed in **Annex I**, have been affixed with the CE mark under the Directive 98/79/EC of the European Parliament and of the Council of 27 October 1998, on *in vitro* diagnostic medical devices and the Royal Decree 1662/2000 of 29 September, transposition to the Spanish legislation, and therefore they can be marketed in Spain and in the rest of the European Union countries, without export restrictions.

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This certificate is issued on the basis of the documentation provided by the company at the time of its issuance and it should not be taken as an AEMPS marketing authorization of the products.

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Signed and issued as evidence thereof to be effective
with the Health Authorities of **India**, in Madrid.

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This Apostille only certifies the authenticity of the signature and the capacity of the person who has signed the public document, and, where appropriate, the identity of the seal or stamp which the public document bears.

This Apostille does not certify the content of the document for which it was issued.

To verify the issuance of this Apostille, see "<https://sede.mjusticia.gob.es/eregister>"

Verification code of the Apostille (*): AD:yeNS-cBZN-egMx-fNrU

This document has been electronically signed in accordance with the provisions of Articles 42 and 43 of Law 40/2015 of October 1st, of Legal Regime of the Public Sector.

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Producto / Product Referencia y modelo / Reference and model	Empresa(s) / Company(s) Domicilio social / Registered office
1 - SARS-CoV-2 Ag Rapid Test Kit (Ref: G-C-0104) (Mod: 1 tests/box)	Hangzhou Jucheng Medical Products Co., Ltd 3rd floor, building 1, No.8 Mingyuan Road, Gaohong Town, Lin'an District, Hangzhou City, Zhejiang Province, REPÚBLICA POPULAR CHINA / PEOPLES REPUBLIC OF CHINA Representante autorizado / Authorised representative CMC Medical Devices & Drugs S.L HORACIO LENGU 18, 29006 - ESPAÑA / SPAIN
2 - SARS-CoV-2 Ag Rapid Test Kit (Ref: G-C-0104) (Mod: 25 tests/box)	Hangzhou Jucheng Medical Products Co., Ltd 3rd floor, building 1, No.8 Mingyuan Road, Gaohong Town, Lin'an District, Hangzhou City, Zhejiang Province, REPÚBLICA POPULAR CHINA / PEOPLES REPUBLIC OF CHINA Representante autorizado / Authorised representative CMC Medical Devices & Drugs S.L HORACIO LENGU 18, 29006 - ESPAÑA / SPAIN
3 - SARS-CoV-2 Ag Rapid Test Kit (Ref: G-C-0104) (Mod: 5 tests/box)	Hangzhou Jucheng Medical Products Co., Ltd 3rd floor, building 1, No.8 Mingyuan Road, Gaohong Town, Lin'an District, Hangzhou City, Zhejiang Province, REPÚBLICA POPULAR CHINA / PEOPLES REPUBLIC OF CHINA Representante autorizado / Authorised representative CMC Medical Devices & Drugs S.L HORACIO LENGU 18, 29006 - ESPAÑA / SPAIN

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Producto / Product Referencia y modelo / Reference and model	Empresa(s) / Company(s) Domicilio social / Registered office
1 - SARS-CoV-2 Ag Rapid Test Kit (Ref: G-C-0104) (Mod: 1 tests/box)	Hangzhou Jucheng Medical Products Co., Ltd 3rd floor, building 1, No.8 Mingyuan Road, Gaohong Town, Lin'an District, Hangzhou City, Zhejiang Province, REPÚBLICA POPULAR CHINA / PEOPLES REPUBLIC OF CHINA Representante autorizado / Authorised representative CMC Medical Devices & Drugs S.L HORACIO LENGÓ 18, 29006 - ESPAÑA / SPAIN
2 - SARS-CoV-2 Ag Rapid Test Kit (Ref: G-C-0104) (Mod: 25 tests/box)	Hangzhou Jucheng Medical Products Co., Ltd 3rd floor, building 1, No.8 Mingyuan Road, Gaohong Town, Lin'an District, Hangzhou City, Zhejiang Province, REPÚBLICA POPULAR CHINA / PEOPLES REPUBLIC OF CHINA Representante autorizado / Authorised representative CMC Medical Devices & Drugs S.L HORACIO LENGÓ 18, 29006 - ESPAÑA / SPAIN
3 - SARS-CoV-2 Ag Rapid Test Kit (Ref: G-C-0104) (Mod: 5 tests/box)	Hangzhou Jucheng Medical Products Co., Ltd 3rd floor, building 1, No.8 Mingyuan Road, Gaohong Town, Lin'an District, Hangzhou City, Zhejiang Province, REPÚBLICA POPULAR CHINA / PEOPLES REPUBLIC OF CHINA Representante autorizado / Authorised representative CMC Medical Devices & Drugs S.L HORACIO LENGÓ 18, 29006 - ESPAÑA / SPAIN

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CERTIFIES: On the basis of the declarations and/or certificates provided:

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Limit of Detection Study Report

Study Summary

This study tested the limit of detection (LOD) of the SARS-CoV-2 Ag Rapid Test Kit on testing SARS-COV-2 nucleocapsid protein.

1.Purpose

This study was to find out the LOD of the SARS-CoV-2 Ag Rapid Test Kit using SARS-COV-2 nucleocapsid protein expressed in vitro and National Standard Reference sample of SARS-CoV-2.

The eligibility standard of the product is: When the SARS-CoV-2 Ag rapid detection kit is used to detect a pre-diluted low-concentration positive sample, and the positive rate is still $\geq 90\%$, it can be determined as the LoD value

2.Reference

The study was performed according to the technical specification of the SARS-CoV-2 Ag Rapid Test Kit.

4.Method

Spike a Sputum/oropharyngeal saliva samples into saline water, respectively. Prepare the supernatant for subsequent use. Spiked a group of serially diluted SARS-COV-2 nucleocapsid protein (500 pg/mL, 100 pg/mL, 50 pg/mL, 10 pg/mL, 0.5 pg/mL, 0 pg/mL) in supernatant described above, respectively. Each supernatant is tested 30 times.

And then, spike a Sputum/oropharyngeal saliva samples into saline water, respectively. Prepare the supernatant for subsequent use. Spiked a group of serially diluted National Standard Reference sample of SARS-CoV-2 (10^6 TCID₅₀/mL, 1×10^5 TCID₅₀/mL, 1×10^4 TCID₅₀/mL, 1×10^3 TCID₅₀/mL, 1×10^2 TCID₅₀/mL, 0TCID₅₀/mL) in supernatant described above, respectively. Each supernatant are tested 30 times.

6.Result

	(Sputum/oropharyngeal saliva) Sample
500 pg/mL	30/30 (100%)
100 pg/mL	30/30 (100%)
50 pg/mL	30/30 (100%)
10 pg/mL	28/30 (93.3%)
0.5 pg/mL	6/30 (20%)
0 pg/mL	0/30 (0%)


	(Sputum/oropharyngeal saliva) Sample
1×10^6 TCID ₅₀ /mL	30/30 (100%)
1×10^5 TCID ₅₀ /mL	30/30 (100%)
1×10^4 TCID ₅₀ /mL	30/30 (100%)
1×10^3 TCID ₅₀ /mL	29/30 (96.7%)
1×10^2 TCID ₅₀ /mL	6/30 (20%)
0 TCID ₅₀ /mL	0/30 (0%)

Test results presented above shows that the LOD of the SARS-CoV-2 Ag Rapid Test Kit is 10 pg/mL SARS-COV-2 nucleocapsid protein. The LOD of the SARS-CoV-2 Ag Rapid Test Kit is 1×10^3 TCID₅₀/mL SARS-COV-2.

7. Conclusion

The SARS-CoV-2 Ag rapid detection kit is superior to the requirements in the product specifications, supporting its feasibility in clinical applications.

Hangzhou Jucheng Medical Products Co., Ltd
Cheng Yinkai

For and on behalf of
HANGZHOU JUCHENG MEDICAL PRODUCTS CO., LTD

Authorized Signature(s)