Striving for excellence in offering
global healthcare solutions

SARS-CoV-2 Ag Rapid Test kit

A rapid, highly reliable and affordable kit, providing an aid in early diagnosis of individuals who are suspected of COVID-19 by their healthcare provider and who are asymptomatic.

Hangzhou Jucheng Medical Products Co., Ltd 3rd floor, building 1, No.8 Mingyuan Road, Ga Lin, an District, Hangzhou City, Zhejiang Province

• Tel: +86-571-63810290

• Web: www.juchengmedical.com

• Email: sales@jcmedical.com

SARS-CoV-2 Antigen Rapid Test

This product is suitable for the qualitative detection of novel coronavirus, or COVID-19, in Saliva. It aids in the diagnosis of infection with novel coronavirus. The novel coronaviruses (SARS-CoV-2) belong to the β genus. COVID-19 is an acute respiratory infectious disease. People are generally susceptible to infection. Currently, the patients infected by the novel coronavirus are the main source of infection; asymptomatic infected people can also be an infectious source. Based on the current epidemiological investigation, the incubation period is 1 to 14 days, particularly 3 to 7 days. The main symptoms include fever, fatigue, and dry cough. Nasal congestion, runny nose, sore throat, myalgia, and diarrhea are also found in some cases.

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Nasal Swab Sample		RT-PCR		Total
		Positive	Negative	
Test reagent	Positive	116	0	116
	Negative	5	105	110
Total		121	105	226

Nasal swab samples: The SARS-CoV-2 Ag Rapid Test Kit showed 95.9% sensitivity and 100% specificity in nasal swab samples.

Clinical sensitivity (%) = $[116 / (116 + 5)] \times 100\% = 95.9\%$

Clinical specificity (%) = $[105 / (0 + 105)] \times 100\% = 100\%$

Total agreement rate (%) = $[(116 + 105) / (116 + 5 + 0 + 105)] \times 100\% = 97.8\%$

Limit of Detection (LoD)

SARS-COV-2 nucleocapsid protein expressed in vitro and National Standard Reference sample of SARS-CoV-2 were used for Limit of Detection (LoD) tests. The LOD of the SARS-CoV-2 Ag Rapid Test Kit is 10 pg/mL SARS-COV-2 nucleocapsid protein. The LOD of the SARS-CoV-2 Ag Rapid Test Kit is 1×10^3 TCID₅₀/mL SARS-COV-2.

N-protein	Saliva	National Standard Reference sample	Saliva
500 pg/mL	30/30 (100%)	1×10^1 TCID ₅₀ /mL	30/30 (100%)
100 pg/mL	30/30 (100%)	1×10^2 TCID ₅₀ /mL	30/30 (100%)
50 pg/mL	30/30 (100%)	1×10^3 TCID ₅₀ /mL	30/30 (100%)
10 pg/mL	28/30 (93.3%)	1×10^4 TCID ₅₀ /mL	29/30 (96.7%)
0.5pg/mL	6/30 (20%)	1×10^5 TCID ₅₀ /mL	6/30 (20%)
0 pg/mL	0/30 (0%)	0 TCID ₅₀ /mL	0/30 (0%)


Recognition performance for mutant viruses

Spike a healthy saliva sample into saline water, respectively. Prepare the supernatant for subsequent use. Spiked different kind of National Standard Reference sample of SARS-CoV-2 mutant virus (1×10^3 TCID₅₀/mL). According to the test results, The detection performance of SARS-CoV-2 Ag Rapid Test Kit is suitable for a variety of SARS-CoV-2 mutant virus strains.

	Saliva		Saliva
B.1.618	50/50 (100%)	B.1.1.7	50/50 (100%)
B.1.617.1	50/50 (100%)	P.1	50/50 (100%)
B.1.617.2	50/50 (100%)	D614G	50/50 (100%)
B.1.1.351	50/50 (100%)	501Y.V2	50/50 (100%)

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1. Cough deeply and collect saliva or mucus which from the deep throat.
2. Continue to gargle with 1 tablespoon of water to further enrich the saliva sample.
3. Put the tampon of the test card into the mouth, keep the horizontal state of the test card, and wait for 1-1.5min.
4. When the wet liquid arrives at the top of the observation window, remove the test card and close the cover. Place it flat on the desktop and wait for 0-15min.
5. Positive results can be judged immediately if it appears no need to spend full 15 minutes.

Ordering Information

Product Name	Catalog NO.	Format	Specimen	Package
SARS-CoV-2 Antigen Rapid Test				25Tests/Kit

EC DECLARATION OF CONFORMITY

Name and address of the manufacturer: **Hangzhou Jucheng Medical Products Co.,Ltd**
3rd floor,building 1,No.8 Mingyuan Road,Gaohong Town,
Lin, an District, Hangzhou City, Zhejiang Province, China

We declare under our sole responsibility that

the medical device: **SARS-CoV-2 Ag Rapid Test Kit**

of class: **Other**
according to article 9 of directive 98/79/EC

meets the provisions of the directive 98/79/EC and its transpositions in national laws which apply to it. The declaration is valid in connection with the "final inspection report" of the device.

Conformity assessment procedure: **Directive 98/79/EC Annex III, excluding Section 6**

Applicable standards: **EN ISO 13485:2016** **EN ISO 15223-1:2016**
EN ISO 18113-1:2011 **EN ISO 18113-2:2011**
EN ISO 23640:2015 **EN ISO 14971:2019**
EN 13641:2002 **EN13612:2002/AC:2002**
EN 13975:2003

Name and address of the authorized representative: **CMC Medical Devices & Drugs S.L**
C/Horacio Lengo N° 18
CP 29006, Málaga-Spain
Tel: +34951214054
Fax: +34952330100

For and on behalf of
HANGZHOU JUCHENG MEDICAL PRODUCTS CO.,LTD

Hangzhou, May 6, 2021.

Place, date

Authorized Signature(s)
Cheng Yinkai Legal Representative

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Kontakt:

E-Mail: sarscov2ivd@pei.de

SARS-CoV-2 Antigen Rapid Detection Kit (Colloidal Gold Method)	Jinan Babio Biotechnology Co., Ltd.
GLINE-2019-nCoV Ag	Shenzhen YHLO Biotech Co. Ltd.
SARS-CoV-2 Rapid Antigen Test (Colloidal Gold Method)	Sansure Biotech Inc.
SARS-CoV-2 Antigen Schnelltestkassette	Hangzhou Sejoy Electronics & Instruments Co.,Ltd.
SARS-CoV-2-Antigen-Schnelltest	Hangzhou Careomedic Tech Co., Ltd.
AS-check COVID-19 Antigen Schnelltest	Asterion Otel Insaat Bilisim Medikal Maden Tic.Ltd.Sti.
Kanzone COVID-19 Antigen Rapid Test	Weihai Kangzhou Biotechnology Engineering Co., Ltd.
COVID-19 Antigen Test Kit	Surge Medical Inc.
GENEDIAN COVID-19 Antigen Test Cassette	Hangzhou DIAN Biotechnology Co., Ltd.
Novel Coronavirus(COVID-19) Antigen Test Kit (Colloidal Gold)	Changzhou Biowin Pharmaceutical Co.,Ltd.
SARS-CoV-2 Antigen Rapid Test (Lateral Flow Assay)	Shenzhen Kang Sheng Bao Bio- Technology Co.,Ltd.
Rapidan Tester COVID-19 Ag Test	Türklab Tibbi Malzemeler San. ve Tic. A.S.
ISIA SARS-CoV-2 Antigen Rapid Test Kit	Chongqing ISIA BIO-Technology Co.,Ltd
Check Up SARS-CoV-2 Nasal Antigen Rapid Test	Cesna Biyoteknoloji Arastirma Gelistirme Laboratuvar Sist. Ins.Müh.Dan.San.Tic-Ltd.Sti (Istanbul, Turkey)
SARS-CoV-2 Antigen Rapid Test (Colloidal Gold Method)	Biohit Healthcare (Hefei) Co., Ltd. (Hefei, China)
Novel Coronavirus (COVID-19) Antigen Detection Kit (Latex Immunochromatography)	Zhejiang Gene Science Co., Ltd. (Hangzhou Bay, China)
SARS-CoV-2 Antigen Rapid Test Cassette	Merlin Biomedical (Xiamen) Co., Ltd. (Xiamen, China)
SARS-CoV-2 Antigen Test Kit (Colloidal Gold)	Shenzhen Dymind Biotechnology Co. Ltd. (Shenzhen, China)
Novel Coronavirus(2019-nCoV) Antigen Rapid Test	Bioscience (Tianjin) Diagnostic Technology Co., Ltd (Tianjin, China)
StrongStep® SARS-COV2 Antigen Rapid Test	Nanjing Liming Bio-Products Co., Ltd. (Nanjing, China)
COVID-19 Rapid Test Cassette Antigen Test Kit	Lifecosc Biotech Limited (Shenzhen, China)
COVID-19 Antigen Detection Kit (Quantum Dots-Based Immunofluorescence Chromatography)	Shenzhen Kingfocus Biomedical Engineering Co., Ltd.
SARS-Cov-2 Antigen Rapid Detection Kit	Shenzhen CAS-Envision Medical Technology Co., Ltd. (Shenzhen, China)
JC+ SARS-CoV-2 Ag Rapid Test Kit	Hangzhou Jucheng Medical Products Co., Ltd
ZandCell COVID-19 Saliva Antigen-Test	ZandCell AB

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Hangzhou Jucheng Medical Products Co., Ltd.
3rd floor, building 1, No. 8 Mingyuan Road, Gaozhong Town, Lin'an District, Hangzhou City, Zhejiang Province, China

EC REP: CMC Medical Devices & Drugs S.L.
C/Horacio Lengua N° 18
CP 29006, Málaga-Spain
Tel: +34951214054
Fax: +34952330100
Email: info@cmcmmedicaldevices.com

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Cat: HY03
Version: 03

Specimens: Saliva
Effective Date: 2021-06

For professional and in vitro diagnostic use only.

【PRODUCT NAME】

SARS-CoV-2 Ag Rapid Test Kit

【PACKING】

Type I (1piece/bag , 1pieces/box ;)
Type II (1piece/bag , 5pieces/box ;)
Type III (1piece/bag , 10pieces/box)
Type IV (1piece/bag , 25pieces/box)
Type V (1piece/bag , 50pieces/box)

【INTENDED USE】

This product is suitable for the qualitative detection of novel coronavirus, or COVID-19, in Saliva. It aids in the diagnosis of infection with novel coronavirus.

【SUMMARY】

The novel coronaviruses (SARS-CoV-2) belong to the β genus. COVID-19 is an acute respiratory infectious disease. People are generally susceptible to infection. Currently, the patients infected by the novel coronavirus are the main source of infection; asymptomatic infected people can also be an infectious source. Based on the current epidemiological investigation, the incubation period is 1 to 14 days, particularly 3 to 7 days. The main symptoms include fever, fatigue, and dry cough. Nasal congestion, runny nose, sore throat, myalgia, and diarrhea are also found in some cases.

【PRINCIPLE】

The SARS-CoV-2 Ag Rapid Test Kit is an immunochromatographic membrane assay that uses highly sensitive monoclonal antibodies to detect nucleocapsid protein and Spiker protein from SARS-CoV-2 in Sputum/oropharyngeal saliva samples. The test strip is composed of the following parts: namely sample pad, reagent pad, reaction membrane, and absorbing pad. The reagent pad contains colloidal-gold conjugated with the monoclonal antibodies against the nucleocapsid protein of SARS-CoV-2 and colloidal-gold conjugated with the monoclonal antibodies against the nucleocapsid protein of SARS-CoV-2; the reaction membrane contains the secondary antibodies for nucleocapsid protein and S protein of SARS-CoV-2. The whole strip is fixed inside a plastic device. When the sample is added into the sample well, conjugates dried in the reagent pad are dissolved and migrate along with the sample. If SARS-CoV-2 antigen presents in the sample, a complex formed between the anti-SARS-2 conjugate and the virus will be captured by the specific anti-SARS-2 monoclonal antibody coated on the test line region (T). Absence of the T line suggests a negative result. To serve as a procedural control a red line will always appear in the control line region (C) indicating that proper volume of sample has been added and membrane wicking has occurred.

【COMPOSITION】

1. Disposable test device.


【STORAGE AND STABILITY】

1. Store as packaged in the hermetically-sealed bag at the temperature (2-30°C or 38-86°F) and avoid direct sunshine. The kit is stable within the expiration date printed on the labeling.
2. Once the sealed bag is opened, the test should be used within one hour. Prolonged exposure to hot and humid environments will cause product deterioration.
3. The lot number and the expiration date are printed on each sealed bag.

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
Allow the test device and specimens to equilibrate to room temperature (15-30°C or 59-86°F) prior to testing.

TEST METHOD 1 :

1. Gargle with 5~10 mL water for about 10 seconds, raise head and cough deeply and collect sputum/oropharyngeal saliva or mucus which from the deep throat at the same time.
2. Put the tampon of the test card into the mouth, gently bite the end of the plastic card housing, and keep the horizontal state of the test card, and wait for 1-3 minutes until the wet liquid reaches the top of the observation window.
3. When the wet liquid arrives at the top of the observation window, remove the test card and close the cover. Place it flat on the desktop and wait for 0-15min.

【INTERPRETATION OF RESULTS (WITHIN 15 MINUTES)】

Positive(+): Both of T and C lines appear within 15 minutes.

Negative(-): C line appears while no T line appeared after 15 minutes.

Invalid: If the C line does not appear, this indicates that the test result is invalid, and you should retest the specimen with another test device.

【NOTES】

1. SARS-CoV-2 Ag Rapid Test Kit is only applicable to Saliva samples. Blood, serum, plasma, urine, and other samples may cause abnormal results. If any sample tests positive, please see your local healthcare authority for further clinical diagnosis and reporting of results.
2. Make sure that the tampon is fully moistened. And wet liquid will arrive at the top of the observation window in 1-1.5 minutes.
3. Keep the horizontal state of the test card during testing.
4. Positive results can be judged immediately if C line and T line appear, and negative results need to spend full 15 minutes.
5. The test device is a disposable product and will contain biohazards after use. Please properly dispose of the test devices, specimens, and all collection materials after use.
6. Must use prior to the expiration date on product labeling.
7. If part of the test membrane containing the reagents is out of the test window, or more than 2 mm of filter paper or latex pad is exposed in the test window, do not use it because the test results will be invalid. Use a new test kit instead.

【PERFORMANCE CHARACTERISTICS】

1. Clinical Performance

A clinical evaluation was carried out to confirm that the sensitivity and specificity of the SARS-CoV-2 Ag Rapid Test Kit for SARS-CoV-2, compare results and RT-PCR. The results are as follows summarized:

(sputum/oropharyngeal saliva) Sample		RT-PCR		Total
		Positive	Negative	
Test reagent	Positive	98	0	98
	Negative	9	105	114
Total		107	105	212

(Sputum/oropharyngeal saliva) samples: The SARS-CoV-2 Ag Rapid Test Kit showed 91.6% sensitivity and 100% specificity in (sputum/oropharyngeal saliva) samples.

Clinical sensitivity (%) = $\frac{98}{98+9} \times 100\% = 91.6\%$, and the 95% confidence interval is 90.01% - 97.53%

Clinical specificity (%) = $\frac{105}{0+105} \times 100\% = 100\%$, and the 95% confidence interval is 97.61% - 100%

Total agreement rate (%) = $\frac{98+105}{98+9+0+105} \times 100\% = 95.8\%$

2. Limit of Detection (LoD)

SARS-CoV-2 nucleocapsid protein expressed in vitro and National Standard Reference sample of SARS-CoV-2 were used for Limit of Detection (LoD) tests. The LOD of the SARS-CoV-2 Ag Rapid Test Kit is 10 pg/mL SARS-CoV-2 nucleocapsid protein. The LOD of the SARS-CoV-2 Ag Rapid Test Kit is 1×10^3 TCID₅₀/mL SARS-CoV-2.

N-protein	Saliva	National Standard Reference sample	Saliva
500 pg/mL	30/30 (100%)	1×10^3 TCID ₅₀ /mL	30/30 (100%)
100 pg/mL	30/30 (100%)	1×10^2 TCID ₅₀ /mL	30/30 (100%)
50 pg/mL	30/30 (100%)	1×10^1 TCID ₅₀ /mL	30/30 (100%)
10 pg/mL	28/30 (93.3%)	1×10^0 TCID ₅₀ /mL	29/30 (96.7%)
0.5pg/mL	6/30 (20%)	1×10^{-1} TCID ₅₀ /mL	6/30 (20%)
0 pg/mL	0/30 (0%)	0 TCID ₅₀ /mL	0/30 (0%)

3. Recognition performance for mutant viruses:

Spiked different kind of National Standard Reference sample of SARS-CoV-2 mutant virus (1×10^3 TCID₅₀/mL) to saliva sample. According to the test results, the detection performance of SARS-CoV-2 Ag Rapid Test Kit is suitable for a variety of SARS-CoV-2 mutant virus strains

	Saliva		Saliva
B.1.618	50/50 (100%)	B.1.1.7	50/50 (100%)
B.1.617.1	50/50 (100%)	P.1	50/50 (100%)
B.1.617.2	50/50 (100%)	D614G	50/50 (100%)
B.1.1.351	50/50 (100%)	501Y.V2	50/50 (100%)

4. Cross-reactivity:

The cross-reactivity with the following organism and virus was examined. The following substances will not produce false positive or false negative reactions when tested with the SARS-CoV-2 Ag Rapid Test Kit for the SARS-CoV-2.

Organism	Concentration (TCID ₅₀ /mL)	Organism	Concentration (TCID ₅₀ /mL)
HKU1	1.5×10^4	Enterovirus D	4×10^3
Oc43	1.5×10^5	Epstein-Barr virus	2.5×10^3
NI63	1.5×10^5	Measles virus	3×10^3
229E	1.5×10^4	Human cytomegalovirus	3×10^3
MERS	1.5×10^3	Rotavirus	5×10^3
Influenza A H1N1	3×10^3	Norovirus	5×10^3
Seasonal Influenza H1N1	2×10^3	Mumps virus	5×10^3
Influenza A H3N2	3×10^3	Rhinovirus C	2.5×10^3
Influenza A H5N1	3×10^3	Adenovirus type 1	5×10^3
Influenza A H7N9	3×10^3	Adenovirus type 2	5×10^3
Influenza B	5×10^3	Adenovirus type 3	5×10^3
Syncytial virus	4×10^3	Adenovirus type 4	3.5×10^3
Rhinovirus A	2.5×10^3	Adenovirus 5	5×10^3
Rhinovirus B	2.5×10^3	Adenovirus type 7	3.5×10^3
Adenovirus 55	4×10^3	Enterovirus B	4×10^3
Enterovirus A	4×10^3	Enterovirus C	4×10^3
Varicella-zoster virus	5×10^3	Chlamydia pneumoniae	4.5×10^3 cells/mL
Human Metapneumovirus (hMPV)	4×10^3	Legionella pneumophila	6×10^3 cells/mL
Parainfluenza virus 1	4×10^3	Staphylococcus aureus	6×10^3 cells/mL
Parainfluenza virus 2	2.5×10^3	Streptococcus pneumoniae	5×10^3 cells/mL
Parainfluenza virus 3	3×10^3	Streptococcus pyogenes	5×10^3 cells/mL
Parainfluenza virus 4	3×10^3	Candida albicans	5×10^3 cells/mL
Respiratory syncytial virus	3.5×10^3	Pooled human sampling site wash	4.5×10^3 cells/mL
Haemophilus influenzae	5×10^3	Bordetella pertussis	4.5×10^3 cells/mL
Mycoplasma pneumoniae	6×10^3 cells/mL		

5. Endogenous/exogenous material interference test

The following substances, which occur naturally in breath samples or which can be artificially introduced into the airways, were evaluated as listed below. The SARS-CoV-2 Ag Rapid Test Kit does not report false positive or false negative.

Substance	Substance	Substance	Substance
Purified	Mucin Total Ig	Mritonavir	Oxymetazoline
Bilirubin	Hematocrit	Abidol	Sodium chloride
Blood lipids	Meropenem	Levofloxacin	Beclomethasone
Hemoglobin	alpha-interferon	Azithromycin	Dexamethasone
Rheumatoid factor	Zanamivir	Ceftriaxone	Flunisolone
Antinuclear antibody	Ribavirin	Fluticasone	Triamcinolone
Antimitochondrial antibody	Oseltamivir	Tobramycin	Budesonide
HAMA	Paramivir	Histamine hydrochloride	Momisson
Total IgG	Lopinavir	Benfurin	

6. Hook effect

The hook effect refers to the false-negative phenomenon caused by the incorrect ratio of antigen to antibody. For SARS-CoV-2 Ag Rapid Test Kit, even if the concentration of SARS-CoV-2 nucleocapsid protein reaches 200 µg/mL, the SARS-CoV-2 Ag Rapid Test Kit still has no hook effect.

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Hangzhou Jucheng Medical Products Co., Ltd
3rd floor, building 1, No. 8 Mingyuan Road, Gaozhong Town,
Lin, an District, Hangzhou City, Zhejiang Province, China

CMC Medical Devices & Drugs S.L
C/Horacio Lengo N° 18
CP 29006, Málaga-Spain
Tel: +34951214054
Fax: +34952330100
Email: info@cmcmedicaldevices.com

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Deutscher Vertreter Glückhouse International GmbH Neubrücker Straße 3, 55768 Hoppstädten-Weiersbach, Deutschland