



## Item 25-326 Used *MBraun™ PURGE 4 port Double Sided Glovebox*

This used *MBraun PURGE* glovebox is ideal for processes using sulfur, chlorine or halides compounds not compatible with the CU catalyst used in recirculating systems. The glovebox includes large & mini antechamber, touch screen controls, H2O sensor and foot switch. The system maintains an inert atmosphere via a constant purge flow. The flow rate is controlled via a manual flow valve. H2O removal can be further enhanced with the included internal VAC *Genesis* purifier. A vacuum pump and gloves are required but not included. 90 Day Warranty for Parts.

Price: **Make offer**



The system is equipped with a PLC operated pressure control system which starts automatically when the system is activated. Upper and lower working pressure parameters are set by the user, provide a comfortable range for the user to work in. Depending upon the application, the system can be operated in the positive or negative pressure state. Operational purge Inlet flow is controlled via the flow meter at a rate between 0-200 liters/min. The user can adjust the inlet flow as required. This purge mode has pressure control. If the pressure is above or below the upper or lower limits, purge will shutdown to regulate box pressure to back within limits. Foot pedals can also be used which will shut down the purge operation when gas or vacuum is called for.



**UsedGloveBox.com**  
SCIENTIFIC MARKET SOLUTIONS



Web: [www.usedglovebox.com](http://www.usedglovebox.com)

Email: [info@usedglovebox.com](mailto:info@usedglovebox.com)

Tel: 702-343-0297






Web: [www.usedglovebox.com](http://www.usedglovebox.com)

Email: [info@usedglovebox.com](mailto:info@usedglovebox.com)

Tel: 702-343-0297



**UsedGloveBox.com**  
SCIENTIFIC MARKET SOLUTIONS




Manual purge flow control valve with indicator and automated flow shut off valves

Web: [www.usedglovebox.com](http://www.usedglovebox.com)

Email: [info@usedglovebox.com](mailto:info@usedglovebox.com)

Tel: 702-343-0297







**UsedGloveBox.com**  
SCIENTIFIC MARKET SOLUTIONS



Web: [www.usedglovebox.com](http://www.usedglovebox.com)

Email: [info@usedglovebox.com](mailto:info@usedglovebox.com)





Tel: 702-343-0297



Top mounted purge exhaust plumbing with pneumatic shut off valve and H2O sensor



H2O sensor



Screen shots of purge touch screen controller

Web: [www.usedglovebox.com](http://www.usedglovebox.com)

Email: [info@usedglovebox.com](mailto:info@usedglovebox.com)

Tel: 702-343-0297




## VAC Genesis Purifier with Included pre-regenerated moisture cartridges



VAC Genesis Purifier blower will sit inside on glovebox floor and enhance H<sub>2</sub>O removal. This can help reduce gas usage and lower the required purge flow rate, especially when H<sub>2</sub>O is the main concern. Replacement cartridges are available from VAC.




## Main Switch





Glovebox leak check conducted with an *Inficon* H2 leak detector using 5% H2 95% N2 gas



Publication: INFICON AB - nina60e1-d (1304) - All information can be modified without prior notice

**Optional On-Site Installation, Start Up Assistance and Training** -We provide the option of on-site installation assistance. Our personnel have combined inert glovebox work experience of over 40 years. We can train your staff on the operation of the glovebox system. We also will provide long term technical support by email and phone at no charge.

**Web:** [www.usedglovebox.com](http://www.usedglovebox.com)

**Email:** [info@usedglovebox.com](mailto:info@usedglovebox.com)

**Tel:** 702-343-0297