



DISCLAIMER:


This content is intended for educational and/or entertainment purposes only. Many of the chemicals used in this research are hazardous and should be handled only by properly trained adults who are familiar with the physical/chemical properties and the techniques for handling/transferring them. We assume no liability for accidents or injuries to individuals or property caused by the use or misuse of our products or information.

Esterification of two cinnamic acids using Dimethylcarbonate

In order to investigate certain synthetic methods for use in natural product synthesis, we needed the methyl esters of two particular aryl-substituted cinnamic acids. Seeking a “green” method and avoiding Fisher esterification (which is incompatible with acid-sensitive substrates), dimethylcarbonate (DMC)¹ was chosen as an alternative to classical (toxic) methylating reagents such as iodomethane and dimethyl sulfate. We referred to a 2013 publication² by Gorin, et al., which used DMSO/catalytic K₂CO₃ in refluxing DMC.


Part 1. Synthesis of Methyl-3,4,5-trimethoxycinnamate

3,4,5-trimethoxycinnamic acid (a natural product itself) and its methyl ester are common precursors to numerous natural products contained in plants used in Ayurvedic medicine. One example is Piperlongumine (*aka* Piplartine), an alkaloid found in Long Pepper (*Piper longum*- which is known as “Pippali” in India) and in Black Pepper³ (*Piper Nigrum*- *aka* “Maricha”). Another example is Colchicine, found in *Colchicum luteum*, *aka* “Suranjana.”



Part 2. Synthesis of Methyl-3-(1,3-benzodioxol-5-yl)-2-propenoate

3,4-(methylenedioxy)cinnamic acid (and its methyl ester) are precursors to fragrances and a multitude of natural products used in homeopathic & Ayurvedic medicine. A well-known example is Piperine, the alkaloid of highest concentration in Black Pepper^{3,4}, which is known to enhance the absorption of supplements or drugs. Another example is Methysticin, a biologically active lactone found in the roots of the Kava (*Piper methysticum*) plant. We were mainly interested in the ester for synthetic method development and for use in the synthesis (Route D) of Berberine.



1. See last page page for list of abbreviations

2. “Catalytic Methyl Transfer from Dimethylcarbonate to Carboxylic Acids.” Yuan Ji, Jessica Sweeney, Jillian Zoglio & David J. Gorin* *J. Org. Chem.*, **2013**, *78*, 11606-11.

3. “Piperine-Type Amides: Review of the Chemical and Biological Characteristics.” *International Journal of Chemistry*; Vol. 5, No. 3; **2013**.

4. “Black Pepper, the “King of Spices”: Chemical composition to applications.” *Arabian Journal of Chemical and Environmental Research* Vol. 06 Issue 1, **2019**, 12–56.

Part 1a. Synthesis of Methyl-3,4,5-trimethoxycinnamate (22 g scale)

3,4,5-Trimethoxycinnamic acid (TMCA, MW 238.26 g/mol) **22 g, 92.3 mmol, 1 eq**

Potassium carbonate K₂CO₃ (MW 138.2 g/mol) **5.1 g, 36.9 mmol, 0.4 eq**

DMSO (MW 78.13 g/mol, d 1.1 g/mL) **92 g, 1.18 mol, 12.75 eq**

DMC (MW 90.08 g/mol, d 1.07 g/mL, BP 90 °C) **92 g, 1.02 mol, 11 eq**

NOTE: Product (TMCA ME) MW = 252.26 g/mol [20329-96-8] or [7560-49-8] Rep MPs: 91-92 °C, 99-100 °C, 108-110 °C

06/19/25			
TIME	EXSBT (°C)	Variac (%)	ACTION
10:25 am	N/A	N/A	Added DMSO followed by TMCA followed by DMC to 350 mL glass PV equipped with stir bar. Stirred. Once TMCA dissolved, added K ₂ CO ₃ . Flushed with N ₂ and sealed.
11 am	"	"	Turned on variac (15%)
12 pm	62	15	
1 pm	73	"	Vented
1:30 pm	75.5	"	↑d variac to 18%
2 pm	85	18	Took photo (2:22 pm)- See PHOTO A
2:30 pm	89	"	
2:45 pm	90	"	Vented
3:15 pm	91	"	
3:30 pm	92	"	
3:45 pm	93	"	
4 pm	94	"	Vented
4:15 pm	"	"	
4:30 pm	"	"	
5 pm	95	"	Vented, backed off variac to 17%
5:30 pm	93	17	
6 pm	"	"	
6:30 pm	"	"	Vented one last time and left for evening
			PHOTO A (06/19/25) 2:22 pm
			PHOTO B (06/20/25) 9:07 am
06/20/25			
8:30 am	91	17	Vented and then carried out TLC (5:1 Hex:EtOAc, PHOTO B : 1_TMCA, 2_TMCA ME from a previous 1 g scale rxn, 3_ALQ from this rxn (diluted ~1:1:1 with IpOAc:H ₂ O)
2:45 pm	92	"	Turned off heat and removed HM; continued stirring
5:15 pm	N/A	N/A	Poured RM into mixture of IpOAc (165 g) + DI water (165 g) in 1 L EFlask and stirred 5 min; Sealed with AL foil and let sit over weekend
06/23/25			
Separated layers & extracted AQ 2 more times (1 x 75 g & 1 x 50 g). Combined IpOAc layers & washed with DI water (2 x 150 g) then brine (1 x 125 g). EOCd in shallow pyrex pan.			
06/24/25			
Mass of (solid) residue (with slight DMS odor): 17 g			
06/25/25			
Added this crude solid to 40 g boiling IPA in a 100 mL beaker then hot-filtered into a clean 250 mL beaker. Washed (extr) filtered solids with 5 g boiling IPA. Sealed beaker with AL foil and let stand.			
06/26/25			
Decanted ML from crystalline mass into clean 100 mL beaker. Dried under stream of air overnight. Sealed beaker containing ML and refrigerated (cooler fridge ~10 °C).			
06/27/25			
Crop1 mass = 13.8 g (MP 100-101 °C)			
06/30/25			
Beaker containing ML was removed from cooler fridge and put on ice bath for 1 h before decanting ML from 2 nd crop of crystals. Crop2 mass = 800 mg (MP 100-102 oC) so total yield = 13.8 g + 0.8 g = 14.6 g (63%).			

Part 1b. Synthesis of Methyl-3,4,5-trimethoxycinnamate (110 g scale)

3,4,5-Trimethoxycinnamic acid (TMCA, MW 238.26 g/mol) **110 g, 462 mmol, 1 eq**

K₂CO₃ (MW 138.2 g/mol) **32 g, 231 mmol, 0.5 eq**

DMSO (MW 78.13 g/mol, d 1.1 g/mL) **541 g, 6.93 mol, 15 eq**

DMC (MW 90.08 g/mol, d 1.07 g/mL, BP 90 °C) **624 g, 6.93 mol, 15 eq**

NOTE: Product (TMCA) MW = 252.26 g/mol; CAS#: 7560-49-8, 20329-96-8; Rep MPs: 91-92, 99-100, 108-110 °C

08/07/25				
TIME	IT (°C)	EXSBT (°C)	VARIAC (%)	ACTION
~2 pm	N/A	N/A	N/A	Weighed DMSO, DMC, K ₂ CO ₃ , & TMCA (PHOTO A) into 3 L RBF equipped with thermowell
~2:30 pm	"	"	"	Assembled glass (Hickman Still under 16" vigreux column with N ₂ inlet) & flushed with N ₂ (used a separate bubbler); Turned on variac & set to 15%
3:30 pm	NC	NC	15	Insulated RBF
4 pm	39	NC	"	↑d variac setting to 16%
4:42 pm	46	54	16	
4:55 pm	49	56	"	↑d variac setting to 18%
6:02 pm	61	68	18	Took photo (PHOTO B)
10:45 pm	88	90	"	Condensation present in neck of Hickman still but nothing (ie. MeOH) collected
08/08/25				
10:18 am	95	96	"	~10 droplets collected in Hickman still; RM has yellowed; questioning if magnetic stirring is adequate
~5 pm	"	"	"	Set timer to cut variac (heat) off ~8:30 am Sunday- so considering as ~60 h at IT ≥88 oC

08/11/25 → 08/14/25

Poured RM into a 4 L FFlask containing (stirring = 3" stir bar) IpOAc & DI water (500 g + 500 g) which was pre-cooled in an in an ice bath to 13 oC (**PHOTO C**). This dilution was exothermic and caused an ↑ in IT (13 oC → 27 oC). The RBF and funnel was rinsed with IpOAc followed by DI water (200 g & 300 g, respectively). Stirred the mixture a few minutes then sealed with AL foil and let sit 24 h. Transferred mixture to 4 L sep funnel and separated layers. Extracted AQ layer again with IpOAc (225 g). Combined IpOAc layers were washed with DI water (2 x 200 g) followed by brine (1 x 250 g). Carried out TLC (**PHOTO D**). Evaporated to a mass of 149 g with a strong smell of DMS.




08/07/25 (~2 pm) PHOTO A
MATERIALS



08/07/25 (6:02 pm)
PHOTO B



08/11/25 (1:21 pm) PHOTO C (ABOVE ↑)
08/13/25 (6:40 pm) PHOTO D (RIGHT →)
TLC (~3:1 Hex:EtOAc): 1_TMCA (SM),
2_TMCA (authentic), 3_DHTMCA (authentic), 4_RM following AQ workup



08/15/25

Further EOCd to mass of 120 g. Still strong smell of DMS. Transferred bulk of crude from glass evap bowl to 600 mL beaker equipped with stir bar. Heated IPA (120 g) in a 600 mL beaker to near boiling and used this hot IPA to extract residue which was adhered to walls of bowl. Did this a 2nd time (120 g IPA) and added both solvent rinses to 600 mL beaker containing the bulk of crude product. Heated to boiling (and measured mass of solvent was 222 g after evap losses) and then hot-filtered (gravity-filter paper) into another clean 600 mL beaker. Sealed with AL foil and let stand. Crystals were observed on the bottom of the beaker in minutes following the hot filtration.

08/18/25

Broke up crystals with spatula, freeing from sides and bottom of beaker. Filtered crystals (aspirator, coarse fritted funnel, no wash); left on aspirator for ~10 min then spread crystals onto shallow glass pan to air dry.

08/19/25

Mass of (dry) Crop1 crystals = 108 g (93%, MP 101-102 oC) **PHOTO E**




PHOTO E

Part 2. Synthesis of Methyl-3-(1,3-benzodioxol-5-yl)-2-propenoate (38.4 g scale)

3,4-Methylenedioxybenzoic acid ("MDCA" MW 192.17 g/mol, MP 242-244 °C) **38.4 g, 0.2 mol, 1 eq**

Potassium carbonate (K₂CO₃) MW 138.2 g/mol **13.8 g, 0.1 mol, 0.5 eq**






DMSO (MW 78.13 g/mol, d 1.1 g/mL) **234 g, 3 mol, 15 eq**

DMC (MW 90.08 g/mol, d 1.07 g/mL, BP 90 °C) **(270 g, 3 mol, 15 eq) + (81 g, 0.9 mol, 4.5 eq) = 19.5 eq**


NOTE: Product (MDCAME) MW = 206.2 g/mol; MP (rep) 133-134 °C

06/10/26			
TIME	IT (°C)	Variac (%)	ACTION
~8:30 pm	N/A-RT	N/A	Took photo of SM, reagents, & solvents (besides extraction solvents) PHOTO A
06/11/26			
~2 pm	"	"	Testing solubility on a small scale prior to rxn: To a 1 oz glass shell vial equipped with stir bar was added DMSO (10 g) followed by MDCA (1.64 g) ← dissolved quickly. DMC (11.5 g) was added ← MDCA remained dissolved. Added K ₂ CO ₃ (590 mg) → remained mostly undissolved (which was expected). Stirred & did not observe any solubility issues (PHOTO B)
2:11 pm	"	"	Added DMSO (224 g) to a 1 L (3-neck) RBF equipped with stir bar & thermowell (PHOTO C)
2:14 pm	"	"	Weighed out MDCA (minus qty in shell vial)- PHOTO D
2:17 pm	"	"	Added MDCA to (stirring) DMSO in 1 L RBF via paper funnel (PHOTO E)
2:23 pm	22	"	MDCA dissolved (PHOTO F)
2:28 pm	"	"	Added DMC (258 g)- PHOTO G & stirred (PHOTO H)
2:30 pm	"	"	Weighed out K ₂ CO ₃ (minus qty in shell vial, PHOTO I) & added to (stirring) solution (PHOTO J)
2:35 pm	"	"	Added contents of shell vial to 1 L RBF (PHOTO K)
2:45 pm	23	"	Added 16" vigreux column/N ₂ inlet/bubbler & flushed system with N ₂ (PHOTO L) & added HM (PHOTO M)
2:51 pm	24	"	Stoppered RBF; adjusted N ₂ with slight positive pressure & powered ON variac/set to 20% (PHOTOS N & O)

space

			
<p>06/10/26 (8:33 pm) PHOTO A</p>	<p>06/11/26 (2:08 pm) PHOTO B</p>	<p>06/11/26 PHOTO C (2:11 pm) Above PHOTO D (2:14 pm) Top Right; PHOTO E (2:17 pm) Bottom Right</p>	

space



06/11/26 (2:23 pm)
PHOTO F


06/11/26 (2:28 pm)
PHOTO G

06/11/26 (2:29 pm)
PHOTO H

06/11/26
PHOTO I (2:30 pm)
TOP; PHOTO J (2:33
pm) BOTTOM

06/11/26 (2:35 pm)
PHOTO K

space



06/11/26 (2:45 pm)
PHOTO L

06/11/26 (2:47 pm)
PHOTO M

06/11/26
PHOTO N (2:51 pm) TOP
PHOTO O (2:52 pm)
BOTTOM

06/11/26
PHOTO P (6:44 pm) ABOVE
PHOTO Q (7:31 pm) TOP
RIGHT
PHOTO R (7:59 pm) TLC
BOTTOM RIGHT (Ignore
circled spots)


space

06/11/26 (continued)

TIME	IT (°C)	Variac (%)	ACTION
3 pm	28	20	
3:30 pm	45	"	

4 pm	58	"	
4:30 pm	67	"	
4:45 pm	70	"	
5:04 pm	74	"	
5:15 pm	75	"	
5:33 pm	76	"	
6:18 pm	79	"	
6:32 pm	"	"	↑d variac to 21%
6:41 pm	"	21	↑d variac to 22%
7 pm	81	22	
7:30 pm	82	"	
8 pm	84	"	
8:15 pm	85	"	
9 pm	88	"	↓d variac to 21%
9:20 pm	87	21	↓d variac to 20%
06/12/26			
7:20 am	85	20	
6:44 p → ~7:30 pm	83	"	Drew a (~1 g) ALQ for TLC using a glass pipet (PHOTO P); Let cool to RT then added ~1 g ice + ~1 g 10% HCl + 1 g Et2O & mixed well; Some insolubles present so transferred to a (larger) 2 nd shell vial & added enough EtOAc to give 2 clear layers (PHOTO Q); ↑d variac to 21%
~7:45 pm	NC	21	Carried out TLC (~2:1 Hex:EtOAc): 1_MDCA (MeOH solution) vs 2_RXN ALQ (PHOTO R): NO RXN YET!!
8:28 pm	87	"	Added Hickman Still between RBF & vigreux column
06/13/26			
9 am	88	"	↑d variac to 24%
5:30 pm	100	24	Took PHOTO S ; 1-2 mL water in Hickman Still
5:35 pm	"	"	Insulated RB & replaced Hickman Still with larger one containing anhydrous Na2SO4 (70 g); SLOWLY added more DMC (81 g) to RM (PHOTO T)
6:10 pm	101	"	Condensate beginning to appear in Hickman Still (PHOTO U)
6:55 pm	103	"	Hickman Still now ~2/3 full of solvent; ↓d variac to 23%
7:02 pm	104	23	Returned the other Hickman Still positioned above the first one (to potentially collect MeOH)- PHOTO V

space



06/13/26 (5:33 pm) PHOTO S

06/13/26 (5:42 pm)
PHOTO T







06/13/26 (6:10 pm) PHOTO U

06/13/26 (7:04 pm) PHOTO V

space

06/14/26			
TIME	IT (°C)	Variac (%)	ACTION
11:33 am	107	23	Maybe 1-2 mL in top Hickman Still
3:41 pm	108	"	~3-4 mL in top Hickman Still (PHOTO W) & distinct color change in RM (PHOTO X)
~3:45 pm	"	"	Drew ALQ (& diluted into EtOAc/5% HCl) & carried out TLC (PHOTO Y)- Looks done!!
4:05 pm	103	"	Turned variac power OFF, removed insulation, & took PHOTO Z ; Left to stir & cool slowly to RT





space

			
	<p>06/14/26 PHOTO Y (4:05 pm) TLC (~2:1 Hex:EtOAc): 1_MDCA; 2_ALQ (Ignore circled spots = permanent stains)</p>	<p>06/14/26 PHOTO Z (4:06 pm) ABOVE 06/15/26 PHOTO ZA (3:50 pm) TOP RIGHT PHOTO ZB (3:55 pm) →</p>	





space

06/15/26			
TIME	IT (°C)	Variac (%)	ACTION
2:50 pm	25	N/A	
3:50 pm	"	"	Decanted most of the supernatant liquid into 4 L beaker via large orange funnel (PHOTO ZA)
3:55 pm	"	"	Added filter paper & filtered remaining liquid + solids (K ₂ CO ₃); Washed (extr) RBF & filtered solids twice with Et ₂ O (200 g + 150 g)- PHOTO ZB
4:36 pm	23	"	Added 3" stir bar & suspended thermometer & stirred mixture at a low setting; Slowly added DI water (~400 g, PHOTOS ZC & ZD); Addition of water was exothermic: IT ↑d to 30 oC (PHOTO ZE)
~4:40 pm	30	"	Stirred vigorously for ~1 min (PHOTO ZF) then let stand; Some solids present at interface (PHOTO ZG)
~5 pm	RT	"	Filtered (thru filter paper in orange funnel) into a 2 L sep funnel (PHOTO ZH)
~5 pm → ~6 pm	"	"	Separated layers & extracted (lower) AQ layer again with IpOAc (250 g); Separated layers & combined Et ₂ O & IpOAc extracts in a 2 L EFlask; Washed combined organic extracts (twice) with DI water (2 x ~300 g) ← PHOTOS ZI & ZJ followed by brine (~300 g) ← PHOTOS ZK & ZL
6:19 pm	"	"	Sealed 2 L EFlask with plastic wrap/AL foil for later processing (PHOTO ZM)






space

			
<p>06/15/26 (4:38 pm) PHOTO ZC</p>	<p>06/15/26 PHOTO ZD (4:39 pm) TOP PHOTO ZE (4:40 pm) BOTTOM</p>	<p>06/15/26 PHOTO ZF (4:41 pm) TOP PHOTO ZG (4:55 pm) BOTTOM</p>	<p>06/15/26 PHOTO ZH (5 pm)</p>

space

			
<p>06/15/26 (5:46 pm) PHOTO ZI</p>	<p>06/15/26 (5:52 pm) PHOTO ZJ</p>	<p>06/15/26 (6:03 pm) PHOTO ZK</p>	<p>06/15/26 (6:08 pm) PHOTO ZL</p>




space

				
06/15/26 (6:19 pm) PHOTO ZM	06/17/26 PHOTO ZN (1:56 pm) TOP PHOTO ZO (2:01 pm) BOTTOM	06/17/26 (2:02 pm) PHOTO ZP	06/17/26 PHOTO ZQ (2:09 pm) TOP PHOTO ZR (2:10 pm) BOTTOM	

space

06/17/26			
TIME	IT (°C)	Variac (%)	ACTION
1:56 pm	RT	N/A	Observed that a small amount of crystals had formed in EFlask (PHOTO ZN)
~2 pm	"	"	Added Hex (195 g, PHOTO ZO) & heated crystals back into solution (PHOTO ZP)
~2:05 pm →2:10 p	"	"	Added ~3 teaspoons anhydrous Na ₂ SO ₄ & stirred warm solution (PHOTO ZQ) before decanting into a 2 L beaker (PHOTO ZR); The EFlask was then rinsed with Et ₂ O (~50 g) & decanted again.
~2:15 pm	"	"	The 2 L beaker was sealed with plastic wrap & AL foil (PHOTO ZS) & let stand to cool slowly to RT
~4:30 pm	"	"	After a couple of hours, the sealed beaker was put on an ice bath in a styrofoam cooler (PHOTO ZT)





space

		
06/17/26 (2:15 pm) PHOTO ZS	06/17/26 (4:29 pm) PHOTO ZT	06/18/26 (9:33 am) PHOTO ZU

space

06/18/26			
TIME	IT (°C)	Variac (%)	ACTION
~9:30 am	~0	N/A	Removed beaker & took PHOTO ZU ; Immediately returned beaker to styrofoam cooler & replenished ice
~4:30 pm → 5 p	"	"	Removed beaker & decanted ML into 2 L EFlask (PHOTO ZV); Left beaker suspended using a chain clamp to let last bit of ML drip into flammable waste bottle (PHOTO ZW).
06/19/26			
1:05 pm	RT	N/A	Mass of crystals (Crop1) = 16 g (PHOTO ZX)
~1:45 pm	"	"	Evaporated volatiles from Crop1 ML to afford a 2 nd crop of 6.9 g (which will require recrystallization) (PHOTO ZY); Yield = 22.9 g (56%)
06/26/26			
~7 pm	"	"	MP (Crop1) = 139-141 °C

space

			
06/18/26 (4:38 pm) PHOTO ZV	06/18/26 (4:49 pm) PHOTO ZW	06/19/26 (1:05 pm) PHOTO ZX Crop1	06/23/26 (6:11 pm) PHOTO ZY Crop2

ABBREVIATIONS:

AQ = aqueous

AL = aluminum

ALQ = aliquot

DI = deionized (ie. distilled)

DMS = dimethyl sulfide

DMSO = dimethyl sulfoxide

EFlask = erlenmeyer flask

EOCd = evaporated outdoors covered

Et₂O = diethyl ether

EtOAc = ethyl acetate

extr = extracted

EXSBT = external sand bath temperature

FFlask = filter flask

g = gram(s)

h = hour(s)

HCl = hydrochloric acid

Hex = hexanes

HM = heating mantle

IPA = isopropyl alcohol

IT = internal temperature

IpOAc = isopropyl acetate

L = Liter(s)

min = minute(s)

mg = milligram(s)

mL = milliliter(s)

ML = mother liquor

MP = melting point

MW = molecular weight

N₂ = nitrogen

PV = pressure vessel

qty = quantity

RB = round bottom

RBF = round bottom flask

RM = reaction mixture

RT = room temperature

RXN = reaction

SEP = separatory

SM = starting material

TLC = thin layer chromatography

↑ = increase

↓ = decrease

↑d = increased

↓d = decreased

