

Summary of Third-Party Testing on Phenotype Model EL-05

Testing results for Phenotype's 5-layer Surgical Masks performed by Third-Party Laboratories are summarized in the following Tables below.

Original Lab Report Documentation is available upon request.

<u>Important Note:</u> Testing results <u>do not represent official certification</u> by NIOSH. No such claim is made by Phenotype Pharmaceuticals. Data is provided as is, for <u>informational purposes only</u>, and is applicable to the mask body only. End user fit determines ultimate effectiveness.

Table 1: Advanced Particle Filtration – Sodium Chloride Aerosol Challenge Performed By: Nelson Labs (Salt Lake City, Utah)

	Test Method	Minimum Passing Standard	Test results
Sodium Chloride Aerosol Challenge – 0.3 μm particles	TEB-APR-STP-0059 /		PASS
	42 CFR Part 84 –	≥95% Filter Efficiency	≥99% Filter Efficiency,
	N95 Respirator		20/20 samples

Table 2: Breathability - Inhalation and Exhalation Resistance Performed By: Nelson Labs (Salt Lake City, Utah)

	Test Method	Minimum Passing	Test results
		Standard	
Inhalation Resistance Pressure	TEB-APR-STP-0007 /		PASS
	42 CFR Part 84 –	≤35mm H2O	≤22mm H2O ,
	N95 Respirator		3/3 samples
Exhalation Resistance Pressure	TEB-APR-STP-0003 /		PASS
	42 CFR Part 84 –	≤25mm H2O	≤19mm H2O,
	N95 Respirator		3/3 samples

Table 3: Fluid Resistance - Resistance to Penetration by Synthetic Blood Performed By: Eurofins Medical Device Testing (Lancaster, Pennsylvania)

	Test Method	Minimum Passing	Test results
		Standard	
Resistance to Synthetic Blood – Level 2 Barrier Rating	ASTM F1862/F1862M-17; ASTM F2100-19 Medical Face Mask	No Fluid Penetration at 120 mmHg Pressure (Level 2)	PASS No Fluid Penetration, Level 2, 32/32 samples



Table 4: Safe for Clothing - Flammability Rating Performed By: Eurofins Medical Device Testing (Lancaster, Pennsylvania)

	Test Method	Minimum Passing Standard	Test results
Flammability Class Rating 1	16 CFR Part 1610	Class 1 – fabric burns ≥ 3.5sec	PASS No Ignition, 9/9 samples