



# Phospho-Flow Cytometric Assays

**--Accurate evaluation of signaling molecules in each single cell**

## ABOUT PHOSPHOFLOW CYTOMETRY

Many proteins involved in cellular signaling cascades are phosphorylated. Fluorescently-labeled antibodies specific to a variety of phosphorylated proteins facilitate the detection of these important proteins at the single cell level with flow cytometry. Additionally, this assay allows for the examination of phosphoproteins in distinct cell subsets simultaneously in samples with multiple cell types. This assay is a more powerful and sensitive approach than traditional assays such as Western Blots which only highlight global protein status in a seemingly homogenous population.

## HOW DOES THE ASSAY WORK?

Phosphoflow cytometric assays can specifically detect the state of protein phosphorylation at the cellular level. Our assays are optimized for multiple cell types including immune cells. A general work flow is as follows:

1. Treated and control cells or tissue/blood samples are processed and stained with

viability dye and cell surface markers antibodies.

2. The cells are fixed and permeabilized.
3. The cells are then stained with the antibodies for intracellular phosphoproteins.
4. The cells are acquired on a flow cytometer.

## ASSAY BENEFITS

- **Superior specificity:** delivers greater specificity than other common technologies like Western blot, potentially revealing heterogeneity in a seemingly homogeneous cell sample.
- **Flexibility:** easily automated for use in routine compound screening
- **Simultaneous detection:** detect multiple targets (phosphoproteins) in the same cell type or in distinct cell subsets in a heterogeneous cell sample.
- **Standardized platform:** 96-well plate format using Cytoflex S flow cytometer.
- **Use of specialized reagents:** Assay conditions are optimized for accurate and reproducible data.
- Measure phosphorylation at a **single cell** level.
- Provides data for **mechanistic insights**.

## SAMPLE REQUIREMENTS


- **Minuscule samples**
- **Compatibility with heterogeneous samples:** works with a variety of sample types such as cultured cells, whole blood, PBMCs, isolated tissue cells (such as skin cells), etc






## SERVICE FEATURES

- ❖ **High throughput:** 96-well format
- ❖ **Simultaneous quantitation:** detect phosphoproteins in multiple cell populations in the same sample and/or multiple phosphoproteins in the same cell type
- ❖ **Highly reproducible**
- ❖ **State-of-art platforms:** Cytoflex S flow cytometer
- ❖ **Large inventory of phosphoprotein specific antibodies available**
- ❖ **Extensive data analysis**
- ❖ **Timely data delivery:** 1-4 weeks or sooner, upon receiving test samples
- ❖ **20+ years of experience:** Expert data analysis and interpretation, high quality scientific and technical support




**Example 1. Multidimensional analyses with phospho-specific flow cytometry.** (A) Gating strategy to identify Tregs and examine STAT5 phosphorylation in a human PBMC samples. Human PBMCs were treated with 60 ng/mL of wildtype (WT), Variant 1 and Variant 2 recombinant IL-2. The cells were then stained with fluorescently-conjugated antibodies against CD3, CD56, CD4, CD8, FoxP3 and phosphorylated STAT5 (pSTAT5). (B) The level of STAT5 phosphorylation in Tregs in response to a fixed concentration of different IL-2 variants was analyzed.





**Example 2. Evaluation of pSTAT5 in multiple cell subsets following IL-2 treatment.** (A) Gating strategy to identify Natural Killer (NK) cells, CD8<sup>+</sup> T cells and Tregs in a PBMC sample. Normal human PBMCs were treated with 60 ng/mL of wildtype (WT) recombinant IL-2. The cells were then stained with fluorescently-conjugated antibodies against CD3, CD56, CD4, CD8, FoxP3 and phosphorylated STAT5 (pSTAT5). (B) The level of STAT5 phosphorylation in Natural Killer (NK) cells, CD8<sup>+</sup> T cells and Tregs in response to WT IL-2 was analyzed.



**Example 3. The effect of IL-2 treatment on STAT5 phosphorylation in Tregs.** Human PBMCs were treated with varying concentrations of wildtype (WT), Variant 1 and Variant 2 recombinant IL-2 for 2 h. The cells were then stained with appropriate fluorophore-conjugated antibodies and acquired on the Cytoteflex S cytometer. The mean fluorescence intensity (MFI) of phosphorylated STAT5 (pSTAT5) was determined at each concentration.

