

Phone: 203-433-4329 Fax: 203-208-0664 Toll Free: 866-702-0260 www.phoenixsongsbio.com

Hepatocyte Culture Medium

Complete Medium for the Maintenance of Mature Human Hepatocytes

Hepatocyte Culture Medium (HCM) is formulated for the maintenance of mature fresh, cryopreserved and even human stem cell-derived human hepatocytes. HCM is a Williams E base medium without phenyl red and contains stabilized L-glutamine and HEPES. After HCM Serum-Free Supplements are added to the basal medium and filtered to ensure sterility, the medium is complete and ready to use for maintenance of human hepatocytes plated on Collagen, Liver Biomatrix or in co-culture with human liver stromal cells. HCM is designed to maintain human hepatocytes for longer than 30 days in co-culture when fed daily. Hepatocyte Plating Medium is complete so you don't need to source additional components. All supplements required for making this medium complete come separate in frozen vials in the supplement pack. The basal medium is shipped at ambient temperature and the supplement packs are shipped on dry ice.

Making Hepatocyte Culture Medium

- 1. Hepatocyte Culture Medium (Cat.# HCM500) includes the following components:
 - a. Hepatocyte Culture basal medium (500ml bottle) is stored at 4°C
 - b. Hepatocyte Culture Supplement Pack is stored at -80°C and contains the following components

| Supplement Pack Components | Amount to Add to 250ml Medium |
|-----------------------------------|-------------------------------|
| Serum-Free Maintenance Supplement | 5.0 ml |
| Dexamethasone | 5.0 µl |
| Penn/Strep | 5.0 ml |

Preparing Hepatocyte Culture Medium

- 1. Warm the basal medium to room temperature
- 2. Thaw the contents of the Serum-Free Supplement Pack
 - a. Aseptically add the supplements to the basal medium
- 3. Filter through 0.22 µm filter unit
- 4. Store the complete medium at 4°C in the dark until ready for use.
- 5. Warm the complete medium to room temperature immediately before use.
 - a. Do not leave the complete medium in a water-bath for longer than 5 minutes as this may result in degradation of the medium components

Note: Degradation of some components is toxic to cells

6. Use complete medium within two weeks

Use Restrictions: The Products and components of the Products are to be used for RESEARCH PURPOSES only and are not for use in humans or animals. No rights expressed or implied are conveyed by the sale of the Products to make, have made, offer to sell or sell the Products alone or in combination with any other products or to use for commercial purposes without prior written consent.