

StemCell Freezing Medium And Instructions

StemCell Freezing Medium contains DMSO and is optimized for the cryopreservation of hepatic and biliary tree stem cells.

Required Supplies

1. Kubota's StemCell Growth medium (Catalog # 11001-250)
2. StemPro Accutase (LifeTech A11105-01)
3. StemCell Freezing medium for hepatic and biliary tree stem cells (Catalog # 11005-050)

Preparing hepatic and biliary tree StemCells for Cryopreservation

1. For all volumes refer to Table 1 below.
2. Working with multiple plates or flasks
 - a. Only work with a stack of 4 x 100mm plates or 4 flasks at a time
 - b. Minimize time out of the incubator during the process
 - c. Work swiftly but carefully to minimize cell loss
3. Remove the media from the plate or flask and add cold Accutase (see table 1 for volumes) (Do not warm Accutase in the waterbath)
4. Take a plate or flask of cells immediately to the microscope so that you can observe the dissociation.

NOTE: Do not place the plates or flasks in an incubator during the dissociation process. Accutase will loose activity if allowed to warm up to 37°C.

5. Tap the plates or flasks against the palm of your hand to dislodge the cells and when the cells are free floating, return the plates to the safety cabinet
6. Gently pipette up and down to ensure single cell suspension, then transfer cell suspension into appropriate size conical tube.
7. Wash the plate or flask with Kubota's StemCell Growth medium and add to the tube, pipetting up and down to ensure a homogenous cell suspension
8. Take an aliquot for cell counting
 - a. Take 20ul cell suspension and add it to 20ul Trypan Blue
 - b. Pipette up and down
 - c. Load 10ul of cell suspension to both counting chambers of a hemacytometer
 - d. Count the cells within the center grid
 - e. Calculate the total cell count

$$\text{Cell count} \times 10,000 \times 2 \times \text{Volume} = \text{Total cell count}$$

Note: 10,000 is a hemacytometer constant and 2 is the Dilution factor

9. Divide the total cell count by 1.5×10^6 to determine the number of cryovials needed for freezing down the cells for long-term storage.
10. Pellet the cells in all the tubes by centrifugation at 500 x g for 5 minutes.
11. While the cells are in the centrifuge print the labels needed for each cryovial
12. Place one label on each cryovial.

13. After the centrifuge stops, resuspend the cell pellet in enough StemCell Freezing Medium to freeze 1.5×10^6 /ml as calculated in step 9 above.
14. Transfer 1ml StemCell Freezing cell suspension into each cryovial
15. Transfer the cryovials into a controlled rate freezer or a styrofoam container then place the container into a -80°C freezer over night
16. Next day transfer the vials of cells into the LN2 freezer for long-term storage

Table 1. Volumes of Accutase and media

Dish/Flask Size	Growth Area (cm ²)	Accutase Volume	Media Volume
100mm	58.1	2	2
60mm	21.3	1	1
35mm	9.6	0.5	0.5
T-75	75	10	10
T-175	175	15	15
T-225	225	20	20

Use Restrictions: The Products and components of the Products are to be used for RESEARCH PURCHASES ONLY and are not for use in humans or animals. No rights expressed or implied are conveyed by the sale of the Products to make, have made, offer to sell or sell the Products alone or in combination with any other products or to use for commercial purposes without PhoenixSongs' prior written consent.