

## Helicobacter ViraStripe® IgG Test Kit

Stripe-Immunoblot for the qualitative detection of **IgG** antibodies against specific **Helicobacter pylori** antigens in human serum resp. for **differentiation** between **high** and **low pathogenic Helicobacter pylori** strains.

The **Helicobacter ViraStripe® IgG Test Kit** is an **immunoblot** based on an enzyme-immunoassay in the line/strip format, carrying the following purified Helicobacter specific antigens: **CagA, VacA, p90, UreA, p26, p24** and **p17**.

### Principle of the assay

During the serum incubation step Helicobacter specific IgG antibodies bind to the immobilised antigens on the test strip. During the conjugate reaction, AP-conjugate binds to the antigen-antibody complex. The alkaline phosphatase converts the chromogen/substrate and thus, stains the antigen-antibody complex on the test strip purple. The washing procedures following serum, conjugate and chromogen/substrate incubation steps remove unbound antibodies and reagents.

The green separation line divides the test strip into a control section and an analytical section. The control section contains the **negative control band**, the **serum control**, **three conjugate controls** (IgG, IgA, IgM) and the **cut off control**.

The test strip code for Helicobacter ViraStripe® IgG test strips is **hG**. Test strips are numbered from **01** to **50**. The analytical section contains the Helicobacter specific antigens.

Order No.:	<b>V-HPSGOK</b>	Order No.:	<b>V-HPSDK (Deca Kit)</b>
Kit size:	<b>1x 50 test strips</b>	Kit size:	<b>10x 50 test strips</b>
Specimen:	<b>20 µl serum</b>	Specimen:	<b>20 µl serum</b>
Time for testing:	<b>approx. 90 minutes</b>	Time for testing:	<b>approx. 90 minutes</b>

### Materials provided

1x or 10x 50 test strips	<b>Helicobacter ViraStripe® Antigen Strips (IgG)</b> Test strips including a control section and Helicobacter specific antigens in the analytical section, ready to use	(Prod. No.: V-HPSGAS)
1x or 10x 9 ml	<b>ViraStripe® / ViraBlot® AP-Anti-Human IgG Conjugate Concentrate</b> , goat	(Order No.: V-UVNGKI)
1x or 10x 100 ml	<b>ViraStripe® / ViraBlot® Diluent / Wash Buffer</b> 10x concentrate	(Order No.: V-UVNUWP)
1x or 10x 5 g	<b>ViraStripe® / ViraBlot® Diluent / Wash Powder</b>	(Order No.: V-UVNUMP)
1x or 10x 90 ml	<b>ViraStripe® / ViraBlot® Chromogen / Substrate Solution</b> Ready to use	(Order No.: V-UVNUCS)
1 or 10 copies	<b>Evaluation Protocol for Helicobacter ViraStripe® IgG Test Kit</b>	

### Additionally available

330 µl	<b>Helicobacter ViraStripe® IgG Positive Control</b> Human, ready to use	(Order No.: V-HPSGPK)
330 µl	<b>Helicobacter ViraStripe® IgG,A,M Negative Control</b> Human, ready to use	(Order No.: V-HPSPNK)
50 copies	<b>Helicobacter ViraStripe® IgG evaluation protocols</b> for automated interpretation with ViraScan® software	(Order No.: V-HPSGEP)

### Preparation of reagents and patient samples

**Bring all reagents to room temperature (20-25°C) prior to use.** Information about stability can be found on page 5.

<b>Diluent / Wash Buffer Working Dilution:</b>	Dilute <b>Diluent / Wash Buffer Concentrate 1:10</b> with distilled or deionised water (100 ml concentrate + 900 ml water). Add Diluent / Wash Powder completely and stir well until all powder is dissolved. If needed, place onto a magnetic stirrer for 10-15 minutes. The pH value should be around pH 7.5 at 20°C.
<b>Antigen Strips:</b>	Carefully separate the required number of test strips by use of <b>forceps</b> at the <b>label</b> and place the test strips in the prepared incubation tray (see assay procedure, step 2). Use test strips directly after removing from packing. Do not touch test strips by hand. Return unused test strips directly into the original packing, seal well and store at 2-8°C.
<b>Patient samples:</b>	Use <b>20 µl patient serum</b> undiluted per test strip.
<b>Controls:</b>	Use <b>100 µl of Positive Control</b> or <b>100 µl of Negative Control</b> undiluted per test strip respectively.
<b>Conjugate Working Dilution:</b>	Prepare <b>Conjugate Concentrate 1:10</b> with Diluent / Wash Buffer Working Dilution (see table 1). Prepare freshly prior to each test run. Do not store for further use.
<b>Chromogen / Substrate Solution:</b>	Ready to use.

## Helicobacter ViraStripe® IgG Test Kit

- 2 -

**Preparation of Conjugate Working Dilution IgG**

Number of strips	Diluent / Wash Buffer Working Dilution	Conjugate Concentrate	Final volume	Number of strips	Diluent / Wash Buffer Working Dilution	Conjugate Concentrate	Final volume		
1	1.35 ml	+	0.15 ml	1.5 ml	26	35.10 ml	+	3.90 ml	39.0 ml
2	2.70 ml	+	0.30 ml	3.0 ml	27	36.45 ml	+	4.05 ml	40.5 ml
3	4.05 ml	+	0.45 ml	4.5 ml	28	37.80 ml	+	4.20 ml	42.0 ml
4	5.40 ml	+	0.60 ml	6.0 ml	29	39.15 ml	+	4.35 ml	43.5 ml
5	6.75 ml	+	0.75 ml	7.5 ml	30	40.50 ml	+	4.50 ml	45.0 ml
6	8.10 ml	+	0.90 ml	9.0 ml	31	41.85 ml	+	4.65 ml	46.5 ml
7	9.45 ml	+	1.05 ml	10.5 ml	32	43.20 ml	+	4.80 ml	48.0 ml
8	10.80 ml	+	1.20 ml	12.0 ml	33	44.55 ml	+	4.95 ml	49.5 ml
9	12.15 ml	+	1.35 ml	13.5 ml	34	45.90 ml	+	5.10 ml	51.0 ml
10	13.50 ml	+	1.50 ml	15.0 ml	35	47.25 ml	+	5.25 ml	52.5 ml
11	14.85 ml	+	1.65 ml	16.5 ml	36	48.60 ml	+	5.40 ml	54.0 ml
12	16.20 ml	+	1.80 ml	18.0 ml	37	49.95 ml	+	5.55 ml	55.5 ml
13	17.55 ml	+	1.95 ml	19.5 ml	38	51.30 ml	+	5.70 ml	57.0 ml
14	18.90 ml	+	2.10 ml	21.0 ml	39	52.65 ml	+	5.85 ml	58.5 ml
15	20.25 ml	+	2.25 ml	22.5 ml	40	54.00 ml	+	6.00 ml	60.0 ml
16	21.60 ml	+	2.40 ml	24.0 ml	41	55.35 ml	+	6.15 ml	61.5 ml
17	22.95 ml	+	2.55 ml	25.5 ml	42	56.70 ml	+	6.30 ml	63.0 ml
18	24.30 ml	+	2.70 ml	27.0 ml	43	58.05 ml	+	6.45 ml	64.5 ml
19	25.65 ml	+	2.85 ml	28.5 ml	44	59.40 ml	+	6.60 ml	66.0 ml
20	27.00 ml	+	3.00 ml	30.0 ml	45	60.75 ml	+	6.75 ml	67.5 ml
21	28.35 ml	+	3.15 ml	31.5 ml	46	62.10 ml	+	6.90 ml	69.0 ml
22	29.70 ml	+	3.30 ml	33.0 ml	47	63.45 ml	+	7.05 ml	70.5 ml
23	31.05 ml	+	3.45 ml	34.5 ml	48	64.80 ml	+	7.20 ml	72.0 ml
24	32.40 ml	+	3.60 ml	36.0 ml	49	66.15 ml	+	7.35 ml	73.5 ml
25	33.75 ml	+	3.75 ml	37.5 ml	50	67.50 ml	+	7.50 ml	75.0 ml

**Table 1: 1:10 dilution of conjugate concentrate with Diluent / Wash Buffer Working Dilution**
**Assay procedure**

- Rinse incubation tray channels once with 1.5 ml Diluent / Wash Buffer Working Dilution, decant the liquid**  
Mark the trays with water-resistant pen. Rinsing removes dust particles.
- Place the needed amount of test strips into the incubation tray - one test strip per channel**  
For each patient serum and each control, carefully separate one test strip by use of forceps at the label and place them into the incubation tray channels. **The side showing the green separation line and the label must face up.**
- Add 1.5 ml Diluent / Wash Buffer Working Dilution and incubate by rocking for 5 minutes at room temperature (RT)**  
Make sure the test strips are completely covered with liquid. Use a 2D rocker with a rocking frequency of approx. 40/min. Avoid spilling of liquid. **Do not decant the liquid after incubation.**
- Add 20 µl of each patient serum or 100 µl of each control**  
Add patient sera and controls directly onto the labelled end of the test strips while the 2D rocker is running or make sure to tilt the incubation tray after adding each serum.
- Incubate by rocking for 30 minutes at RT**  
Make sure the test strips are completely covered with liquid. Use a 2D rocker with a rocking frequency of approx. 40/min.
- Decant the liquid**  
Carefully tap the incubation tray on absorbent paper to remove the remaining liquid. **Test strips adhere to the incubation tray when liquid is decanted.**
- 3 x washing:**  
- add 1.5 ml Diluent / Wash Buffer Working Dilution  
- incubate by rocking for 5 minutes at RT  
- decant the liquid  
Wash on the 2D rocker. Carefully tap the incubation tray on absorbent paper to remove the remaining liquid.
- Add 1.5 ml fresh Conjugate Working Dilution**  
Make sure the test strips are completely covered with Conjugate Working Dilution
- Incubate by rocking for 15 minutes at RT**  
Make sure the test strips are completely covered with liquid. Use a 2D rocker with a rocking frequency of approx. 40/min.
- Decant the liquid**  
Carefully tap the incubation tray on absorbent paper.
- 3 x washing as in step 7**  
Wash on the 2D rocker.
- Add 1.5 ml distilled or deionised water and incubate by rocking for 1 minute at RT**  
Make sure the test strips are completely covered with liquid. Use a 2D rocker with a rocking frequency of approx. 40/min.
- Decant the liquid**  
Carefully tap the incubation tray on absorbent paper.
- Add 1.5 ml Chromogen / Substrate Solution**  
Make sure the test strips are completely covered with liquid.
- Incubate by rocking at RT**  
Stop the reaction as soon as the cut off control becomes visible. **The Cut off control is located in the test strip control section. Caution:** Prolonged incubation causes background staining.
- Stop the reaction by decanting the liquid**  
Carefully tap the incubation tray on absorbent paper.
- Wash 3 x with 1.5 ml distilled or deionised water**  
Wash without incubation time.
- Dry test strips for interpretation**  
Carefully tap the incubation tray on absorbent paper to remove the remaining liquid. Place wet test strips with forceps on unbleached absorbent paper and allow to air dry before interpretation.

## Assay interpretation

- Evaluation protocol:** Record data on the evaluation protocol. Glue the test strips on the evaluation protocol. Place the green separation line of the test strips exactly onto the separation line printed on the evaluation protocol.
- Validity of test strips:** A test strip is considered as valid if the following bands are visible:
  - The **serum control**.
  - The **conjugate control** of the conjugate class being used. If more than one of the three conjugate controls becomes visible, the strongest band must indicate the appropriate conjugate class.
  - The **cut off control**.**and** if the following band is **not** visible:
  - The **negative control band**.
 Do not assess invalid test strips!
- Assignment of antigen bands:** The green separation line of the test strips indicates position and orientation for the assignment of bands with the bandlocator on the evaluation protocol. Assign bands and record results according to 4.
- Assessment of bands:** According to quality laboratory guidelines, the use of a cut off control for each run is recommended (7). **The cut off control of the Helicobacter ViraStripe® IgG is located in the control section on each test strip.** The intensity of the cut off control indicates the threshold of which bands are being assessed:
 

A band is considered as **distinct** if its intensity is **equal** to or **higher** than the intensity of the cut off control. Mark bands with **X** in the evaluation protocol appropriately.

A band is not assessed if it is **not present** or if its intensity is **lower** than the intensity of the cut off control.
- Interpretation of patient bands:** Patient bands have to be considered as symptoms of the disease. A final clinical diagnosis should always be made considering anamnesis, clinical manifestations and laboratory data. Bands of the following antigens are considered as **highly specific** for Helicobacter pylori: **CagA, VacA, p90, p30, UreA (specific), p26, p24 and p17**


## IgG Interpretation criteria

General note: Distinct bands must have a minimum intensity ( $\geq$  cut off), which has to be determined by the cut off control. The cut off control is located in the control section of each test strip.

Identified bands	Result	Interpretation
At least <b>one distinct</b> band out of: <b>CagA, VacA</b>	<b>Positive for type I</b>	Specific antibodies against <b>Helicobacter pylori type I</b> detectable. An infection with a high pathogenic Helicobacter pylori strain is probable.
At least <b>two distinct</b> bands out of: <b>p90, p30, UreA, p26, p24, p17</b>	<b>Positive for type II</b>	Specific antibodies against <b>Helicobacter pylori type II</b> detectable. An infection with a low pathogenic Helicobacter pylori strain is probable.
<b>One distinct</b> band out of: <b>p90, p30, p26, p24, p17</b>	<b>Equivocal</b>	Specific antibodies against <b>Helicobacter pylori</b> detectable. Evidence of a Helicobacter pylori infection. As control, check a second sample for IgG and IgA specific antibodies after 3-4 weeks.
<b>No band(s)</b> or singular distinct UreA band	<b>Negative</b>	No specific antibodies against Helicobacter pylori detectable.

If bands for both type I and type II are present, the result is considered positive for type I.

## IgG test strip



**Figure 1:** Schematic drawing of a Helicobacter ViraStripe® IgG test strip in full scale.

## Nomenclature and description of Helicobacter species bands from literature

Band nomenclature	Antigen	Comments
CagA	Highly specific Cytotoxin-associated-gene A	Highly specific for H. pylori strains type I (1); known as a marker for ulcer, detected in 80-100% of patients with gastric or duodenal ulcers, less frequent among patients showing gastritis without ulcer (2, 3, 4).
VacA	Highly specific Vacuolating-Cytotoxin A	Highly specific for H. pylori strains type I (1); known as a marker for ulcer and frequently co-expressed with CagA (3). Different variants of this protein are described (5).
p90	Highly specific	Protein with a molecular weight of 90 kD.
p30	Highly specific	Protein with a molecular weight of 30 kD.
UreA	Specific Urease-subunit A	UreA shows some similarity with urease-subunit A of other organisms, specific for H. pylori (6).
p26	Highly specific	Protein with a molecular weight of 26 kD.
p24	Highly specific	Protein with a molecular weight of 24 kD.
p17	Highly specific	Protein with a molecular weight of 17 kD.

## Diagnostic significance of antibodies against Helicobacter pylori

**1. High pathogenic (type I strains) and low pathogenic strains (type II strains)** are described for Helicobacter pylori (1). According to the WHO, high pathogenic strains are classified as human carcinogen class I. High and low pathogenic strains differ in two highly specific proteins which are mainly described for high pathogenic strains: **CagA** (Cytotoxin associated gene A) and **VacA** (Vacuolating cytotoxin A). In case of H. pylori infection with type II strains, antibodies against these proteins are generally missing. Both strains have in common the proteins Urease subunit A and five not further characterised proteins with a molecular weight of 90 kD, 30 kD, 26 kD, 24 kD and 17 kD.

**2.** Approximately 50% of the world population are carriers of the gram negative Helicobacter pylori bacterium, with chronic gastritis as major symptom. Severe infections with symptoms like peptic ulcers and gastric cancer are only seen in rare cases. The virulence of different Helicobacter pylori strains is mainly regulated by expression of several mediators, themselves being present in multiple genetic variations (8). Helicobacter pylori has developed active mechanisms to avoid

opsonisation and is able to retard immunogenic phagocytosis, resulting in permanent bacterial persistence (9).

**3. IgG antibodies** usually appear several weeks to months after an infection and might not be detectable in the early stage of the infection. IgG and IgA antibodies should be checked in case of a suspected recent infection / colonisation and a second sample should be analysed at a later time point. Patients with a chronic disease are usually positive for IgG antibodies. Titers steadily decrease after eradication therapy (11).

**4. IgM antibodies** usually appear 1-8 weeks after beginning of the disease.

**5.** In some patients **IgA antibodies** indicate an inflammation of the gastric mucosa.

**6.** An early antibiotic therapy can suppress the development of antibodies.

**7.** Medication and immunoglobulin therapy can cause unspecific antibody reactions (11).

## IgG performance data

### Sensitivity:

The **clinical sensitivity** of the Helicobacter ViraStripe® IgG Test Kit was determined using 15 clinical defined Helicobacter pylori positive sera:

	Helicobacter ViraStripe® IgG Test Kit, % (n)
H. pylori - positive, (n= 15)	<b>100% (15)</b>

80 sera found positive by the reference test Helicobacter Virablot® IgG Test Kit, were used to determine the sensitivity of the Helicobacter ViraStripe® IgG Test Kit:

	Helicobacter ViraStripe® IgG Test Kit, % (n)
H. pylori - positive with reference test, (n= 80)	<b>100% (80)</b>

### Specificity:

94 sera found negative by the reference test Helicobacter Virablot® IgG Test Kit, were used to determine the specificity of the Helicobacter ViraStripe® IgG Test Kit:

	Helicobacter ViraStripe® IgG Test Kit, % (n)
H. pylori - negative with reference test, (n= 94)	<b>95% (89)</b>

## Helicobacter ViraStripe® IgG Test Kit

- 5 -

**Warnings and precautions**

- All human serum components were tested for HCV, HIV1,2 antibodies and HBs antigens and found to be negative. Nevertheless, all human kit components as well as the patient samples should be considered as potentially infectious and handled according to safety precautions. While working with potentially infectious/hazardous materials, all national and international rules, regulations, guidelines and laws must be taken into account. This also applies to storage and disposal of chemicals and reagents being used.
- While working with hazardous or toxic substances/ biological agents precautions have to be applied following national biosafety regulations. In general, biological and chemical agents should be handled according to „Good Laboratory Practice (GLP)“ guidelines. Precautions among others are:
  - Do not pipette by mouth.
  - Wear disposable gloves and safety glasses while working.
  - Do not eat, drink or smoke in the working area.

**Storage and stability of reagents**

- Test strips:** In closed bags stable until the expiration date if stored at 2-8°C.
- Conjugate Concentrate:** Stable until the expiration date if stored at 2-8°C.
- Conjugate Working Dilution:** Prepare freshly prior to each run. Do not store for further use.
- Diluent / Wash Buffer Concentrate, 10x:** Stable until the expiration date if stored at 2-8°C.

**Specimen indications**

- The **Helicobacter ViraStripe® IgG Test Kit** must be used with human serum.
- Only clear, non-hemolysed, non-microbially contaminated specimens must be used.
- Using icteric, lipemic, hemolytic and/or heat-inactivated serum may lead to false results.

**Limitation of the procedure**

- To ensure reliable results, follow carefully the Instruction for Use and “Good Laboratory Practice”.
- A positive result is based on elevated specific antibody titers and should be considered as a symptom. The correlation to a disease is only conditionally possible.
- A negative result does not exclude a contact with the pathogen or the presence of a disease.
- Adequately trained personnel only should perform the assay procedure.

**Literature**

- XIANG, Z. et al.: Analysis of expression of CagA and VacA virulence factors in 43 strains of Helicobacter pylori reveals that clinical isolates can be divided into two major types and that CagA is not necessary for expression of the vacuolating cytotoxin. *Infect. Immun.* 1995, 63 (1): 94-98
- COVER, T. L. et al.: Helicobacter pylori: A bacterial cause of gastritis, peptic ulcer disease and gastric cancer. *ASM News*, 1995, 61 (1): 21-26
- COVACCI, A. S.: Molecular characterisation of the 128 kDa immunodominant antigen of Helicobacter pylori associated with cytotoxicity and duodenal ulcer. *PNAS*, 1993, 90: 5791
- COVER, T. L. et al.: Serologic detection of infection with cagA+ helicobacter pylori strains. *J Clin. Microbiol.*, 1995, 33 (6): 1496-1500
- COVER, T. L. et al.: Divergence of genetic sequences for the vacuolating cytotoxin among Helicobacter pylori strains. *J Biol. Chem.* 1994, 269 (14): 10556-10573
- MALFERTHEINER, P., Helicobacter pylori – Von der Grundlage zur Therapie, Thieme, 1996, 2. Aufl.
- RILI-BÄK: Bäk-Richtlinie zur Qualitätssicherung quantitativer laboratoriumsmedizinischer Untersuchungen, (01.04.2008), [www.bundesaeztekammer.de](http://www.bundesaeztekammer.de)
- OLIVARES, D. et al.: Factors involved in the pathogenesis of Helicobacter pylori infection. *Rev Esp Enferm Dig* 2006; 98: 374-386
- ALLEN, L. A.: Phagocytosis and persistence of Helicobacter pylori. *Cellular Microbiology* (2007) 9(4), 817–828
- DIN 58967-40: Serodiagnostik von Infektions- und Immunkrankheiten, Immunoblot (IB), Begriffe, Allgemeine methodenspezifische Anforderungen
- THOMAS, L.: 2008, Labor und Diagnose, Med. Verlagsgesellschaft Marburg

- The chromogen/substrate solution contains BCIP and NBT. Avoid contact with skin and mucous membranes. In case of contact with skin and eyes wash immediately with large quantities of water.
- Samples and all potentially contaminated materials must be decontaminated using validated laboratory techniques, e.g. by autoclaving 20 minutes at 121°C under humid conditions. Liquid disposals can be mixed with sodium hypochlorite to a final concentration of 1% sodium hypochlorite. Incubate 30 min for complete disinfection.
- Please refer to material safety data sheets for detailed information on potential risks, first aid guidelines, accidental release measures, handling and storage recommendations, personal protective equipment, directions for disposal and indications to toxicology.

- Diluent / Wash Buffer Working Dilution:** Stable for 2 weeks if stored at 2-8°C. For longer storage, aliquot and freeze at -20°C.
- Diluent / Wash Powder:** Stable until the expiration date if stored at 2-8°C.
- Chromogen / Substrate Solution:** Stable until the expiration date if stored at 2-8°C. Avoid exposure to light!









- Normally, human serum can be stored up to 5 days at 2-8°C. Specimens may be stored at -20°C (or below) for long term storage.
- Prior test processing, specimens should have reached room temperature. Mix specimens carefully after thawing. Precipitates in specimens can be removed by centrifugation.
- Avoid multiple freeze and thaw cycles.

- The detection of specific antibodies can vary within different assays from different manufacturers and can lead to different results due to different sensitivity, specificity and assay methodologies.
- Test strips showing a high background level should not be interpreted, especially if band intensities are lighter than the background level.
- In vitro* diagnostics must not be used beyond expiration date as reliable results may not be possible.
- Efficient washing after each incubation step is essential for consistent results; insufficient washing may lead to false results.

Helicobacter ViraStripe® IgG Test Kit

- 6-

**Symbols used**

	Manufacturer	<b>REF</b>	Order Number
	Refer to Instructions for Use		Use by / Expiration Date
<b>IVD</b>	<i>In-Vitro</i> Diagnostic Medical Device		Temperature Limitation (Storage)
<b>LOT</b>	Test Kit Lot Number	<b>CONTROL +</b>	Positive Serum Control
	Sufficient for 50 Tests	<b>CONTROL -</b>	Negative Serum Control
	Room Temperature in °C	<b>CONTROL</b>	Control
	User	<b>DATE</b>	Date
<b>#</b>	Serum Number	 <b>SUBSTRATE</b>	Chromogen/Substrate Incubation Time in Minutes
<b>PROTOCOL</b>	Evaluation Protocol	<b>No</b>	Protocol Number