



Hospicare™ 70% v/v Isopropyl Alcohol Disinfectant Wipes

Hospicare™ 70% v/v Isopropyl Alcohol Disinfectant Wipes contains 70% Isopropyl alcohol solution. Its fast-acting solution is able to effectively kill 99.99% harmful germs. Provides an effective disinfection of environmental and clinical surfaces. Tested according to EN 1040, EN 1276, EN 13624, EN 13727, EN 14348, EN 14476, EN 14561, AATCC 100 Test Method.

- Biodegradable solution
- Kills 99.99% bacteria, fungi and viruses
- Antimicrobial effectiveness tested.
- Tested according to EN 1040, EN 1276, EN 13624, EN 13727, EN 14348, EN 14476, EN 14561, EN 1650.
- EN 14476 Tested – Effective against enveloped viruses including Coronaviruses and Monkeypox virus in 30 seconds

PRODUCT SPECIFICATIONS

Product Code: HOSPICARE-IPA70-0003

Product Description:

Hospicare™ 70% v/v Isopropyl Alcohol Disinfectant Wipes With Cap

HS Code: 48182000

Product Barcode: 8887245225496

Carton Barcode: 8887245235785

ORDERING INFORMATION

Product Packing: 200 Sheets/ Canister

Carton Packing: 200 Sheets x 12 Canisters

PATENTED PRIMARY PACKAGING

Canister Material: HDPE with Cap

Canister Dimension: H250 x Ø115 mm

INGREDIENT

70% v/v Isopropyl Alcohol, Deionized Water

OUTER PACKAGING

Carton Material: Thick Double Wall

Carton Dimension: 440(L) x 330(W) x 255(H) mm

NONWOVEN

Paper material: Nonwoven 23gsm, Paper Dimension: 260(L) x 210(W) mm



High Absorbency



No Fluorescer



Soft And Delicate



Both for Wet and Dry

APPLICATION PROPERTIES

- For general cleaning and disinfection of hard surfaces e.g. tables, counter tops, equipment.

DIRECTIONS FOR USE:

- 1) Test on small corners to check for colour fades or runs. Clean the surface in an S shape.
Do not use on multiple surface. Let the surface air dry.
- 2) Use a new wipe for a new surface.
- 3) Discard used wipes in a waste bin.

USABLE LIFE AND STORAGE

This product has a usable life of 2 years from the date of production. Store in a cool, dry and well-ventilated area to prevent drying. Keep away from open flame and heat.

PRECAUTIONS

This is for surfaces use only. Isopropyl alcohol will attack some forms of plastic, rubber, and coatings. Test a small portion of the surface for some colour change and reaction prior to use. Allow the tools and surfaces to dry before use. Flammable, keep away from heat and flame. Keep out of reach of children

EFFICACY DETAILS

TEST ORGANISM(S)	TEST METHOD	CONTACT TIME
BACTERIA		
Methicillin resistant Staphylococcus aureus (MRSA)	EN 13727	30 Seconds
	EN 1276	2 Minutes
Vancomycin resistant Enterococcus faecalis (VRE)	EN 13727	30 Seconds
	EN 1276	2 Minutes
Staphylococcus aureus	EN 13727, EN 14561, EN 1040	5 Minutes
Pseudomonas aeruginosa	EN 1276	30 Seconds
	EN 13727, EN 14561, EN 1040	5 Minutes
Enterococcus hirae	EN 1276	30 Seconds
	EN 13727, EN 14561	5 Minutes
Escherichia coli	EN1276	30 Seconds
Klebsiella pneumoniae	EN1276	30 Seconds
Salmonella typhimurium	EN1276	30 Seconds
Campylobacter jejuni	EN1276	30 Seconds
MYCOBACTERIA		
Mycobacterium terrae (TB bacterium)	EN 14348	30 Seconds
FUNGI		
Candida albicans	EN 1650	30 Seconds
	EN 13624	60 Minutes
Saccharomyces cerevisiae	EN 1650	30 Seconds
Aspergillus niger	EN 13624	60 Minutes
VIRUSES		
Vaccinia virus, Ankara strain (covers all Enveloped viruses including Coronaviruses)	EN 14476	30 Seconds

Test organism(s)	Method	Results (% Reduction)		
		30 secs	1 min	2 mins
Bacteria				
Staphylococcus aureus (MRSA)	EN13727	≥99.999	≥99.999	≥99.999
Enterococcus faecalis (VRE)	EN13727	≥99.999	≥99.999	≥99.999
Pseudomonas aeruginosa	EN1276	≥99.999	≥99.999	≥99.999
Escherichia coli	EN1276	≥99.999	≥99.999	≥99.999
Enterococcus hirae	EN1276	≥99.999	≥99.999	≥99.999
Klebsiella pneumoniae	EN1276	≥99.999	≥99.999	≥99.999
Salmonella typhimurium	EN1276	≥99.999	≥99.999	≥99.999
Campylobacter jejuni	EN1276	≥99.999	≥99.999	≥99.999
Mycobacteria				
Mycobacterium terrae	EN 14348	≥99.99	≥99.999	≥99.999
Fungicidal				
Candida albicans	EN1650	≥99.999	≥99.999	≥99.999
Saccharomyces cerevisiae	EN1650	≥99.999	≥99.999	≥99.999
Virucidal				
Vaccinia virus, Ankara strain (covers all Enveloped viruses including Coronaviruses)	EN14476	≥99.999	≥99.999	≥99.999

PRODUCTION INFORMATION

Our state-of-the-art machinery allows us to produce a wide range of products in various material, size, weight, fold configurations and packaging. This includes single sachet wipes, resealable wipes, wipes in canisters or jumbo bags, as well as bottling solutions for liquid products

RESEARCH & DEVELOPMENT

Our in-house R&D and Quality Control teams, close partnership with experts and consultants, as well as quality suppliers from around the world ensures that we are not only able to consistently develop new and innovative products, improved packaging and revolutionary systems, but also assist our clients in manufacturing a product based on their contract manufacturing requirements and provide support in obtaining necessary regulatory and tests documentations.

QUALITY MANAGEMENT SYSTEM CERTIFICATION

