



1. Location: There are three possible locations that will give a successful stimulation of the common peroneal nerve (dorsiflexion). The **A** position aligns with the fibula head, while **B** and **C** are aligned with a prominent lateral tendon. The reasons for choosing an alternative fitting position include: skin sensitivity or lesions at the other position(s), better stimulation of the nerve because of anatomical variations, better comfort or ease of fitting, or the position of dressings.




2. Cleaning: Wash and dry the skin where the device will be fitted.

3^A. Fitting: Remove the film from the geko™ device and place the marker line >>●<<< over the fibula head.



3^B. Fitting: Place the marker line >>●<<< on the lateral/outer/biceps femoris tendon below the crease of the knee.



3^C. Fitting: Place the marker line >>●<<< on the lateral/outer tendon just above the crease of the knee.



4. Turning On: To turn on, use a short press of the ⊕ button.



5. Settings: Use the ⊕ button to increase the setting and ⊖ button to decrease. Increase the setting until you get a rhythmic upwards and outwards movement of the foot.



6. Switching Off: To turn it off, hold ⊖ button down for 3 seconds.



7. Removing: Remove carefully in one piece, to avoid damaging the skin.

