


# OhmniClean™

Autonomous UV-C Disinfection Solution

## UV-C Efficacy and Microbiology Brief

# How does OhmniClean Disinfect?

OhmniClean delivers **UV-C light** to high-touch surfaces through **8 powerful UV-C bulbs**.



Cover multiple surfaces in a room through **autonomous disinfection**.




# What is UV-C disinfection?

UV-C is ultraviolet light with a wavelength from 200-280 nanometers

Destroys or inactivates bacteria and viruses by deactivating the molecular bonds of their DNA and RNA.

Used to disinfect water, air and surfaces for over a century.

Valuable technique used to reduce Hospital Acquired Infections (HAIs).



# How does UV-C disinfection work?

1


The UV-C device emits UV-C light and hits the surface being disinfected.

2

Continued electromagnetic radiation from the UV-C light modifies the microorganism's genetic material and renders DNA/RNA inactive.

3

Once enough UV-C light is delivered, the ability of the microorganism to reproduce is limited, and prevents its continued spread on the surface.



# How Much UV-C Light is Needed to Disinfect Microorganisms?

Every microorganism requires certain levels of **UV-C dose** to be significantly disinfected (measured through Log Reductions).

Received  
UV-C Dose

=

UV-C Light Power x


UV-C Exposure  
Time

( UV-C Bulb Power , Distance from Surface )

Resistant pathogens require **more UV-C Dose** and require **more Time** to be disinfected.

Getting **closer** to surfaces **increases Power** and allows disinfection to be done in **less Time**.

# How Effective is OhmniClean?



OhmniClean provides a more powerful UV-C dose by moving close to all target surfaces when compared to other emitters.

The strength of disinfection by UV-C light follows the **Inverse Square Law** (a surface **twice** as far from the light source receives **25%** of the original **UV-C** strength).

# How Effective is OhmniClean?



OhmniClean has been independently tested with microbiology lab EMSL Analytical Inc. to prove disinfection efficacy:

- E.Coli: Log6 reduction
- C.Diff: > Log4 reduction
- MRSA: Log6 reduction
- VRE: > Log6 reduction
- P. aeruginosa: Log7 reduction

**\*Results achieved by OhmniClean in under 3 min**



**EMSL ANALYTICAL, INC.**

# What other Microorganisms is UV-C Effective Against?

- UV-C has been used for hospital-based disinfection for over 15 years.
- A summary of UV-C thresholds for relevant microorganisms is shown in the table below:

| Pathogen  | Disinfection Thresholds<br>(254nm UV-C)            | Sources   | Pathogen  | Disinfection Thresholds<br>(254nm UV-C)  | Sources  |
|---|--|---|---|--|--|
| Vaccina Virus<br>(Poxvirus)                               | Log 2 reduction (99%) 0.3mJ/cm <sup>2</sup>        | <a href="#">McDevitt et al. 2007</a> - Applied and Environmental Microbiology   | Clostridioides<br>difficile (C.<br>Diff)        | Log 3 reduction (99.9%): 46 mJ/cm <sup>2</sup>   | <a href="#">Liscynsky et al. 2017</a> - Infection Control & Hospital Epidemiology    |
| SARS-CoV-2 (all<br>variants)                              | Log 3 reduction (99.9%): 4-7 mJ/cm <sup>2</sup>    | <a href="#">Ma et al. 2021</a> - Journal of Applied and Environmental Microbiology  |   | Log 5 reduction (99.999%): ~100 mJ/cm <sup>2</sup>   | <a href="#">Cadnum et al. 2019</a> - Open Forum Infectious Diseases                  |
|   | Log 4 reduction (99.99%): 16-32 mJ/cm <sup>2</sup> | <a href="#">Storm et al. 2020</a> - Nature, Scientific Reports<br><a href="#">Biasin et al. 2020</a> - Nature, Scientific Reports<br><a href="#">Sabino et al. 2020</a> - Elsevier Journal of Photodiagnosis and Photodynamic Therapy, September 2020 |   |  | <a href="#">Yang et al. 2019</a> - Journal of Microbiology, Immunology and Infection |
| Staphylococcus<br>Aureus                                  | Log 4 reduction (99.99%): 10mJ/cm <sup>2</sup>     | <a href="#">McKinney et al. 2012</a> - Environmental Science and Technology   | Candida Auris<br>(C. Auris)                     | Log 3 reduction (99.99%) ~45 mJ/cm <sup>2</sup>  | <a href="#">Marijita et al. 2022</a> - MicrobiologyOpen                              |
| Methicillin-resista<br>nt Staphylococcus<br>Aureus (MRSA) | Log 4 reduction (99.99%): 31 mJ/cm <sup>2</sup>    | <a href="#">Ploydaeng et al. 2020</a> - Photodermatology, Photoimmunology, Photomedicine  | Log 4 reduction (99.99%) ~60 mJ/cm <sup>2</sup> | <a href="#">Rutala et al. 2022</a> - Infection Control & Hospital Epidemiology)  |  |
|   | Log 5 reduction (99.999%): 50 mJ/cm <sup>2</sup>   | <a href="#">Wong et al. 2015</a> - American Journal of Infection Control<br><a href="#">Yang et al. 2019</a> - Journal of Microbiology, Immunology and Infection  |   | <a href="#">Lemons et al. 2019</a> - American Journal of Infection Control<br><a href="#">Groot et al. 2019</a> - Mycoses, Diagnosis, Therapy and Prophylaxis of Fungal Diseases |  |

# What other Microorganisms is UV Effective Against?

- Log Reduction thresholds for other Microorganisms can be referenced with the comprehensive peer-reviewed paper:

## Fluence (UV Dose) Required to Achieve Incremental Log Inactivation of Bacteria, Protozoa, Viruses and Algae

Mayaleri et al. 2016, International Ultraviolet Association

| Spore   | Lamp Type      | Fluence (UV dose) (mJ/cm <sup>2</sup> ) for a given log reduction without photoreactivation |     |                   |     |    | Protocol? | Notes | Reference                  |
|---|----------------|---|-----|-------------------|-----|----|-----------|-------|----------------------------|
|   |                | 1   | 2   | 3                 | 4   | 5  |           |       |                            |
| <i>Aspergillus brasiliensis</i> (previously known as <i>Aspergillus niger</i> ) ATCC 16404 (dark culture) | LP             | 122   | 226 | 293               |     |    | yes       |       | Taylor-Edmonds et al. 2015 |
| <i>Aspergillus niger</i> ATCC 32625   | LP             | 116   | 245 | 370               | 560 |    | yes       |       | Clauß 2006                 |
| ATCC 32625  | Excimer 222 nm | 90  | 220 | 325               | 430 |    | yes       |       | Clauß 2006                 |
| <i>Bacillus anthracis</i>   |                |   |     |                   |     |    |           |       |                            |
| Sterne  | LP             | 28  | 37  | 52                |     |    | yes       |       | Nicholson & Galeano 2003   |
| Sterne  | LP             | 23  | 30  |                   |     |    | yes       |       | Blatchley III et al. 2005  |
| Arnes   | LP             | 25  | ~40 | >120 with tailing |     |    | yes       |       | Rose & O'Connell 2009      |
| 34F2 (Sterne) method: soil extract-peptone-beef extract agar  | LP             | 23  | ~40 | >120 with tailing |     |    | yes       |       | Rose & O'Connell 2009      |
| 34F2 (Sterne) method: Schaeffer's sporulation medium  | LP             | 23  | 36  | 80                |     |    | yes       |       | Rose & O'Connell 2009      |
| <i>Bacillus atrophaeus</i> ATCC 9372  | LP             | 22  | 38  | 55                | 71  |    | yes       |       | Zhang et al. 2014          |
|   | LP             | 10  | 16  | 26                | 39  |    | yes       |       | Sholtes et al. 2016        |
|   | UV-LED 260 nm  | 6   | 10  | 14                | 19  | 31 | yes       |       | Sholtes et al. 2016        |

| Bacterium  | Lamp Type      | Fluence (UV dose) (mJ/cm <sup>2</sup> ) for a given log reduction without photoreactivation |     |     |      |     |     | Protocol? | Notes | Reference                 |
|--|----------------|---|-----|-----|------|-----|-----|-----------|-------|---------------------------|
|  |                | 1   | 2   | 3   | 4    | 5   | 6   |           |       |                           |
| <i>Aeromonas hydrophila</i> ATCC 7966            | LP             | 1.1   | 2.5 | 4.0 | 5.5  | 6.9 | 8.4 | yes       |       | Wilson et al. 1992        |
| <i>Aeromonas salmonicida</i> AL 2017             | LP             | 1.5   | 2.7 | 3.1 | 5.9  |     |     | yes       |       | Litved & Landfald 1996    |
| <i>Arthrobacter nicotinovorus</i>                |                |   |     |     |      |     |     |           |       |                           |
| ATCC 49919                                       | LP             | 8   | 10  | 12  | 14   |     |     | yes       |       | Clauß 2006                |
| ATCC 49919                                       | Excimer 222 nm | 10  | 15  | 18  | 20   |     |     | yes       |       | Clauß 2006                |
| <i>Bacillus cereus</i> (veg. bacteria)           |                |   |     |     |      |     |     |           |       |                           |
| ATCC 11778                                       | LP             | 6   | 7   | 9   | 12   |     |     | yes       |       | Clauß 2006                |
| ATCC 11778                                       | Excimer 222 nm | 9   | 11  | 14  | 18   |     |     | yes       |       | Clauß 2006                |
| <i>Bacillus megaterium</i> (veg. cells) QMB 1551 | 265 nm         | 4.6   |     |     |      |     |     | no        |       | Donnellan & Stafford 1968 |
| <i>Burkholderia mallei</i>                       |                |   |     |     |      |     |     |           |       |                           |
| M3   | LP             | 1.0   | 2.4 | 3.8 | 5.2  |     |     | yes       |       | Rose & O'Connell 2009     |
| M13  | LP             | 1.2   | 2.7 | 4.1 | 5.5  |     |     | yes       |       | Rose & O'Connell 2009     |
| <i>Bruceella melitensis</i>                      |                |   |     |     |      |     |     |           |       |                           |
| ATCC 23456                                       | LP             | 2.8   | 5.3 | 7.8 | 10.3 |     |     | yes       |       | Rose & O'Connell 2009     |
| IL195  | LP             | 3.7   | 5.8 | 7.8 | 9.9  |     |     | yes       |       | Rose & O'Connell 2009     |
| <i>Burkholderia pseudomallei</i>                 |                |   |     |     |      |     |     |           |       |                           |
| ATCC 11668                                       | LP             | 1.7   | 3.5 | 5.5 | 7.4  |     |     | yes       |       | Rose & O'Connell 2009     |
| CA650  | LP             | 1.4   | 2.8 | 4.3 | 5.7  |     |     | yes       |       | Rose & O'Connell 2009     |
| <i>Brucella suis</i>                             |                |   |     |     |      |     |     |           |       |                           |

| Virus                       | Host                                | Lamp Type | Fluence (UV dose) (mJ/cm <sup>2</sup> ) for a given log reduction without photoreactivation |     |     |     |     |     | Protocol? | Notes | Reference                |
|-----------------------------|-------------------------------------|-----------|---|-----|-----|-----|-----|-----|-----------|-------|--------------------------|
|                             |                                     |           | 1   | 2   | 3   | 4   | 5   | 6   |           |       |                          |
| <i>Adenovirus</i>           |                                     |           |   |     |     |     |     |     |           |       |                          |
| Type 1 method: MPN          | PLC/ PRF/5 and HeLa cell line       | LP        | 35  | 69  | 103 | 138 |     |     | yes       |       | Nwachuku et al. 2005     |
| Type 2                      | PLC/ PRF/5                          | LP        | 40  | 78  | 119 | 160 | 195 | 235 | yes       |       | Gerbe et al. 2002        |
| Type 2                      | Human lung cell line                | LP        | 35  | 55  | 75  | 100 |     |     | yes       |       | Balister & Malley 2004   |
| Type 2                      | A549 cell line                      | LP        | 20  | 45  | 80  | 110 |     |     | yes       |       | Shin et al. 2005         |
| Type 2                      | A549 cell line                      | LP        | -30   | -60 |     |     |     |     | yes       |       | Linden et al. 2007       |
| Type 2                      | A549 cell line                      | MP        | -10   | -20 | -30 | -40 | -50 |     | yes       |       | Linden et al. 2007       |
| Type 2                      | A549 cell line                      | MP        | <240 nm blocked   |     |     |     |     |     |           |       |                          |
| Type 2                      | A549 cell line                      | LP        | -15   | -30 | -45 | -60 |     |     | yes       |       | Linden et al. 2007       |
| Type 2                      | A549 cell line                      | LP        | 8   | 31  | 50  | 80  | 117 |     | yes       |       | Eischeid et al. 2009     |
| Type 2 method: TCID50       | A549 cell line                      | LP        | 35  | 78  | 126 | 168 |     |     | yes       |       | Linden et al. 2009       |
| Type 2 method: TCID50       | A549 cell line                      | MP        | 14  | 29  | 44  | 80  | 120 |     | yes       | (3)   | Linden et al. 2009       |
| Type 2 method: cell culture | HEK293 cells human embryonic kidney | LP        | 37  | 88  | 120 |     |     |     | yes       |       | Baxter et al. 2007       |
| Type 2 adenoc 6 (VR-846)    | A-549 cell line (CCL-165)           | LP        | 42  | 83  | 124 | 166 |     |     | yes       |       | Sirikanchana et al. 2008 |

# Which Microorganisms are more Susceptible to UV-C?

- **Viruses** are generally the most susceptible to UV-C disinfection and require the least UV-C dose.
- **Spores** are generally the least susceptible to UV-C disinfection and require the most UV-C dose.

| <b>Viruses</b>  | <b>Vegetative Bacteria</b>  | <b>Mycobacteria</b>  | <b>Vegetative Fungi</b>  | <b>Bacterial Spores</b>   |
|---|---|--|--|---|
| <b>EASIER TO DISINFECT</b>  |   | <b>HARDER TO DISINFECT</b>   |  |   |
| <b>Coronavirus</b><br>Influenza viruses<br>Measles<br>Norovirus<br>Smallpox | <b><i>Staphylococcus aureus</i></b><br><i>Escherichia coli</i><br><i>Pseudomonas aeruginosa</i><br><i>Serratia marcescens</i> | <b><i>Mycobacterium tuberculosis</i></b><br><i>Mycobacterium Abscessus</i><br><i>Mycobacterium bovis</i> | <b><i>Candida Auris</i></b><br><i>Fusarium Species</i><br><i>Aspergillus Species</i> | <b><i>Clostridioides difficile</i></b><br><i>Bacillus species</i> |

# How does OhmniClean target different Microorganisms?

- Staff can easily adjust a **Disinfection Intensity** for any target pathogen
- **Disinfection intensities** can be set hospital-wide for easier staff operation

| Viruses   | Vegetative Bacteria   | Mycobacteria  | Vegetative Fungi  | Bacterial Spores                                    |
|---|---|---|---|---|
| EASIER TO DISINFECT   |   | HARDER TO DISINFECT   |   |   |
| <b>Coronavirus</b><br>Influenza viruses<br>Measles<br>Norovirus<br>Smallpox | <b>Staphylococcus aureus</b><br>Escherichia coli<br>Pseudomonas aeruginosa<br>Serratia marcescens | <b>Mycobacterium tuberculosis</b><br>Mycobacterium Abscessus<br>Mycobacterium bovis | <b>Candida Auris</b><br>Fusarium Species<br>Aspergillus Species | <b>Clostridioides difficile</b><br>Bacillus species |



# How fast does OhmniClean disinfect different Microorganisms?

- Disinfection Time depends on the: **Size of the room, Disinfection path, Disinfection intensity**
- OhmniClean provides an accurate **estimated disinfection time** before every cycle



## Example Disinfection Times for ~200sqft Patient Rooms

| Viruses             | Vegetative Bacteria          | Mycobacteria                      | Vegetative Fungi     | Bacterial Spores                |
|---------------------|------------------------------|-----------------------------------|----------------------|---------------------------------|
| Coronavirus         | <i>Staphylococcus aureus</i> | <i>Mycobacterium tuberculosis</i> | <i>Candida Auris</i> | <i>Clostridioides difficile</i> |
| EASIER TO DISINFECT |                              |                                   | HARDER TO DISINFECT  |                                 |
| <b>1-3min</b>       | <b>4-6min</b>                | <b>7-9min</b>                     | <b>8-12min</b>       | <b>9-15min</b>                  |



OhmniClean™
Operator Master Admin


PREPARE FOR DISINFECTION

Disinfection Intensity

Estimated disinfection time

# 3m

## Example Disinfections in Other Hospital Areas



# OhmniClean - Additional Academic References

[Haag et al. 2022](#) - Comparing UV-C dosages of emitter placement strategies in a community hospital setting.

[Raggi et al. 2018](#) - Clinical, operational, and financial impact of an ultraviolet-C terminal disinfection intervention at a community hospital

[Wong et al. 2015](#) - Postdischarge decontamination of MRSA, VRE, and Clostridium difficile isolation rooms using ... ultraviolet-C-emitting devices.

[Anderson et al. 2017](#) - Enhanced terminal room disinfection ... caused by multidrug-resistant organisms and Clostridium difficile ...

**OhmniClean's Autonomous UV-C applied a much larger UV-C dose throughout an entire room compared to Manual UV-C.**

**Routine UV-C resulted in 19% lower HAI, and a direct cost savings of \$1,219,878 over a 12-month period.**

**UV-C use reduced MRSA, VRE, and C.Diff bioburden in Patient Rooms 8-10x more than regular Manual Disinfection.**

**UV-C use in high-risk rooms led to a hospital-wide decrease in C. Diff and VRE patient exposure over 28 months.**