












## **LM Periodontics**

Online Category: [http://issuu.com/lminstruments/docs/lm\\_catalog](http://issuu.com/lminstruments/docs/lm_catalog)

 A yellow-handled dental scaler with a curved blade and a hook-like tip.	<p><b>Laplander</b></p> <ul style="list-style-type: none"><li>● Similar to H6-H7 but with a lateral bending</li><li>● Easy access to molar area and proximal surfaces</li></ul> <p>For calculus removal, place the tip third of the scaler's cutting edge against the tooth.</p> <p>Tilt the instrument toward the tooth to achieve 70° to 80° angle between the tooth and the blade.</p> <p>Apply lateral pressure and activate the scaler by using vertical, diagonal or horizontal pull strokes.</p> <p>To maintain control, use short 2-3 mm long strokes.</p>
 A silver-handled dental scaler with a curved blade and a hook-like tip.	<p><b>McCall 11A-12A</b></p> <p>For calculus removal, place the tip third of the scales cutting edge against the tooth.</p> <p>Tilt the instrument toward the tooth to achieve 70° to 80° angle between the tooth and the blade.</p> <p>Apply lateral pressure and activate the scaler by using vertical, diagonal or horizontal pull strokes.</p> <p>To maintain control, use short 2-3 mm long strokes.</p>
 A silver-handled dental scaler with a curved blade and a hook-like tip.	<p><b>Micro Sickle</b></p> <ul style="list-style-type: none"><li>● Especially for tight interdental spaces.</li><li>● 204SD design</li><li>● Delicate working end</li><li>● Elongated shank</li></ul>

	<p><b>Mini Kaplan</b> Same angle but with a more delicate blade than in the Crane-Kaplan.</p>
	<p><b>Mini Sickle</b></p> <ul style="list-style-type: none"> <li>• For all tooth surfaces, especially proximal surfaces.</li> <li>• Fine, slightly angled blade</li> <li>• For all tooth surfaces</li> </ul>
	<p><b>Push Scaler H4-H5</b> The chisel-like working end is used by pushing. The jacquette working end used by pulling.</p> <ul style="list-style-type: none"> <li>• 90° working surface</li> </ul>
	<p><b>Scaler H6-H7</b> The scaler H6/H7 is used for the removal of calculus and plaque in the anterior and premolar regions.</p> <p>Scalers are suitable for removal of calculus up to a pocket depth of 3mm. In deeper pockets, the chances that the sharp tip of the instrument can cause damage to the surrounding tissue is too great.</p>

 A dental scaler with a blue, double-lobed handle and a single curved working end.	<p><b>Scaler U15 Towner</b> Single-ended.</p>
 A dental sickle with a blue, double-lobed handle and two curved working ends.	<p><b>Sickle LM204S</b> Sharply angulated shank and a very fine working end.</p>
 A dental sickle with a blue, double-lobed handle and two curved working ends.	<p><b>Sickle LM23</b> Has a shorter blade and a longer, lower shank than in the Mini Sickle.</p>
 A dental sickle with an orange, double-lobed handle and a single curved working end.	<p><b>Sickle Scalers</b> Single-ended.</p>