

LM™ Extraction

LM

feel the
difference



Extract with confidence

LM LiftOut™

LM LiftOut™


Luxating instruments for atraumatic extraction

LM-LiftOut™ instruments allow tooth extraction to be done as atraumatically as possible. This is important in order to enable rapid healing and future implant placement. The thin tip of the instrument is introduced into the periodontal space and slowly advanced toward the apex of the root while twisting the instrument gently back and forth. LM-LiftOut™ is not suitable for use as an elevator.




Angulations and sizes

Straight blade



Curved blade



Extract with confidence & luxate atraumatically

Familiar LM™ ergonomics and tactile sensitivity also in extracting instruments.

LM™ Feature	Benefit
Ergonomic design	– comfortable to hold and easy to rotate
Secure LM-ErgoTouch™ surface	– light, non-slip grip
Well-balanced	– controlled grip
Lightweight	– less strain on the hands
Clear color-coding	– easy identification
LM-DuraGradeMAX® supersteel	– strong and sharp blades
Tight seam at blade juncture	– hygienic



LM™ extraction instruments rated as Best Product 2014.




Scan the QR for the complete evaluation!

LM SlimLift™

Luxating instruments designed for extremely narrow places

The LM-SlimLift™ instruments are especially designed to reach places you cannot go with traditional luxating instruments. It allows tooth extraction to be done as atraumatically as possible. This is important in order to enable rapid healing. The extremely fine tip of the instrument is introduced into the very narrow periodontal space. The instrument is slowly advanced towards the apex of the root while turning the instrument gently back and forth. The LM-SlimLift™ is not suitable for use as an elevator.

The extremely slim design of the blade allows for the most atraumatic extractions.



Thin LM-SlimLift™ instruments – a genuinely more atraumatic solution

MD, DDS, D.Med.Sc, Oral Surgeon
Peter Jungell, Helsinki, Finland



"I use the LM-SlimLift™ instruments currently for nearly all extractions. I also use my older preferred luxating instruments while the use of the traditional ones has become less frequent.

LM™'s thin luxating instrument can enter even particularly narrow periodontal spaces. When wide instruments are used, the luxating force is spread over a wider area to the sides and deeper towards the tip of the roots so that the fracturing risk of the fragile roots, in particular, is limited considerably. The majority of extractions can be performed with luxating instruments without any need for any other operative measures.

The design and material of the LM-SlimLift™ instrument make it possible to get an excellent grip

on the instrument. At first touch, the handles may even feel too big but in clinical work I have noticed they fit the hand perfectly. The instrument can be held either on the palm of your hand or held only with your fingers. It is easy to pick a suitable instrument for any clinical need from the broad size selection. The wide blades are appropriate for the extraction of the lower molars while the narrow tips are perfect for instance for the extraction of the narrow roots of the upper teeth.

I can recommend the LM-SlimLift™ instruments for all dentists: it makes sense to abandon the old and poorly functioning instruments altogether. Too often you hear of cases where the treatment of the root has taken close to two hours."

Angulations and sizes

Straight blade



Curved blade



LM TwistOut™

LM TwistOut™

Elevators for forceful extraction

LM-TwistOut™ instruments are suitable for tooth extraction in situations where strong force or torque must be applied.



“This elevator from LM™ is specifically designed for torque applications”

Bite Magazine, August 2014
by Dr. Davis Lakatos,
Complete Dentistry, Kilcoy, QLD
Australia

“I came across this elevator at the Australian Dental Congress in Melbourne last year. As we’re a rural practice and probably see a few more extractions than in a city area it’s a very handy tool.

I find that luxators are often quite fine and if a little too much pressure is applied, they tend to bend or break. This elevator from LM™ is specifically designed for torque applications.

I use this elevator for all extractions and have the full range of straight and angled models in all sizes. It’s a very sturdy, simple tool that makes extractions much easier to perform.


I had often wondered why there wasn’t a more sturdy luxator on the market but, as it turns out, this elevator was exactly what I was looking for.”

Angulations and sizes

Straight blade



Curved blade





Easy to maintain

Cassette LM 810001

LM-LiftOut™ and LM-TwistOut™ instruments are supplied in a convenient cassette that protects both the instrument and the handler during the maintenance cycle. The cassette keeps the instrument from puncturing the sterilization pouch, so sterility is more likely to be assured. The cassette can be color-coded.

Sharpening stone LM 818002



The blades, hand-finished and made from LM-DuraGradeMAX® supersteel, will stay sharp and last a long time when properly used and maintained. Sharpen the extracting instruments occasionally using this round sharpening stone as a simple and quick way to keep them in top condition and always ready to use.



Root tip instruments

Root Tip Pick LM 769-771 XSI/ES

- for the removal of root tip
- suits well in all difficult to reach areas
- available with LM-ErgoSense® (ES) and LM-ErgoMax™ (XSI) handles



LM-RootOut™ LM 812210

- designed for delicate root teasing procedures
- not recommended for more substantial root tip removal
- similar to Heidbrink root tip picks
- ergonomically designed handle for a good grip

LM Dental™

 LM-Instruments Oy

Norrbyn rantatie 8 | FI-21600 Parainen | Finland
info@lm-dental.com | www.lm-dental.com

