

Safety first!

Pro-Nox Scavanger Options







Agenda

- Why Nitrous?
- Pro-Nox highlights
- NIOSH / OSHA requirements
- Nitrous Oxyde Properties
- Side effects / Physiological effects
- Methods to minimize environmental exposure
- Who needs a scavenger system?
- Active Scavenger
- Passive Scavanger

Why Nitrous?

- Quick Onset enables a more natural alternative
 - Perfect first defense (non-invasive) tool in your tool box for pain & anxiety relief
- Patient Administered allows patients a sense of control
- Safe for Patients— many studies done worldwide
- Proven Effective proven over decades of use worldwide
- Quick Recovery allows patient to actively participate in their procedure while being relaxed
- Doesn't change/add to your patient monitoring levels





Pro-Nox System Highlights

- Super Safe, Easy to Use
- Audible & Visual Alarms
- Superior Infection Control
- Portable, lightweight, can move from room to room
- Customizable Wall Mount or Roll Stand, 3 sizes of masks, ready-made kits
- Internal On-Demand Valve
- 2 scavenger options Active with wall suction or Passive with a pump









NIOSH / OSHA Requirements

- The U.S. has the most stringent limits on acceptable exposure in the world – 25 ppm over an 8-hr shift. European standards are 100 ppm.
- Long-term exposure to levels above 500 ppm have been associated with decreased fertility.
- Use of a properly maintained scavenging system can reduce exposure more than 70%.

Controlling Exposures to Nitrous Oxide During Anesthetic Administration; DHHS (NIOSH) Publication Number 94-100 - NIOSH (National Institute for Occupational Safety & Health) OSHA - Occupational Safety and Health Administration







Nitrous Oxide Gas Properties

- Clear, colorless vapor with a slightly sweet odor and taste
- Under pressure, N2O can be stored as a liquid
- An e-cylinder can hold 1590 liters (O2 = 625 liters)
- Like oxygen, N2O supports combustion and caution must be exercised near an ignition source
- N2O Vapor is 1.5 times as heavy as air at 77° F








Nitrous Oxide – Physiological Effects Side Effects

- Effects response varies between patients
 - Decreased sensation of pain
 - A disconnection from the pain they are feeling
 - Less scared, anxious, or stressed
- Side Effects (from Likis FE et al; 2012 Vanderbilt review of the literature – mostly labor & delivery)
 - Nausea (0-28%)
 - Vomiting (0-14%) No aspiration
 - Dizziness, Light-headedness (3-23%)
 - Drowsiness (0-67%)
 - Unconsciousness (<1%)
 - Dry Mouth, Pins and needles, numbness





Methods to minimize environmental exposure


- Scavenger
- Room ventilation
- Patient Coaching
- Demand valve for intermittent applications



Methods to minimize environmental exposure

- Safety Features:
 - Internal Demand Valve

- Protection of healthcare workers
 - 2 scavenger options available
 - Nitrous Badges to test exposure





Active Scavenger – Wall Suction

Passive Scavenger & Pump-No-Wall Suction



Who needs a Scavenger? ASK QUESTIONS!

- Duration of procedures?
- Number of procedures per day?
- Size of treatment room?
- Offer measurement when in doubt
 - Badge to be mounted on stand for about 8h (1 working day)
 - Send to a laboratory
 - Get the result after 3 days




Active Scavenger

- Wall Suction
 - Part # GD 1925
 - Mount at the Rolling Stand
 - Relieve Valves / Ports (over- and underpressure)
 - Gas Reservoire









Passive Scavenger

- When no Wall Suction available
 - PRO-VAK pump
 - Part # GD-14610A
 - PRO-VAK-WK Wall Install kit

or

PRO-VAK-USK Under the Sink Install





PRO-VAK pump (Ventilator)



GD-14610A (Safety valve)

PRO-VAK-WK Wall Install kit





PRO-VAK-USK Under the Sink Install







PRO-VAK-USK Under the Sink Install



Pearls of Wisdom

- Win-Win for everyone
- Low Cost to implement
- Low Cost to run program
- Low Cost to maintain
- Patient Satisfier
- Patient Centered Pain Relief Option
- Safety
- Effective
- Positive patient experience









The BEST Part About The TREATMENT Is The



Laughing Gas!



