AVID IF₂TM

Patient Manual

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Avid IF2 Interferential Stimulator / Estimulador interferencial Avid IF2

Avid App Login Information

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ENGLISH

VQ OrthoCare is not liable for any misuse or misunderstanding of the Avid IF2™ product or operating manual. Please call your local representative if any additional assistance is required regarding this product and its operating instructions.

This patient manual is proprietary and confidential. No part of this document may be reproduced or transmitted in any form or by any means, electronic or mechanical, including photocopy, recording or information storage and retrieval system, without permission in writing from VQ OrthoCare.

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Table of Contents

Letter to the Patient	2
System Components and Accessories.....	3
Precautions and Prescription Information.....	4
Controls and Features.....	7
Starting Your Treatment Sessions	8
Normal Running Screens.....	10
Completing Your Treatment Sessions.....	12
Treatment Outcome Questions	12
Avid IF2 App.....	13
Electrode Use and Care	16
Care and Maintenance.....	17
Technical Data and Specifications	18
Equipment Classification and Device Symbol Descriptions..	20
EMC Compliance Table.....	23
U.S. and International Safety Standards.....	25
Waveform Diagrams	26
Troubleshooting.....	27
Warranty.....	28
Notes	57

Letter to the Patient

Your physician has prescribed the Avid IF2™ as part of your treatment plan. The Avid IF2 is an interferential stimulator. It will help reduce pain and other symptoms of your condition.

The Avid IF2 sends small pulses of current through electrodes* on your skin. These pulses pass through the skin and activate nerves. The Avid IF2 activates sensory nerves at low levels. This causes a tingling sensation and helps with pain. Motor nerves are activated at higher levels. This helps increase circulation.

Always follow the instructions of your physician to get the most out of your treatment. He or she knows your specific condition and the Avid IF2 settings. He or she will instruct you on the proper settings for your individual condition.

If you have any questions about the Avid IF2 or need any assistance in setting up or maintaining this equipment, or to report unexpected operations or events, contact VQ OrthoCare at 800.452.7993.

NOTE: The Avid IF2 is not a substitute for full medical evaluation and treatment. Talk to your physician or therapist if you have any questions about the Avid IF2 or your treatment plan. The Avid IF2 should only be used by the person for whom it is prescribed.

* For use only with stimulating electrodes, commonly referred to as "electrodes" throughout this manual

System Components & Accessories

The Avid IF₂™ Interferential Stimulation System is operated with:

- Avid IF₂ Interferential device
- Non-Sterile, reusable electrodes (See “Treatment Sessions” and “Electrode Use and Care”)
- Lead wire that connects the device to the electrodes

Accessories

- AC Adapter/Charger
- Adhesive Remover Towelettes
- Additional packages of non-sterile, reusable electrodes (depending on your doctor’s protocol)
- Patient Manual (this booklet)
- Protocol

Precautions & Prescription Information

CAUTION: Federal law restricts this device to sale by, or on the order of, a practitioner licensed by the law of the State in which he/she practices to use or order the use of this device.

NOTE: Please read the following prescription information carefully before using the Avid IF2™. If you have any questions regarding this information, consult with your physician or therapist before proceeding.

Interferential Stimulation can be used in the following applications:

- Symptomatic relief of post surgical and/or post traumatic acute pain
- Symptomatic relief of chronic intractable pain
- Relaxation of muscle spasms
- Maintaining and increasing range of motion
- Increasing local blood circulation

CONTRAINDICATION:

Stimulators should not be used on patients with cardiac demand pacemakers.

WARNINGS:

- The long-term effects of chronic electrical stimulation are unknown.
- Stimulation should not be applied over the carotid sinus nerves, particularly in patients with a known sensitivity to the carotid sinus reflex.
- Stimulation should not be applied over the neck or mouth; severe spasm of the laryngeal and pharyngeal muscles may occur and the contractions may be strong enough to close the airway or cause difficulty in breathing
- Stimulation should not be applied directly on the eyes.
- Electrodes should not be placed on the chest in such a configuration that would allow electrical current delivered by the device to flow through the chest because it may cause a cardiac arrhythmia.
- Stimulation should not be applied transcerebrally.
- Stimulation should not be applied over swollen, infected, or inflamed areas or skin eruptions, e.g. phlebitis, thrombophlebitis, varicose veins, etc.
- Stimulation should not be applied over, or in proximity to, cancerous lesions.
- Due to the potential for causing skin burns, this device should only be used with electrodes which limit the current density to 0.25W/cm²

WARNINGS (CONT.):

- Stimulation should not be applied transthoracically in that the introduction of electrical current may cause cardiac arrhythmias.
- Simultaneous connection of a patient to a high frequency surgical medical electronic equipment may result in burns at the site of the stimulator electrodes and possibly damage to the Avid IF2.
- Application of electrodes near the thorax may increase the risk of cardiac fibrillation.
- Any electrodes that have current densities exceeding 2mA/cm² may require the special attention of the operator.

PRECAUTIONS:

- Safety of stimulators for use during pregnancy has not been established.
- Caution should be used for patients with suspected or diagnosed heart problems.
- Stimulation should not be applied across or through the head, directly on the eyes, covering the mouth, on the front of the neck (especially the carotid sinus), or from electrodes placed directly on the chest and upper back or crossing over the heart.
- Caution should be used for patients with suspected or diagnosed epilepsy.
- Caution should be used in the presence of the following:
 - When there is a tendency to hemorrhage following acute trauma or fracture
 - Following recent surgical procedures when muscle contraction may disrupt the healing process
 - Over the menstruating or pregnant uterus
 - Over areas of the skin which lack normal sensation
- Some patients may experience skin irritation or hypersensitivity due to electrical stimulation or electrical conductive medium. The irritation can usually be reduced by using an alternative conductive medium or alternate electrode placement.
- Electrode placement and stimulation setting should be based on the guidance of the prescribing practitioner.
- Keep out of the reach of children.
- Keep lead wires away from children; the length of the lead wire could cause entanglement of a baby or small child, potentially causing strangulation.

Precautions & Prescription Information (Cont.)

PRECAUTIONS (CONT.):

- Use only with the lead wires and electrodes recommended for use by the manufacturer.
- Portable stimulators should not be used while driving, operating machinery, or during any activity in which involuntary muscle contractions may put the user at undue risk of injury.
- Do not trim electrodes as cut edges may affect even distribution of stimulation.
- Do not immerse electrodes in water.
- Electronic monitoring equipment (such as ECG monitors and ECG alarms) may not operate properly when the Avid IF2 is in use.
- Wireless communications equipment such as wireless home network devices, mobile phones, cordless telephones and their base stations, walkie-talkies can affect the Avid IF2 and should be kept at least 1 (one) meter away from the equipment. In addition, the presence of RFID and security systems (e.g., electromagnetic anti-theft systems, and metal detectors) might interfere with the operation of the Avid IF2 stimulator.

ADVERSE REACTIONS:

Skin irritation and burns beneath the electrodes have occasionally been reported with the use of electrical stimulators.


If you experience any symptoms of skin irritation or any unusual reactions while using the Avid IF2, discontinue use immediately and contact VQ OrthoCare's Patient Care Department at 800.452.7993 or consult your physician.

Controls & Features

The Avid IF2™ is a touch screen device, allowing the user to control the device on the screen and providing information regarding its operation and status. Virtually all of the functions of the device are viewed and controlled via a finger touch on the screen.

The “” button at the bottom will turn the device on. Pressing and holding this button will turn the device **OFF**.

Main Display Screen

 The battery symbol in the upper left corner will display the battery status. Press on the battery and the screen will show you the battery life available and whether or not the device is using the AC adapter.

 This button accesses the settings for the device.

The Time is displayed in the upper right corner of the screen.

The number in the middle of the screen is the **protocol**. You may have one or more to choose from, depending on your physician's prescription.

 Increases the protocol number on the screen.


 Decreases the protocol number on the screen.

 Selects the setting shown on the screen.

Settings Menu

- 1. Display:** Language, Brightness, Audible beep, Night Mode, Compliance Display.
- 2. Bluetooth:** Enable connection (Avid IF2 is bluetooth ready) and Last Connection Date. During stimulation, the bluetooth function is automatically disabled.
- 3. Time/Date:** Set clock to 12 or 24 hr, set time and date
- 4. Tools:** This is used by clinicians and technicians to program your device.

The settings function will only display when treatment is **not** running.

Starting Your Treatment Sessions

1. Prepare the skin. Electrode application sites must be clean, dry, unbroken skin surfaces. To reduce the risk of skin irritation, follow these instructions for proper skin care before treatment.
 - ✓ Wash site of electrode placement with mild soap and water before applying electrodes.
 - ✓ Rinse all soap off the skin before continuing.
 - ✓ Dry the skin thoroughly.
 - ✓ Trim excess body hair from electrode sites with scissors. Be careful not to cut (or break) the skin. **DO NOT SHAVE THE SKIN IMMEDIATELY PRIOR TO STIMULATION.** Shaving will create small cuts in the surface layer of skin.

CAUTION:

DO NOT apply lotions, oils or other ointments to the skin prior to the application of electrodes. Topical agents may increase the risk of skin irritation.

DO NOT clean the skin with anything other than soap and water. Other cleaning agents may interfere with proper conduction.

USE OF electrodes smaller than 2" round or 2x2" square is not recommended as increased power density could cause skin irritation or burn.

USE OF electrodes for long periods of time without removal may cause skin irritation.

In the event of skin irritation, discontinue use immediately and contact VQ OrthoCare's Patient Care Department at 800.452.7993 or consult your physician or therapist.

2. Inspect each electrode and lead wire before application. Do not proceed if there is any doubt about the integrity or proper function of any electrode or lead wire. Connect the lead wire to the electrodes. Make sure there is no bare metal exposed at the connection points. Do not use unnecessary force in connecting the electrodes to the lead wires. If they do not fit, use a different set of electrodes and call VQ OrthoCare to report the situation and receive a replacement set of electrodes. Save the non-working electrodes for evaluation by VQ OrthoCare's Quality Assurance Department; **DO NOT THROW THEM AWAY!**
3. Gently remove the electrodes from the plastic liner. Place the electrodes on the sites prescribed by your physician or therapist and press them firmly onto the skin. **Electrodes should be placed at least one inch apart from each other. Skin irritation or burn could result.**




4. Insert the lead wire into the Avid IF2™ with the strain relief pointing towards the bottom of the Avid IF2.
5. Turn the Avid IF2 on by pressing the on/off (⏻) button. The device will run through a self-test, displaying the VQ OrthoCare logo and then show the protocol screen.
6. Scroll through the presets by pressing the ▶ or ◀ buttons until the number indicated for the protocol is shown. Press the to select a preset, as directed by your physician or therapist. Most devices will have two default presets to choose from unless additional presets were prescribed by your physician.
7. Increase the stimulation intensity level by repeatedly pressing and releasing the + button until you feel a comfortable tingling sensation or as directed by your physician or therapist.

CAUTION: If using pain medication and/or other pain management devices, exercise caution when increasing the amplitude of the Avid IF2.

Normal Running Screens


1. The top of the screen displays the battery status, the “home” button and the time. While treatment is running, the **Home** button (🏠) will return you to the Select Protocol Screen. When treatment is off, the **Settings** icon (⚙️) will be in the top middle of the screen.
2. The middle of the screen displays the intensity of the stimulation.
3. Below the intensity is the timer. For programs with a set time, it will count down. For continuous programs, the timer will count up.
4. The bottom right corner is the increase **+** stimulation button.
5. The bottom left corner is the decrease **-** stimulation button.
6. The bottom middle is the pause button **||**.
7. These are the screens that appear during most of your treatment time. Press the **+** or **-** button to change stimulation intensity. For your safety, the keypad locks after 30 seconds of inactivity. If you see the lock in the middle of the lower line, pressing the display anywhere below the top row will show the unlock icon (🔓), which can be pressed to unlock the device) You can then increase the intensity.


8. When the intensity is equal to or greater than 10 mA, a “**Vector**” symbol (⌘) will appear under the clock and the stimulation can be **vectored**. Press the ⌘ symbol, then adjust the stimulation to the right or left until you feel it concentrated where you have the most pain. The intensity of one channel will increase while the other channel will decrease. This targets the stimulation so that it concentrates where you have the most pain. The ⌘ symbol rotates to indicate channels are vectored.

9. At any time, you can press the pause button to pause stimulation. The screen will show a pause sign at the bottom and a “**Play**” button (▶) in the middle to resume treatment. The intensity will slowly increase from zero up to its last setting. If locked, you need to touch screen once to wake the device and second time to make a change.
10. If at any time you wish to discontinue treatment, press the 🏠 button at the top. The screen will show a stop sign and a play button. Press the stop sign to stop treatment and return to the main treatment screen. Press the play button to resume treatment. Alternatively, you can press and hold the on/off button and the device will turn off completely.


Completing Your Treatment Sessions

After the treatment is complete, device will beep and the screen will show **“TREATMENT COMPLETE.”** Press the  button to return to the Select Protocol Screen or the  button at the bottom of the device to turn it off. After 10 minutes of inactivity, the device will shut off to conserve power.

Disconnect the lead wires from the electrodes. Remove the electrodes from your skin. Place them on the plastic liner and place the plastic liner in the resealable pouch. See also Electrode Use and Care on page 16.

Grip the lead wire close to the socket and pull straight out to disconnect the lead wire from the device.

Treatment Outcome Questions

Your device is programmed to ask outcome questions related to your treatment to track your progress and benefit from your treatment. Use the **+** and **-** to choose the correct answer and to select your answer. The range of answers is listed under the question.

Avid App

The Avid IF2 is Bluetooth enabled to provide important treatment and outcome related information to your prescribing physician via an app. Use of the app is voluntary and will not interfere with your treatment in any way. However, providing this information to your physician will help him or her facilitate changes to your treatment plan to promote optional outcomes for your injury as well as providing valuable information to ensure authorization for continued supplies.

To set up your Avid IF2 device with the app:

1. Download from Google Play or the Apple App Store. Search for **“Avid IF2.”** Open the app and click on **“Sign Up”** at the bottom of the screen.
2. Fill out the requested information. A user name and password is required as well as your e-mail address. An additional e-mail can be entered for your physician to receive reports. Enter your name as you want it to appear on your reports. If no name is entered, your user name will appear on reports to your physician. Enter your VQ OrthoCare account number (It is a 5-digit alpha-numeric code located on your delivery contract at the top of the page and on your patient orientation checklist, if you received one).
3. At **“Scan Serial Number,”** turn on your device, go to Settings  and then press **“Bluetooth.”** Make sure the Bluetooth is enabled. At the bottom of the screen it says **“Serial Number”** and has a bar code. Press **“Scan”** on your phone and scan the bar code on the device (Figure A).
4. Press **“Sign Up”** at the bottom. You should see a pop up that says **“Registered Successfully.”** Press **“OK”** (Figure B).
5. The app will open to the connection page and you should see your device listed. Go to the Bluetooth screen on the device and open it. The app should connect automatically or press **“Connect”** on the app.


Figure A

Figure B

- The app will send you an e-mail to confirm your sign-up. Open the e-mail, click on the link to verify your e-mail and your set up is complete (Figure C).


Figure C

If you do not verify your e-mail address, the app will not let you sign in again (Figure D).

- Keep the app open and it will connect to the device overnight. The app will **NOT** disrupt your treatment and will **NOT** shut the device off to upload information.
- If your device has not uploaded information in a week or more, you will get a message reminding you to open you app again to update information. To upload your data, press “**Send Data**.” This will take you to the Bluetooth screen. Follow the instructions below. To start a treatment instead, press “**Skip**” (Figure E).

Figure D

Figure E

To manually connect to the the app:

- Turn on your Avid IF2 Device and go to Settings . Choose Bluetooth. Make sure Bluetooth is enabled.
- Open the app on your phone. It should connect automatically. If not, press “**Connect**” next to your device. After upload, press the two plugs at the top to disconnect (Figure F).
- If you are already connected to the app and want to refresh your data, press “**Get Data**” at the bottom of the Device Information Screen.

Figure F

To manually connect to the the app:

1. The app can be used to adjust your display settings. It can only adjust the device when it is connected.
2. From the Device Information Screen (Figure G), choose “**Settings**” at the bottom of the screen. You can change the language, brightness, audible alarm, and sleep mode. These settings are also adjustable on the device.

Figure G

To view your data on the app:

1. From the device Information Screen, press Usage Data at the bottom of the screen. You will see each preset that has been used in treatment as well as the questions that are asked and answered on the device.
2. To only see treatments, press “**Usage**” at the top of the page. Press any individual treatment line to see the details of that treatment.
3. To only see outcome questions, press “**Questions**” at the top of the page. Press any individual question line to see specific answers. If a question was not asked, you will see “**UA**” for “**Unasked.**”

To disconnect your device from the app:

From the Device: Turn the Avid IF2 off by pressing and holding the .

From the App: Press the two plugs in the upper right corner of the screen.

 Connected

 Not connected

If you have any problems with the app, please contact **Customer Service** at **(800) 452-7993**.

Electrode Use and Care

WARNING: This product should be used only under the direct supervision of a registered healthcare practitioner. Product should be used according to the instructions in this manual. Improper use of this or any electrode may result in a patient burn. Due to the potential for causing skin burns, devices must be limited to 0.25W per cm². Electrodes should not be placed on the body so that the Avid IF₂'s electrical current flows across the chest. This may cause cardiac arrhythmia.

NOTE: The electrodes supplied with this Avid IF₂ have been designed specifically for use with this product. They must not be manipulated, cut, or altered in any way. This may affect the even distribution of stimulation and/or its intensity beyond the recommended 0.25W per cm² maximum.

NOTE: Do not place electrodes on cut, broken or irritated skin. If irritation occurs, consult with your physician or therapist.

To apply self-adhesive, reusable electrodes:

- ✓ Attach the electrodes to the lead wires provided with the Avid IF₂.
- ✓ Grasp an electrode by the edge and gently peel away the electrode from its plastic liner.
- ✓ Apply the electrodes to the treatment area by firmly pressing the electrode against the skin. Check for proper adhesion. Your physician or therapist will show you where to apply the electrodes.
- ✓ Always make sure electrodes are placed **at least one inch apart from each other**.

To remove the electrodes:

- ✓ Grasp the electrode by the edge and gently remove it from the skin. Return it to the plastic liner.
- ✓ Replace the electrodes on the plastic liner in the resealable pouch provided.

CAUTION: Never remove the electrode by pulling on the wire. This may damage your electrode.

CAUTION: Be careful when using the Avid IF₂ while sleeping. Wires may tangle and become loose.

Replacing Electrodes:

- ☑ Electrode life varies due to skin condition, skin preparation, storage, climate and care. For most patients, electrodes will last for approximately 7-10 applications and removals.
- ☑ If your electrodes start to dry out, gently rub a few drops of water into the gel on the electrode to remoisten.
- ☑ After repeated use, reusable electrodes lose their adhesion, deliver less stimulation and shorten battery life. **Replace electrodes that do not stick uniformly or cause discomfort.**

Care and Maintenance

CAUTION: DO NOT attempt to clean the Avid IF2 while in use.

Maintenance, care and cleaning of the Avid IF2 is quick and easy when these instructions are followed:

Case: Clean the outside of the Avid IF2 case as needed with a damp cloth and mild soap. Never immerse the Avid IF2 in any liquid. This could seriously damage the internal electronics. Avoid using cleaning fluids or solvents to remove stains and dirt. This may damage the casing. **DO NOT** use abrasive materials on the LCD viewing area.

Lead wire: Do not bend or crimp the lead wire. Do not use the lead wire if it is moist or wet – allow to air dry for at least one hour. Always remove the lead wire by firmly grasping the connector. If possible, tuck in excess length of the lead wires to avoid interference with normal activity. These wires only need to be cleaned if they become dirty. Wipe with a damp cloth and dry thoroughly.

AC Adapter: The Avid IF2 is powered by a non-removable, rechargeable, internal battery or an AC Adapter. When the AC Adapter is plugged in, it will charge the internal battery and supply power for operating the device. Use only the AC Adapter provided (Globtek Model # GTM96180-1811-2.0). Use of other AC Adapters may damage the Avid IF2. The AC adapter serves as the only means to disconnect the device from the mains supply.

Technical Data & Specifications

Physical:

Power source:	Internal, rechargeable, 7.4V. Lithium Ion battery. Optional VQ Orthocare AC Adapter.	
Size:	2.90" (73 mm) x 4.95" (125 mm) x 1.00" (25 mm) without belt clip	
Weight:	6.8 oz. (193 g)	
Lead wire length:	22.0" (55.88 cm) unstretched	
Electrical tolerances:	All electrical specifications are $\pm 10\%$ into a 500 ohm load for each channel and using the AC Adapter, unless stated otherwise	
Treatment timer:	1-30 minutes in one-minute increments 35 minutes - 8 hours in five-minute increments or continuous	
Compliance timer:	0 min - 999 days, in 1-min increments	
Number of preset programs:	1 - 20, according to physician orders	
Cycles:	1 - 99	
Next preset:	1 - 20 or OFF	
Expected Service Life	2.5 years	
Operating Conditions	Temperature:	5° C to 40° C
	Relative Humidity:	15% - 93%
	Atmospheric Pressure:	700 to 1060 hPa
Storage and Transport	Temperature:	-25° C to +70° C
	Relative Humidity:	30% - 93%
	Atmospheric Pressure:	700 to 1060 hPa

* At the end of the product service life, you may dispose of the device according to your local governing ordinances or recycling plans. Contact your local authorities to determine the proper method for disposal since this is electronic equipment contains a lithium-ion battery and potentially biohazardous parts and accessories.

Interferential Stimulation:

Pulse amplitude:	0 - 50mA
Pulse shape:	With four (4) electrodes – symmetric, biphasic square wave With two (2) electrodes – symmetric, biphasic rectangular wave
Pulse width:	With four (4) electrodes – 125 μ s per phase With two (2) electrodes – 5 μ s to 125 μ s per phase

Output Configuration:

Quadrapolar: Ch. 1 delivers a fixed-frequency carrier pulse train. Ch. 2 delivers a user-adjustable frequency that consists of the carrier frequency, plus a beat frequency. Four (4) electrodes are required.

Bipolar: Ch. 1 delivers a pulse train that is pre-modulated inside the unit. Ch. 2 is not used and has no stimulation pulses. Two (2) electrodes are required.

Carrier frequency: 4000 Hz

Adjustable frequency:

Continuous mode range: 4001 Hz - 4150 Hz, in 10 Hz increments (beat frequency of 1 - 150 Hz)

Sweep modes	1	1 - 10 Hz over a 30-sec period
	2	1 - 150 Hz over a 30-sec period
	3	80 - 150 Hz over a 30-sec period
	4	1 - 10 Hz over a 15-sec period
	5	1 - 150 Hz over a 15-sec period
	6	80 - 150 Hz over a 15-sec period

Frequency shift mode range: 4000.6 Hz - 4240.0 Hz

- All displayed values are $\pm 10\%$ of maximum setting
- Load impedance of less than 200 Ohms or greater than 2000 Ohms will result in lower total power output from the device

Technical Data & Specifications (Cont.)

Pulse coordination modes:

Continuous: Stimulation is provided at the user-selected beat frequency.

- 11:** After one second of stimulation, the beat frequency abruptly changes from 60% of user setting to 160% of user setting. One second later, the beat frequency abruptly changes back to 60% of user setting. The cycle then starts over.
- 616:** After six seconds of stimulation, the beat frequency abruptly changes from 60% of user setting to 160% of user setting. Six seconds later, the beat frequency abruptly changes back to 60% of user setting. The cycle then starts over.
- 6/6:** The beat frequency gradually changes from 60% of user setting to 160% of user setting over a six-second period. Over the next six seconds, the beat frequency gradually changes back to 60% of user setting. The cycle then starts over.

Equipment Classification and Device Symbol Descriptions

Symbol	Meaning
	Stimulator, lead wires and electrodes are type BF Applied parts
	Indicates the manufacturer
	Indicates the manufacturer's serial number
	Refer to instructional manual
	Warning symbol

Symbol	Meaning
	Storage humidity range
	Storage temperature range
	Storage pressure range
<p data-bbox="181 610 380 683">IP 22</p>	The Avid IF2 is rated IP22, indicating that the device is protected from harmful effects against ingress of fingers or similar objects and that vertically dripping water shall have no harmful effect when the enclosure is tilted at an angle up to 15° from its normal position
	Do not throw in trash
	Requires prescription
	Does not contain dry natural rubber latex
	Class II device
	Indoor use (AC adapter)
	Keep dry

Equipment Classification and Device Symbol Descriptions (cont.)

Equipment Classification

- IEC 60529 enclosure rating IP22
- Equipment not suitable for use in an oxygen rich environment
- Suitable for continuous use

⚠ WARNING: No modification of this equipment is allowed. There are no user-serviceable parts in the Avid IF2. Attempting to modify or repair this device by anyone other than the manufacturer could damage the device and will void the warranty.

- The Avid IF2 is not affected by lint, dust or light (including sunlight).
- The Avid IF2 should be kept away from pets and children
- The Avid IF2 should be kept away from areas with pests

⚠ WARNING: This product can expose you to chemicals including acrylamide, which are known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to [www. P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).

Symbols not from Standards

Symbol	Meaning
	On/off symbol
	Decrease value or stimulation
	Increase value or stimulation
	Home button
	Settings
	Display
	Tools
	Bluetooth
	Battery
	Select
	Vector

Symbol	Meaning
	Indicates another program will follow the current one when timer finishes counting down.
	Timer
	Date/time set
	Continue stimulation/ play/start
	Pause
	Stop
	Lock/unlock
	Delete

EMC Compliance Table

Medical electrical equipment needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided below.

Portable and mobile RF communications equipment can affect medical electrical equipment. Shortwave or microwave therapy equipment may produce instability in the stimulator output if operated in close proximity to the Avid IF2 (e.g. within 1 meter).

- Simultaneous connection of a patient to a high frequency surgical medical electronic equipment may result in burns at the site of the Avid IF2 electrodes and possible damage to the stimulator.
- Operation in close proximity (e.g. 1 m) to a shortwave or microwave therapy medical electronic equipment may produce instability in the Avid IF2 output.

Use of lead wires and other accessories other than those supplied by VQ OrthoCare may result in increased emissions or decreased immunity of the Avid IF2.

Guidance and manufacturer's declaration - electromagnetic emissions		
The Avid IF2™ is intended for use in the electromagnetic environment specified below. The customer or user of the Avid IF2 should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic Environment - Guidance
RF Emissions	Group1	The Avid IF2 uses RF energy only for its internal function. RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF Emissions	Class B	The Avid IF2 is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Complies	
Voltage fluctuations / flicker emissions IEC 61000-3-3	Complies	

Guidance and manufacturer's declaration - electromagnetic immunity			
The Avid IF2™ is intended for use in the electromagnetic environment specified below. The customer or user of the Avid IF2 should assure that it is used in such an environment.			
Immunity test	IEC 60601 Test Level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact ± 15 kV air	± 8 kV contact ± 15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.

Guidance and manufacturer's declaration - electromagnetic immunity			
Immunity test	IEC 60601 Test Level	Compliance level	Electromagnetic environment - guidance
Electrical fast transient / burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input / output lines	Not Applicable	
Surge IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	Not Applicable	Not Applicable
Voltage dips, short	< 5 % U _T	Not Applicable	Not Applicable
Interruptions and voltage variations on power supply input lines IEC 61000-4-11	(>95 % dip in U _T) for 0,5 cycle 40 % U _T (60 % dip in U _T) for 5 cycles 70 % U _T (30 % dip in U _T) for 25 cycles < 5 % U _T (>95 % dip in U _T) for 5 s		
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	Not Applicable	Not Applicable

Note: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a. Field strengths from fixed transmitters such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Avid IF2 is used exceeds the applicable RF compliance level above, the Avid IF2 should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the Avid IF2.

b. Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

NOTE: U_T is the a.c. mains voltage prior to application of the test level.

U.S. and International Safety Standards

When traveling with the Avid IF2, users may be requested to provide documentation identifying it as a medical device and substantiating its electrical safety. The following facts are provided for this purpose.

The Avid IF2, is a FDA Class II Medical Device prescribed by a licensed physician or other medical professional. This Device has been cleared by the U.S. Food and Drug Administration under 510(k)K183692.

The Avid IF2 meets or exceeds federal and international standards for electrical safety, including:

- IEC 60601-1 : 1988 Standards for Medical Electrical Equipment; Part 1: General Requirements for Safety, Amendment A1 : 1991, Amendment A2 : 1995
- IEC 60601-2-10 : 1987 Standards for Medical Electrical Equipment; Part 2: Particular Requirements for the Safety of Nerve and Muscle Stimulators, Amendment A1 : 1991
- ANSI/AAMI NS4-1985 American National Standards for Transcutaneous Electrical Nerve Stimulator; Environmental Testing
- IEC60601-1-2:2007 Medical electrical equipment - General requirements for basic safety and essential performance - Collateral standard: Electromagnetic compatibility - Requirements and tests
- IEC60601-1-11:2010 General requirements for basic safety and essential performance - Collateral Standard: Requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment

Further information regarding this medical device may be obtained directly from VQ OrthoCare at 800.444.1456 (inside the US) or 949.261.3000 (outside of the US).

Waveform Diagrams


Figure 1 – Premodulated IF stimulation

This drawing shows the output waveform in the single-channel premodulated IF stimulation mode at maximum amplitude (50mA). The waveform is shown with a purely resistive load of 500 ohms. The two signals created in two-channel mode are internally added together and output on a single channel.



Figure 2 – Two-channel IF stimulation

This drawing shows the output waveform in the two-channel IF stimulation mode at maximum amplitude (50mA). The waveform is shown with a purely resistive load of 500 ohms. Ch. 1 is shown; Ch. 2 is identical except that the period is slightly less since the rate is slightly greater.

Troubleshooting

Problem	Solution
<p>Display shows “warning sign”</p> 	<p>In most cases, this screen will appear when something becomes disconnected (e.g. lead wires or electrodes). Check the connections between the electrodes and your skin, the electrodes and the lead wire and the lead wire and the device. This error may also appear if the electrodes are too dry to adhere properly to your skin. In this case, replace the electrodes with a new set. When finished, press (▶) to proceed as directed on the screen. If none of these corrects the problem, contact Patient Care or your healthcare provider.</p>
<p>The device will not increase the intensity</p>	<p>This is most likely caused by the “lock out” feature that prevents accidental increases to your intensity. Press the (🔒) button to unlock the keypad, then (+) to increase the intensity to the desired level.</p>
<p>Stimulation is weak, even at maximum intensity</p>	<p>This is most likely due to low battery charge. Plug in the AC Adapter. This can also be caused by a decrease in sensation over certain areas of the skin (e.g., those without good nerve conduction, immediately after ice application, after pain medications are taken or injected). Use caution when increasing intensity in these situations.</p>
<p>An audible beep and flashing battery with the percentage (< 10%)</p>	<p>The battery is low. Plug in the AC Adapter.</p>

Warranty

VQ OrthoCare warrants the Avid IF2™ to be free from electrical circuit defects in workmanship and materials for a period of one (1) year upon delivery to the customer. VQ OrthoCare will repair or replace, at its facilities, any Avid IF2 unit found to have become defective within the warranty period. Do not attempt to repair or modify the device in any way as there are no user-serviceable parts and it may void the warranty.

This warranty does not apply to accessories; specifically lead wires, electrodes, tapes, gel or AC Adapter, nor does it apply to Avid IF2 units that have been damaged due to misuse, or repaired or altered by parties other than VQ OrthoCare at its facilities. This warranty is in lieu of any other warranties expressed or implied. No person or entity is authorized to bind VQ OrthoCare to any representation of warranty other than those specifically set forth herein.

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© 2019 VisionQuest Industries, Inc.