


Rezzimax® Tuner Pro II Rezzipes: Dogs, Cats & Horses



TUNER PRO
II





Introduction

This Rezzipe Book is a step-by-step guide intended to assist you in using the Rezzimax Tuner Pro II on your pet. Follow each step as you would a recipe. It will provide you with basic techniques for supporting your pet's natural healing. It is divided nose-to-tail based on where pain or tension is currently felt. After following these, feel free to develop your own Rezzipes and share with us what works best for your pet! For additional support, [view our online training videos](#).

**Do not leave pets unsupervised with the Tuner Pro II or any of its accessories. Do not allow them to ingest any part of the Tuner Pro II or its accessories.

Disclaimer: These techniques are not intended to diagnose, cure, or treat diseases.

Table of Contents

The Basics	4
Dog	
Whole Body	6
Head & Neck	10
Back, Chest, & Abdomen	14
Gastrointestinal System	16
Legs	16
Hind End, Pelvis, & Tail	18
Appendix	41
Cat	
Whole Body	20
Head & Neck	22
Back, Chest, & Abdomen	24
Gastrointestinal System	26
Legs	27
Hind End, Pelvis, & Tail	29
Appendix	41
Horse	
Whole Body	31
Head & Neck	33
Back, Chest, & Abdomen	35
Gastrointestinal System	36
Legs	37
Hind End, Pelvis, & Tail	39
Appendix	42

The Basics:

It is important to start introducing your pet to the device slowly and with rewards to create a positive experience. Show them that it is okay for them by using the Tuner Pro II on yourself first. They will map your response to the device and will be more receptive to it. Follow the instructions below or watch the [Rezzimax Pro Pain Tuner II orientation video](#) to learn how to use your device.

1. Use the arrow buttons to change between manual levels. To switch from extra gentle to normal levels and vice versa, with Tuner Pro II turned off, press the two circle buttons, and let go (it flashes to signify it has changed).
2. For each step begin with lowest levels and increase intensity as tolerated. The 4 modes are blue, green, red, and white that change to the next by pressing the button between the two arrows. Completion of a mode is signaled by two distinct pulses. Avoid use of modes until your pet can tolerate high intensity.
3. As you work with your pet, watch for the following:
 - Licking the lips
 - Yawning
 - Chewing
 - Taking a deep breath
 - Blowing out of the mouth (horses)

When your animal shows these responses but are not showing signs of anxiety or stress, then these responses signify that their nervous system is becoming more in tune with their body. This primes their body for healing.



Dogs



Cats




Horses



Dogs: Whole Body [Click to Watch Video](#)

This Rezzipe supports your dog's natural healing response.



Step 1:

- Tuner: Accessory Tongs
- Mode/Level: Low Level
- Duration: 2 min
- Pressure: Light
- Location: Neck and Shoulders

Attach tongs to Tuner and place tongs on both sides of the neck with the Tuner resting against the chest for one minute. Then place tongs on shoulders with Tuner on the spine for one minute.



Step 2:

- Tuner: Accessory Tongs
- Mode/Level: Low Level
- Duration: 1-2 min
- Pressure: Light
- Location: Hips

Place tongs on hips with Tuner on the spine, pointed toward the head.

Step 3:

- Tuner: Wing
- Mode/Level: Low Level
- Duration: 1-2 min per leg
- Pressure: Light-Moderate
- Location: Paw and Leg

Remove Tongs, place Clear Tip, and scrape along the top of a paw with wing for 15 seconds. Then press the tip of the wing in between each toe for a few seconds. Scrape along all remaining sides of the paw and leg moving up to the body. Repeat on all legs.



Step 4:

- Tuner: Accessory Clear Tip
- Mode/Level: Low Level
- Duration: 1-2 min
- Pressure: Light-Moderate
- Location: Tail

Place Clear Tip on areas around and on the tail.



Step 5:

- Tuner: Wing
- Mode/Level: Low Level
- Duration: 2 min
- Pressure: Light-Moderate
- Location: Back of Head and Spine

Scrape along each side of the spine from neck to tail for one minute. Then scrape from behind the ears and back of head to base of the neck for one minute.



Step 6:

- Tuner: Wing
- Mode/Level: Low Level/Blue
- Duration: 2 min
- Pressure: Light-Moderate
- Location: Shoulders


Scrape all sides of the shoulders



Step 7:

- Tuner: Tuner Base
- Mode/Level: Lowest Level
- Duration: 2-3 min
- Pressure: Light
- Location: Head

Place the base of the Tuner under the chin, then move it along the cheeks, jaw, and under the ear to the back of the head.



Head and Neck [Click to Watch Video](#)

Additional Items Needed



- Rubber/latex glove or other cover to keep the Clear Tip Attachment clean
- Recommended Treat: Peanut butter, wet food, pumpkin puree, non-fat plain Greek yogurt, or other similar consistency
- (Optional) May substitute the Blue Brain with a soft toy or pillow if needed



Step 1:

- Tuner: Accessory Blue Brain
- Mode/Level: Low Level
- Duration: 3-4 min
- Pressure: Light
- Location: Head

Hold the Blue Brain on your dog's back for 2 minutes or until they relax. Then hold the Blue Brain around the face and back of the head for 10-20 seconds on each area.

Step 2:

- Tuner: Accessory Clear Tip
- Mode/Level: Low Level
- Duration: 1-2 min
- Pressure: Light
- Location: Tongue

Remove the Blue Brain and place the Clear Tip. Cover the tip with a rubber or latex glove. Apply a soft treat to the covering. Allow your dog to lick at the treat while the Tuner is running.



Step 3:

- Tuner: Wing
- Mode/Level: Low Level
- Duration: 2-3 min
- Pressure: Light on Face, Light-Moderate on Neck and Shoulders
- Location: Face, Neck, and Shoulders

Scrape along the cheeks, top and back of the head, then the neck for 15-30 seconds each. Then scrape around the shoulder area for one minute.



Step 4:

- Tuner: Tuner Base
- Mode/Level: Lowest Level
- Duration: 2-3 min
- Pressure: Light
- Location: Head

Place base of the Tuner under the chin, then move it along the cheeks, jaw, and under the ear to the back of the head.





Step 5:

- Tuner: Accessory Tongs
- Mode/Level: Low Level
- Duration: 2 min
- Pressure: Light
- Location: Jaw and Neck


Replace the clear tip with the tongs. Place tongs on both sides of the jaw for one minute. Then place tongs on both sides of the base of the neck with Tuner resting down the middle of the chest for one minute.



Step 6:

- Tuner: Accessory Tongs
- Mode/Level: Low Level
- Duration: 2-3 min
- Pressure: Light
- Location: Head and Neck

Move tongs on top of the head from front to back 7-8 times. Then place the tongs on back of the head behind the ears with Tuner resting down the spine for 1-2 minutes.



Step 7:

- Tuner: Accessory Tongs
- Mode/Level: Low Level/Blue
- Duration: 2-3 min
- Pressure: Light-Moderate
- Location: Neck and Shoulders

Place tongs on base of neck with Tuner on the spine for one minute. Then place tongs on shoulders with Tuner on the spine for one to two minutes.




Back, Chest, and Abdomen [Click to Watch Video](#)



Step 1:

- Tuner: Accessory Tongs
- Mode/Level: Low Level
- Duration: 1-2 min
- Pressure: Light-Moderate
- Location: Base of the Neck

Place tongs on base of neck with Tuner on the spine.



Step 2:

- Tuner: Accessory Tongs
- Mode/Level: Lowest Level
- Duration: 2 min
- Pressure: Light
- Location: Base of Neck

Place Tuner with one prong at the base of the neck and the other below the throat.



Step 3:

- Tuner: Accessory Tongs
- Mode/Level: Blue
- Duration: 2 min
- Pressure: Moderate
- Location: Shoulders and Spine


Place tongs behind shoulders with Tuner on the spine for one minute, then move Tuner from head to tail 4-5 times.



Step 4:

- Tuner: Accessory Tongs
- Mode/Level: Low Level
- Duration: 2 min
- Pressure: Moderate
- Location: Pectoral Muscles

Place tongs on pectoral muscles with Tuner resting on the chest.



Step 5:

- Tuner: Wing
 - Mode/Level: Blue
 - Duration: 3-5 min
 - Pressure: Moderate-Firm
 - Location: Shoulders, Back, Chest, and Abdomen
- Remove tongs and scrape around the shoulders, back, sides, chest, and abdomen with the wing.

Gastrointestinal System [Click to Watch Video](#)



Step 1:

- Tuner: Accessory Tongs
- Mode/Level: Blue
- Duration: 3-5 min
- Pressure: Light-Moderate
- Location: Abdomen

Place tongs at the base of the ribcage below the diaphragm.

Legs [Click to Watch Video](#)

Note: Use caution with dogs that are sensitive to their paws and tail being touched.

Step 1:

- Tuner: Wing
- Mode/Level: Low Level
- Duration: 1-2 min per leg
- Pressure: Light-Moderate
- Location: Paw and Leg

Scrape along the top of a paw with wing for 15 seconds. Then press the tip of the wing in between each toe for a few seconds. Scrape along all remaining sides of the paw and leg moving up to the body. Repeat on all legs.



Step 2:

- Tuner: Accessory Tongs
- Mode/Level: Low Level
- Duration: 1 minute per paw
- Pressure: Light-Moderate
- Location: Paw

Place tongs on the Tuner and put one tong on the top and bottom of the paw for 30 seconds. Then place tongs on sides of the paw for 30 seconds. Repeat on remaining paws.



Step 3:

- Tuner: Accessory Tongs
- Mode/Level: Low Level
- Duration: 1-2 min
- Pressure: Light-Moderate
- Location: Shoulders

Place tongs on shoulders with Tuner on spine.



Step 4:

- Tuner: Accessory Tongs
- Mode/Level: Low Level
- Duration: 2 min
- Pressure: Light-Moderate
- Location: Hips and Tail

Place tongs on hips with Tuner on spine pointed toward the head for one minute. Rotate Tuner so that tongs remain on hips, but the Tuner is now against the tail. Hold the tail against the Tuner for one minute.



Hind End, Pelvis, and Tail [Click to Watch Video](#)

Note: Use caution with dogs who are sensitive to their tail being touched.



Step 1:

- Tuner: Accessory Tongs
- Mode/Level: Low Level
- Duration: 2 min
- Pressure: Light-Moderate
- Location: Hips and Tail

Place tongs on hips with Tuner on spine pointed toward the head for one minute. Rotate Tuner so that tongs remain on hips, but the Tuner is now against the tail. Hold the tail against the Tuner for one minute.



Step 2:

- Tuner: Wing
- Mode/Level: Low Level
- Duration: 2-3 min
- Pressure: Moderate-Deep
- Location: Hips

Remove tongs and replace with Clear Tip. Scrape the underside and around the hip using the wing. Hold in tense areas. Repeat with the other hip.



Step 3:

- Tuner: Accessory Clear Tip
- Mode/Level: Low Level
- Duration: 1-2 min
- Pressure: Light
- Location: Tail

Place Clear Tip around the base of the tail hitting points all around it for 10-15 seconds each, then hold on the tip of the tail for 15-30 seconds.



Cats: Whole Body [Click to Watch Video](#)


This Rezzipe supports your cat's natural healing response.



Step 1:

- Tuner: Accessory Tongs
- Mode/Level: Low Level
- Duration: 2 min
- Pressure: Light
- Location: Neck and Shoulders


Attach tongs to Tuner and place tongs on both sides of the neck with the Tuner resting against the chest for one minute. Then place tongs on shoulders with Tuner on the spine for one minute.



Step 2:

- Tuner: Accessory Tongs
- Mode/Level: Low Level
- Duration: 1-2 min
- Pressure: Light
- Location: Hips

Place tongs on hips with Tuner on the spine, pointed toward the head.



Step 3:

- Tuner: Wing
- Mode/Level: Low Level
- Duration: 1-2 min per leg
- Pressure: Light-Moderate
- Location: Paw and Leg

Remove Tongs, place Clear Tip, and scrape along the top of a

paw with wing for 15 seconds. Then press the tip of the wing in between each toe for a few seconds. Scrape along all remaining sides of the paw and leg moving up to the body. Repeat on all legs.

Step 4:

- Tuner: Accessory Clear Tip
- Mode/Level: Low Level
- Duration: 1-2 min
- Pressure: Light-Moderate
- Location: Tail

Place Clear Tip on areas around and on the tail.



Step 5:

- Tuner: Wing
- Mode/Level: Low Level
- Duration: 2 min
- Pressure: Light-Moderate
- Location: Back of Head and Spine

Scrape along each side of the spine from neck to tail for one minute. Then scrape from behind the ears and back of head to base of the neck for one minute



Step 6:

- Tuner: Wing
- Mode/Level: Low Level/Blue
- Duration: 2 min
- Pressure: Light-Moderate
- Location: Shoulders

Scrape all sides of the shoulders.



Step 7:

- Tuner: Tuner Base
- Mode/Level: Lowest Level
- Duration: 1-2 min
- Pressure: Light
- Location: Head

Place base of Tuner under the chin for 10-15 seconds, then move it along the cheeks, jaw, and under the ear to back of the head.



Note: You may find a spot where your cat becomes jumpy. In this case, hold the Tuner Pro II lightly to that area until they calm down and the area moves freely. Use caution for cats that are sensitive to their face, feet, or tail being touched.

Head and Neck [Click to Watch Video](#)

Additional Items Needed

- Rubber/latex glove or other cover to keep the Clear Tip Attachment clean
- Recommended Treat: Wet food, Churu®, or other similar consistency



Step 1:

- Tuner: Accessory Blue Brain
- Mode/Level: Low Level
- Duration: 3-4 min
- Pressure: Light
- Location: Head

Hold the Blue Brain on your cat's back for 2 minutes or until they relax. Then hold the Blue Brain around the face and back of the head for 10-20 seconds on each area.



Step 2:

- Tuner: Accessory Clear Tip
- Mode/Level: Low Level
- Duration: 1-2 min
- Pressure: Light
- Location: Tongue

Remove Blue Brain and place Clear Tip. Cover the tip with a rubber or latex glove. Apply a soft treat to the covering. Allow your cat to lick at the treat while the Tuner is running.



Step 3:

- Tuner: Wing
- Mode/Level: Low Level
- Duration: 2-3 min
- Pressure: Light on Face, Light-Moderate on Neck and Shoulders
- Location: Face, Neck, and Shoulders

Scrape along the cheeks, top and back of the head, then the neck for 15-30 seconds each. Then scrape around the shoulder area for one minute.

Step 4:

- Tuner: Tuner Base
- Mode/Level: Lowest Level
- Duration: 2-3 min
- Pressure: Light
- Location: Head


Place base of the Tuner under the chin, then move it along the cheeks, jaw, and under the ear to the back of the head.



Step 5:

- Tuner: Accessory Tongs
- Mode/Level: Low Level
- Duration: 2 min
- Pressure: Light
- Location: Jaw and Neck

Replace clear tip with tongs. Place tongs on both sides of the jaw for one minute. Then place tongs on both sides of base of the neck with Tuner resting down the middle of the chest for one minute.



Step 6:

- Tuner: Accessory Tongs
- Mode/Level: Low Level
- Duration: 2-3 min
- Pressure: Light
- Location: Head

Move tongs on top of head from front to back 7-8 times. Then place tongs on both sides of back of head behind the ears with Tuner resting down the spine for 1-2 minutes.



Step 7:

- Tuner: Accessory Tongs
- Mode/Level: Low Level/Blue
- Duration: 2-3 min
- Pressure: Light-Moderate
- Location: Neck and Shoulders


Place the tongs on the base of the neck with the Tuner on the spine for about one minute. Then place the tongs on the shoulders with the Tuner on the spine for one minute.



Variation to Step 1:

- Tuner: Accessory Tongs
- Mode/Level: Low Level
- Duration: 1-2 min
- Pressure: Light
- Location: Head

Place the Tongs Attachment. Hold Tuner and allow your cat to rub their muzzle, cheeks, and other areas around the face on tongs



Back, Chest, and Abdomen [Click to Watch Video](#)



Step 1:

- Tuner: Accessory Tongs
- Mode/Level: Low Level
- Duration: 1-2 min
- Pressure: Light
- Location: Base of the Neck

Place tongs on the base of the neck with the Tuner on the spine.



Step 2:

- Tuner: Accessory Tongs
- Mode/Level: Lowest Level
- Duration: 2 min
- Pressure: Light
- Location: Base of Neck

Place Tuner with one prong at back of neck and other below the throat.

Step 3:

- Tuner: Accessory Tongs
- Mode/Level: Low/Blue
- Duration: 2 min
- Pressure: Moderate
- Location: Shoulders and Spine

Place tongs behind shoulders with Tuner on the spine for one minute, then move the Tuner from head to tail 4-5 times



Step 4:

- Tuner: Accessory Tongs
- Mode/Level: Low Level
- Duration: 2 min
- Pressure: Moderate
- Location: Pectoral Muscles

Place tongs on pectoral muscles with Tuner resting on the chest.



Step 5:

- Tuner: Wing
- Mode/Level: Blue
- Duration: 3-5 min
- Pressure: Moderate-Deep
- Location: Shoulders, Back, Chest, and Abdomen

Remove Tongs and scrape around the shoulders, back, sides, chest, and abdomen with the wing.



Step 6:

- Tuner: Body of the Tuner
- Mode/Level: Low/Blue
- Duration: 3-5 min
- Pressure: Moderate-Deep
- Location: Sides and Abdomen

Hold Tuner flat along your cat's side. Repeat on the other side, and then on the abdomen.



Gastrointestinal System [Click to Watch Video](#)



Step 1:

- Tuner: Accessory Tongs
- Mode/Level: Low/Blue
- Duration: 3-5 min
- Pressure: Light-Moderate
- Location: Abdomen

Place tongs on both sides of base of ribcage below the diaphragm and base of Tuner on the lower abdomen.



Step 2:

- Tuner: Accessory Blue Brain
- Mode/Level: Low/Blue
- Duration: 2-3 min
- Pressure: Light
- Location: Abdomen

Remove tongs and replace with Blue Brain. Hold Blue Brain against your cat's abdomen.

Legs [Click to Watch Video](#)

This Rezzipe supports your cat's natural healing response in the legs.

Note: Use caution with cats that are sensitive to their paws and tail being touched.

Step 1:

- Tuner: Wing
- Mode/Level: Low Level
- Duration: 1-2 min per leg
- Pressure: Light-Moderate
- Location: Paw and Leg

Scrape along the top of a

paw with wing for 15 seconds. Then press the tip of the wing in between each toe for a few seconds. Scrape along all remaining sides of the paw and leg moving up to the body. Repeat on all legs.



Step 2:

- Tuner: Accessory Tongs
- Mode/Level: Low Level
- Duration: 1 minute per paw
- Pressure: Light-Moderate
- Location: Paw

Place tongs on the Tuner and

put one tong on the top and bottom of the paw for 30 seconds. Then place tongs on sides of the paw for 30 seconds. Repeat on remaining paws.



Step 3:

- Tuner: Accessory Tongs
- Mode/Level: Low Level
- Duration: 1-2 min
- Pressure: Light-Moderate
- Location: Shoulders

Place tongs on shoulders with Tuner on the spine.



Step 4:

- Tuner: Accessory Tongs
- Mode/Level: Low Level
- Duration: 2 min
- Pressure: Light-Moderate
- Location: Hips and Tail

Place tongs on hips with Tuner on the spine pointed toward the head for one minute. Rotate Tuner so that tongs remain on hips, but Tuner is now against the tail. Hold tail against the Tuner for one minute.



Hind End, Pelvis, and Tail [Click to Watch Video](#)

Note: Use caution with cats who are sensitive to their tail being touched.

Step 1:

- Tuner: Accessory Tongs
- Mode/Level: Low Level
- Duration: 2 min
- Pressure: Light-Moderate
- Location: Hips and Tail

Place tongs on hips with Tuner on the spine pointed toward the head for one minute. Rotate Tuner so that tongs remain on hips, but Tuner is now against the tail. Hold tail against the Tuner for one minute.



Step 2:

- Tuner: Wing
- Mode/Level: Low Level
- Duration: 2-3 min
- Pressure: Moderate-Firm
- Location: Hips

Remove tongs and place Clear Tip. Scrape underside and around the hip using the wing. Hold in tense areas. Repeat with other hip.



Step 3:

- Tuner: Accessory Clear Tip
- Mode/Level: Low Level
- Duration: 1-2 min
- Pressure: Light
- Location: Tail

Place the Clear Tip around base of the tail hitting points all around it, then hold on the tip of the tail.




Horses: Whole Body [Click to Watch Video](#)

This Rezzipe supports your horse's natural healing response.

Note:

- Remember to start by showing the horse it is okay by using the Tuner on yourself prior to using it on the horse.
- You may find a spot where your horse becomes jumpy. In this case, hold the Tuner Pro II lightly to that area until they calm down and the area moves freely.
- Use caution for horses that are sensitive to their face, hooves, or tail being touched.



Step 1:

- Tuner: Accessory Tongs
- Mode/Level: Low Level
- Duration: 2 min
- Pressure: Light-Moderate
- Location: Neck and Shoulders

Attach tongs to Tuner and place tongs on both sides of neck with Tuner resting against the chest for one minute. Then place tongs on shoulders with Tuner on the spine for one minute.

Step 2:

- Tuner: Accessory Tongs
- Mode/Level: Low Level
- Duration: 1-2 min
- Pressure: Light-Moderate
- Location: Hips

Place tongs on hips with Tuner on the spine, pointed toward the head.



Step 3:

- Tuner: Wing
- Mode/Level: Low Level
- Duration: 2-3 min per leg
- Pressure: Moderate-Deep
- Location: Leg

Remove tongs, place Clear Tip, and scrape along all sides of the leg starting at the coronet band and moving up to the body using the wing. Repeat on all legs



Step 4:

- Tuner: Accessory Clear Tip
- Mode/Level: Low Level
- Duration: 1-2 min
- Pressure: Light-Moderate on tail, Firm-Deep around base of tail
- Location: Tail

Place Clear Tip on areas around and on the tail.



Step 5:

- Tuner: Wing
- Mode/Level: Low Level
- Duration: 3-5 min
- Pressure: Deepest Tolerated
- Location: Back of Head and Spine

Scrape along each side of the spine from neck to tail for about two to three minutes. Then scrape from behind the ears and back of the head to the base of the neck for one to two minutes.



Step 6:

- Tuner: Wing
- Mode/Level: Low Level/Blue
- Duration: 2-3 min
- Pressure: Firm-Deep
- Location: Shoulders

Scrape all sides of the shoulders



Step 7:

- Tuner: Tuner Base
- Mode/Level: Lowest Level
- Duration: 1-2 min
- Pressure: Light
- Location: Head

Place base of Tuner under chin, then move it along the cheeks, jaw, and under the ear to back of the head.



Head and Neck [Click to Watch Video](#)

Note: Use caution with horses sensitive to their ears being touched.

Additional Items Needed

- Rubber/latex glove or other cover to keep the Clear Tip Attachment clean
- Recommended Treat: Carrot Puree, Apple Puree, or other similar consistency
- Dish Towel/Small Rope - See Variation on Page 34

Step 1:

- Tuner: Wing/Accessory Clear Tip
- Mode/Level: Low Level
- Duration: 3-4 min
- Pressure: Light
- Location: Head

Scrape your horse's back for 2 minutes or until they relax. Then hold the Clear Tip around the face and back of the head for 10-20 seconds on each area.



Step 2:

- Tuner: Accessory Clear Tip
- Mode/Level: Low Level
- Duration: 1-2 min
- Pressure: Light
- Location: Tongue

Cover the Clear Tip with a rubber or latex glove. Apply a soft treat to the covering. Allow your horse to lick at the treat while the Tuner is running.



Step 3:

- Tuner: Wing
- Mode/Level: Low Level
- Duration: 2-3 min
- Pressure: Light-Moderate on Face, Firm-Deep on Neck and Shoulders
- Location: Face, Neck, and Shoulders

Scrape along the cheeks, top and back of the head, then the neck for 15-30 seconds each. Then scrape around the shoulder area for one minute.



Step 4:

- Tuner: Accessory Tongs
- Mode/Level: Low Level
- Duration: 2 min
- Pressure: Light
- Location: Jaw and Pectorals


Replace the clear tip with the tongs. Place tongs under jaw for one minute. Then place tongs on the pectorals with Tuner resting down the middle of the chest for one minute.



Step 5:

- Tuner: Accessory Tongs
- Mode/Level: Low Level
- Duration: 2-3 min
- Pressure: Light
- Location: Head and Spine

Move tongs on top of head from front to back 7-8 times. Then place tongs on back of head behind the ears with Tuner resting down the spine for one to two minutes.



Step 6:

- Tuner: Accessory Tongs
- Mode/Level: Low Level/Blue
- Duration: 2-3 min
- Pressure: Light-Moderate
- Location: Neck and Shoulders

Place tongs on base of the neck with Tuner on the spine for about one minute. Then place the tongs on the shoulders with the Tuner on the spine for about one minute.



Step 7:

- Tuner: Tuner & Belt Pack
- Mode/Level: Comfortable Level
- Duration: 1 min per side
- Pressure: N/A
- Location: Neck

Remove tongs attachment. Place Tuner in the Belt Pack with buttons facing the mesh back. Clip the Belt Pack so Tuner is laying against the side of the neck near the shoulder. Repeat on the other side.



Optional Variation for Step 2:

- Tuner: Base of Tuner
- Mode/Level: Low Level
- Duration: 1-2 min
- Pressure: Light
- Location: Tongue

Being careful not to hurt the mouth, place a dish towel/rope in the mouth behind the front teeth. Tie Tuner to ends of the towel outside of the mouth and let the vibration run through the tongue and mouth.



Back, Chest, and Abdomen [Click to Watch Video](#)



Step 1:

- Tuner: Accessory Tongs
- Mode/Level: Low Level
- Duration: 1-2 min
- Pressure: Light
- Location: Base of the Neck

Place tongs on base of the neck with Tuner on the spine



Step 2:

- Tuner: Accessory Tongs
- Mode/Level: Low Level
- Duration: 2 min
- Pressure: Moderate
- Location: Pectorals

Place tongs on upper pectoral muscles with Tuner resting on the chest for one minute. Then move tongs down to lower pectoral muscles for one minute.

Step 3:

- Tuner: Wing
- Mode/Level: Blue
- Duration: 3-5 min
- Pressure: Moderate-Deep
- Location: Back, Shoulders, Chest, and Abdomen

Remove tongs and scrape around the shoulders, back, sides, chest, and abdomen with the wing.

- Note: If your horse's back seems particularly painful or if your horse suffers from Swayback/Lordosis, then scrape back and forth between chest and lower abdominal area for 1-2 minutes.



Gastrointestinal System [Click to Watch Video](#)

Step 1:

- Tuner: Accessory Tongs
- Mode/Level: Low/Blue
- Duration: 3-5 min
- Pressure: Light-Moderate
- Location: Abdomen

Place tongs on both sides of base of ribcage below the diaphragm and base of Tuner on the abdomen.



Step 2:

- Tuner: Wings
- Mode/Level: Med-High
- Duration: 2-3 min
- Pressure: Moderate-Firm
- Location: Abdomen

Scrape from chest to lower abdomen with the wing.



Legs [Click to Watch Video](#)

Note:

- Use caution with horses that are sensitive to their hooves and tail being touched.
- If there is a particularly tight or painful area and your horse does not tolerate it being touched, try scraping the opposite side of the leg for a minute and re-evaluate. Repeat as needed.

Step 1:

- Tuner: Wing
- Mode/Level: Low Level
- Duration: 2-3 min per leg
- Pressure: Moderate-Firm
- Location: Leg

Scrape along all sides of the leg starting at the coronet band and moving upward using the wing. Repeat on all legs



Step 2:

- Tuner: Accessory Tongs
- Mode/Level: Low Level
- Duration: 1-2 min
- Pressure: Light-Moderate
- Location: Shoulders

Place tongs on shoulders with Tuner on the spine.




Step 3:

- Tuner: Accessory Tongs
- Mode/Level: Low Level
- Duration: 3 min
- Pressure: Light-Moderate
- Location: Hips and Tail

Place tongs on hips with Tuner on the spine pointed toward the head for one minute. Rotate Tuner so tongs remain on hips, but Tuner is now toward tail for one minute. Then hold tail against Tuner for one minute.

- Note: Can strap Tuner to tail using provided strap.



Step 4 (Optional):

- Tuner: Strap & Tuner
- Mode/Level: Low Level/Blue
- Duration: 3-5 min
- Pressure: Light
- Location: Leg

Strap Tuner to leg on desired location



Hind End, Pelvis, and Tail [Click to Watch Video](#)


Note: Use caution with horses who are sensitive to their tail being touched.

Step 1:

- Tuner: Accessory Tongs
- Mode/Level: Low Level
- Duration: 3 min
- Pressure: Light-Moderate
- Location: Hips and Tail

Place tongs on hips with Tuner on the spine pointed toward the head for one minute. Rotate Tuner so tongs remain on hips, but Tuner is now toward tail for one minute. Then hold tail against Tuner for one minute.

- Note: Can strap Tuner to tail using provided strap.



Step 2:

- Tuner: Wing
- Mode/Level: Low Level
- Duration: 2-3 min
- Pressure: Light
- Location: Hips

Remove tongs and place Clear Tip Attachment. Scrape underside and around hip using the wing. Hold in tense areas. Repeat with other hip.



Step 3:

- Tuner: Accessory Clear Tip
- Mode/Level: Low Level
- Duration: 1-2 min
- Pressure: Light
- Location: Tail

Place Clear Tip around base of tail hitting points all around it for 30 seconds to one minute, then hold on the tip of the tail for 30 seconds to one minute.



Appendix

Note: The techniques included may be helpful for the following conditions. Page numbers are where the Rezzipe starts. Complete the whole Rezzipe for best results.

Appendix - Dog

Anxiety	6
Arthritis	
• Neck	10
• Spine	14
• Legs	16
• Osteochondritis Dissecans (OCD)	
Behavioral Issues	6, 10
Brain Diseases/Issues	6, 10
Degenerative Myelopathy	6, 16, 18
Diarrhea	16
Digestive Problems	16
Emotional Issues	6
Head Tilt	10
Hunching	14
Hydrocephalus	10
Irritable Bowel Syndrome (IBS)	16
Immune Mediated Conditions	6, 10
• Lupus	10
• Pemphigus	10
Inflammation	
• Abdominal	16
• Inflammatory Bowel Disease (IBD)	
• Back or Chest	14
• General	6
• Head or Neck	10
• Salivary Glands	10
• Leg	16
• Tail	18
Intervertebral Disc Disease (IVDD)	14
• Cervical IVDD	10
Ligament Disease	16

Limber Tail Syndrome	18
Limping	16
Masticatory Myositis	10
Muscular Dystrophy	14
Nasal Tumor	10
Neurological Issues	6, 10
Pain	
• Abdominal	16
• Back or Chest	14
• General	6
• Head or Neck	10
• Leg	16
• Nerve	6
• Tail	18
Pelvic Issues	18
Pharyngeal Paralysis	10
Reverse Sneezing	10
Rhinitis	10
Sialadenosis	10
Seizures	10, 16
Soft Tissue Injury	16
Stomatitis	10
Stress	6
Temporomandibular Joint Disorder (TMD)	10
Training Help	6, 10
Trauma	
• Abdominal	16
• Back or Chest	14
• Emotional	6
• General	6
• Head or Neck	10

• Leg	16
• Tail	18
Vestibular Syndrome	10
Weakness	
• Back or Chest	14
• General	6
• Head or Neck	10
• Leg	16
• Tail	18

Appendix - Cat

Anxiety	20
Arthritis	
• Neck	22
• Spine	24
• Legs	27
• Osteochondritis Dissecans (OCD)	
Behavioral Issues	20, 22
Cystitis	24, 26
Diarrhea	26
Digestive Problems	26
Dystonia	24
Emotional Issues	20
Epaxial Muscle Disease	22
Feline Audiogenic Reflex	22
Feline Orofacial Pain Syndrome	22
Irritable Bowel Syndrome (IBS)	26
Immune-Mediated Disease	22
Inflammation	
• Abdominal	26

• Inflammatory Bowel Disease (IBD)		• Osteochondritis Dissecans (OCD)		Swelling	
• Back or Chest	24	Cervical Dysfunction	33	• Abdominal	36
• General	20	Colic	36	• Back or Chest	35
• Head or Neck	22	Colitis	36	• Head or Neck	33
• Leg	27	Degenerative Joint Disease	37	• Leg	37
• Tail	29	Degenerative Suspensory Ligament		Tension	31
• Intervertebral Disc Disease (IVDD)	24	Desmitis (DSLDD)	37	Temporomandibular Joint Disorder	
• Cervical IVDD	22	Diarrhea	36	(TMD)	33
Ligament Disease	27	Digestive Problems	36	Tongue Paralysis	33
Myoclonus/Muscle Contraction		Emotional Issues	31	Trauma	
Disease	24	Enteritis	36	• Abdominal	36
Pain		Facial Paralysis	33	• Back or Chest	35
• Abdominal	26	Foundering	37	• Emotional	31
• Back or Chest	24	Head Shaking Caused by		• General	31
• Head or Neck	22	Neuropathic Pain Response	33	• Head or Neck	33
• Leg	27	Immune-Mediated Disease	33	• Leg	37
• Nerve	20	Inflammation		• Tail	39
• Tail	29	• Abdominal	36	Weakness	
Pelvic Issues	29	• Inflammatory Bowel Disease (IBD)		• Back or Chest	35
Rhinitis/Sinusitis	22	• Back or Chest	35	• General	31
Seizures	22, 27	• General	31	• Head or Neck	33
Soft Tissue Injury	27	• Head or Neck	33	• Leg	37
Stomatitis/Resorptive Lesions	22	• Leg	37	• Tail	39
Stress	20	• Tail	39		
Temporomandibular Joint Disorder		Laminitis	37		
(TMD)	22	Ligament Disease	37		
Trauma		Muscle Strain	35, 37		
• Abdominal	26	Neck Soreness	33		
• Back or Chest	24	Neurological Issues	31		
• Emotional	20	Pain			
• General	20	• Abdominal	36		
• Head or Neck	22	• Back or Chest	35		
• Leg	27	• Head or Neck	33		
• Tail	29	• Neck Joint Pain	33		
Vestibular Disease/Head Tilt	22	• Leg	37		
Weakness		• Nerve	31		
• Back or Chest	24	• Tail	39		
• General	20	Pelvic Issues	39		
• Head or Neck	22	Seizures	33, 37		
• Leg	27	Sinusitis	33		
• Tail	29	Soft Tissue Injury	37		
		Soreness			
		• Neck	33		
		• Back or Chest	35		
		• Leg	37		
Appendix - Horse		Stenotic Myelopathy (Wobbler			
Anxiety	31	Syndrome)	33		
Arthritis		Stress	31		
• Neck	33	Swayback/Lordosis	35		
• Spine	35				
• Leg	37				

