

Determining Fatty Acid – Membrane Partition Coefficients

Synopsis

ADIFAB2 can be used to determine the partition coefficient of a fatty acid between a membrane phase and aqueous solution. Simply add fatty acid to a cuvette containing ADIFAB2 and a membrane and measure the fluorescence ratio (550/457 upon excitation at 375 nm).

Procedure

For details on measuring the ADIFAB2 ratio and calculating [FFA] see [Determining the ADIFAB2 Ratio](#). To determine R_0 , add ADIFAB2 and a membrane of known lipid concentration to a cuvette containing measuring buffer (20 mM HEPES, 140 nM NaCl, 5 mM KCl, 1 mM Na_2HPO_4 , at pH 7.4), and measure the fluorescence ratio (550/457 nm). Titrate the solution with fatty acid aliquots of known concentration and measure the R value after each addition —be sure to wait at least 5-10 minutes for equilibrium before measuring R . For each R measured, calculate the partition coefficient, K_p , using Eq. (1):

$$K_p = \frac{V_a}{V_m} \cdot \frac{([\text{FA}]_{\text{total}} - [\text{FFA}])}{[\text{FFA}]} \quad (1)$$

where $[\text{FA}]_{\text{total}}$ is the *total* fatty acid concentration in the cuvette after each addition, $[\text{FFA}]$ is the *free* fatty acid concentration, and V_a and V_m are the volumes of the aqueous and membrane phases, respectively. (Note: to accurately report $[\text{FA}]_{\text{total}}$, measure the concentration of the fatty acid stock according to [Determining the Concentration of Fatty Acid in an Aqueous Solution](#).) With sufficient membrane present, no correction for wall binding to the cuvette walls is necessary because for typical conditions >95% of the fatty acid will be bound to the membrane; very little will be free, bound to ADIFAB2 or bound to the walls.

Example

For an example experiment using *ADIFAB* see [Determining Fatty Acid – Membrane Partition Coefficients:ADIFAB](#).