

# IMC EMERGENCY DEPARTMENT - X-RAY #2 - UPGRADE

INTERMOUNTAIN HEALTHCARE  
5125 SOUTH COTTONWOOD STREET, MURRAY UT 84107



524 SOUTH 600 EAST  
SALT LAKE CITY, UT 84102  
801.575.8800 | VCBO.COM



**owner**

STEVE BROWN  
INTERMOUNTAIN HEALTHCARE  
38 S. STATE ST. # 2100  
S.L.C. UT 84111  
Steve.Brown@imail.org  
801.718.2411

**architect**

JEFF PINEGAR  
VCBO ARCHITECTURE  
524 SOUTH 600 EAST  
SALT LAKE CITY, UT 84102  
jpinegar@vcbo.com  
801.575.8800

**structural engineer**


JESSICA CHAPPELL  
BEAWELEY ENGINEERS  
675 EAST 500 SOUTH #400  
S.L.C. UT 84102  
jchappell@beaweley.com  
801.486.3883

**electrical engineer**

PETER JOHANSEN  
SPECTRUM ENGINEERS  
324 S STATE ST., Suite 400  
S.L.C. UT 84111  
pj@\_spectrum-engineers.com  
801.328.5151

PERMIT SET / BID SET  
06/08/2020






REV	DATE	DESCRIPTION



VCBO NUMBER: 20370  
DATE: 06/08/2020

**IMC EMERGENCY DEPARTMENT - X-RAY #2 - UPGRADE**  
INTERMOUNTAIN HEALTHCARE  
5125 SOUTH COTTONWOOD STREET, MURRAY UT 84107  
PERMIT SET / BID SET



**A2 BUILDING 5 - LEVEL LL-1 - OVERALL PLAN**  
SCALE: 1/32" = 1'-0"

**A4 PLAN - SITE - OVERALL**  
SCALE: 1" = 80'-0"



PROJECT AREA @ LEVEL LL-1:  
REFER TO SHEETS AD110.0 &  
A110. PLANS FOR ADDITIONAL  
INFORMATION IN THIS AREA

KEYED NOTES	
KEY VALUE	KEYNOTE TEXT
215.1	EXISTING SINK, REMOVE AND DISPOSE
215.2	EXISTING FAUCET, REMOVE AND DISPOSE
215.3	EXISTING FOOT PEDALS, REMOVE AND DISPOSE
219.2	REMOVE AND DISPOSE EXISTING PLASTIC LAMINATE COUNTERTOP, PREPARE FOR NEW SOLID SURFACE COUNTERTOP (PER PLANS)
219.3	REPLACE CABINET FINISHES (BASE AND UPPER CABINETS), PROTECT BODY BOX IN PLACE, REPAIR AS NEEDED, SALVAGE AND PROTECT HARDWARE FOR REUSE
220.2	REMOVE AND DISPOSE EXISTING CEILING TILES AND CEILING GRIDS, PREPARE FOR NEW FINISHES
220.3	REMOVE, SALVAGE AND REINSTALL DIFFUSERS AND FIRE SPRINKLER HEADS AS NEEDED, PROTECT AS NEEDED
220.4	REMOVE AND DISPOSE EXISTING METAL STRUT SYSTEM ABOVE CEILING, PREPARE STRUCTURE FOR NEW METAL STRUT SYSTEM
225.0	REMOVE AND DISPOSE EXISTING WALL PROTECTION AND CORNER GUARDS AROUND THIS ROOM, PROTECT WALLS FOR FUTURE FINISHES, REPAIR AS NEEDED PRIOR TO NEW INSTALLATION
240.0	DEMO FLOOR COVERING
250.0	SLAB LEVELING WORK IN THIS AREA, REFER TO STRUCTURAL AND CARESTREAM DRAWINGS FOR ADDITIONAL INFORMATION
260.0	EXISTING TABLE ANCHOR SYSTEM TO BE REMOVED, PATCH AND FILL AS NEEDED FOR FUTURE FINISHES
261.0	EXISTING FLOOR ELECTRICAL BOX TO BE REMOVED, PATCH AND FILL AS NEEDED FOR FUTURE FINISHES
270.0	EXISTING SOAP DISPENSER, REMOVE AND SALVAGE FOR REUSE, REPAIR AND REPAINT WALL AS REQUIRED
270.1	EXISTING PAPER TOWEL DISPENSER, REMOVE AND SALVAGE FOR REUSE, REPAIR AND REPAINT WALL AS REQUIRED

**GENERAL DEMOLITION NOTES**


- FIELD VERIFY DIMENSIONS AND CONDITIONS INCLUDING EXISTING UTILITIES PRIOR TO BIDDING. BRING DIFFERING DIMENSIONS AND CONDITIONS TO ARCHITECT'S ATTENTION PRIOR TO BIDDING.
- A HAZARDOUS MATERIAL SURVEY IS AVAILABLE FROM THE OWNER. ABATEMENT MUST BE COMPLETED PRIOR TO DEMOLITION OF BUILDINGS OR CONSTRUCTION ELEMENTS.
- PROVIDE DUSTPROOF ENCLOSURES AT PERIMETER OF CONSTRUCTION & DEMOLITION FOR PROTECTION OF ADJACENT SPACES.
- COORDINATE MAINTENANCE OF FIRE EGRESS FOR OCCUPANTS IN EXISTING BUILDING WITH THE OWNER AND FIRE MARSHAL. PROVIDE NECESSARY TEMPORARY WALLS OR ENCLOSURES, EMERGENCY LIGHTS, ETC., FOR THE DURATION OF CONSTRUCTION.
- BRING TO ARCHITECT'S ATTENTION EXISTING CONDITIONS THAT PRESENT ANY CODE VIOLATIONS, INCORRECT CONSTRUCTION OR SAFETY PROBLEMS.
- MAINTAIN EXISTING FIRE RATINGS, AND ASSOCIATED FIRE PROTECTION SYSTEMS (I.E. FIRE SPRINKLERS AND FIRE ALARM SYSTEMS) THROUGHOUT CONSTRUCTION. COORDINATE ANY INTERRUPTION TO THESE SYSTEMS WITH THE OWNER AND FIRE MARSHAL. PROVIDE FIRE WATCH REQUIREMENTS ASSOCIATED WITH INTERRUPTIONS TO THESE SYSTEMS.
- PROTECT EXISTING STRUCTURE, FINISHES, AND SITE ELEMENTS NOT SCHEDULED FOR DEMOLITION. RESTORE DAMAGED ITEMS TO THEIR ORIGINAL CONDITION OR REPLACE AT CONTRACTOR'S EXPENSE.
- REMOVE AND DISPOSE SELECTIVE DEMOLITION MATERIAL PER CITY REQUIREMENTS.
- SALVAGE MATERIAL WHERE INDICATED. REMOVE ITEMS FROM CURRENT LOCATIONS & PREPARE FOR TRANSPORT BY THE OWNER.

**GENERAL PLAN DEMOLITION NOTES**

- REFER TO ELECTRICAL AND MECHANICAL PLANS FOR REQUIRED ADDITIONAL DEMOLITION
- MAINTAIN EXISTING FIRE RATINGS THROUGHOUT CONSTRUCTION
- DO NOT DISTURB EXISTING FIRE RATED ELEMENTS INCLUDING FIREPROOFING, PATCH/REPAIR DAMAGED OR DISTURBED ITEMS.
- AFTER DEMOLITION, PRIOR TO FINISH, PATCH AND REPAIR EXISTING WALLS TO PROVIDE SMOOTH SURFACE SUITABLE FOR PAINTING OR WALL COVERING.
- PATCH & LEVEL EXISTING CONCRETE SLABS FOR NEW FINISHES WITH FLOOR LEVELING COMPOUND.
- FIELD VERIFY AND COORDINATE SAW CUTTING OF THE CONCRETE FLOOR SLAB WITH PLUMBING AND ELECTRICAL.
- REPLACE SLAB AND TRENCH BY COMPACTING CLEAN GRAVEL IN 8 INCH LIFTS. DRILL #4 EPOXY-COATED REBAR INTO EXISTING SLAB @ 12 INCHES OC. POUR SLAB TO PROVIDE A SMOOTH EVEN FLOOR.
- WHERE ELECTRICAL CIRCUIT CONTINUITY IS INTERRUPTED, BUT MUST BE MAINTAINED, MAKE NECESSARY MODIFICATIONS TO MAINTAIN CIRCUIT INTEGRITY
- REMOVE ELECTRICAL BOXES BEHIND RELOCATED MILLWORK AND CAP AS REQUIRED.
- CAP EXISTING DUCT WORK FOR DUST CONTROL.

**DEMOLITION LEGEND**

- HALF-TONE LINE DENOTES ITEMS TO REMAIN
- - - DASHED LINE DENOTES ITEMS TO BE DEMOLISHED
- [ ] AREA TO REMAIN UNDISTURBED DURING CONSTRUCTION
- [X] AREA OF DISTURBANCE DURING CONSTRUCTION (COORDINATE W. OWNER)
- [ ] AREA OUT OF SCOPE OF WORK



REV	DATE	DESCRIPTION


VCBO NUMBER: 20370  
DATE: 06/08/2020

**IMC EMERGENCY DEPARTMENT - X-RAY #2 - UPGRADE**  
INTERMOUNTAIN HEALTHCARE  
5125 SOUTH COTTONWOOD STREET, MURRAY UT 84107  
PERMIT SET / BID SET




X-RAY ROOM ENTRANCE  
(Existing signage to remain)

DEMOLITION PLAN & DEMO REFLECTED CEILING PLAN - ENLARGED  
**AD110.1**



**A1 DEMO PLAN - ENLARGED**  
SCALE: 3/8" = 1'-0"  
NORTH



**A4 DEMO REFLECTED CEILING PLAN - ENLARGED**  
SCALE: 3/8" = 1'-0"  
NORTH

REFER TO STRUCTURAL AND ELECTRICAL DEMOLITION DRAWINGS FOR ADDITIONAL INFORMATION


KEYED NOTES	
KEY VALUE	KEYNOTE TEXT
306.2	REINFORCED CONCRETE SLAB ON GRADE, 15" THICK, TO MATCH EXISTING
330.0	LEVEL WITH A TOTAL TOLERANCE OF 1/8" (SHARED AREA)
945.0	REPAIR AND PREPARE WALL FOR NEW PAINT, AS NEEDED (AT EXISTING TALL CABINET)
950.0	PAINT EXISTING DOOR/WINDOW FRAME

### GENERAL NOTES

- IT IS THE CONTRACTOR'S RESPONSIBILITY TO REVIEW AND COORDINATE THE WORK OF ALL SUB-CONTRACTORS, TRADES AND SUPPLIERS WITH THE REQUIREMENTS OF THE CONTRACT DOCUMENTS BEFORE COMMENCING CONSTRUCTION, AND TO ASSURE THAT ALL PARTIES ARE AWARE OF ALL REQUIREMENTS, REGARDLESS OF WHERE THE REQUIREMENTS OCCUR IN THE CONTRACT DOCUMENTS, WHICH MIGHT AFFECT THE WORK OF THAT PARTY.
- AS PART OF THE CONTRACTOR'S RESPONSIBILITY TO COORDINATE THE WORK OF ALL SUB-CONTRACTORS, TRADES AND SUPPLIERS, THE CONTRACTOR SHALL ENDEAVOR TO IDENTIFY AND NOTIFY THE ARCHITECT OF ANY CONFLICTS BETWEEN THE WORK OF DIFFERENT PARTIES AT THE EARLIEST POSSIBLE DATE SO AS TO ALLOW REASONABLE AND ADEQUATE TIME FOR THE CONFLICT TO BE RESOLVED WITHOUT DELAYING THE WORK. ALL DEVIATIONS FROM THAT WHICH IS REQUIRED BY THE CONTRACT DOCUMENTS MUST BE APPROVED IN ADVANCE BY THE ARCHITECT.
- THE ARCHITECTURAL DRAWINGS ESTABLISH AND COORDINATE THE FINISHED APPEARANCE AND EXACT LOCATION OF ALL EXPOSED ELEMENTS OF THE WORK OF ALL THE TRADES, INCLUDING THAT WORK WHICH IS ILLUSTRATED PRIMARILY ON DRAWINGS OF OTHER DISCIPLINES. QUANTITIES ARE TO BE PROVIDED AS SHOWN ON DRAWINGS OF OTHER DISCIPLINES BUT LOCATIONS SHOWN ON OTHER DRAWINGS ARE SCHEMATIC UNLESS OTHERWISE NOTED ON THE ARCHITECTURAL DRAWINGS. THE ARCHITECTURAL DRAWINGS TAKE PRECEDENCE FOR THE FINISHED APPEARANCE AND EXACT LOCATION OF ALL PARTS OF THE WORK.
- EXCEPT WHERE DIRECTED TO PLACE ITEMS OF WORK AT THE APPROXIMATE LOCATION SHOWN, DO NOT SCALE DRAWINGS FOR DIMENSIONAL INFORMATION. ALL ELEMENTS OF THE DRAWINGS MAY NOT BE DRAWN TO EXACT SCALE. ALL DIMENSIONS REQUIRED ARE SHOWN OR MAY BE DERIVED FROM THOSE SHOWN ON THE FLOOR PLANS, DETAIL PLANS, ELEVATIONS, SECTIONS, DETAILS, SCHEDULES AND SPECIFICATIONS. IF DIMENSIONS ARE NOT PRESENT, THE ARCHITECT IS TO BE NOTIFIED SO THAT A CLARIFICATION CAN BE ISSUED.
- CONTRACTOR TO FOLLOW CURRENT ANSI 117-1 STANDARDS AS REPRESENTED ON SHEET C301, GENERAL ACCESSIBILITY GUIDELINES. NOTIFY ARCHITECT IF THE DESIGN DRAWINGS CONFLICT WITH THIS SHEET.

### DIMENSION NOTES

- ALL DIMENSIONS ARE TO CENTER OF STUD WALL OR FACE OF CONCRETE, MASONRY OR ROUGH OPENING UNLESS NOTED OTHERWISE. WHERE THE END OF A WALL IS INDICATED THE DIMENSION IS TO THE FINISH SURFACE OF THE WALL END.
- UNLESS DIMENSIONED OTHERWISE, THE DIMENSION FROM THE BUCK OF A DOOR FRAME IS TO BE 4" TO THE WALL CORNER.
- EXCEPT WHERE DIRECTED TO PLACE ITEMS OF WORK AT THE APPROXIMATE LOCATION SHOWN, DO NOT SCALE DRAWINGS FOR DIMENSIONAL INFORMATION. ALL ELEMENTS OF THE DRAWINGS MAY NOT BE DRAWN TO EXACT SCALE. ALL DIMENSIONS REQUIRED ARE SHOWN OR MAY BE DERIVED FROM THOSE SHOWN ON THE FLOOR PLANS, DETAIL PLANS, ELEVATIONS, SECTIONS, DETAILS, SCHEDULES AND SPECIFICATIONS. IF DIMENSIONS ARE NOT PRESENT, THE ARCHITECT IS TO BE NOTIFIED SO THAT A CLARIFICATION CAN BE ISSUED.



REV	DATE	DESCRIPTION

VCBO NUMBER: 20370  
DATE: 06/08/2020


## IMC EMERGENCY DEPARTMENT - X-RAY #2 - UPGRADE

INTERMOUNTAIN HEALTHCARE  
5125 SOUTH COTTONWOOD STREET, MURRAY UT 84107


PERMIT SET / BID SET

LEVEL LL-1 SLAB PLAN +  
ANNOTATED & DIMENSION  
PLAN - ENLARGED


**A110.1**



**A2 SLAB PLAN - ENLARGED**  
SCALE: 1/2" = 1'-0"



**A5 LL-1 ANNOTATED & DIMENSION PLAN - ENLARGED**  
SCALE: 1/2" = 1'-0"




LEGEND - FINISH					
MARK	PRODUCT DESCRIPTION	MANUFACTURER	NAME	COLOR	COMMENTS
<b>FLOOR</b>					
F1	LINOLEUM	MANNINGTON	BIOSPEC SR 67361	SANDRIFT	MAIN COLOR
F2	LINOLEUM	MANNINGTON	BIOSPEC SR 67203	FLAX	ACCENT COLOR (AT X-RAY AREA)
F3	LINOLEUM	MANNINGTON	BIOSPEC SR 67369	BEDROCK	ACCENT COLOR (AT DOOR & AROUND X-RAY)
<b>BASE</b>					
B1	LINOLEUM	MANNINGTON		SANDRIFT / FLAX / BEDROCK	4" COVED - CARRY OVER SAME COLOR AS FLOORING WHERE IT OCCURS
<b>PAINT</b>					
P1	LATEX PAINT - EGGSHELL	SHERWIN WILLIAMS		PURE WHITE SW 7005	BASE COLOR
P2	LATEX PAINT - EGGSHELL	SHERWIN WILLIAMS		AUSTERE GRAY SW6184	ACCENT COLOR
P3	LATEX PAINT - EGGSHELL	SHERWIN WILLIAMS		PORPOISE SW7047	DOOR FRAME & WINDOW FRAME
<b>SURFACE</b>					
PL1	PLASTIC LAMINATE	WILSONART	PREMIUM LAMINATE	8213K-28 PHANTOM COCOA	GLOSS LINE FINISH
SS1	SOLID SURFACE - CORIAN	DUPONT	CORIAN	CONCRETE	USE NON-FERROUS FASTENERS / ATTACHMENTS

### GENERAL FINISH NOTES

- ALL FLOOR TRANSITIONS TO BE LOCATED AT CENTER OF DOOR, U.N.O.
- FIELD VERIFY ALL DIMENSIONS BEFORE FABRICATION OF MILLWORK.
- COORDINATE ALL MILLWORK WITH APPLIANCES BEFORE FABRICATION.
- AT SOFFITS RECEIVING COLOR. PAINT ALL SIDES OF SOFFIT.
- ALL WOOD TRIM TO BE STAINED TO MATCH DOOR STAIN.
- ALL COUNTERTOP, BACKSPLASHES, AND EDGE BANDING TO HAVE COORDINATING FINISHES.
- PROVIDE A SMOOTH TRANSITION AT ALL FLOOR MATERIALS - CONTRACTOR TO INSTALL ALL FLOOR FINISHES AT SAME LEVEL, DESPITE DIFFERENT THICKNESS. PROVIDE FLOOR TRANSITION WHERE OCCURS.
- GYPSUM BOARD SOFFITS TO BE PAINTED WHITE.
- ALL EXPOSED CEILING TO BE PAINTED (WHERE OCCURS). REFER TO REFLECTED CEILING PLANS. COORDINATE WITH ARCHITECT FOR PAINT COLOR.

### FINISH PLAN SYMBOLS





### GENERAL CEILING NOTES

- REFER TO DETAIL A2/A500 FOR TYPICAL CEILING SUSPENSION & SEISMIC BRACING
- REFER TO DETAIL A3/A500 FOR TYPICAL SUSPENDED GYP. BOARD CEILING
- ALL UNIDENTIFIED CEILING TYPES ON THE PLANS SHALL BE TYPE "A" AT 9'-0" A.F.F.
- GRID SUSPENSION SYSTEMS SHALL BE CENTERED WITHIN AREAS INDICATED, UNLESS NOTED OTHERWISE
- PAINT BLACK ALL EXPOSED STRUCTURE, MECHANICAL, DUCTS, ELECTRICAL WORK, PIPING, ETC.
- REFER TO ARCHITECTURAL DRAWINGS FOR LOCATION OF MECHANICAL GRILLES, AND TO MECHANICAL DRAWINGS FOR QUANTITIES AND TYPES
- REFER TO ARCHITECTURAL DRAWINGS FOR LOCATIONS OF LIGHT FIXTURES AND TO ELECTRICAL DRAWINGS FOR QUANTITY AND TYPES WITH SPRINKLER CONTRACTOR TO AVOID CONFLICTS IN FIELD
- MECHANICAL AND ELECTRICAL CONTRACTORS TO COORDINATE WORK WITH SPRINKLER CONTRACTOR TO AVOID CONFLICTS IN FIELD
- ALL CEILING HEIGHTS ARE ELEVATION ABOVE TOP OF CONCRETE FLOOR SLAB



### KEYED NOTES

KEY VALUE	KEYNOTE TEXT
900.1	NON-STRUCTURAL METAL STUD FRAMING WALL WITH (1) LAYER OF GYPSUM BOARD ON EACH SIDE. PREPARE FOR NEW FINISHES
914.2	SCHEDULED PAINT. EXISTING CEILING / SOFFIT. REPAIR, PATCH AND PREPARE EXISTING CEILING / SOFFIT FOR NEW PAINT



### CEILING LEGEND



**A2 REFLECTED CEILING PLAN - ENLARGED**



**A5 FINISH PLAN - ENLARGED**




REV	DATE	DESCRIPTION

VCBO NUMBER: 20370  
DATE: 06/08/2020

**IMC EMERGENCY DEPARTMENT - X-RAY #2 - UPGRADE**  
INTERMOUNTAIN HEALTHCARE  
5125 SOUTH COTTONWOOD STREET, MURRAY UT 84107  
PERMIT SET / BID SET


LEVEL LL-1 REFLECTED CEILING PLAN + FINISH PLAN - ENLARGED

**A110.2**




**D2 X-RAY ROOM - WEST ELEVATION**

SCALE: 1/2" = 1'-0"




**C2 X-RAY ROOM - EAST ELEVATION**

SCALE: 1/2" = 1'-0"




**D4 X-RAY ROOM - SOUTH ELEVATION**

SCALE: 1/2" = 1'-0"




**C4 X-RAY ROOM - NORTH ELEVATION**

SCALE: 1/2" = 1'-0"



**A2 X-RAY ROOM - 3D VIEW 1**

SCALE:



**A4 X-RAY ROOM - 3D VIEW 2**

SCALE:

**ARCHITECTURAL MILLWORK KEY**

EXAMPLE:

B 1 1 1 \*\* ACCESSORY (IF NOTED)  
ONE ADJUSTABLE SHELF  
ONE DOOR  
ONE DRAWER  
BASE CABINET \*

W 0 2 1 \*\* ACCESSORY (IF NOTED)  
ONE ADJUSTABLE SHELF  
TWO DOORS  
ZERO DRAWERS  
WALL CABINET \*

T 0 2 5 L ACCESSORY (IF NOTED)  
FIVE ADJUSTABLE SHELVES  
TWO DOORS  
ZERO DRAWERS  
TALL CABINET \*

\* NOTE:  
SB DENOTES SINK BASE CABINET  
TW DENOTES TALL WARDROBE CABINET  
FMD DENOTES FRAMED PERFORATED METAL PANEL DOOR  
FGD DENOTES FRAMED GLASS PANEL DOOR

\*\* NOTE:  
L DENOTES LOCKS AT DOOR (S) / DRAWER (S)  
F DENOTES FILE DRAWER (S)

CABINET MEASUREMENTS SHOWN ARE ACTUAL SIZES. BASE CABINET HEIGHTS ALLOW FOR A COUNTERTOP 1 1/2" THICK. CABINET DEPTHS ARE MEASURED FROM THE BACK TO THE FACE OF THE DOOR OR DRAWER FRONT (WHERE APPLICABLE)

**MILLWORK LEGEND**


- MILLWORK DIMENSION NUMBERS ARE WIDTH X HEIGHT X DEPTH.
- ALL MILLWORK DIMENSIONED FROM BASE TO TOP OF IDENTIFIED COUNTERTOP, TYP.
- CABINET DEPTHS ARE MEASURED FROM THE WALL TO THE FACE OF THE DOOR OR DRAWER FRONT (WHERE APPLICABLE).
- PROVIDE BASE AT ALL CABINET TOE SPACE, UNLESS NOTED OTHERWISE.
- PROVIDE GROMMET WHERE "G" IS LABELED ON PLANS OR ELEVATIONS.
- ALL COUNTERTOPS TO HAVE A 4" BACKSLASH, UNLESS NOTED OTHERWISE, TO MATCH COUNTERTOP, ON BACK AND SIDE WALLS.
- PROVIDE FILLER PANELS TO SEAL SIDES AND TOPS OF ALL CABINETS PLACED AT AN ANGLE TO ADJACENT WALL(S).
- ALL MILLWORK TO FINISHED ON ENDS, TYP.
- CONTRACTOR TO PROVIDE BLOCKING BEHIND ALL CABINETS, COAT RACKS, PENCIL SHARPENER BLOCKS, T.V. BRACKETS AND PROJECTION SCREENS AS WELL AS ALL WALL MOUNTED ACCESSORIES, INCLUDING WHITE BOARDS, TACKBOARDS, TOILET AND URINAL PARTITIONS AND TOILET ROOM ACCESSORIES, ETC. NOTE: ONLY 2X WOOD BLOCKING IS ACCESSIBLE BEHIND MILLWORK AND TOILET ROOM PARTITIONS.
- REFER TO SHEET A1103 FOR FINISH COLORS ON ALL MILLWORK AND CASEWORK.

**KEYED NOTES**

KEY VALUE	KEYNOTE TEXT
904.0	GYPSUM BOARD PAINTED WHERE EXPOSED
912.0	SCHEDULED BASE
1022.0	STAINLESS STEEL GROMMET, GLOVE BOX DISPENSER ACCESS HOLE
A10101	


REV	DATE	DESCRIPTION

VCBO NUMBER: 20370  
DATE: 06/08/2020




**D2 CONTROL ROOM - WEST ELEVATION**

SCALE: 1/2" = 1'-0"




**D4 CONTROL ROOM - SOUTH ELEVATION**

SCALE: 1/2" = 1'-0"



**C2 CONTROL ROOM - EAST ELEVATION**


SCALE: 1/2" = 1'-0"



**C4 CONTROL ROOM - NORTH ELEVATION**

SCALE: 1/2" = 1'-0"

**ARCHITECTURAL MILLWORK KEY**



CABINET MEASUREMENTS SHOWN ARE ACTUAL SIZES. BASE CABINET HEIGHTS ALLOW FOR A COUNTERTOP 1 1/2" THICK. CABINET DEPTHS ARE MEASURED FROM THE BACK TO THE FACE OF THE DOOR OR DRAWER FRONT (WHERE APPLICABLE)

**MILLWORK LEGEND**

- MILLWORK DIMENSION NUMBERS ARE WIDTH X HEIGHT X DEPTH.
- ALL MILLWORK DIMENSIONED FROM BASE TO TOP OF IDENTIFIED COUNTERTOP, TYP.
- CABINET DEPTHS ARE MEASURED FROM THE WALL TO THE FACE OF THE DOOR OR DRAWER FRONT (WHERE APPLICABLE).
- PROVIDE BASE AT ALL CABINET TOE SPACE, UNLESS NOTED OTHERWISE.
- PROVIDE GROMMET WHERE "G" IS LABELED ON PLANS OR ELEVATIONS.
- ALL COUNTERTOPS TO HAVE A 4" BACKSPLASH, UNLESS NOTED OTHERWISE, TO MATCH COUNTERTOP, ON BACK AND SIDE WALLS.
- PROVIDE FILLER PANELS TO SEAL SIDES AND TOPS OF ALL CABINETS PLACED AT AN ANGLE TO ADJACENT WALL(S).
- ALL MILLWORK TO FINISHED ON ENDS, TYP.
- CONTRACTOR TO PROVIDE BLOCKING BEHIND ALL CABINETS, COAT RACKS, PENCIL SHARPENER BLOCKS, T.V. BRACKETS AND PROJECTION SCREENS AS WELL AS ALL WALL MOUNTED ACCESSORIES, INCLUDING WHITE BOARDS, TACKBOARDS, TOILET AND URINAL PARTITIONS AND TOILET ROOM ACCESSORIES, ETC. NOTE: ONLY 2X WOOD BLOCKING IS ACCESSIBLE BEHIND MILLWORK AND TOILET ROOM PARTITIONS.
- REFER TO SHEET A110.3 FOR FINISH COLORS ON ALL MILLWORK AND CASEWORK.

**KEYED NOTES**

KEY VALUE	KEYNOTE TEXT
900.1	NON-STRUCTURAL METAL STUD FRAMING WALL WITH (1) LAYER OF GYPSUM BOARD ON EACH SIDE. PREPARE FOR NEW FINISHES
902.0	1 1/2" METAL STUD FRAMING
904.0	GYPSUM BOARD PAINTED WHERE EXPOSED
912.0	SCHEDULED BASE
A10101	



REV	DATE	DESCRIPTION

VCBO NUMBER: 20370  
 DATE: 06/08/2020

**IMC EMERGENCY DEPARTMENT - X-RAY #2 - UPGRADE**

INTERMOUNTAIN HEALTHCARE  
 5125 SOUTH COTTONWOOD STREET, MURRAY UT 84107  
 PERMIT SET / BID SET

X-RAY ROOM - INTERIOR ELEVATIONS - 3D VIEWS





1. Design Criteria

- 1.1. Governing Building Code ..... 2018 International Building Code (IBC)
A. Risk Category ..... IV
1.2. Earthquake
A. Seismic Design Category ..... D
1. Note that all Unistrut connections must meet code requirements for a high seismic application.
B. Spectral Response Accelerations
Ss = 1.488 g
S1 = 0.992 g
C. Soil Site Class ..... D
Fv = 1.0
D. Importance Factor, Ie ..... 1.5

2. Concrete

- 2.1. Materials shall comply with the Standards specified in American Concrete Institute (ACI) 318-14, "Building Code Requirements for Structural Concrete."
A. Concrete mix design requirements shall be as follows:
Table with columns: Location, 28 days (psi), Max WC Ratio, Air Content (%), Max Aggregate Size, Exposure Classes (F, S, C)
B. Cementitious Materials:
1. Portland Cement (ASTM C150):
a. Type I or II for exposure class S0.
b. Type I or V for exposure class S1.
c. Type V for exposure class S2 and S3.
2. Fly Ash (ASTM C618, Class C or F): maximum fly ash content as a percentage of total weight of cementitious materials shall be 25 percent.
C. Concrete Density (Maximum Air Dry Weight):
1. Normal weight concrete shall be approximately 145 to 155 pounds per cubic foot. Aggregate shall be ASTM C33.
2. Lightweight concrete shall not exceed 110 pounds per cubic foot and shall be made of lightweight coarse aggregates and either lightweight and/or normal weight fines meeting ASTM C330.
D. Admixtures:
1. Air-entraining admixtures, comply with ASTM C 260 (when used).
a. Tolerance on air content as delivered shall be +/- 1.5%.
b. When air content of a trowel finished floor slab exceeds 3%, there is an increased risk for delaminations and blistering to occur. When this situation is present, the Contractor shall pay special attention to the finishing procedures to help minimize such risks. Refer to ACI 302.1R-15 "Guide for Concrete Floor and Slab Construction" for proper finishing guidelines.
2. Corrosion Inhibiting admixture, comply with ASTM C1582 (when used):
a. Corrosion inhibiting additive containing a minimum of 30 percent calcium nitrite dosed at 3 gallons per cubic yard shall be added to all reinforced concrete with exposure class C2.
3. The use of super plasticizers and water reducers is allowed, but not required.
4. Calcium chloride or admixtures containing calcium chloride shall not be added to the concrete mix.
E. Chloride Ion: Maximum water soluble chloride ion concentrations in hardened concrete at age between 28 and 42 days contributed from the ingredients including water, aggregates, cementitious materials, and admixtures shall not exceed a maximum, by weight of cement, of 1.00% for concrete with exposure class C0, 0.30% for concrete with exposure class C1, 0.15% for concrete with exposure class C2, and 0.06% for all prestressed concrete.
F. Slump Limit: 4 inches, maximum for all concrete prior to the addition of plasticizers and water reducing admixtures. The concrete supplier shall indicate the final slump of each concrete mix in the submitted mix design.

3. Structural Steel

- 3.1. Material:
A. W-Shapes: ASTM A992 (Fy = 50 ksi), except as noted otherwise
B. All Other Shapes and Plates: ASTM A36 (Fy = 36 ksi), except as noted otherwise
3.2. Fabrication and construction shall comply with the following Codes and Standards:
A. American Institute of Steel Construction (AISC) 360-16, "Specification for Structural Steel Buildings"
B. AISC 903-16, "Code of Standard Practice for Steel Buildings and Bridges" excluding the following: Section 3.3 (last sentence of first paragraph), Section 4.4, Section 4.4.1, Section 4.4.2, Section 4.5, and Section 7.13.3
1. The architectural drawings are the prime contract drawings. Consultants' drawings by other disciplines are supplementary to the architectural drawings. The structural drawings shall be used in conjunction with the architectural drawings. Detailing and shop drawing production for structural elements (including dimensions) contained in architectural, structural, and/or other consultants' drawings. Refer to the Special Instructions section of the general notes, below.
C. AISC/RISC 2014, "Specification for Structural Joints Using High-Strength Bolts"
D. American Welding Society (AWS) D1.1:2015, "Structural Welding Code - Steel" (specific items do not apply when they conflict with the AISC requirements)
E. American Welding Society (AWS) D1.8:2016, "Structural Welding Code - Seismic Supplement" (specific items do not apply when they conflict with the AISC requirements)
3.3. Structural shapes and plates shall be fabricated from newly rolled (milled) one-piece sections without splices, unless specifically noted otherwise on the structural drawings. Connections for structural steel shall comply with the structural drawings, unless written approval is given by the Structural Engineer.
3.4. Welding:
A. It is recommended the steel erector contractor and steel fabricator contact the Quality Assurance Agency prior to beginning any welds. A program of joint preparation and welding procedures should be worked out between the two parties before the welding is started so that correct welds will be made from the beginning.
B. Certification of Welders: All shop and field welding shall be executed by AWS certified welders who have been specifically certified for the process of welding being performed. The welder's certification will be considered as being current unless the welder is not engaged in the process of welding being performed for a period exceeding six months or there is a specific reason to question a welder's ability as required by AWS. Certification and records must comply with AWS Standards. Certification and appropriate records must be provided to the Architect prior to beginning work.
C. Electrodes: E-70 XX or as noted otherwise. E60 XX may be used for welding steel floor and roof decks.
D. Minimum Welds: All intersecting steel shapes that are not bolted shall be connected by a fillet weld all around, unless noted otherwise. Fillet weld sizes that are not shown shall be 1/16" less than the thinnest of the connected parts for thicknesses 1/4" and larger. Fillet welds on plates less than 1/4" shall be of the same size as the thinnest of the connected parts.
3.5. Bolted Connections:
A. Provide snug tightened joints with Group A (threads not excluded) bolts for steel to steel connections, unless noted otherwise. Snug tightened joints shall be used in connections for simple span framing and beam (or girder) to bearing plate connections. Snug tight is the condition that exists when all of the plates in a connection have been pulled into firm contact by the bolts in the joint and all of the bolts in the joint have been tightened sufficiently to prevent the removal of the nuts without the use of a wrench. The snug tightened condition is typically achieved with a few impacts of an impact wrench, application of an electric torque wrench until the wrench begins to slow, or the full effort of a worker on an ordinary spud wrench.
B. Provide hardened washers beneath the turned element of all bolts or nuts. Provide hardened beveled washers, to compensate for the lack of parallelism, where the outer face of the bolted parts has a slope greater than one in twenty with respect to the plane normal to the bolt axis. Hardened washers or plates installed over oversized holes or slotted holes shall be at least 5/16" thick and shall conform to ASTM F436. Plates or bars installed at slotted holes shall have a size sufficient to completely cover the slot after installation.
C. Bolts, nuts and washers shall not be reused.

4. Miscellaneous

- 4.1. Post-Installed Anchors in Concrete
A. Anchorage to hardened concrete shall include all mechanical and adhesive anchors and epoxy dowelled reinforcing bars of size, quantity, spacing, and embedment as shown on the drawings. Additional anchors shall not be used without approval from the Engineer prior to installation.
B. Special inspection is required during the installation of all post-installed anchors. Refer to applicable code evaluation reports and the Quality Assurance and Statement of Special Inspections sections of the General Structural Notes.
C. Anchorage to Concrete:
1. All post-installed anchors into hardened concrete shall be selected from the following pre-approved products, unless noted otherwise:
Table with columns: Steel Screw Anchor, Evaluation Report
Table with columns: Steel Expansion/Wedge Anchor, Evaluation Report
Table with columns: Adhesive Anchor System, Evaluation Report
2. Adhesive anchors shall be installed into concrete having a minimum age of 21 days. For installations sooner than 21 days, consult the adhesive manufacturer.
D. Alternate anchors or adhesives are permitted with approval of the Engineer. The Contractor shall submit the proposed anchor product data and code evaluation report demonstrating the anchor is equivalent or exceeds the capacity of the specified anchor.
E. Installation of adhesive anchors horizontally or upwards inclined to support sustained tension loads shall be performed by personnel certified by an applicable certification program. Certification shall include written and performance tests in accordance with the ACI/CRSI Adhesive Anchor Installer Certification program, or equivalent. Proof of current certification shall be submitted to the Engineer for approval prior to commencement of installation.
F. Anchors shall be installed according to the Manufacturer's Printed Installation Instructions and applicable code evaluation reports including:
1. Hole diameter, depth, and cleaning procedure
2. Adhesive mixing, preparation, and placement
3. Installation torque

- G. Locate all existing reinforcement and embedded items prior to drilling into concrete or masonry elements. Do not damage rebar or embeds while drilling or installing anchors.
H. Grout all defective or abandoned holes with non-shrink grout or an injectable epoxy adhesive matching the surrounding concrete compressive strength. Consult the Architect for additional requirements at architecturally exposed concrete.
I. Drilled anchors are not allowed in post-tensioned concrete without approval of the Architect and Engineer.
J. Carbon steel anchors are limited to use in dry, interior locations.

5. Special Instructions

- 5.1. The project specifications are not superseded by the General Structural Notes but are intended to be complementary to them. Consult the specifications for additional requirements in each section. Notes and specific details on the drawings shall take precedence over General Structural Notes and typical details.
5.2. The architectural drawings are the prime contract drawings. Consultant drawings by other disciplines are supplementary to the architectural drawings. All omissions or conflicts, including dimensions between the various elements of the consultants' drawings and/or specifications shall be brought to the attention of the Architect before proceeding with any work involved. In case of conflict, follow the most stringent requirement as directed by the Architect without additional cost to the Owner. Any work done by the Contractor after discovery of such discrepancy shall be done at the Contractor's risk.
5.3. The structural drawings shall be used in conjunction with the architectural drawings. Primary structural elements and overall structural layout are indicated within the structural plans and details. Some secondary elements, architectural layouts, alcoves, elevations, slopes, depressions, curbs, mechanical equipment and electrical equipment, are not indicated within the structural drawings. Detailing and shop drawing production for structural elements will require information (including dimensions) contained in the architectural, structural and/or other consultants' drawings.
5.4. Shoring and Bracing Requirements:
A. Floor and Roof Structures - The General Contractor is responsible for the method and sequence of all structural erection. The Contractor shall provide temporary shoring and bracing as the method of erection requires to provide adequate vertical and lateral support. Shoring and bracing shall remain in place as the chosen method requires until all permanent members are in place and all final connections are completed, including all roof and floor attachments. The building shall not be considered stable until all connections are complete.
B. Foundation walls must be braced until the complete floor or roof systems is completed. Do not backfill until floor or roof systems are in place.
C. Walls above grade shall be braced until the structural system is complete. Walls shall not be considered to be self-supporting.
5.5. Submittals: A copy of all shop drawings that have been submitted for review must be kept at the construction site for reference. These drawings must bear the appropriate review stamps. The shop drawing review shall relieve the Contractor of the responsibility of completing the project according to the contract documents. The General Contractor shall review and mark all shop drawings prior to submitting them to the Architect for review. Shop Drawings made from reproductions of (these) contract drawings will be rejected.
5.6. Project Coordination: It shall be the responsibility of the General Contractor to coordinate with all trades any and all items that are to be integrated into the structural system. Openings or penetrations through, or attachments to the structural system that are not indicated on these drawings shall be the responsibility of the General Contractor and shall be coordinated with the Architect/Engineers. The order of construction is the responsibility of the General Contractor. It is the Contractor's obligation to provide all items necessary for the chosen procedure.
5.7. Contractor shall field verify all dimensions, and conditions. If the contract drawings do not represent actual conditions, Contractor shall notify Architect/Engineer prior to fabrication or construction within that area.
5.8. Notice of Copyright: The structural drawings, plans, schedules, notes and details are hereby copyrighted by Reaveley Engineers. Submission or distribution of documents to meet official regulatory requirements or for similar purposes in connection with the project is not to be construed as publication in derogation of Reaveley Engineers' reserved rights. The documents defining the structure are instruments of service prepared by Reaveley Engineers for one use only. Furthermore, these documents shall not be reproduced, or copied, in whole or in part by the Contractor or subcontractors for preparation of shop drawings or other submittals.

6. Quality Assurance

- 6.1. Quality Assurance Agency Requirements:
A. The Owner shall engage a qualified Quality Assurance Agency (QAA) to provide all special inspection and quality assurance testing for the project. The QAA shall provide all information necessary for the building official to determine that the agency meets the applicable requirements.
1. The QAA shall be objective, competent and independent from the Contractor responsible for the work being inspected. The agency shall disclose to the building official and the registered design professional in responsible charge possible conflicts of interest so that objectivity can be confirmed.
2. The QAA shall have adequate equipment to perform required tests. The equipment shall be periodically calibrated.
3. The QAA shall employ experienced personnel educated in conducting, supervising and evaluating tests and special inspections. Experience or training shall be considered relevant where the documented experience or training is related in complexity to the same type of special inspection or testing activities for projects of similar complexity and material qualities.
4. The QAA shall send copies of all inspection and testing reports to the building official, Owner, Architect, Engineers shall indicate that the work inspected was or was not completed in conformance to the approved construction documents. Discrepancies shall be brought to the immediate attention of the Contractor for correction. If they are not corrected, the discrepancies shall be brought to the attention of the Architect and Engineer.
5. The QAA shall submit a final report documenting required special inspections and tests, and correction of any discrepancies noted in the inspections or tests. The final report shall be distributed to the building official, Owner, Architect and Engineer in a timely manner prior to the completion of the project.
6.2. Contractor Responsibilities:
A. The Contractor shall submit a written statement of responsibility to the building official and the Owner or the owner's authorized agent prior to the commencement of work on the systems or components listed in the statement of special inspections. The Contractor's statement of responsibility shall contain acknowledgement or awareness of the special requirements contained in the statement of special inspections.
B. Notification of QAA: The Contractor shall notify the QAA in a timely manner so that inspection and testing may be performed as outlined in the statement of special inspections.
6.3. Structural Observations by the Engineer of Record:
A. The Engineer of Record will perform structural observations at critical phases of the project. Observations will be made on a periodic basis throughout the construction of the structural system. Copies of the Engineer's report will be distributed to the Architect, Contractor, Owner, and building official.
B. Observation visits to the site by the Engineer's field representatives shall not be construed as inspection or approval of construction.

7. Statement of Special Inspections

- 7.1. The following materials, systems and components require special inspection or testing per Chapter 17 of the International Building Code (IBC).
7.2. For items requiring continuous inspection, a special inspector must be present onsite during the performance of that task. In most cases, periodic inspections/tests shall be performed prior to commencing the task, intermittently during the task, and at the completion of the task. Frequency marked with (E) designates periodic inspections that must be performed prior to or upon completion of every task.
7.3. Re-purposed anchor points require special inspection.

Concrete Construction per IBC Sections 1705.3 & 1705.12

Table with columns: Item, Frequency, Detailed Instructions. Lists post-installed adhesive anchors and mechanical anchors with their respective inspection frequencies and instructions.

PLAN LEGEND table showing symbols for STEEL BRACED FRAME, STEEL BEAM OR GIRDER, STEEL JOIST OR PURLIN, HORIZONTAL BRIDGING, EXISTING CONCRETE SHEAR WALL, FOUNDATION WALL OR RETAINING WALL, EXISTING STEEL COLUMN - WIDE FLANGE, EXISTING STEEL BEAM OR GIRDER, and EXISTING STEEL JOIST OR PURLIN.

ABBREVIATIONS table listing various abbreviations and their full names, such as @ AT, AB ANCHOR BOLT (S), BLW BELOW, etc.

ABBREVIATIONS table listing various abbreviations and their full names, such as REQD REQUIRED, SFRS SEISMIC FORCE RESISTING SYSTEM, SHT SHEET, etc.

PLAN MARKS table listing various plan marks and their descriptions, such as BF# BRACED FRAME, CB# CONCRETE BEAM, CC# CONCRETE COLUMN, etc.

STRUCTURAL DRAWING LIST table with columns SHT NO. and SHT NAME, listing drawing sheets like S-001 GENERAL STRUCTURE NOTES & LEGENDS & ABBREVIATIONS.





Table with columns REV, DATE, DESCRIPTION.

VCBO NUMBER: 20370
DATE: 06/08/2020

IMC EMERGENCY DEPARTMENT - X-RAY #2 - UPGRADE
INTERMOUNTAIN HEALTHCARE
5125 SOUTH COTTONWOOD STREET, MURRAY UT 84107
PERMIT SET / BID SET






REV	DATE	DESCRIPTION


VCBO NUMBER: 20370  
DATE: 06/08/2020

**IMC EMERGENCY DEPARTMENT - X-RAY #2 - UPGRADE**  
INTERMOUNTAIN HEALTHCARE  
5125 SOUTH COTTONWOOD STREET, MURRAY UT 84107  
PERMIT SET / BID SET


LEVEL 01 FRAMING PLAN



**A1**  
**S-101** Level LL-1 FLOOR SLAB PLAN  
SCALE: 1/4" = 1'-0"




**A5**  
**S-101** LEVEL 01 FLOOR FRAMING PLAN  
SCALE: 1/4" = 1'-0"




REV	DATE	DESCRIPTION


VCBO NUMBER: 20370  
DATE: 06/08/2020




**B3** UNISTRUT CONNECTION TO STEEL BEAM  
S-501 NO SCALE




**B4** UNISTRUT CONNECTION DETAIL  
S-501 NO SCALE



**B5** TYPICAL UNISTRUT CONNECTION TO CONC SLAB OVER FLOOR DECK  
S-501 NO SCALE



**A2** RADIOLOGY EQUIPMENT SUPPORT DETAIL  
S-501 NO SCALE



**A4** RADIOLOGY EQUIPMENT SUPPORT DETAIL  
S-501 NO SCALE

**IMC EMERGENCY DEPARTMENT - X-RAY #2 - UPGRADE**

INTERMOUNTAIN HEALTHCARE  
5125 SOUTH COTTONWOOD STREET, MURRAY UT 84107  
PERMIT SET / BID SET

SYMBOLS LEGEND	
SYMBOL	DESCRIPTION
<b>REFERENCE AND LINE SYMBOLS</b>	
01	DETAIL INDICATOR: A5 INDICATES DETAIL NUMBER, E-501 INDICATES DRAWING SHEET WHERE DETAIL IS SHOWN.
02	ELEVATION OR SECTION INDICATOR, EXTERIOR: A5 INDICATES ELEVATION OR SECTION NUMBER, E-201 INDICATES DRAWING SHEET WHERE ELEVATION OR SECTION IS SHOWN.
03	ELEVATION OR SECTION INDICATOR, INTERIOR: A5 INDICATES ELEVATION OR SECTION NUMBER, E-201 INDICATES DRAWING SHEET WHERE ELEVATION OR SECTION IS SHOWN.
04	ROOM IDENTIFIER WITH ROOM NAME AND NUMBER.
05	KEYNOTE INDICATOR.
06	REVISION INDICATOR.
07	EQUIPMENT INDICATOR.
08	BREAK, STRAIGHT: TO BREAK PARTS OF DRAWING
09	BREAK, ROUND
10	NEW LINE: MEDIUM LINE.
11	HIDDEN FEATURES LINE: HIDDEN, THIN LINE
12	EXISTING TO REMAIN LINE: THIN LINE.
13	DEMOLITION LINE: DASHED, MEDIUM LINE
<b>WIRING METHODS</b>	
01	WIRING.
02	BRANCH CIRCUIT HOME RUN TO PANELBOARD: NUMBER OF ARROWS INDICATES NUMBER OF CIRCUITS. LETTER AND NUMBER NOTATIONS IDENTIFY PANEL AND CIRCUIT NUMBERS. USE #12 CONDUCTORS. EXCEPT #10 CONDUCTORS SHALL BE INSTALLED IF DISTANCES EXCEED THOSE SPECIFIED IN THE ELECTRICAL SPECIFICATIONS.
03	BRANCH CIRCUIT HOME RUN TO PANELBOARD: NUMBER OF ARROWS INDICATES NUMBER OF CIRCUITS. LETTER AND NUMBER NOTATIONS IDENTIFY PANEL AND CIRCUIT NUMBERS. NUMBER IN BOX REFERS TO THE CONDUCTOR AND CONDUIT SCHEDULE. FOR BRANCH WIRING USE #12 CONDUCTORS. EXCEPT #10 CONDUCTORS SHALL BE INSTALLED IF DISTANCES EXCEED THOSE SPECIFIED IN THE ELECTRICAL SPECIFICATIONS.
04	WIRING AND/OR RACEWAY: THIN LINE. WHERE "X" = : CATV = CABLE TELEVISION NC = NURSE CALL CCTV = CLOSED CIRCUIT TV P = POWER TELEVISION RC = RIGID CONDUIT FA = FIRE ALARM S = SOUND FO = FIBER OPTICS T = TELEPHONE I = INTERCOM TV = TELEVISION
OTHERS AS NOTED IN OTHER SCHEDULES. RACEWAYS AND WIRING SHALL BE SIZED AS SHOWN AND/OR SPECIFIED.	
05	LOW VOLTAGE WIRING: DIVIDE, MEDIUM LINE.
06	CONDUIT STUB. DIMENSION RECORD DRAWINGS AND MARK.
07	CONDUCTOR & CONDUIT ("CC") SCHEDULE INDICATOR. REFER TO ONE-LINE DIAGRAM.
08	ADA ACCESS PUSH PLATE
09	JUNCTION BOX.
10	JUNCTION BOX, SYSTEMS FURNITURE COMMUNICATION CONNECTION.
11	JUNCTION BOX, SECURITY SYSTEM. PROVIDE CONDUIT AND ROUGH-IN PER SECURITY DRAWINGS.
12	CABLE TRAY ABOVE ACCESSIBLE CEILING.
13	EARTH GROUND (ONE-LINE DIAGRAM).
14	JUNCTION BOX, CEILING.
15	LADDER RACK.
16	MECHANICAL EQUIPMENT CONNECTION. REFER TO EQUIPMENT SCHEDULE FOR REQUIREMENTS.
<b>LIGHTING (REFER TO FIXTURE SCHEDULE FOR SYMBOLS)</b>	
01	FIXTURE IDENTIFICATION: (W-3) INDICATES FIXTURE TYPE AS SCHEDULED.
02	FIXTURE IDENTIFICATION, EMERGENCY WITH BATTERY PACK, CONNECTED TO GENERATOR AS INDICATED: (W-3) INDICATES FIXTURE TYPE AS SCHEDULED.
03	EGRESS DIRECTION ARROW (EXIT SIGNS).
04	EXIT SIGN: SINGLE FACE; CEILING MOUNTED
05	EXIT SIGN: SINGLE FACE; WALL MOUNTED
06	EXIT SIGN: DOUBLE FACE; CEILING MOUNTED
07	EXIT SIGN: DOUBLE FACE; WALL MOUNTED
<b>LIGHTING CONTROL</b>	
01	OCCUPANCY SENSOR, DUAL TECHNOLOGY, OMNI-DIRECTIONAL, CEILING.
02	OCCUPANCY SENSOR, DUAL TECHNOLOGY, WALL.
03	VACANCY SENSOR, DUAL TECHNOLOGY, OMNI-DIRECTIONAL, CEILING.
04	VACANCY SENSOR, DUAL TECHNOLOGY, WALL.
05	PHOTOCELL.
06	LOW VOLTAGE DIGITAL LIGHTING CONTROL SWITCH: LETTER "a,b" INDICATES ZONING WHERE SHOWN (REFER TO PLANS, SCHEDULES, AND DETAILS FOR EXACT BUTTON CONFIGURATION AND PROGRAMMING REQUIREMENTS)
07	DIGITAL LIGHTING DIMMING CONTROLLER
08	DIGITAL PLUG LOAD CONTROLLER
09	DIGITAL LIGHTING ROOM CONTROLLER
10	LIGHTING SPACE CONTROL TYPE. X INDICATES TYPE. SEE SCHEDULE / DIAGRAM.
<b>CCTV</b>	
01	CCTV CABLE, POWER.
02	CCTV CABLE, VIDEO SIGNAL.
03	CCTV HEADEND EQUIPMENT.
04	CCTV MONITOR.
05	CCTV CAMERA/ENCLOSURE WITH LENS, TYPICAL. SEE SCHEDULE.

SYMBOLS LEGEND	
SYMBOL	DESCRIPTION
<b>WIRING DEVICES</b>	
01	RECEPTACLE, DUPLEX, NEMA 5-20R.
02	RECEPTACLE, DUPLEX, ABOVE COUNTER. NEMA 5-20R.
03	RECEPTACLE, DUPLEX, CEILING. NEMA 5-20R.
04	RECEPTACLE, DUPLEX WITH GROUND FAULT CIRCUIT INTERRUPTER, DRINKING FOUNTAIN: CONCEAL WATER COOLER RECEPTACLE BEHIND WATER COOLER. SEE MECHANICAL PLUMBING SHOP DRAWINGS FOR INSTALLATION REQUIREMENTS.
05	RECEPTACLE, DUPLEX, HOSPITAL GRADE. NEMA 5-20R.
06	RECEPTACLE, DUPLEX ON EMERGENCY POWER. NEMA 5-20R.
07	RECEPTACLE, DUPLEX, HOSPITAL GRADE ON EMERGENCY POWER. NEMA 5-20R.
08	RECEPTACLE, DUPLEX WITH GROUND FAULT CIRCUIT INTERRUPTER. NEMA 5-20R.
09	RECEPTACLE, DUPLEX WITH GROUND FAULT CIRCUIT INTERRUPTER, HOSPITAL GRADE. NEMA 5-20R.
10	RECEPTACLE, DUPLEX WITH GROUND FAULT CIRCUIT INTERRUPTER, HOSPITAL GRADE ON EMERGENCY POWER. NEMA 5-20R.
11	RECEPTACLE, QUADRAPLEX, HOSPITAL GRADE. NEMA 5-20R.
12	RECEPTACLE, QUADRAPLEX ON EMERGENCY POWER. NEMA 5-20R.
13	RECEPTACLE, QUADRAPLEX, HOSPITAL GRADE ON EMERGENCY POWER. NEMA 5-20R.
14	RECEPTACLE, QUADRAPLEX WITH GROUND FAULT CIRCUIT INTERRUPTER. NEMA 5-20R.
15	RECEPTACLE, SPECIAL PURPOSE. PROVIDE RECEPTACLE TO MATCH EQUIPMENT PLUG.
16	RECEPTACLE, SPECIAL PURPOSE ON EMERGENCY POWER. PROVIDE RECEPTACLE TO MATCH EQUIPMENT PLUG.
17	RECEPTACLE, DRYER. NEMA 14-30R.
18	RECEPTACLE, RANGE. NEMA 14-50R.
19	MULTI-OUTLET ASSEMBLY. NEMA 5-20R.
20	DROP CORD. SEE DETAIL.
21	FLUSH FLOOR BOX. "F" SHOWN ON DRAWINGS. REFER TO WIRING DEVICE SCHEDULE IN THE ELECTRICAL SPECIFICATIONS FOR CONFIGURATION AND DEVICES.
22	POWER POLE. "P" SHOWN ON DRAWINGS. REFER TO WIRING DEVICE SCHEDULE IN THE ELECTRICAL SPECIFICATIONS FOR CONFIGURATION AND DEVICES.
23	FLUSH FIRE RATED POKE THRU. "F" SHOWN ON DRAWINGS. REFER TO WIRING DEVICE SCHEDULE IN THE ELECTRICAL SPECIFICATIONS FOR CONFIGURATION AND DEVICES.
24	SWITCH, DIMMER.
25	SWITCH, SINGLE POLE ("X" INDICATES FIXTURES CONTROLLED).
26	SWITCH, DOUBLE POLE ("X" INDICATES FIXTURES CONTROLLED).
27	SWITCH, THREE-WAY ("X" INDICATES FIXTURES CONTROLLED).
28	SWITCH, FOUR-WAY ("X" INDICATES FIXTURES CONTROLLED).
29	SWITCH, MOMENTARY.
30	RECEPTACLE, QUADRAPLEX WITH GROUND FAULT CIRCUIT INTERRUPTER, HOSPITAL GRADE. NEMA 5-20R.
31	RECEPTACLE, QUADRAPLEX WITH GROUND FAULT CIRCUIT INTERRUPTER, HOSPITAL GRADE ON EMERGENCY POWER. NEMA 5-20R.
32	RECEPTACLE, SINGLE PLEX, WITH USB OUTLET
33	RECEPTACLE, DUPLEX, RECESSED, NEMA 5-20R, AUTOMATICALLY CONTROLLED THROUGH TIME OR OCCUPANCY BASED CONTROLS (REFER TO PLANS FOR CONTROL METHOD)
34	RECEPTACLE, QUADRAPLEX, RECESSED, NEMA 5-20R, AUTOMATICALLY CONTROLLED THROUGH TIME OR OCCUPANCY BASED CONTROLS (REFER TO PLANS FOR CONTROL METHOD)
35	INDICATES A RECEPTACLE IS AUTOMATICALLY CONTROLLED THROUGH TIME OR OCCUPANCY BASED CONTROLS (REFER TO PLANS FOR CONTROL METHOD)
<b>STRUCTURED CABLING IHC</b>	
01	IHC COMMUNICATIONS DEVICE (1 DATA).
02	IHC COMMUNICATIONS DEVICE (1 DATA / 1 ANALOG).
03	IHC COMMUNICATIONS DEVICE (1 DATA WALL PHONE).
04	IHC COMMUNICATIONS DEVICE (2 DATA).
05	IHC COMMUNICATIONS DEVICE (3 DATA).
06	IHC COMMUNICATIONS DEVICE (4 DATA).
07	IHC COMMUNICATIONS DEVICE (6 DATA).
08	IHC COMMUNICATIONS DEVICE PHYSIOLOGICAL MONITOR (1 DATA).
09	IHC COMMUNICATIONS DEVICE WIRELESS ACCESS POINT (2 DATA).
<b>NURSE CALL</b>	
01	CORRIDOR LIGHT.
02	BATHROOM PULL CORD STATION.
03	DUTY STATION.
04	EMERGENCY ASSISTANCE CALL STATION.
05	EMERGENCY ASSISTANCE CODE BLUE CALL STATION.
06	PATIENT STATION.
07	STAFF STATION.
08	TOUCH SCREEN NURSE CALL MASTER STATION.
<b>TV DISTRIBUTION</b>	
01	TV DISTRIBUTION CABLE, INDIVIDUAL DROPS.
02	TV DISTRIBUTION CABLE, TRUNK.
03	COMBINER.
04	DIRECTIONAL COUPLER.
05	TV OUTLET.
06	TERMINATOR, 75 OHM (TV DISTRIBUTION).

SYMBOLS LEGEND	
SYMBOL	DESCRIPTION
<b>ELECTRICAL POWER AND DISTRIBUTION</b>	
01	FUSE WITH RATING (ONE-LINE DIAGRAM).
02	DISCONNECT, FUSED (ONE-LINE DIAGRAM).
03	DISCONNECT, NONFUSED (ONE-LINE DIAGRAM).
04	DISCONNECT WITH FUSE AND MOTOR STARTER COMBINATION (ONE-LINE DIAGRAM).
05	OVERLOAD RELAY (ONE-LINE DIAGRAM).
06	STARTER (ONE-LINE DIAGRAM).
07	CIRCUIT BREAKER, MOLDED CASE (ONE-LINE DIAGRAM).
08	CIRCUIT BREAKER, MOLDED CASE WITH SHUNT TRIP (ONE-LINE DIAGRAM).
09	CIRCUIT BREAKER, SOLID STATE (ONE-LINE DIAGRAM).
10	CIRCUIT BREAKER, SOLID STATE WITH GROUND FAULT PROTECTION (ONE-LINE DIAGRAM).
11	CIRCUIT BREAKER, SOLID STATE WITH GROUND FAULT PROTECTION (ONE-LINE DIAGRAM).
12	MOTOR.
13	TRANSFORMER (ONE-LINE DIAGRAM).
14	PANELBOARD WITH MAIN LUGS ONLY. BUS SIZE AND PHASE AS SHOWN (ONE-LINE DIAGRAM).
15	PANELBOARD WITH MAIN CIRCUIT BREAKER. SIZE AND PHASE AS SHOWN (ONE-LINE DIAGRAM).
16	PANELBOARD WITH MAIN AND SUB FEED CIRCUIT BREAKER (ONE-LINE DIAGRAM).
17	PANELBOARD WITH SUB FEED LUGS (ONE-LINE DIAGRAM).
18	PANELBOARD WITH CIRCUIT BREAKER AND SUB FEED LUGS (ONE-LINE DIAGRAM).
19	CT CABINET PER UTILITY'S REQUIREMENTS (ONE-LINE DIAGRAM).
20	TRANSFER SWITCH (ONE-LINE DIAGRAM).
21	DIGITAL MULTIMETER (ONE-LINE DIAGRAM).
22	SERVICE ENTRANCE SURGE PROTECTION (ONE-LINE DIAGRAM).
23	GENERATOR, POWER (ONE-LINE DIAGRAM).
24	METER.
25	VARIABLE FREQUENCY MOTOR CONTROLLER (ONE-LINE DIAGRAM).
26	DISCONNECT SWITCH, FUSED.
27	STARTER, COMBINATION WITH DISCONNECT SWITCH.
28	STARTER OR MOTOR CONTROLLER.
29	PUSHBUTTON.
30	PANELBOARD CABINET, MOTOR CONTROL.
31	PANELBOARD CABINET, FLUSH MOUNTED.
32	PANELBOARD CABINET, SURFACE MOUNTED, 1 SECTION.
33	PANELBOARD CABINET, SURFACE MOUNTED, 2 SECTION.
34	DISTRIBUTION PANEL OR SWITCHBOARD.
35	LIGHTING RELAY, CONTACTOR PANEL, OR DIMMING ENCLOSURE.
36	LIGHTING CONTROL STATION.
37	SWITCH, TOGGLE MOTOR STARTER WITH OVERLOAD PROTECTION.
38	TRANSFORMER: NUMBER INDICATES KVA.
<b>TECHNOLOGY SYSTEMS</b>	
01	TECHNOLOGY SYSTEM CABLE. SEE SPECIFIC JOB EQUIPMENT LIST FOR APPLICABLE DESIGNATIONS. EXAMPLES: C = CONTROL CABLE G = GROUND CABLE, 10 AWG, 1 CONDUCTOR, GREEN INSULATED M = MICROPHONE CABLE S = SPEAKER CABLE, 18 OHM SYSTEM Z = SPEAKER CABLE, 8 OHM SYSTEM
02	SPEAKER, CEILING MOUNTED.
03	SPEAKER, WALL MOUNTED.
04	EQUIPMENT CABINET.
05	CONNECTION PANEL.




SYMBOLS LEGEND	
SYMBOL	DESCRIPTION
<b>FIRE ALARM</b>	
01	FIRE SYSTEM ANNUNCIATOR.
02	FIRE ALARM CONTROL PANEL, SEMI-RECESSED.
03	FIRE ALARM NOTIFICATION POWER SUPPLY.
04	FIRE ALARM TRANSDUCER OR TRANSMITTER.
05	SMOKE CONTROL PANEL.
06	AUTOMATIC DOOR CLOSERS: DOOR CLOSERS SHALL BE FURNISHED WITH DOOR HARDWARE AND CONNECTED TO BY FIRE ALARM INSTALLERS.
07	CONTROL MODULE.
08	MONITOR MODULE.
09	FIRE ALARM MANUAL PULL STATION.
10	SHUT DOWN RELAY. INSTALL RELAY IN CONTROL CIRCUIT OF EQUIPMENT TO BE CONTROLLED IN THE EVENT OF A FIRE.
11	MAGNETIC DOOR HOLDER.
12	FIRE SERVICE OR EMERGENCY TELEPHONE STATION, ACCESSIBLE.
13	FIRE SERVICE OR EMERGENCY TELEPHONE STATION, HANDSET.
14	FIRE SERVICE OR EMERGENCY TELEPHONE STATION, JACK.
15	DETECTOR, SMOKE.
16	DETECTOR, SMOKE WITH AUXILIARY CONTACT.
17	DETECTOR, SMOKE, BEAM RECEIVER.
18	DETECTOR, SMOKE, BEAM TRANSMITTER.
19	DETECTOR, SMOKE, ELEVATOR RECALL DESIGNATION.
20	DETECTOR, SMOKE WITH GUARD.
21	DETECTOR, SMOKE, RESIDENTIAL.
22	DETECTOR, SMOKE, DUCT WITH HOUSING AND SAMPLING TUBE.
23	INDICATOR LAMP.
24	STROBE.
25	STROBE. SUBSCRIPT INDICATES CANDELA RATING.
26	ALARM, HORN/SPEAKER, WEATHERPROOF.
27	ALARM, HORN/SPEAKER, ONE ASSEMBLY.
28	ALARM, HORN/SPEAKER, ONE ASSEMBLY. SUBSCRIPT INDICATES CANDELA RATING.
29	ALARM, CHIME/STROBE, ONE ASSEMBLY.
30	ALARM, HORN/SPEAKER WITH GUARD, ONE ASSEMBLY.
31	ALARM, MINI HORN/SPEAKER, ONE ASSEMBLY.
32	SPEAKER, EVACUATION.
33	SPEAKER, EVACUATION, COMBINATION STROBE.
34	DETECTOR, FLOW SWITCH: FLOW SWITCHES SHALL BE PROVIDED AND INSTALLED WITH FIRE SPRINKLER SYSTEM AND SHALL BE CONNECTED TO LOCATIONS SHOWN ON THE FIRE SPRINKLER SHOP DRAWINGS.
35	DETECTOR, TAMPER SWITCH WITH VALVE: TAMPER SWITCHES SHALL BE PROVIDED AND INSTALLED TO LOCATIONS SHOWN ON THE FIRE SPRINKLER SHOP DRAWINGS.
36	SMOKE DAMPER.
37	FIRE AND SMOKE DAMPER.
38	BELL (GONG).
39	DETECTOR, CARBON MONOXIDE.
40	DETECTOR, SMOKE/STROBE, RESIDENTIAL.
41	ALARM, HORN/SPEAKER, ONE ASSEMBLY, CEILING MOUNTED. SUBSCRIPT INDICATES CANDELA RATING.
42	ALARM, STROBE, CEILING MOUNTED. SUBSCRIPT INDICATES CANDELA RATING.
<b>SECURITY</b>	
01	SECURITY CABLE. SEE EQUIPMENT SCHEDULE FOR CABLE TYPE.
02	ACCESS CONTROL HEADEND EQUIPMENT.
03	SECURITY CONTROL PANEL.
04	INTRUSION DETECTION HEADEND EQUIPMENT.
05	CARD ACCESS DOOR TYPE #1 OR AS NOTED. SEE SCHEDULE.
06	CARD READER.
07	KEYPAD/CARD READER COMBINATION.
08	DOOR SWITCH, BALANCED MAGNETIC CONTROL.
09	EXIT REQUEST.
10	REMOTE DOOR RELEASE BUTTON.
11	BELL.
12	BUZZER.
13	BUZZER, COMBINATION BELL.
14	PANIC DURESS SWITCH.

ABBREVIATIONS			
NOTE: ALL ABBREVIATIONS MAY NOT BE USED.			
1P	SINGLE POLE	KV	KILOVOLT
1PH	SINGLE-PHASE	KVA	KILOVOLT AMPERE
1WAY	ONE-WAY	KVAR	KILOVOLT AMPERE REACTIVE
2C	TWO-CONDUCTOR	KW	KILOWATT
2WAY	TWO-WAY	KWH	KILOWATT HOUR
3/C	THREE-CONDUCTOR	KWD	LIQUID TIGHT FLEXIBLE DIODE
3WAY	THREE-WAY	LFD	LIQUID TIGHT FLEXIBLE METAL CONDUIT
4OUT	QUADRUPLE RECEPTACLE OUTLET	LFLNC	LIQUID TIGHT FLEXIBLE NONMETALLIC CONDUIT
4PDT	FOUR-POLE DOUBLE THROW	LPS	LOW PRESSURE SODIUM
4PST	FOUR-POLE SINGLE THROW	LRA	LOCKED ROTOR AMPS
4W	FOUR-WIRE	LTV	LIGHTING
4WAY	FOUR-WAY	LV	LOW VOLTAGE
A	ABOVE COUNTER	MA TV	MASTER ANTENNA TELEVISION SYSTEM
AC	ARMORED CABLE	MAX	MAXIMUM
ADA	AMERICANS WITH DISABILITIES ACT	MC	METAL CLAD
ADJ	ADJACENT	MCC	MINIMUM CIRCUIT AMPS
AFF	ABOVE FINISHED FLOOR	MCB	MAIN CIRCUIT BREAKER
AFG	ABOVE FINISHED GRADE	MCC	MOTOR CONTROL CENTER
AIC	AMPERE INTERRUPTING CAPACITY	MCP	MOTOR CIRCUIT PROTECTION
ALUM	ALUMINUM	MDP	MAIN DISTRIBUTION PANEL
AMP	AMPERE	MG	MOTOR GENERATOR
ANN	ANNUNCIATOR	MH	MANHOLE
AP	ACCESS POINT (WIRELESS DATA)	MIN	MINIMUM
AR	AS REQUIRED	ML	MAIN LUGS ONLY
AS	AS REQUIRED	MOC	MAXIMUM OVERCURRENT PROTECTION
ASC	AMPS SHORT CIRCUIT AUTOMATIC TRANSFER SWITCH	MTS	MANUAL TRANSFER SWITCH
ATC	AUDIO VISUAL	NA	NOT APPLICABLE
AV	AMERICAN WIRE GAGE	NC	NORMALLY CLOSED
BB	BUCK-BOOST TRANSFORMER	NEC	NATIONAL ELECTRICAL CODE
XB	CEILING MOUNTED COMMUNITY ANTENNA TELEVISION	NEMA	NATIONAL ELECTRICAL MANUFACTURERS ASSOCIATION
CB	CIRCUIT BREAKER	NFC	NATIONAL FIRE CODE
CBCA	CUSTOM COLOR AS SELECTED BY ARCHITECT	NFPA	NATIONAL FIRE PROTECTION ASSOCIATION
CCTV	CLOSED CIRCUIT TELEVISION CONTRACTOR FURNISHED/ CONTRACTOR INSTALLED	NIC	NOT IN CONTRACT
CF/CI	CONTRACTOR FURNISHED/ CONTRACTOR INSTALLED	NO	NORMALLY OPEN
CF/CI	CONTRACTOR FURNISHED/ OWNER INSTALLED	NTS	NOT TO SCALE
CFBA	CUSTOM FINISH AS SELECTED BY ARCHITECT	OC	ON CENTER
CKT	CIRCUIT	OC	OVER CURRENT PROTECTION
CM	CONSTRUCTION MANAGER	OF/CI	OWNER FURNISHED/ CONTRACTOR INSTALLED
CND	CONDUIT	OF/CI	OWNER FURNISHED/ OWNER INSTALLED
CO	CONVENIENCE OUTLET	OP	OBTAIN FROM PLANS
COR	CONTRACTING OFFICER'S REPRESENTATIVE	OH DR	OVERHEAD (COILING) DOOR
CP	CURRENT TRANSFORMER	OL	OVERLOAD
CTV	CABLE TELEVISION	OS	OVERSIGHT
CU	COPPER	PF	POWER FACTOR
DBA	UNIT OF SOUND LEVEL	PH	PHASE
DPTD	DOUBLE POLE, DOUBLE THROW	PNL	PANEL
DS	DISCONNECT SWITCH	PT	POTENTIAL TRANSFORMER
EA	EACH	PTZ	PANTILT/ZOOM
EM	EMERGENCY	QTY	QUANTITY
EMT	ELECTRICAL METALLIC TUBING	R	REMOVE
ENT	ELECTRIC NONMETALLIC TUBING	RCP	REFLECTED CEILING PLAN
EPO	EMERGENCY POWER OFF EQUIPMENT	RMC	RIGID METAL CONDUIT
EA	EXISTING	RNC	RIGID NONMETAL CONDUIT
EX	FURNITURE MOUNTED	RPM	REVOLUTIONS PER MINUTE
FA	FIRE ALARM	RR	REMOVE AND RELOCATE
FCA	FIRE ALARM CONTROL PANEL	S	START/STOP
FCL	FULL LOAD AMPS	SCA	SHORT CIRCUIT AMPS
FMC	FLEXIBLE METAL CONDUIT	SCBA	STANDARD COLOR AS SELECTED BY ARCHITECT
FOB	FREIGHT ON BOARD	SF	SQUARE FOOT (FEET)
FVNR	FULL VOLTAGE NON-REVERSING	SFBA	STANDARD FINISH AS SELECTED BY ARCHITECT
FVR	GENERATOR	SPD	SURGE PROTECTIVE DEVICE
GEN	GROUND FAULT INTERRUPTER	SPDT	SINGLE POLE, DOUBLE THROW
GFCI	GROUND FAULT PROTECTION	SPEC	SPECIFICATION
GND	GROUND	SPST	SINGLE POLE, SINGLE THROW
HD	HEAVY DUTY	SWBD	SWITCHBOARD
HD	HIGH INTENSITY DISCHARGE	SWGR	SWITCHGEAR
HDA	HAND-OFF-AUTOMATIC	TL	TWIST LOCK
HPS	HORSE POWER	TP	TELEPHONE POLE
HPE	HIGH POWER FACTOR	TP	TWISTED PAIR
HPS	HIGH PRESSURE SODIUM	TTB	TELEPHONE TERMINAL BOARD
HV	HIGH VOLTAGE	TVSS	TRANSIENT VOLTAGE SURGE SUPPRESSER
HZ	HERTZ	TYP	TYPICAL
IG	INPUT/OUTPUT	UF	UNDERFLOOR
IG	ISOLATED GROUND	UGND	UNDERGROUND
IMC	INTERMEDIATE METAL CONDUIT	UNPS	UNINTERRUPTIBLE POWER SUPPLY
INS	INSULATED/ISOLATED	V	VOLTS
IR	INFRARED	VFCVFD	VARIABLE FREQUENCY MOTOR CONTROLLER
J-BOX	JUNCTION BOX	W	WITH
		WO	WITHOUT
		WP	WEATHERPROOF
		XFMR	TRANSFORMER

GENERAL ELECTRICAL NOTES	
1.	CLARIFICATION METHODS: AT THE TIME OF BIDDING, BIDDERS SHALL FAMILIARIZE THEMSELVES WITH THE DRAWINGS AND SPECIFICATIONS. ANY QUESTIONS, MISUNDERSTANDINGS, CONFLICTS, DELETIONS, DISCONTINUED PRODUCTS, CATALOG NUMBER DISCREPANCIES, DISCREPANCIES BETWEEN THE EQUIPMENT SUPPLIED AND THE INTENT OR FUNCTION OF THE EQUIPMENT, ETC. SHALL BE SUBMITTED TO THE ARCHITECT/ENGINEER IN WRITING FOR CLARIFICATION PRIOR TO ISSUANCE OF THE FINAL ADDENDUM AND BIDDING OF THE PROJECT. WHERE DISCREPANCIES OR MULTIPLE INTERPRETATIONS OCCUR, THE MOST STRINGENT (WHICH IS GENERALLY RECOGNIZED AS THE MOST COSTLY) THAT MEETS THE INTENT OF THE DOCUMENTS SHALL BE ENFORCED.
2.	OWNER FURNISHED ITEMS: THE OWNER WILL FURNISH MATERIAL AND EQUIPMENT AS INDICATED IN THE CONTRACT DOCUMENTS TO BE INCORPORATED INTO THE WORK. THESE ITEMS ARE ASSIGNED TO THE INSTALLER AND COSTS FOR RECEIVING, HANDLING, STORAGE, IF REQUIRED, AND INSTALLATION ARE INCLUDED IN THE CONTRACT SUM. A. THE INSTALLER'S RESPONSIBILITIES ARE THE SAME AS IF THE INSTALLER FURNISHED THE MATERIALS OR EQUIPMENT. B. THE OWNER WILL ARRANGE AND PAY FOR DELIVERY OF OWNER FURNISHED ITEMS FREIGHT ON BOARD JOB SITE AND THE INSTALLER WILL INSPECT DELIVERIES FOR DAMAGE. IF OWNER FURNISHED ITEMS ARE DAMAGED, DEFECTIVE OR MISSING, DOCUMENT DAMAGED ITEMS WITH THE



### E2 RECEPTACLE MOUNTING DETAILS

SCALE: NTS




### D2 ADA DETAIL

SCALE: NTS




### D3 SWITCH MOUNTING DETAILS

SCALE: NTS




### C1 CMU DEVICE MOUNTING ALIGNMENT DETAIL

SCALE: NTS




### C2 LIGHTING MOUNTING DETAILS

SCALE: NTS




### C3 COMMUNICATIONS MOUNTING DETAILS

SCALE: NTS




### B1 FIRE ALARM MOUNTING DETAILS

SCALE: NTS




### A1 BOX MOUNTING DETAILS

SCALE: NTS



### A2 TYPICAL WALL MOUNTED DEVICES ALIGNMENT DETAIL

SCALE: NTS



### GENERAL SHEET NOTES

1. DETERMINE MOUNTING HEIGHTS OF ELECTRICAL AND ELECTRONIC EQUIPMENT IN THE FOLLOWING ORDER OF PRIORITY:
  - 1 - ELEVATIONS (ARCHITECTURAL, ELECTRICAL, MECHANICAL, ETC).
  - 2 - EQUIPMENT SHOP DRAWINGS.
  - 3 - FIELD INSTRUCTIONS.
2. LOCATE RECEPTACLES SERVING THE SAME TYPE OF USE AT A UNIFORM HEIGHT UNLESS DIRECTED OTHERWISE.
3. MECHANICAL, ELECTRICAL, AND COMMUNICATION ROOMS: COORDINATE LOCATION OF LIGHTING AND POWER RECEPTACLES WITH EQUIPMENT, PIPING, AND DUCTWORK. DO NOT INSTALL RECEPTACLES BEHIND EQUIPMENT OR WHERE OTHERWISE INACCESSIBLE. POSITION LIGHTING REGARDLESS OF WHERE SHOWN ON DRAWING TO PROVIDE PROPER ILLUMINATION.
4. MOUNT RECEPTACLE BOXES FOR SWITCHES AND RECEPTACLES WITH LONG AXIS OF THE DEVICE VERTICAL UNLESS OTHERWISE INDICATED.
5. SET BOXES WITH PLASTER RINGS FLUSH WITH FINISHED SURFACE.
6. LOCATE BOX COVERS OR DEVICE PLATES SO THEY WILL NOT SPAN DIFFERENT TYPES OF BUILDING FINISHES EITHER VERTICALLY OR HORIZONTALLY.
7. VERIFY ALL DOOR CONDITIONS ON ARCHITECTURAL DRAWINGS PRIOR TO INSTALLING SWITCHES.
8. LOCATE WIRING DEVICES WHICH ARE ADJACENT AND ARE COMPATIBLE VOLTAGES IN ONE PLATE.
9. WHERE DEVICES ARE LOCATED IN CLOSE PROXIMITY OF THE SAME VERTICAL PLANE, ALIGN DEVICES VERTICALLY PER THE TYPICAL WALL MOUNTED DEVICES ALIGNMENT DETAIL, UNLESS OTHERWISE INDICATED.

### SHEET KEYNOTES


1. LOCATE RECEPTACLES BEHIND DRINKING FOUNTAINS.
2. REFER TO ARCHITECTURAL ELEVATIONS FOR PLACEMENT OF OUTLETS.
3. LOCATE AT BOTTOM OF BEAMS (OR JOISTS) OR AT CEILING. (REDUCE SPACING BY 3 PERPENDICULAR TO BEAM OR JOIST DIRECTION.) FOR OTHER CONDITIONS, REFER TO NFPA 72.
4. LOCATE DETECTOR ANYWHERE IN SHADED AREA BUT NOT IN TOP 4\"/>

REV	DATE	DESCRIPTION
-----	------	-------------

VCBO NUMBER: 20370  
DATE: 06/08/2020


REV	DATE	DESCRIPTION

VCBO NUMBER: 20370  
DATE: 06/08/2020




- NOTES:
1. TYPICAL FOR WOOD AND METAL STUD ROUGH-IN.
  2. PLASTER RINGS NOT SHOWN.
  3. LOCATE ALL OUTLET BOXES IN ACCORDANCE WITH ARCHITECTURAL AND MECHANICAL DRAWINGS AND WITH ALL APPLICABLE SHOP DRAWINGS.
  4. IN ACCORDANCE WITH IBC 714.3.2 EXCEPTION 1, OUTLETS ON OPPOSITE SIDES OF WALLS OR PARTITIONS IN THE SAME STUD SPACE IN A RATED FIRE SEPARATION WALL MUST BE SEPARATED BY A MINIMUM OF 24" HORIZONTAL DISTANCE OR LISTED, SOUND AND FIRE RATED PUTTY PADS SHALL BE USED ON THE OUTLET BOXES.
  5. IN NON-RATED WALLS, OUTLETS ON OPPOSITE SIDES OF WALLS OR PARTITIONS MUST BE SEPARATED BY 16" FOR SOUND ATTENUATION.


**1** TYPICAL ROUGH-IN REQUIREMENTS DETAIL  
SCALE: NTS




**2** TYPICAL RACEWAY SUPPORT METHODS DETAIL  
SCALE: NTS




**4** RECESSED FIXTURE MOUNTING DETAIL  
SCALE: NTS




**3** TYPICAL CONDUIT RACK DETAIL  
SCALE: NTS



**7** FIRE STOP FOR METAL CONDUIT THROUGH GYPSUM WALL BOARD  
SCALE: NTS



**6** TYPICAL FIRE STOP FOR CABLES/CONDUIT THROUGH CONCRETE FLOORING  
SCALE: NTS



**5** TYPICAL FIRE STOP FOR CABLES/CONDUIT THROUGH CONCRETE WALLS  
SCALE: NTS

**CONNECTIVITY REQUIREMENTS**

Broadband Connections are necessary during the installation process and going forward to ensure full support from the Engineering Teams for the customers system. Maximum performance and availability for the customers system is maintained and closely monitored during the lifetime of the system. Proactive and reactive maintenance is available utilizing the wide range of digital tools using the connectivity solutions listed below:


- Site-to-Site VPN/GE Solution
- Site-to-Site VPN/Custom Solution
- Connection through Dedicated Service Network
- Internet Access - connectivity for InSite 2.0

The requirements for these connectivity solutions are explained in the broadband solutions catalogue (separate document).

**ELECTRICAL NOTES**

- All wires specified shall be copper stranded, flexible, thermo-plastic, color coded, cut 10 foot long at outlet boxes, duct termination points or stubbed conduit ends. All conductors, power, signal and ground, must be run in a conduit or duct system. Electrical contractor shall ring out and tag all wires at both ends. Wire runs must be continuous copper stranded and free from splices.
- Aluminum or solid wires are not allowed.
- Wire sizes given are for use of equipment. Larger sizes may be required by local codes.
- It is recommended that all wires be color coded, as required in accordance with national and local electrical codes.
- Conduit sizes shall be verified by the architect, electrical engineer or contractor, in accordance with local or national codes.
- Convenience outlets are not illustrated. Their number and location are to be specified by others. Locate at least one convenience outlet close to the system control, the power distribution unit and one on each wall of the procedure room. Use hospital approved outlet or equivalent.
- General room illumination is not illustrated. Caution should be taken to avoid excessive heat from overhead spotlights. Damage can occur to ceiling mounting components and wiring if high wattage bulbs are used. Recommend low wattage bulbs no higher than 75 watts and use dimmer controls (except MR). Do not mount lights directly above areas where ceiling mounted accessories will be parked.
- Routing of cable ductwork, conduits, etc., must run direct as possible otherwise may result in the need for greater than standard cable lengths (refer to the interconnection diagram for maximum usable lengths point to point).
- Conduit turns to have large, sweeping bends with minimum radius in accordance with national and local electrical codes.
- A special grounding system is required in all procedure rooms by some national and local codes. It is recommended in areas where patients might be examined or treated under present, future, or emergency conditions. Consult the governing electrical code and confer with appropriate customer administrative personnel to determine the areas requiring this type of grounding system.
- The maximum point to point distances illustrated on this drawing must not be exceeded.
- Physical connection of primary power to GE equipment is to be made by customers electrical contractor with the supervision of a GE representative. The GE representative would be required to identify the physical connection location, and insure proper handling of GE equipment.
- GEHC conducts power audits to verify quality of power being delivered to the system. The customer's electrical contractor is required to be available to support this activity.

- All junction boxes, conduit, duct, duct dividers, switches, circuit breakers, cable tray, etc., are to be supplied and installed by customers electrical contractor.
- Conduit and duct runs shall have sweep radius bends
- Conduits and duct above ceiling or below finished floor must be installed as near to ceiling or floor as possible to reduce run length.
- Ceiling mounted junction boxes illustrated on this plan must be installed flush with finished ceiling.
- All ductwork must meet the following requirements:
  1. Ductwork shall be metal with dividers and have removable covers.
  2. Ductwork shall be certified/rated for electrical power purposes.
  3. Ductwork shall be electrically and mechanically bonded together in an approved manner.
  4. PVC as a substitute must be used in accordance with all local and national codes.
- All openings in access flooring are to be cut out and finished off with grommet material by the customers contractor.
- General contractor to insert pull cords for all cable run conduits between the equipment room and the operators control room.
- 10 foot pigtails at all junction points.
- Grounding is critical to equipment function and patient safety. Site must conform to wiring specifications shown on this plan.



ELECTRICAL LAYOUT ITEM LIST	
1	One 3" stubbed conduit
2	Flush 12"x4"x4" box for Control
3	Flush pull box for Chest Unit - size per local code
4	Flush pull box for Access Unit - size per local code
5	Flush pull box for TB - size per local code
6	Flush pull box for Generator - size per local code
7	18" x 3 1/2" Flush vertical wall duct with minimum 2 dividers
8	Box above ceiling size per local code
9	Box below floor size per local code
10	One 1" cnd
11	One 1 1/2" cnd
12	One 2" cnd
13	One 2 1/2" cnd
14	Main Disconnect Panel

ITEM		Outlet Legend for GE Equipment	
☑	System emergency off (SEO), (recommended height 1.2m [48"] above floor)		
⊗	X-Ray room warning light control panel		
⊕	X-Ray ON lamp (L1) - 24V		
⊗	Door interlock switch (needed only if required by state/local codes)		
⊕	Duplex hospital grade, dedicated wall outlet 120-v, single phase power		
⊕	Dedicated telephone line(s)		
⊕	Network outlet		


Conduit Legend	
----	Above Ceiling
---	Below Floor

Additional Conduit Runs (Contractor Supplied and Installed)				
From	To	Qty	In.	mm
3 phase power	Main disconnect	1	AS REQ'D	AS REQ'D
Main disconnect	Emergency off	1	1/2"	16
Main disconnect	Systems Cabinet	1	AS REQ'D	AS REQ'D
Warning light	Warning light control	1	1/2"	16
1 phase power		1	AS REQ'D	AS REQ'D
		1	1/2"	16
Systems Cabinet	Door Switch	1	1/2"	16
	Tether Interface Box	2	2"	53
	Access Point	1	1"	27
Operators Console	Tether Interface Box	1	1"	27
	Access Point	1	2"	53

Intermountain Medical Center | OPTIMA XR646 HD (G3) | RAD-M210912-FIN-00-A.DWG | Rev A | Date 13/May/2020 | E1 - Electrical Notes | 12/16

Intermountain Medical Center | OPTIMA XR646 HD (G3) | RAD-M210912-FIN-00-A.DWG | 1/4"-1"0" Rev A | Date 13/May/2020 | E2 - Electrical Layout | 13/16



REV	DATE	DESCRIPTION

VCBO NUMBER: 20370  
DATE: 06/08/2020

**IMC EMERGENCY DEPARTMENT - X-RAY #2 - UPGRADE**  
INTERMOUNTAIN HEALTHCARE  
5125 SOUTH COTTONWOOD STREET, MURRAY UT 84107  
PERMIT SET / BID SET

GE DRAWINGS  
**EE701**

**CONNECTIVITY REQUIREMENTS**

Broadband Connections are necessary during the installation process and going forward to ensure full support from the Engineering Teams for the customers system. Maximum performance and availability for the customers system is maintained and closely monitored during the lifetime of the system. Proactive and reactive maintenance is available utilizing the wide range of digital tools using the connectivity solutions listed below:


- Site-to-Site VPN/GE Solution
- Site-to-Site VPN/Custom Solution
- Connection through Dedicated Service Network
- Internet Access - connectivity for InSite 2.0

The requirements for these connectivity solutions are explained in the broadband solutions catalogue (separate document).

**ELECTRICAL NOTES**

- All wires specified shall be copper stranded, flexible, thermo-plastic, color coded, cut 10 foot long at outlet boxes, duct termination points or stubbed conduit ends. All conductors, power, signal and ground, must be run in a conduit or duct system. Electrical contractor shall ring out and tag all wires at both ends. Wire runs must be continuous copper stranded and free from splices.
- Aluminum or solid wires are not allowed.
- Wire sizes given are for use of equipment. Larger sizes may be required by local codes.
- It is recommended that all wires be color coded, as required in accordance with national and local electrical codes.
- Conduit sizes shall be verified by the architect, electrical engineer or contractor, in accordance with local or national codes.
- Convenience outlets are not illustrated. Their number and location are to be specified by others. Locate at least one convenience outlet close to the system control, the power distribution unit and one on each wall of the procedure room. Use hospital approved outlet or equivalent.
- General room illumination is not illustrated. Caution should be taken to avoid excessive heat from overhead spotlights. Damage can occur to ceiling mounting components and wiring if high wattage bulbs are used. Recommend low wattage bulbs no higher than 75 watts and use dimmer controls (except MR). Do not mount lights directly above areas where ceiling mounted accessories will be parked.
- Routing of cable ductwork, conduits, etc., must run direct as possible otherwise may result in the need for greater than standard cable lengths (refer to the interconnection diagram for maximum usable lengths point to point).
- Conduit turns to have large, sweeping bends with minimum radius in accordance with national and local electrical codes.
- A special grounding system is required in all procedure rooms by some national and local codes. It is recommended in areas where patients might be examined or treated under present, future, or emergency conditions. Consult the governing electrical code and confer with appropriate customer administrative personnel to determine the areas requiring this type of grounding system.
- The maximum point to point distances illustrated on this drawing must not be exceeded.
- Physical connection of primary power to GE equipment is to be made by customers electrical contractor with the supervision of a GE representative. The GE representative would be required to identify the physical connection location, and insure proper handling of GE equipment.
- GEHC conducts power audits to verify quality of power being delivered to the system. The customer's electrical contractor is required to be available to support this activity.

- All junction boxes, conduit, duct, duct dividers, switches, circuit breakers, cable tray, etc., are to be supplied and installed by customers electrical contractor.
- Conduit and duct runs shall have sweep radius bends
- Conduits and duct above ceiling or below finished floor must be installed as near to ceiling or floor as possible to reduce run length.
- Ceiling mounted junction boxes illustrated on this plan must be installed flush with finished ceiling.
- All ductwork must meet the following requirements:
  1. Ductwork shall be metal with dividers and have removable covers.
  2. Ductwork shall be certified/rated for electrical power purposes.
  3. Ductwork shall be electrically and mechanically bonded together in an approved manner.
  4. PVC as a substitute must be used in accordance with all local and national codes.
- All openings in access flooring are to be cut out and finished off with grommet material by the customers contractor.
- General contractor to insert pull cords for all cable run conduits between the equipment room and the operators control room.
- 10 foot pigtails at all junction points.
- Grounding is critical to equipment function and patient safety. Site must conform to wiring specifications shown on this plan.



ELECTRICAL LAYOUT ITEM LIST	
1	One 3" stubbed conduit
2	Flush 12"x4"x4" box for Control
3	Flush pull box for Chest Unit - size per local code
4	Flush pull box for Access Unit - size per local code
5	Flush pull box for TB - size per local code
6	Flush pull box for Generator - size per local code
7	18" x 3 1/2" Flush vertical wall duct with minimum 2 dividers
8	Box above ceiling size per local code
9	Box below floor size per local code
10	One 1" cnd
11	One 1 1/2" cnd
12	One 2" cnd
13	One 2 1/2" cnd
14	Main Disconnect Panel

ITEM		Outlet Legend for GE Equipment	
☑	System emergency off (SEO), (recommended height 1.2m [48"] above floor)		
⊗	X-Ray room warning light control panel		
⊕	X-Ray ON lamp (L1) - 24V		
⊗	Door interlock switch (needed only if required by state/local codes)		
⊕	Duplex hospital grade, dedicated wall outlet 120-v, single phase power		
⊕	Dedicated telephone line(s)		
⊕	Network outlet		


Conduit Legend	
----	Above Ceiling
---	Below Floor

Additional Conduit Runs (Contractor Supplied and Installed)				
From	To	Qty	In.	mm
3 phase power	Main disconnect	1	AS REQ'D	AS REQ'D
Main disconnect	Emergency off	1	1/2"	16
Main disconnect	Systems Cabinet	1	AS REQ'D	AS REQ'D
Warning light	Warning light control	1	1/2"	16
1 phase power		1	AS REQ'D	AS REQ'D
		1	1/2"	16
Systems Cabinet	Door Switch	1	1/2"	16
	Tether Interface Box	2	2"	53
	Access Point	1	1"	27
Operators Console	Tether Interface Box	1	1"	27
	Access Point	1	2"	53

Intermountain Medical Center | OPTIMA XR646 HD (G3) | RAD-M210912-FIN-00-A.DWG | Rev A | Date 13/May/2020 | E1 - Electrical Notes | 12/16

Intermountain Medical Center | OPTIMA XR646 HD (G3) | RAD-M210912-FIN-00-A.DWG | 1/4"-1"0" Rev A | Date 13/May/2020 | E2 - Electrical Layout | 13/16



REV	DATE	DESCRIPTION

VCBO NUMBER: 20370  
DATE: 06/08/2020

**IMC EMERGENCY DEPARTMENT - X-RAY #2 - UPGRADE**  
INTERMOUNTAIN HEALTHCARE  
5125 SOUTH COTTONWOOD STREET, MURRAY UT 84107  
PERMIT SET / BID SET

GE DRAWINGS  
**EE701**

**CONNECTIVITY REQUIREMENTS**

Broadband Connections are necessary during the installation process and going forward to ensure full support from the Engineering Teams for the customers system. Maximum performance and availability for the customers system is maintained and closely monitored during the lifetime of the system. Proactive and reactive maintenance is available utilizing the wide range of digital tools using the connectivity solutions listed below:


- Site-to-Site VPN/GE Solution
- Site-to-Site VPN/Custom Solution
- Connection through Dedicated Service Network
- Internet Access - connectivity for InSite 2.0

The requirements for these connectivity solutions are explained in the broadband solutions catalogue (separate document).

**ELECTRICAL NOTES**

- All wires specified shall be copper stranded, flexible, thermo-plastic, color coded, cut 10 foot long at outlet boxes, duct termination points or stubbed conduit ends. All conductors, power, signal and ground, must be run in a conduit or duct system. Electrical contractor shall ring out and tag all wires at both ends. Wire runs must be continuous copper stranded and free from splices.
- Aluminum or solid wires are not allowed.
- Wire sizes given are for use of equipment. Larger sizes may be required by local codes.
- It is recommended that all wires be color coded, as required in accordance with national and local electrical codes.
- Conduit sizes shall be verified by the architect, electrical engineer or contractor, in accordance with local or national codes.
- Convenience outlets are not illustrated. Their number and location are to be specified by others. Locate at least one convenience outlet close to the system control, the power distribution unit and one on each wall of the procedure room. Use hospital approved outlet or equivalent.
- General room illumination is not illustrated. Caution should be taken to avoid excessive heat from overhead spotlights. Damage can occur to ceiling mounting components and wiring if high wattage bulbs are used. Recommend low wattage bulbs no higher than 75 watts and use dimmer controls (except MR). Do not mount lights directly above areas where ceiling mounted accessories will be parked.
- Routing of cable ductwork, conduits, etc., must run direct as possible otherwise may result in the need for greater than standard cable lengths (refer to the interconnection diagram for maximum usable lengths point to point).
- Conduit turns to have large, sweeping bends with minimum radius in accordance with national and local electrical codes.
- A special grounding system is required in all procedure rooms by some national and local codes. It is recommended in areas where patients might be examined or treated under present, future, or emergency conditions. Consult the governing electrical code and confer with appropriate customer administrative personnel to determine the areas requiring this type of grounding system.
- The maximum point to point distances illustrated on this drawing must not be exceeded.
- Physical connection of primary power to GE equipment is to be made by customers electrical contractor with the supervision of a GE representative. The GE representative would be required to identify the physical connection location, and insure proper handling of GE equipment.
- GEHC conducts power audits to verify quality of power being delivered to the system. The customer's electrical contractor is required to be available to support this activity.

- All junction boxes, conduit, duct, duct dividers, switches, circuit breakers, cable tray, etc., are to be supplied and installed by customers electrical contractor.
- Conduit and duct runs shall have sweep radius bends
- Conduits and duct above ceiling or below finished floor must be installed as near to ceiling or floor as possible to reduce run length.
- Ceiling mounted junction boxes illustrated on this plan must be installed flush with finished ceiling.
- All ductwork must meet the following requirements:
  1. Ductwork shall be metal with dividers and have removable covers.
  2. Ductwork shall be certified/rated for electrical power purposes.
  3. Ductwork shall be electrically and mechanically bonded together in an approved manner.
  4. PVC as a substitute must be used in accordance with all local and national codes.
- All openings in access flooring are to be cut out and finished off with grommet material by the customers contractor.
- General contractor to insert pull cords for all cable run conduits between the equipment room and the operators control room.
- 10 foot pigtails at all junction points.
- Grounding is critical to equipment function and patient safety. Site must conform to wiring specifications shown on this plan.



ELECTRICAL LAYOUT ITEM LIST	
1	One 3" stubbed conduit
2	Flush 12"x4"x4" box for Control
3	Flush pull box for Chest Unit - size per local code
4	Flush pull box for Access Unit - size per local code
5	Flush pull box for TB - size per local code
6	Flush pull box for Generator - size per local code
7	18" x 3 1/2" Flush vertical wall duct with minimum 2 dividers
8	Box above ceiling size per local code
9	Box below floor size per local code
10	One 1" cnd
11	One 1 1/2" cnd
12	One 2" cnd
13	One 2 1/2" cnd
14	Main Disconnect Panel

ITEM		Outlet Legend for GE Equipment	
☑	System emergency off (SEO), (recommended height 1.2m [48"] above floor)		
⊗	X-Ray room warning light control panel		
⊕	X-Ray ON lamp (L1) - 24V		
⊗	Door interlock switch (needed only if required by state/local codes)		
⊕	Duplex hospital grade, dedicated wall outlet 120-v, single phase power		
⊕	Dedicated telephone line(s)		
⊕	Network outlet		


Conduit Legend	
----	Above Ceiling
---	Below Floor

Additional Conduit Runs (Contractor Supplied and Installed)				
From	To	Qty	In.	mm
3 phase power	Main disconnect	1	AS REQ'D	AS REQ'D
Main disconnect	Emergency off	1	1/2"	16
Main disconnect	Systems Cabinet	1	AS REQ'D	AS REQ'D
Warning light	Warning light control	1	1/2"	16
1 phase power		1	AS REQ'D	AS REQ'D
		1	1/2"	16
Systems Cabinet	Door Switch	1	1/2"	16
	Tether Interface Box	2	2"	53
	Access Point	1	1"	27
Operators Console	Tether Interface Box	1	1"	27
	Access Point	1	2"	53

Intermountain Medical Center | OPTIMA XR646 HD (G3) | RAD-M210912-FIN-00-A.DWG | Rev A | Date 13/May/2020 | E1 - Electrical Notes | 12/16

Intermountain Medical Center | OPTIMA XR646 HD (G3) | RAD-M210912-FIN-00-A.DWG | 1/4"-1"0" Rev A | Date 13/May/2020 | E2 - Electrical Layout | 13/16



REV	DATE	DESCRIPTION

VCBO NUMBER: 20370  
DATE: 06/08/2020

**IMC EMERGENCY DEPARTMENT - X-RAY #2 - UPGRADE**  
INTERMOUNTAIN HEALTHCARE  
5125 SOUTH COTTONWOOD STREET, MURRAY UT 84107  
PERMIT SET / BID SET

GE DRAWINGS  
**EE701**

**CONNECTIVITY REQUIREMENTS**

Broadband Connections are necessary during the installation process and going forward to ensure full support from the Engineering Teams for the customers system. Maximum performance and availability for the customers system is maintained and closely monitored during the lifetime of the system. Proactive and reactive maintenance is available utilizing the wide range of digital tools using the connectivity solutions listed below:


- Site-to-Site VPN/GE Solution
- Site-to-Site VPN/Custom Solution
- Connection through Dedicated Service Network
- Internet Access - connectivity for InSite 2.0

The requirements for these connectivity solutions are explained in the broadband solutions catalogue (separate document).

**ELECTRICAL NOTES**

- All wires specified shall be copper stranded, flexible, thermo-plastic, color coded, cut 10 foot long at outlet boxes, duct termination points or stubbed conduit ends. All conductors, power, signal and ground, must be run in a conduit or duct system. Electrical contractor shall ring out and tag all wires at both ends. Wire runs must be continuous copper stranded and free from splices.
- Aluminum or solid wires are not allowed.
- Wire sizes given are for use of equipment. Larger sizes may be required by local codes.
- It is recommended that all wires be color coded, as required in accordance with national and local electrical codes.
- Conduit sizes shall be verified by the architect, electrical engineer or contractor, in accordance with local or national codes.
- Convenience outlets are not illustrated. Their number and location are to be specified by others. Locate at least one convenience outlet close to the system control, the power distribution unit and one on each wall of the procedure room. Use hospital approved outlet or equivalent.
- General room illumination is not illustrated. Caution should be taken to avoid excessive heat from overhead spotlights. Damage can occur to ceiling mounting components and wiring if high wattage bulbs are used. Recommend low wattage bulbs no higher than 75 watts and use dimmer controls (except MR). Do not mount lights directly above areas where ceiling mounted accessories will be parked.
- Routing of cable ductwork, conduits, etc., must run direct as possible otherwise may result in the need for greater than standard cable lengths (refer to the interconnection diagram for maximum usable lengths point to point).
- Conduit turns to have large, sweeping bends with minimum radius in accordance with national and local electrical codes.
- A special grounding system is required in all procedure rooms by some national and local codes. It is recommended in areas where patients might be examined or treated under present, future, or emergency conditions. Consult the governing electrical code and confer with appropriate customer administrative personnel to determine the areas requiring this type of grounding system.
- The maximum point to point distances illustrated on this drawing must not be exceeded.
- Physical connection of primary power to GE equipment is to be made by customers electrical contractor with the supervision of a GE representative. The GE representative would be required to identify the physical connection location, and insure proper handling of GE equipment.
- GEHC conducts power audits to verify quality of power being delivered to the system. The customer's electrical contractor is required to be available to support this activity.

- All junction boxes, conduit, duct, duct dividers, switches, circuit breakers, cable tray, etc., are to be supplied and installed by customers electrical contractor.
- Conduit and duct runs shall have sweep radius bends
- Conduits and duct above ceiling or below finished floor must be installed as near to ceiling or floor as possible to reduce run length.
- Ceiling mounted junction boxes illustrated on this plan must be installed flush with finished ceiling.
- All ductwork must meet the following requirements:
  1. Ductwork shall be metal with dividers and have removable covers.
  2. Ductwork shall be certified/rated for electrical power purposes.
  3. Ductwork shall be electrically and mechanically bonded together in an approved manner.
  4. PVC as a substitute must be used in accordance with all local and national codes.
- All openings in access flooring are to be cut out and finished off with grommet material by the customers contractor.
- General contractor to insert pull cords for all cable run conduits between the equipment room and the operators control room.
- 10 foot pigtails at all junction points.
- Grounding is critical to equipment function and patient safety. Site must conform to wiring specifications shown on this plan.



ELECTRICAL LAYOUT ITEM LIST	
1	One 3" stubbed conduit
2	Flush 12"x4"x4" box for Control
3	Flush pull box for Chest Unit - size per local code
4	Flush pull box for Access Unit - size per local code
5	Flush pull box for TB - size per local code
6	Flush pull box for Generator - size per local code
7	18" x 3 1/2" Flush vertical wall duct with minimum 2 dividers
8	Box above ceiling size per local code
9	Box below floor size per local code
10	One 1



**POWER REQUIREMENTS**

POWER SUPPLY	380/400/415/440/460/480V ±10%, THREE-PHASE + G
FREQUENCIES	50/60Hz ± 3Hz
POWER DEMAND	97KVA
MAXIMUM LINE RESISTANCE PER 2 PHASES (Ohm)	380V : 0.118 / 400V : 0.131 / 415V : 0.138 440V : 0.154 / 480V : 0.185

- Power supply should come into a power distribution box (MDP) containing the protective units and controls.
- The section of the supply cable should be calculated in accordance with its length and the maximum permissible voltage drops.
- There must be discrimination between supply cable protective material at the beginning of the installation (main low-voltage transformer side) and the protective devices in the MDP.

**SUPPLY CHARACTERISTICS**

- Power input must be separated from any others which may generate transients (elevators, air conditioning, radiology rooms equipped with high speed film changers...)
- All equipment (lighting, power outlets, etc...) installed with GE system components must be powered separately.

**GROUND SYSTEM**

- Equipotential : the equipotential link will be by means of an equipotential bar. This equipotential bar should be connected to the protective earth conductors in the ducts of the non GE cableways and to additional equipotential connections linking up all the conducting units in the rooms where GE units are located.

**CABLES**

- Power and cable installation must comply with the distribution diagram below.

- All cables must be isolated and flexible.
- Cable color codes must comply with standards for electrical installation.

Case PDB furnished by GE : The cables for signals and remote control (DLK1, SEO, XRL1...) will go to MDP with a pigtail length of 1.5m [4.9 ft.], and will be connected during installation. Each conductor will be identified and isolated (screw connector).

**CABLEWAYS**

The general rules for laying cableways should meet the conditions laid down in current standards and regulations, with regard to:

- Protecting cables against water (cableways should be waterproof)
- Protecting cables against abnormal temperatures (proximity to heating pipes or ducts)
- Protecting cables against temperature shocks
- Replacing cables (cableways should be large enough for cables to be replaced) metal cableways should be grounded.

**FEEDER TABLE**

MIN. FEEDER WIRE SIZE, AWG OR MCM (sq. mm)/VAC	MINIMUM FEEDER WIRE LENGTH - ft. (m)								
	50 (15)	100 (30)	150 (46)	200 (61)	250 (76)	300 (91)	350 (107)	400 (122)	450 (137)
480 VAC	4 (21)	4 (21)	4 (21)	2 (14)	1 (45)	1/0 (54)	1/0 (54)	2/0 (68)	3/0 (85)



**GENERAL NOTES**

In all cases qualified personnel must verify that the feeder (at the point of take-off) and the run to the GE system meet all the requirements stated in the PIM

For a single unit installation, the minimum transformer size is 112.5 kVA. Synthesized power feed is not acceptable

Grounding conductor will be of the same size as the feeder. This ground will run from the equipment back to the power source/main grounding point and always travel in the same conduit with the feeders.

**POWER DISTRIBUTION**



REV	DATE	DESCRIPTION

VCBO NUMBER: 20370  
DATE: 06/08/2020

**IMC EMERGENCY DEPARTMENT - X-RAY #2 - UPGRADE**  
INTERMOUNTAIN HEALTHCARE  
5125 SOUTH COTTONWOOD STREET, MURRAY UT 84107  
PERMIT SET / BID SET


GE DRAWINGS

### GENERAL SHEET NOTES


- 1 UNLESS NOTED OTHERWISE REMOVE ALL LIGHTING FIXTURES DEVICES AND EQUIPMENT SHOWN DASHED. REMOVE CONDUIT AND WIRING BACK TO PANELBOARD OF ORIGIN OR TO FIRST ACTIVE DEVICE THAT REMAINS.
- 2 SALVAGE ALL LIGHT FIXTURES, TWIST-LOCK RECEPTACLES AND WALLPLATES, CEILING SPEAKERS AND SECURITY AND FIRE ALARM DEVICES TO OWNER. PROTECT SALVAGED EQUIPMENT FROM DAMAGE.
- 3 PRIOR TO SUBMITTING BID, VISIT THE SITE AND FIELD VERIFY THE EXTENT OF ELECTRICAL DEMOLITION WORK TO MEET THE INTENT OF THE BID DOCUMENTS AND INCLUDE ALL COSTS IN BID.
- 4 PRIOR TO REMOVAL OF ANY ELECTRICAL EQUIPMENT OR WIRING, FIELD VERIFY THAT THE EQUIPMENT OR WIRING IS INACTIVE OR NO LONGER IN USE.
- 5 REMOVE ALL DEVICES, RACEWAYS AND WIRING FROM WALLS TO BE REMOVED. WHERE ACTIVE RACEWAYS OCCUR IN WALLS TO BE REMOVED, RE-ROUTE THE RACEWAY WITH ASSOCIATED WIRING TO KEEP THE CIRCUIT OPERATIONAL.
- 6 REMOVE ALL FIRE ALARM DEVICES WHERE EXISTING WALLS AND CEILINGS ARE BEING REMOVED. WITH ASSOCIATED CONDUIT AND WIRING. EXISTING FIRE ALARM DEVICES AND SYSTEM NOT INDICATED FOR REMOVAL SHALL REMAIN ACTIVE THROUGHOUT DEMOLITION AND CONSTRUCTION UNTIL THE NEW SYSTEM IS TESTED AND OPERATIONAL. MAINTAIN ALL CLASS A FIRE ALARM INITIATING AND INDICATING LOOPS WHERE EXISTING DEVICES ARE REMOVED.
- 7 REMOVE ALL ABANDONED RACEWAY, CONDUIT, WIRING AND CABLING WHETHER ABANDONED PREVIOUS TO THIS PROJECT OR AS A RESULT OF THIS PROJECT. NOT ALL ABANDONED ITEMS ARE SHOWN ON THESE PLANS AND FIELD VERIFICATION OF DEMOLITION SCOPE EXTENT IS REQUIRED.
- 8 DEVICES MARKED "RR" ARE TO BE REMOVED AND RELOCATED PER NEW PLANS. EXTEND CIRCUITING AS REQUIRED FOR RELOCATION.
- 9 REFER TO ARCHITECTURAL DRAWINGS FOR REMOVAL OF MOTORS, CONDUIT, CONDUCTOR AND CONTROL WIRING ASSOCIATED WITH EXISTING MOTORIZED DOORS, PARTITIONS AND LIGHTING.
- 10 PROVIDE DEDICATED NEUTRALS FOR ALL BRANCH CIRCUITS.
- 11 PROVIDE NEW TYPED PANEL SCHEDULE FOR ALL PANELS AFFECTED BY PROJECT.
- 12 REFER TO THE GE X-RAY DRAWINGS FOR ADDITIONAL CONTRACTOR RESPONSIBILITIES.

### SHEET KEYNOTES


- 1 CONNECT TO THE EXISTING LIGHTING CIRCUIT THAT PREVIOUSLY FED THE X-RAY IN USE SIGN.
- 2 CONNECT TO THE LIGHTING CIRCUIT THAT PREVIOUSLY FED THIS SPACE.
- 3 PROVIDE A MAG HOLD OPEN FOR THE DOOR, MOUNTED HORIZONTALLY IN THE TOEKICK OF THE CABINETS.




**3 LL-1 CEILING DEMOLITION PLAN**  
SCALE: 1/4" = 1'-0"



**4 LL-1 ELECTRICAL DEMOLITION PLAN**  
SCALE: 1/4" = 1'-0"




**1 LL-1 LIGHTING PLAN**  
SCALE: 1/4" = 1'-0"



**2 LL-1 POWER PLAN**  
SCALE: 1/4" = 1'-0"

**SHEET KEYNOTES**

- GROUNDING CONDUCTOR TO BE THE SAME SIZE AS CURRENT CARRYING CONDUCTORS.



**COPPER CONDUCTOR AND CONDUIT SCHEDULE**


SCHEDULE NUMBER		(E.G.) 5 IG				
SUBSCRIPT (NOTE 5)		IG				
SYM	AMP	CONDUIT SIZE	CONDUCTOR(NOTE 1) QTY SIZE G	IG/HH	SBJ	NOTES
1	20	75	2 12	12 8	2	
2	20	75	3 12	12 8	2,3	
3	20	75	4 12	12 8	2,3	
4	30	75	2 10	10 8	2	
5	30	75	3 10	10 8	2	
6	30	75	4 10	10 8	2	
7	40	1	2 8	8 6	2	
8	40	1	3 8	8 6	2	
9	40	1	4 8	8 6	2	
10	55	1	2 6	6 4	2	
11	55	1	3 6	6 4	2	
12	55	1.25	4 6	6 4	2	
13	70	1	2 4	4 2	2	
14	70	1.25	3 4	4 2	2	
15	70	1.25	4 4	4 2	2	
16	85	1.25	2 3	3 2	2	
17	85	1.25	3 3	3 2	2	
18	85	1.25	4 3	3 2	2	
19	95	1.25	3 2	2 2	2	
20	95	1.50	4 2	2 2	2	
21	130	1.50	3 1	1 2	2	
22	130	1.50	4 1	1 2	2	
23	150	2	3 1/0	6 2 1/0	2	
24	150	2	4 1/0	6 2 1/0	2	
25	175	2	3 2/0	6 2 2/0	2	
26	175	2	4 2/0	6 2 2/0	2	
27	200	2	3 3/0	6 2 2/0	2	
28	200	2.50	4 3/0	6 2 2/0	2	
29	230	2.50	3 4/0	4 2 2/0	2	
30	230	2.50	4 4/0	4 2 2/0	2	
31	255	2.50	3 250	4 1 2/0	2	
32	255	2.50	4 250	4 1 2/0	2	
33	310	3	3 350	3 1/0 3/0	2	
34	310	3	4 350	3 1/0 3/0	2	
35	380	3.50	3 500	3 3/0 3/0	2	
36	380	4	4 500	3 3/0 3/0	2	
37	400	2 EA 2	3 3/0	3 3/0 3/0	2	
38	400	2 EA 2.50	4 3/0	3 3/0 3/0	2	
39	510	2 EA 2.50	3 250	1 4/0 3/0	2	
40	510	2 EA 3	4 250	1 4/0 3/0	2	
41	620	2 EA 3	3 350	1/0 4/0 3/0	2,4	
42	620	2 EA 3	4 350	1/0 4/0 3/0	2,4	
43	760	2 EA 3.50	3 500	1/0 4/0 3/0	2,4	
44	760	2 EA 4	4 500	1/0 4/0 3/0	2,4	
45	855	3 EA 3	3 300	2/0 4/0 3/0	2,4	
46	855	3 EA 3	4 300	2/0 4/0 3/0	2,4	
47	1000	3 EA 3.50	3 400	2/0 4/0 3/0	4	
48	1000	3 EA 3.50	4 400	2/0 4/0 3/0	4	
49	1140	3 EA 4	3 500	3/0 4/0 3/0	4	
50	1140	3 EA 4	4 500	3/0 4/0 3/0	4	
51	1240	4 EA 3	3 350	3/0 4/0 3/0	4	
52	1240	4 EA 3	4 350	3/0 4/0 3/0	4	
53	1675	5 EA 4	4 400	4/0 4/0 4/0	4	
54	2010	6 EA 4	4 400	250 250 250	4	
55	2660	7 EA 4	4 500	350 350 350	4	
56	3040	8 EA 4	4 500	500 500 500	4	
57	4180	11 EA 4	4 500	500 500 500	4	
58		5 EA 4			6	
59		5			6	
60		10 EA 4			6	

REV	DATE	DESCRIPTION


VCBO NUMBER: 20370  
DATE: 06/08/2020

**IMC EMERGENCY DEPARTMENT - X-RAY #2 - UPGRADE**  
 INTERMOUNTAIN HEALTHCARE  
 5125 SOUTH COTTONWOOD STREET, MURRAY UT 84107  
 PERMIT SET / BID SET

- CONDUCTOR AND CONDUIT SCHEDULE NOTES**
- CONDUCTORS SHOWN ARE SHOWN FOR EACH CONDUIT WITH MODIFICATIONS AS NOTED IN NOTE 5. ALL CONDUCTORS SHOWN ARE THWN UNLESS OTHERWISE NOTED.
  - PROVIDE EQUIPMENT GROUND CONDUCTORS PER TABLE 250-122 WHEN CIRCUIT BREAKERS ARE SIZED GREATER THAN AMPERE RATING SHOWN IN TABLE.
  - PROVIDE #10 NEUTRALS FOR MULTIWIRE BRANCH CIRCUITS SERVING COMPUTERS.
  - GROUND CONDUCTOR SHALL BE OMITTED BETWEEN THE UTILITY TRANSFORMER AND THE FIRST OVERCURRENT PROTECTIVE DEVICE.
  - SYMBOL SUBSCRIPTS:
    - "2N": INCLUDE TWO NEUTRAL CONDUCTORS, SIZED AS SCHEDULED FOR PHASED AND NEUTRAL CONDUCTORS.
    - "FG": FULL SIZE GROUND, SIZE EQUIPMENT GROUNDING CONDUCTOR TO BE THE SAME SIZE AS THE PHASE CONDUCTOR.
    - "HH": NEUTRAL CURRENTS EXIST DUE TO HIGH HARMONIC "NONLINEAR" LOADS. CURRENT CARRYING CONDUCTORS DERATED ACCORDINGLY. PROVIDE THE IG/HH SIZE FOR THE EQUIPMENT GROUNDING CONDUCTOR.
    - "IG": INCLUDE IG (INSULATED/ISOLATED GROUND CONDUCTOR) SCHEDULED ALONG WITH GROUND OF EQUIPMENT GROUND CONDUCTOR.
    - "SBJ": SUBSTITUTE "SBJ" CONDUCTOR FOR "G" CONDUCTOR SHOWN, WHICH IS SIZED FOR THE SYSTEM BONDING JUMPER OF THE SEPARATELY DERIVED SYSTEM.
  - RACEWAY ONLY. CONDUCTORS PROVIDED BY UTILITY.



**1 DEMOLITION ONE-LINE**  
SCALE: NTS





**2 NEW WORK ONE-LINE**  
SCALE: NTS

## INTERIOR LIGHTING FIXTURE SCHEDULE

### ABBREVIATIONS

#### MOUNTING

B - BASE  
 C - CEILING  
 F - FLANGE  
 G - GRID  
 P - PENDANT  
 PL - POLE  
 R - RECESSED  
 S - SURFACE  
 W - WALL



#### LUMINAIRE OPTIONS

ARHR - AIR RETURN AND HEAT REJECTION  
 DL - DAMP LOCATION  
 EQC - EARTHQUAKE CLIPS  
 F - FUSING  
 HLD - HINGED AND LATCHED DOOR  
 HS - HOUSE SIDE SHIELD  
 PS - PHOTOCELL SWITCH  
 QRS - QUARTZ RESTRIKE  
 ST - STATIC  
 WG - WIRE GUARD  
 WL - WET LOCATION

#### FINISH

MW - MATTE WHITE  
 BL - BLACK  
 SL - SILVER  
 CL - GOLD  
 CL - CLEAR  
 PW - PAINTED WHITE  
 EA - EXTRUDED ALUMINUM  
 S - STEEL  
 GS - GALVANIZED STEEL  
 C - CAST  
 CBA - COLOR BY ARCHITECT  
 SCBA - STANDARD COLOR BY ARCHITECT  
 CCA - CUSTOM COLOR BY ARCHITECT  
 FS - MEETS FEDERAL STANDARD 209D  
 TP - THERMALLY PROTECTED  
 FL - FLUSH  
 R - REGRESS  
 M - MITERED

#### DIFFUSER/LENS

#A - ACRYLIC #THICK  
 #DA - ACRYLIC #THICK (OPAL)  
 GC - GLASS (CLEAR)  
 GO - GLASS (OPAL)  
 GF - GLASS (FROSTED)  
 SOL - SOFT GLOW LENS  
 HPL - HIGH PERFORMANCE LENS  
 DO - DROP OPAL  
 CGL - CONVEX GLASS LENS  
 S - SATIN LENS

#### NOTES


#### REFLECTOR

OP - NONE/OPEN  
 SP - SPECULAR  
 SS - SEMI-SPECULAR  
 D - DIFFUSE (WHITE ENAMEL)  
 SC - SPECULAR (COLORED)  
 PR - PRISMATIC  
 FDR - FULL DEPTH REFLECTOR  
 DS - DIFFUSE (SEMI SPECULAR) SILVER  
 LI - LOW REDUCENT  
 IR - IRRIDISCENT  
 SL - SILVER  
 GL - GOLD  
 CA - CLEAR ALZAK

### GENERAL NOTES

- PROVIDE UNIT PRICES AND FIXTURE BRAND SELECTED FOR ADD/DELETE CHANGES FOR EACH FIXTURE TYPES SHOWN WITHIN 48 BUSINESS HOURS OF THE BID DATE. FAILURE TO COMPLY WITH THIS REQUIREMENT MAY DISQUALIFY THE PRODUCTS AND EMPOWER THE ENGINEER TO DETERMINE FAIR VALUE FOR FIXTURE AND INSTALLATION CHANGES, WITHOUT FURTHER INPUT FROM THE CONTRACTOR OR INSTALLER.
- CONTRACTOR ALLOWANCE PRICES ARE ACCURATE WHEN THIS JOB WAS SPECIFIED, CONTRACTOR AND ELECTRICAL DISTRIBUTOR SHALL VERIFY THIS ALLOWANCE AND REPORT ANY PROBLEMS TO THE ENGINEER BEFORE THE BID. ALLOWANCE PRICE MAY OR MAY NOT INCLUDE LAMP(S) OR FREIGHT AS NOTED, AND DO NOT INCLUDE ANY TAXES.
- SUBSTITUTIONS AND/OR EQUAL FIXTURES MUST RECEIVE APPROVAL PRIOR TO BIDDING. THEY MUST BE SUBMITTED TO THE ENGINEER NO LESS THAN 2 WEEKS PRIOR TO BID OPENING.
- SAMPLES MUST BE PROVIDED FOR ANY AND ALL FIXTURES UPON A/E REQUEST PRIOR TO RELEASING FIXTURES.
- ALL FIXTURES SHALL BE LISTED AND APPROVED FOR THEIR INTENDED USE AND LOCATION.
- VERIFY THE PROPER MOUNTING KITS OR ACCESSORIES TO FACILITATE INSTALLATION AS SHOWN AT EACH LOCATION ON THE DRAWINGS.
- COMPLY WITH THE "INTERIOR LIGHTING" SECTION OF THE SPECIFICATIONS.
- REFER TO SPECIFICATIONS FOR IMPORTANT TECHNICAL REQUIREMENTS FOR LIGHTING FIXTURES, DRIVERS, AND LAMPS.
- ALL LIGHT FIXTURES TO BE EITHER "DLC" OR "LIGHTING FACTS" LISTED OR TO BE APPROVED BY ARCHITECT/ENGINEER AND OWNER.

ID	DESCRIPTION	NOMINAL SIZE				MOUNTING	TYPE	COLOR TEMP	CRI	DRIVER CONFIGURATION	VOLTAGE	WATTS	FINISH	FIGURE LUMENS	DIFFUSER/LENS	REFLECTOR	OPTIONS	NOTES	MANUFACTURER (CATALOG SERIES)			
		LENGTH	DEPTH	HEIGHT	DIAMETER/APERTURE														OPTION 1	OPTION 2	OPTION 3	
(DX-2)	6" ROUND, RECESSED LED DOWNLIGHT, SEMI-SPECULAR REFLECTOR, WHITE TRIM FINISH	-	-	-	6" - 6"	CR	LED	3500K		0-10V DIMMING (10%)	120/277	23	-	2000						GOTHAM (EVO-3520-6AR-WO-LSS-MVOLT-EZ10-TWR)	LITON (LHALD625COT1-D10/LRALD6SWF151-860-T35)	PORTFOLIO (LD6220D10/EU/68-10208 0356LBW2H HB26)
(UC-5)	24" LED UNDERCABINET LIGHT	2' - 0"	-	-	-		LED	3500K		ELV DIMMING	120/277	8		600						DAY-BRITE (LINGS10E1-28-885-UNV-WHG-DIM)	KENALL (AUCLED-1-MW-1-1L35K-2-4277)	AIREY-THOMPSON (19HC-N-39K-24-2-3-D11)



REV	DATE	DESCRIPTION

VCBO NUMBER: 20370  
 DATE: 06/08/2020

**IMC EMERGENCY DEPARTMENT - X-RAY #2 - UPGRADE**  
 INTERMOUNTAIN HEALTHCARE  
 5125 SOUTH COTTONWOOD STREET, MURRAY UT 84107  
 PERMIT SET / BID SET