

- ### # KEYED NOTES
1. STACK PIPING AND ROUTE CHILLED WATER PIPES ALONG SIDE EXISTING STACKED PIPING.
 2. ROUTE CHILLED WATER PIPING ALONG EXISTING CONCRETE WALL.
 3. MAKE CONNECTIONS TO CHILLER PER MANUFACTURER'S RECOMMENDATIONS. PIPING TO BE COPPER, INSULATED WITH ALUMINUM JACKET RATED FOR EXTERIOR USE.
 4. NEW CHILLER PROVIDED BY MRI SUPPLIER AND INSTALLED BY THIS CONTRACTOR. INSTALL ON 6" CONCRETE PAD AND ANCHOR CHILLER TO MEET SEISMIC REQUIREMENTS FOR 1P = 1.5.
 5. ROUTE CHILLED WATER ABOVE EXISTING. REMOVE, UPDATE, AND REPLACE EXISTING METAL PIPE COVER AS NEEDED.
 6. CHILLED WATER PIPING RISE ALONG WITH EXISTING CHILLED WATER PIPE. SHAFT TO BE ADAPTED TO FIT NEW PIPING, SEE ARCHITECTURAL SHEETS.
 7. FIELD VERIFY AND COORDINATE PIPE ROUTING THROUGH THE CEILING SPACE.
 8. INSTALL CHILLER WITH FRONT, OR AIRFLOW SIDE, TOWARD THE BUILDING AS SHOWN.
 9. MAINTAIN CHILLER CLEARANCES.
 10. MOUNT GLYCOL FEEDS ON AVAILABLE WALL SPACE. FIELD VERIFY.
 11. MOUNT EXPANSION TANK ON AVAILABLE WALL SPACE. FIELD VERIFY.
 12. PUMPS STACKED. VERIFY INSTALL LOCATION ALLOWS ACCESS.
 13. HEAT EXCHANGER TO BE INSTALLED AT THE SAME ELEVATION AS THE CHILLER. SNOW SHIELD SHALL BE EXTENDED TO COVER NEW PIPING AND HEAT EXCHANGER.

1 LL-2 Mechanical Piping Plan
SCALE: 1/4" = 1'-0"




REV	DATE	DESCRIPTION

VCBO NUMBER: 19480
CLIENT NUMBER:
DATE: 07/15/2019

KEYED NOTES

1. FIELD VERIFY AND COORDINATE PIPE ROUTING THROUGH THE CEILING SPACE.
2. CHILLED WATER PIPES TO SERVE MRI EQUIPMENT. COORDINATE CONNECTION AND ROUTING WITH SITE SPECIFIC MRI DRAWING. SEE DETAIL 5M501 FOR ALL COMPONENTS.
3. PIPING DOWN TO LL2, SEE SHEET MP101 FOR CONTINUATION.
4. COORDINATE VFD LOCATION WITH AVAILABLE WALL SPACE WITH NEW MRI EQUIPMENT LAYOUT.




1 - LL1 Mechanical Piping Plan
SCALE: 1/4" = 1'-0"

REV	DATE	DESCRIPTION


VCBO NUMBER: 19480
CLIENT NUMBER:
DATE: 07/15/2019

KEYED NOTES

1. REMOVE EXISTING 3/4" RPPB AND ASSOCIATED PIPING SERVING PREVIOUS MRI.
2. CAP REMOVED DOMESTIC COLD WATER LINE.



1 LL1 PLUMBING DEMOLITION PLAN
SCALE: 1/4" = 1'-0"




REV	DATE	DESCRIPTION


VCBO NUMBER: 19480
CLIENT NUMBER:
DATE: 07/15/2019

KEYED NOTES

1. DISCHARGE FROM EQUIPMENT ROOM ABOVE. COORDINATE EXACT LOCATION WITH EQUIPMENT LAYOUT IN EQUIPMENT ROOM.
2. DROP LINE TO EXISTING FLOOR SINK AND TERMINATE WITH AIR GAP AT FLOOR SINK.
3. MECHANICAL EQUIPMENT. SEE SHEET MP101.



1 LL-2 PLUMBING PLAN
SCALE: 1/4" = 1'-0"




REV	DATE	DESCRIPTION

VCBO NUMBER: 19480
CLIENT NUMBER:
DATE: 07/15/2019

KEYED NOTES

1. EXTEND WATER LINE TO HEAT EXCHANGER PER MRI SITE SPECIFIC DRAWINGS. SEE DETAIL 5/M501.
2. DISCHARGE LINE FROM HEAT EXCHANGER. DROP LINE DOWN TO FLOOR BELOW AND RUN TO FLOOR SINK. SEE SHEET P101 FOR CONTINUATION.
3. CONNECT TO EXISTING 2" DOMESTIC COLD WATER LINE. FIELD VERIFY EXISTING CONDITIONS AND LOCATION.




1 LL1 PLUMBING PLAN
SCALE: 1/4" = 1'-0"

REV	DATE	DESCRIPTION



VCBO NUMBER: 19480
CLIENT NUMBER:
DATE: 07/15/2019

KEYED NOTES

1. REMOVE EXISTING MV, MA, AND O2 BRANCH LINES SERVING MRI ROOM BACK TO MAIN.



1 LL1 MEDICAL GAS DEMOLITION PLAN
SCALE: 1/4" = 1'-0"




REV	DATE	DESCRIPTION



VCBO NUMBER: 19480
CLIENT NUMBER:
DATE: 07/15/2019

KEYED NOTES

- ROUTE NEW MV, MA, AND O2 BRANCH LINES SERVING MRI ROOM THROUGH NEW PENETRATION PANEL. TERMINATE AT NEW MEDICAL GAS OUTLET, SEE ARCHITECTURAL ELEVATIONS.



1 LL1 MEDICAL GAS PLAN
SCALE: 1/4" = 1'-0"




REV	DATE	DESCRIPTION



VCBO NUMBER: 19480
CLIENT NUMBER:
DATE: 07/15/2019

KEYED NOTES

1. EXISTING SINGLE INTERLOCK PREACTION SYSTEM VALVE LOCATION.
2. DEMOLISH ALL FIRE PROTECTION PIPING AND HEADS WITHIN THE MRI ROOM.




1 LL1 FIRE PROTECTION DEMOLITION PLAN
SCALE: 1/4" = 1'-0"




REV	DATE	DESCRIPTION

VCBO NUMBER: 19480
CLIENT NUMBER:
DATE: 07/15/2019



- # KEYED NOTES
1. PREACTION PIPING SHALL BE GALVANIZED, TYPICAL. PROVIDE COPPER PIPING AND HARDWARE FOR MRI.
 2. AREA PROTECTED BY EXISTING SINGLE INTERLOCK. PREACTION SYSTEM WITH ELECTRIC ACTIVATION. COORDINATE HEAD PLACEMENT WITH ARCHITECTURAL REFLECTED CEILING PLANS AND MRI SPECIFIC SITE DRAWINGS. COORDINATE PIPE ROUTING WITH OTHER CEILING ELEMENTS.
 3. CONNECT TO EXISTING FIRE PIPING AND ROUTE THROUGH PENETRATION PANEL.

1 LL1 FIRE PROTECTION PLAN
SCALE: 1/4" = 1'-0"




VAN BOERUM & FRANK ASSOCIATES, INC.
CONSULTING ENGINEERS
WWW.VBFA.COM
Salt Lake City, Logan - St. George - Temple
330 South 300 East 801.530.3148 T
Salt Lake City, UT 84111 801.530.3150 F
VBFA Project Number: 19303

REV	DATE	DESCRIPTION


VCBO NUMBER: 19480
CLIENT NUMBER:
DATE: 07/15/2019

IMED BUILDING 5 - MRI CARING SUITE
INTERMOUNTAIN HEALTHCARE
5125 SOUTH COTTONWOOD STREET, MURRAY UT 84107
BID SET / CONSTRUCTION DOCUMENTS


E2 RECEPTACLE MOUNTING DETAILS
SCALE: NTS




D2 ADA DETAIL
SCALE: NTS




D3 SWITCH MOUNTING DETAILS
SCALE: NTS




C1 CMU DEVICE MOUNTING ALIGNMENT DETAIL
SCALE: NTS




C2 LIGHTING MOUNTING DETAILS
SCALE: NTS



C3 COMMUNICATIONS MOUNTING DETAILS
SCALE: NTS




B1 FIRE ALARM MOUNTING DETAILS
SCALE: NTS



A1 BOX MOUNTING DETAILS
SCALE: NTS



A2 TYPICAL WALL MOUNTED DEVICES ALIGNMENT DETAIL
SCALE: NTS



GENERAL SHEET NOTES

1. DETERMINE MOUNTING HEIGHTS OF ELECTRICAL AND ELECTRONIC EQUIPMENT IN THE FOLLOWING ORDER OF PRIORITY:
1 - ELEVATIONS (ARCHITECTURAL, ELECTRICAL, MECHANICAL, ETC.)
2 - EQUIPMENT SHOP DRAWINGS.
3 - FIELD INSTRUCTIONS.
2. LOCATE RECEPTACLES SERVING THE SAME TYPE OF USE AT A UNIFORM HEIGHT UNLESS DIRECTED OTHERWISE.
3. MECHANICAL, ELECTRICAL, AND COMMUNICATION ROOMS: COORDINATE LOCATION OF LIGHTING AND POWER RECEPTACLES WITH EQUIPMENT, PIPING, AND DUCTWORK. DO NOT INSTALL RECEPTACLES BEHIND EQUIPMENT OR WHERE OTHERWISE INACCESSIBLE. POSITION LIGHTING REGARDLESS OF WHERE SHOWN ON DRAWING TO PROVIDE PROPER ILLUMINATION.
4. MOUNT RECEPTACLE BOXES FOR SWITCHES AND RECEPTACLES WITH LONG AXIS OF THE DEVICE VERTICAL UNLESS OTHERWISE INDICATED.
5. SET BOXES WITH PLASTER RINGS FLUSH WITH FINISHED SURFACE.
6. LOCATE BOX COVERS OR DEVICE PLATES SO THEY WILL NOT SPAN DIFFERENT TYPES OF BUILDING FINISHES EITHER VERTICALLY OR HORIZONTALLY.
7. VERIFY ALL DOOR CONDITIONS ON ARCHITECTURAL DRAWINGS PRIOR TO INSTALLING SWITCHES.
8. LOCATE WIRING DEVICES WHICH ARE ADJACENT AND ARE COMPATIBLE VOLTAGES IN ONE PLATE.
9. WHERE DEVICES ARE LOCATED IN CLOSE PROXIMITY OF THE SAME VERTICAL PLANE, ALIGN DEVICES VERTICALLY PER THE TYPICAL WALL MOUNTED DEVICES ALIGNMENT DETAIL, UNLESS OTHERWISE INDICATED.

SHEET KEYNOTES


1. LOCATE RECEPTACLES BEHIND DRINKING FOUNTAINS.
2. REFER TO ARCHITECTURAL ELEVATIONS FOR PLACEMENT OF OUTLETS.
3. LOCATE AT BOTTOM OF BEAMS (OR JOISTS) OR AT CEILING. (REDUCE SPACING BY 5 PERPENDICULAR TO BEAM OR JOIST DIRECTION.) FOR OTHER CONDITIONS, REFER TO NFPA 72.
4. LOCATE DETECTOR ANYWHERE IN SHADED AREA BUT NOT IN TOP 4" OF PEAK.
5. LOCATE AT BOTTOM OF BEAMS IF DIM < 1 OR W/H < 4. OTHERWISE, LOCATE IN BEAM POCKET. FOR D > 4 REDUCE SPACING .33 PERPENDICULAR TO BEAMS.

REV	DATE	DESCRIPTION


VCBO NUMBER:	19480
CLIENT NUMBER:	00000
DATE:	07/15/2019

REV	DATE	DESCRIPTION


VCBO NUMBER: 19480
CLIENT NUMBER: 00000
DATE: 07/15/2019




C1 TYPICAL FIRE STOP FOR CABLES/CONDUIT THROUGH CONCRETE FLOORING
SCALE: NTS



C2 TYPICAL FIRE STOP FOR CABLES/CONDUIT THROUGH CONCRETE WALLS
SCALE: NTS




C3 FIRE STOP FOR METAL CONDUIT THROUGH GYPSUM WALL BOARD
SCALE: NTS




- NOTES:
1. TYPICAL FOR WOOD AND METAL STUD ROUGH-IN.
 2. PLASTER RINGS NOT SHOWN.
 3. LOCATE ALL OUTLET BOXES IN ACCORDANCE WITH ARCHITECTURAL AND MECHANICAL DRAWINGS AND WITH ALL APPLICABLE SHOP DRAWINGS.
 4. IN ACCORDANCE WITH IBC 714.3.2 EXCEPTION 1, OUTLETS ON OPPOSITE SIDES OF WALLS OR PARTITIONS IN THE SAME STUD SPACE IN A RATED FIRE SEPARATION WALL MUST BE SEPARATED BY A MINIMUM OF 24" HORIZONTAL DISTANCE OR LISTED. SOUND AND FIRE RATED PUTTY PADS SHALL BE USED ON THE OUTLET BOXES.
 5. IN NON-RATED WALLS, OUTLETS ON OPPOSITE SIDES OF WALLS OR PARTITIONS MUST BE SEPARATED BY 16" FOR SOUND ATTENUATION.


A1 TYPICAL ROUGH-IN REQUIREMENTS DETAIL
SCALE: NTS



NOTE:
TIE WIRE SHALL NOT BE USED AS A COMPONENT OF ANY RACEWAY HANGER SYSTEM.



A2 TYPICAL RACEWAY SUPPORT METHODS DETAIL
SCALE: NTS



A4 RECESSED FIXTURE MOUNTING DETAIL
SCALE: 1/8" = 1'-0"

LIGHTING REQUIREMENTS

- All lighting fixtures and associated components must meet all RF shielded room and RF grounding requirements (e.g., track lighting is not recommended due to possible RF noise).
- All lighting must use direct current (the DC must have less than 5% ripple).
- 300 lux must be provided at the front of the magnet for patient access and above the magnet for servicing.
- Fluorescent lighting must not be used in the magnet room.
- Lighting must be adjusted using a discrete switch or a variable DC lighting controller.
- Sr dimmers or rheostats must not be used.
- DC led lighting may be used if the power source is located outside the magnet room RF.
- Battery chargers (e.g., used for emergency lighting) must be located outside the magnet RF room.
- Short filament length bulbs are recommended.
- Linear lamps are not recommended due to the high burnout rate.

CONNECTIVITY REQUIREMENTS

Broadband Connections are necessary during the installation process and going forward to ensure full support from the Engineering Teams for the customers system. Maximum performance and availability for the customers system is maintained and closely monitored during the lifetime of the system. Proactive and reactive maintenance is available utilising the wide range of digital tools using the connectivity solutions listed below:


- Site-to-Site VPN/GE Solution
- Site-to-Site VPN/Customer Solution
- Connection through Dedicated Service Network
- Internet Access - connectivity for InSite 2.0

The requirements for these connectivity solutions are explained in the broadband solutions catalogue (separate document).

ELECTRICAL NOTES

- All wires specified shall be copper stranded, flexible, thermo-plastic, color coded, cut 10 foot long at outlet boxes, duct termination points or stubbed conduit ends. All conductors, power, signal and ground, must be run in a conduit or duct system. Electrical contractor shall ring out and tag all wires at both ends. Wire runs must be continuous copper stranded and free from splices.
 - Aluminum or solid wires are not allowed.
 - Wire sizes given are for use of equipment. Larger sizes may be required by local codes.
 - It is recommended that all wires be color coded, as required in accordance with national and local electrical codes.
- Conduit sizes shall be verified by the architect, electrical engineer or contractor, in accordance with local or national codes.
- Convenience outlets are not illustrated. Their number and location are to be specified by others. Locate at least one convenience outlet close to the system control, the power distribution unit and one on each wall of the procedure room. Use hospital approved outlet or equivalent.
- General room illumination is not illustrated. Caution should be taken to avoid excessive heat from overhead spotlights. Damage can occur to ceiling mounting components and wiring if high wattage bulbs are used. Recommend low wattage bulbs no higher than 75 watts and use dimmer controls (except mr). Do not mount lights directly above areas where ceiling mounted accessories will be parked.
- Routing of cable ductwork, conduits, etc., must run direct as possible otherwise may result in the need for greater than standard cable lengths (refer to the interconnection diagram for maximum usable lengths point to point).
- Conduit turns to have large, sweeping bends with minimum radius in accordance with national and local electrical codes.
- A special grounding system is required in all procedure rooms by some national and local codes. It is recommended in areas where patients might be examined or treated under present, future, or emergency conditions. Consult the governing electrical code and confer with appropriate customer administrative personnel to determine the areas requiring this type of grounding system.
- The maximum point to point distances illustrated on this drawing must not be exceeded.
- Physical connection of primary power to GE equipment is to be made by customers electrical contractor with the supervision of a GE representative. The GE representative would be required to identify the physical connection location, and insure proper handling of GE equipment.
- GENC conducts power audits to verify quality of power being delivered to the system. The customer's electrical contractor is required to be available to support this activity.

- All junction boxes, conduit, duct, duct dividers, switches, circuit breakers, cable tray, etc., are to be supplied and installed by customers electrical contractor.
- Conduit and duct runs shall have sweep radius bends.
- Conduits and duct above ceiling or below finished floor must be installed as near to ceiling or floor as possible to reduce run length.
- Ceiling mounted junction boxes illustrated on this plan must be installed flush with finished ceiling.
- All ductwork must meet the following requirements:
 - Ductwork shall be metal with dividers and have removable, accessible covers.
 - Ductwork shall be certified/rated for electrical power purposes.
 - Ductwork shall be electrically and mechanically bonded together in an approved manner.
 - PVC as a substitute must be used in accordance with all local and national codes.
- All openings in access flooring are to be cut out and finished off with grommet material by the customers contractor.
- General contractor to insert pull cords for all cable run conduits between the equipment room and the operators control room.
- 10 foot signals at all junction points.
- Grounding is critical to equipment function and patient safety. Site must conform to wiring specifications shown on this plan.







ITEM	DESCRIPTION (CONTRACTOR SUPPLIED & INSTALLED)
1	Cable ladder 18" x 6"
2	Cable ladder 18" x 6" for gradient cables
3	Non-ferrous cable ladder 18" x 6"
4	Non-ferrous cable ladder 18" x 6" for gradient cables
5	Non-ferrous unistrut cable support
6	Existing 12" x 8" x 6" Junction box
7	Existing 4" x 4" x 4" Junction box
8	Main disconnect panel (reuse existing location)
9	Existing 2" conduit above RF screen
10	Existing 3" conduit above RF screen

ITEM	Outlet Legend for GE Equipment
⚡	System emergency off (SEO), (recommended height 48" above floor)
⚡	Door interlock switch (needed only if required by state/local codes)
⚡	Emergency exhaust fan switch, (recommended height 48" above floor)
⚡	Duplex hospital grade, dedicated wall outlet 120-v, single phase power
⚡	Network outlet
⚡	Dedicated telephone lines/network connection
⚡	Duplex hospital grade, dedicated outlet 120-v emergency, single phase power, 15a
⚡	Duplex hospital grade, dedicated outlet 120-v, single phase outlet routed through RF filter

EXISTING ELECTRICAL NOTE:
USE EXISTING DUCT/CONDUITS WHERE POSSIBLE. ADDITIONAL DUCT/CONDUIT RUNS MAY BE NECESSARY IF EXISTING SYSTEM IS INADEQUATE IN SIZE AND/OR LOCATION FOR THIS INSTALLATION. VERIFY EXISTING SIZE AND LOCATION.

Additional Conduit Runs (Contractor Supplied and Installed)				
From	To	Qty	Size (in)	Size (mm)
Main Disconnect Panel	Facility power	1	as Req'd	
	Integrated Systems Cabinet	1	as Req'd	
Main Disconnect Panel	Integrated Cooling Cabinet	1	as Req'd	
	System emergency off	1	1/2	16
System emergency off	Cooling Penetration Panel	1	1/2	16
Door Switch	Integrated Systems Cabinet	1	3/4	20
System emergency off	Cooling Penetration Panel	1	3/4	20
Magnet Run Down Unit	Magnet	1	1	25
	RF filter	1	as Req'd	
	RF filter	1	as Req'd	
Room Light	RF filter	1	as Req'd	
	RF filter	1	as Req'd	
	Facility emergency power	1	as Req'd	



- Cable Tray Requirements (Side-By-Side)**
- Ceiling Floor
 - Finished Floor
 - Magnet Isocenter
 - Minimum cable tray height required at back of Magnet: 2581 mm [101.5 in].
Tray height may be lower at other points to avoid obstructions.
 - Maximum height from floor to top of tray (anywhere in Magnet room): 3251 mm [128 in].
 - Minimum distance from top of cable tray to ceiling: 254 mm [10 in].
Minimum distance on either side of an obstruction: 254 mm [10 in], minimum distance from top of cable tray to Obstruction: 178 mm [7 in].
 - Tray end to iso-center: 1245 ±12 mm [49 ±0.5 in].
 - Other cable termination to iso-center: 864 ±12 mm [34 ±0.5 in].
 - Minimum distance between trays: 12 mm [0.5 in].
 - Non-ferrous cable support
 - The center of the gradient cable group is 89 mm [3.5 in] from the inside edge of the tray, in line with the magnet center.

POWER REQUIREMENTS

INTRODUCTION
The system requires two independent power inputs:
 • main power supply
 • uninterrupted power supply

SPECIFICATIONS OF MAIN POWER INPUT

POWER SUPPLY	380/400/415/480V +/-10%, THREE-PHASE + G
FREQUENCIES	50/60Hz ± 3Hz
MAXIMUM INPUT POWER (50 msec MAX)	349kVA
INSTALLED LOAD	183kVA
STAND-BY POWER	< 17 kVA

- Power input must be separated from any others which may generate transients (elevators, air conditioning, radiology rooms equipped with high speed film changers...).
- Total harmonic distortion less than 2.5%.
- Phase imbalance must not exceed 2%.

SPECIFICATIONS OF EMERGENCY POWER REQUIREMENT

Magnet Monitor	
POWER INPUT	EMERGENCY LIFE-SAFETY POWER, SINGLE PHASE + G
POWER DEMAND	2.0 A
VOLTAGE	110/220
FREQUENCY	50/60Hz
POWER INPUT	Shield Cryocooler Compressor
POWER DEMAND	AC 380, 400, 415, 460, 480 3Ø
POWER DEMAND	Minimum 9kVA / Recommended 12kVA
FREQUENCY	50/60Hz ± 3Hz


CABLES

- Power and cable installation must comply with the distribution diagram.
- Size of the MDP power input cable is determined by the customer, taking its length and admissible voltage drops into consideration.
- All cables must be isolated and flexible, cable color codes must comply with standards for electrical installation.
- The cables from signaling and remote control (Y,EO#,L...) will go to MDP with a pigtail length of 1.5m, and will be connected during installation.
- Each conductor will be identified and isolated (screw connector).

GROUND SYSTEM

- The equipotential link will be by means of an equipotential bar.
- The grounding point of MDP is directly connected to the building's ground by an isolated copper cable.

POWER DISTRIBUTION



NOTE:

- Cryocooler Compressor (CRY) must operate 24 hours per day, 7 days per week.
- Runs E3030 and M3030 are GE supplied cables, all other wiring is customer supplied and installed.
- Two remote Emergency Off Buttons are supplied with the MDP. Emergency Off removes power from all outputs when activated.
- All MDP output circuits drop out on loss of power, the Cryocooler (CRY) circuit will automatically restart upon restoration of power.
- GE MDP Short circuit current rating is 25,000 amperes at 480 VAC.
- GE MDP is NRTL labeled.
- All feeder circuits require ground wires.


Accepts following range of stranded conductors. All wire types, color and stg to be selected in accordance with governing electrical code(s).

Item	Phase		Ground	
	AWG/kcmil	sq mm	AWG/kcmil	sq mm
A	8-350	6-150	6-250	16-120
B	8-350	6-150	6-250	16-120
C	12-3/0	2.5-70		
D	10-14	4-1.5		

MAIN FEEDER CALCULATION

Irms =	Continuous Power
FLA =	Vph-ph * √3
Ipk =	Irms * 1.25
Vloss =	Peak Power
Rtotal =	Vph-ph * √3
	Vloss
	Ipk * √3
Cable Resistance =	Rtotal
	Cable length

INTERCONNECTIONS



CABLES ROUTING

Configuration	Equipment Room	Magnet Room
A	Short	Short
B	Long	Short
C	Short	Long

CABLES ROUTING FOR OPTIONS


OPTION	FROM	TO	CABLE LENGTH m (ft)
MRE	MRE	Magnet Isocenter	Nominal: 7.31 (24) Maximum: 10.06 (33)
	MRE	ISC cabinet	15.24 (50)
	MRE	Ethernet Hub in ISC	15.24 (50)
	MRE	Customer Supplied Outlet	60Hz: 6.10 (20) 50Hz: 7.62 (25)

GENERAL SHEET NOTES


- 1 UNLESS OTHERWISE INDICATED, REMOVE ALL LIGHTING FIXTURES, OUTLETS, DEVICES AND EQUIPMENT IN HATCHED AREAS. REMOVE ASSOCIATED CONDUIT AND WIRING BACK TO THE PANELBOARD OF ORIGIN. SYSTEMATICALLY CHECK EACH BRANCH PANELBOARD CIRCUIT TO VERIFY EACH THAT CIRCUIT BREAKER NO LONGER HAS ANY ACTIVE LOADS. DISCONNECT THE WIRING AND TURN THE CIRCUIT BREAKER OFF. ANY REMAINING ACTIVE LOADS SHALL BE LABELED AT THE PANELBOARD AS TO WHAT LOAD IS SERVED.
- 2 UNLESS NOTED OTHERWISE REMOVE ALL LIGHTING FIXTURES DEVICES AND EQUIPMENT SHOWN DASHED. REMOVE CONDUIT AND WIRING BACK TO PANELBOARD OF ORIGIN OR TO FIRST ACTIVE DEVICE THAT REMAINS.
- 3 SALVAGE ALL POWER POLES, LIGHT FIXTURES, TWIST-LOCK RECEPTACLES AND WALLPLATES, CEILING SPEAKERS AND SECURITY AND FIRE ALARM DEVICES TO OWNER. PROTECT SALVAGED EQUIPMENT FROM DAMAGE.
- 4 PRIOR TO SUBMITTING BID, VISIT THE SITE AND FIELD VERIFY THE EXTENT OF ELECTRICAL DEMOLITION WORK TO MEET THE INTENT OF THE BID DOCUMENTS AND INCLUDE ALL COSTS IN BID.
- 5 PRIOR TO REMOVAL OF ANY ELECTRICAL EQUIPMENT OR WIRING, FIELD VERIFY THAT THE EQUIPMENT OR WIRING IS INACTIVE OR NO LONGER IN USE.
- 6 REMOVE ALL DEVICES, RACEWAYS AND WIRING FROM WALLS TO BE REMOVED. WHERE ACTIVE RACEWAYS OCCUR IN WALLS TO BE REMOVED, REROUTE THE RACEWAY WITH ASSOCIATED WIRING TO KEEP THE CIRCUIT OPERATIONAL.
- 7 REMOVE ALL FIRE ALARM DEVICES WHERE EXISTING WALLS AND CEILINGS ARE BEING REMOVED, WITH ASSOCIATED CONDUIT AND WIRING. EXISTING FIRE ALARM DEVICES AND SYSTEM NOT INDICATED FOR REMOVAL SHALL REMAIN ACTIVE THROUGHOUT DEMOLITION AND CONSTRUCTION UNTIL THE NEW SYSTEM IS TESTED AND OPERATIONAL. MAINTAIN ALL CLASS A FIRE ALARM INITIATING AND INDICATING LOOPS WHERE EXISTING DEVICES ARE REMOVED.
- 8 REMOVE ALL ABANDONED RACEWAY, CONDUIT, WIRING AND CABLING WHETHER ABANDONED PREVIOUS TO THIS PROJECT OR AS A RESULT OF THIS PROJECT. NOT ALL ABANDONED ITEMS ARE SHOWN ON THESE PLANS AND FIELD VERIFICATION OF DEMOLITION SCOPE EXTENT IS REQUIRED.
- 9 DEVICES MARKED "RR" ARE TO BE REMOVED AND RELOCATED PER NEW PLANS. EXTEND CIRCUITING AS REQUIRED FOR RELOCATION.

SHEET KEYNOTES

- 1 REMOVE EXISTING LIGHTING CONTROLS FOR MRI 3T EXAM ROOM #2.
- 2 REMOVE EXISTING DISCONNECT PANEL FOR MRI.



A2 LOWER LEVEL 1 ELECTRICAL DEMOLITION PLAN
SCALE: 1/4" = 1'-0"



A4 LOWER LEVEL 1 CEILING DEMOLITION PLAN
SCALE: 1/4" = 1'-0"



GENERAL SHEET NOTES

- 1 ALL RECEPTACLES LOCATED WITHIN 6" OF THE EDGE OF A SINK MUST BE GFCI PROTECTED.
- 2 PROVIDE A DEDICATED NEUTRAL FOR ALL BRANCH CIRCUITS.
- 3 PROVIDE NEW TYPED PANEL SCHEDULES FOR ALL PANELS AFFECTED DURING CONSTRUCTION.
- 4 REFER TO GE AND PDC DRAWINGS FOR ADDITIONAL CONTRACTOR RESPONSIBILITIES.

SHEET KEYNOTES




REV	DATE	DESCRIPTION

VCBO NUMBER: 19480
 CLIENT NUMBER: 00000
 DATE: 07/15/2019

A1 LOWER LEVEL 1 OVERALL POWER PLAN
 SCALE: 1/8" = 1'-0"

IMED BUILDING 5 - MRI CARING SUITE
 INTERMOUNTAIN HEALTHCARE
 5125 SOUTH COTTONWOOD STREET, MURRAY UT 84107
 BID SET / CONSTRUCTION DOCUMENTS




GENERAL SHEET NOTES

- 1 ALL RECEPTACLES LOCATED WITHIN 6" OF THE EDGE OF A SINK MUST BE GFCI PROTECTED.
- 2 PROVIDE A DEDICATED NEUTRAL FOR ALL BRANCH CIRCUITS.
- 3 PROVIDE NEW TYPED PANEL SCHEDULES FOR ALL PANELS AFFECTED DURING CONSTRUCTION.
- 4 REFER TO GE AND PDC DRAWINGS FOR ADDITIONAL CONTRACTOR RESPONSIBILITIES.

SHEET KEYNOTES

- 1 PROVIDE 200/3 FUSED DISCONNECT AND A FRS-125 FUSE WITH A NEMA 3R ENCLOSURE MOUNTED ADJACENT TO EQUIPMENT.



19480
00000
07/15/2019


A1 LOWER LEVEL 2 OVERALL POWER PLAN
SCALE: 3/32" = 1'-0"

GENERAL SHEET NOTES


- 1 ALL RECEPTACLES LOCATED WITHIN 6' OF THE EDGE OF A SINK MUST BE GFCI PROTECTED.
- 2 PROVIDE A DEDICATED NEUTRAL FOR ALL BRANCH CIRCUITS.
- 3 PROVIDE NEW TYPED PANEL SCHEDULES FOR ALL PANELS AFFECTED DURING CONSTRUCTION.
- 4 REFER TO GE AND PDC DRAWINGS FOR ADDITIONAL CONTRACTOR RESPONSIBILITIES.

SHEET KEYNOTES

- 1 CIRCUIT NEW FOUR-PLEX TO THE SAME CIRCUIT THAT FED THE FOUR-PLEX THAT WAS DEMOLISHED.
- 2 PROVIDE 120V CIRCUIT TO RF DOOR CONTROL BOX.
- 3 PROVIDE UBS8200XX WALL OUTLETS.
- 4 PROVIDE LEVITON 1081 OR EQUIVALENT WALL SWITCHES.



A1 LOWER LEVEL 1 POWER PLAN
SCALE: 1/4" = 1'-0"



19480
00000
07/15/2019

IMED BUILDING 5 - MRI CARING SUITE

INTERMOUNTAIN HEALTHCARE
5125 SOUTH COTTONWOOD STREET, MURRAY UT 84107

BID SET / CONSTRUCTION DOCUMENTS

LOWER LEVEL 1 POWER PLAN

EP111

SHEET KEYNOTES

- GROUNDING CONDUCTOR TO BE THE SAME SIZE AS CURRENT CARRYING CONDUCTORS.
- MRI MDP WILL BE FURNISHED WITH MRI EQUIPMENT, AND INSTALLED BY ELECTRICAL CONTRACTOR.
- PROVIDE NEW CIRCUIT BREAKER IN EXISTING GE DISTRIBUTION PANEL.

BRANCH CIRCUIT CONDUCTOR AND CONDUIT SIZING TABLE

CIRCUIT AMPACITY/VOLTAGE	CIRCUIT LENGTH	CONDUCTOR SIZE (PHASE, NEUTRAL AND GR)	CONDUIT SIZE
20A/120V	0' - 60'	#12 AWG	0.75" Ø
20A/120V	60' - 95'	#10 AWG	0.75" Ø
20A/120V	95' - 150'	#8 AWG	1" Ø
20A/120V	150' - 240'	#6 AWG	1.25" Ø
20A/277V	0' - 140'	#12 AWG	0.75" Ø
20A/277V	140' - 220'	#10 AWG	0.75" Ø
20A/277V	220' - 350'	#8 AWG	1" Ø
20A/277V	350' - 550'	#6 AWG	1.25" Ø

NOTES:

- WIRE SIZING IS BASED ON COPPER CONDUCTORS SUPPLYING A 20A, 120V CIRCUIT AT THE INDICATED VOLTAGE, ASSUMED TO BE 80% LOADED (16A), WITH MAXIMUM VOLTAGE DROP OF 3% AT THE LOAD.
- DOWN-SIZED WIRE AT DEVICE/LOAD AS REQUIRED AND TERMINATE CONDUCTORS IN A SAFE AND CODE COMPLIANT MANNER.
- CONDUIT SIZE IS BASED ON A MAXIMUM OF 3 CIRCUITS PER CONDUIT, EACH WITH A SEPARATE NEUTRAL CONDUCTOR.



COPPER CONDUCTOR AND CONDUIT SCHEDULE

SCHEDULE NUMBER: (E.G.) 5 IG
SUBSCRIPT (NOTE 5)

SYM	AMP	CONDUIT SIZE	CONDUCTOR(NOTE 1)	QTY	SIZE	IG/HH	SBJ	NOTES
1	20	.75	2	12	12	12	8	2
2	20	.75	3	12	12	12	8	2,3
3	20	.75	4	12	12	12	8	2,3
4	30	.75	2	10	10	10	8	2
5	30	.75	3	10	10	10	8	2
6	30	.75	4	10	10	10	8	2
7	40	1	2	8	10	8	6	2
8	40	1	3	8	10	8	6	2
9	40	1	4	8	10	8	6	2
10	55	1	2	6	10	8	4	2
11	55	1	3	6	10	8	4	2
12	55	1.25	4	6	10	8	4	2
13	70	1	2	4	8	4	2	2
14	70	1.25	3	4	8	4	2	2
15	70	1.25	4	4	8	4	2	2
16	85	1.25	2	3	8	3	2	2
17	85	1.25	3	3	8	3	2	2
18	85	1.25	4	3	8	3	2	2
19	95	1.25	3	2	8	2	2	2
20	95	1.50	4	2	8	2	2	2
21	130	1.50	3	1	6	2	2	2
22	130	1.50	4	1	6	2	2	2
23	150	2	3	1/0	6	2	1/0	2
24	150	2	4	1/0	6	2	1/0	2
25	175	2	3	2/0	6	2	2/0	2
26	175	2	4	2/0	6	2	2/0	2
27	200	2	3	3/0	6	2	2/0	2
28	200	2.50	4	3/0	6	2	2/0	2
29	230	2.50	3	4/0	4	2	2/0	2
30	230	2.50	4	4/0	4	2	2/0	2
31	255	2.50	3	250	4	1	2/0	2
32	255	2.50	4	250	4	1	2/0	2
33	310	3	3	350	3	1/0	3/0	2
34	310	3	4	350	3	1/0	3/0	2
35	380	3.50	3	500	3	3/0	3/0	2
36	380	4	4	500	3	3/0	3/0	2
37	400	2 EA 2	3	3/0	3	3/0	3/0	2
38	400	2 EA 2.50	4	3/0	3	3/0	3/0	2
39	510	2 EA 2.50	3	250	1	4/0	3/0	2
40	510	2 EA 3	4	250	1	4/0	3/0	2
41	620	2 EA 3	3	350	1/0	4/0	3/0	2,4
42	620	2 EA 3	4	350	1/0	4/0	3/0	2,4
43	760	2 EA 3.50	3	500	1/0	4/0	3/0	2,4
44	760	2 EA 4	4	500	1/0	4/0	3/0	2,4
45	855	3 EA 3	3	300	2/0	4/0	3/0	2,4
46	855	3 EA 3	4	300	2/0	4/0	3/0	2,4
47	1000	3 EA 3.50	3	400	2/0	4/0	3/0	4
48	1000	3 EA 3.50	4	400	2/0	4/0	3/0	4
49	1140	3 EA 4	3	500	3/0	4/0	3/0	4
50	1140	3 EA 4	4	500	3/0	4/0	3/0	4
51	1240	4 EA 3	3	350	3/0	4/0	3/0	4
52	1240	4 EA 3	4	350	3/0	4/0	3/0	4
53	1675	5 EA 4	4	400	4/0	4/0	4/0	4
54	2010	6 EA 4	4	400	250	250	250	4
55	2660	7 EA 4	4	500	350	350	350	4
56	3040	8 EA 4	4	500	500	500	500	4
57	4180	11 EA 4	4	500	500	500	500	4
58		5 EA 4						6
59		5						6
60		10 EA 4						6

CONDUCTOR AND CONDUIT SCHEDULE NOTES

- CONDUCTORS SHOWN ARE SHOWN FOR EACH CONDUIT WITH MODIFICATIONS AS NOTED IN NOTE 5. ALL CONDUCTORS SHOWN ARE THWN UNLESS OTHERWISE NOTED.
- PROVIDE EQUIPMENT GROUND CONDUCTORS PER TABLE 250-122 WHEN CIRCUIT BREAKERS ARE SIZED GREATER THAN AMPERE RATING SHOWN IN TABLE.
- PROVIDE #10 NEUTRALS FOR MULTI-WIRE BRANCH CIRCUITS SERVING COMPUTERS.
- GROUND CONDUCTOR SHALL BE OMITTED BETWEEN THE UTILITY TRANSFORMER AND THE FIRST OVERCURRENT PROTECTIVE DEVICE.
- SYMBOL SUBSCRIPTS:
 - "2N": INCLUDE TWO NEUTRAL CONDUCTORS, SIZED AS SCHEDULED FOR PHASE AND NEUTRAL CONDUCTORS.
 - "FG": FULL SIZE GROUND, SIZE EQUIPMENT GROUNDING CONDUCTOR TO BE THE SAME SIZE AS THE PHASE CONDUCTORS.
 - "HH": NEUTRAL CURRENTS EXIST DUE TO HIGH HARMONIC "NON-LINEAR" LOADS. CURRENT CARRYING CONDUCTORS DERATED ACCORDINGLY. PROVIDE THE IG/HH SIZE FOR THE EQUIPMENT GROUNDING CONDUCTOR.
 - "IG": INCLUDE IG (INSULATED/ISOLATED GROUND CONDUCTOR) SCHEDULED ALONG WITH GROUND OF EQUIPMENT GROUND CONDUCTOR.
 - "SBJ": SUBSTITUTE "SBJ" CONDUCTOR FOR "G" CONDUCTOR SHOWN, WHICH IS SIZED FOR THE SYSTEM BONDING JUMPER OF THE SEPARATELY DERIVED SYSTEM.
- RACEWAY ONLY. CONDUCTORS PROVIDED BY UTILITY.



1 DEMOLITION ONE-LINE
SCALE: NTS

2 NEW WORK ONE-LINE
SCALE: NTS



REV	DATE	DESCRIPTION

VCBO NUMBER: 19480
CLIENT NUMBER: 00000
DATE: 07/15/2019

IMED BUILDING 5 - MRI CARING SUITE
INTERMOUNTAIN HEALTHCARE
5125 SOUTH COTTONWOOD STREET, MURRAY UT 84107
BID SET / CONSTRUCTION DOCUMENTS




GENERAL SHEET NOTES

SHEET KEYNOTES

1 LIGHTING TO BE PROVIDED BY PDC AND INSTALLED BY CONTRACTOR. REFER TO PDC DRAWINGS.




19480
00000
07/15/2019



A4 LOWER LEVEL 1 LIGHTING PLAN
SCALE: 1/4" = 1'-0"

IMED BUILDING 5 - MRI CARING SUITE
INTERMOUNTAIN HEALTHCARE
5125 SOUTH COTTONWOOD STREET, MURRAY UT 84107
BID SET / CONSTRUCTION DOCUMENTS

LOWER LEVEL 1 LIGHTING PLAN



A4 LOWER LEVEL 1 AUXILIARY PLAN
SCALE: 1/4" = 1'-0"

GENERAL SHEET NOTES

SHEET KEYNOTES



REV	DATE	DESCRIPTION

VCBO NUMBER: 19480
CLIENT NUMBER: 00000
DATE: 07/15/2019