

Asept.2X Gold – Mobile UV Rapid Disinfectors



The ASEPT.2X Mobile UV is a mobile rapid disinfectant that provides operator initiated closed room UVC disinfection.

Placement of 1-2 units within the room to be disinfected, allows for simultaneous sterilization of a room from all angles, minimizing the problematic shadow areas left by the conventional single unit UV sterilizers.

User-friendly remote Wi-Fi controls which can be used on any smart device.

The ASEPT.2X units used in pairs simultaneously will sterilize up to 99.9999% of a 22' x 22' room (7mx7m) in 10 minutes. Shorter cycle of disinfection (5) minutes will apply for smaller room dimensions.

View product video at: <https://www.youtube.com/watch?v=VafgVMJa0YY&t=27s>

Specifications

Disinfection Time	5-10 min per room
Safety System	4 passive infrared sensors, E-stop
Power Requirements	120-240V AC, 50-60Hz, 7.5A per tower
Overall Dimensions	712 mm W x 712 mm L x 1645 mm H (28.0" W x 28.0" L x 64.8" H)
Weight	45.5 kg (100lbs)
Communication	Wi-Fi
UV Wavelength	254 nm

Features/Deliverables

- UVC lamps (8): T6 quartz 40 in (101.6 cm) Teflon encapsulated, mercury vapour
- Aluminum construction and medical grade stainless steel
- Multiple infrared motion detectors
- Fan-based lamp temperature management
- 360° handle for easy maneuverability
- Heavy-duty lockable swivel casters
- Retractable medical grade cord
- Scalable software with updates
- Chronological data collection
- Integrated WEB server
- User-friendly WEB interface (user/administrator), upgradable and customizable

Benefits

- Automated disinfection with minimal human intervention
- Simultaneous operation of a single unit for quicker disinfection
- 99.99% disinfection in 5 minutes
- Adjustable disinfection time
- Paired units to maximize disinfection dose and time


Pricing

Item	Description	Qty	Price (USD)
1	ASEPT.2X Includes: 1 tower, Android 7in tablet configured for use, Caution floor sign/coupler and (2) door caution panels	1	
2	UVC LED Floor Disinfection Module	1	Upon Request
Optional Accessories			
3	uv cense Dosimeter (UVC254)	1	
4	Intelligo UVC Indicators - pkg of 100		
	0-100 mJ/cm ²	1	
	0-1000 mJ/cm ²	1	
	Safety Indicators (0-6 mJ/cm ²)	1	
5	Aranet4 Pro Air Quality Monitor	1	

Terms and Conditions

Terms:	50% with Purchase Order 30% at Shipping 20% Net 30
Warranty:	All parts are warranted against manufacturer's defects for a period of 12 months from date of installation, not including wear and tear parts.
Taxes:	All taxes extra
F.O.B.	Prescient ^x Warehouse – Cambridge, Ontario, Canada
Returns:	Stock Items: RGA Required; Restocking Charge Applies; Non-Stock Items: No Returns
Manufactured Items:	Cancellation of orders will result in a 10% penalty through submittal stage. After approved submittals 20% penalty until manufacturing begins. After manufacturing begins, cancellations will not be permitted.
Prices:	Quotation valid for 90 days.

Sanuvox EPA

United States				OMB Control No. 2070-0078 Expires on 1/31/2020	
 ENVIRONMENTAL PROTECTION AGENCY Washington, DC 20460 Pesticide Report for Pesticide-Producing and Device-Producing Establishments Office of Enforcement and Compliance Assurance <small>Section 7, Federal Insecticide, Fungicide, and Rodenticide Act (7 U.S.C. sec. 136e)</small> http://www.epa.gov/compliance/pesticide-establishment-registration-and-reporting					
1. Establishment Name SANUVOX TECHNOLOGIES INC.			2. EPA Est. No. 91660-CAN-1		
3. Establishment Site Address 146 BARR					
4. City SAINT LAURENT				5. State	
6. Province/Region QUEBEC PROVINCE		7. Country CAN		8. Zip/Postal Code H4T 1Y4	
9. Establishment Mailing Address ATTN: KEITH JORDAN, 12725 Outlook Ave					
10. City FORT WORTH				11. State TX	
12. Province/Region		13. Country USA		14. Zip/Postal Code 76244	
15. Telephone Number 8887268869		16. Email mec@sanuvox.com			
17. Check if your establishment's name or address changed <input checked="" type="checkbox"/>			18. Did you produce or distribute in 2019? <input checked="" type="checkbox"/> Yes <input type="checkbox"/> No		
19. Company Name SANUVOX TECHNOLOGIES INC.					
20. Company Site Address 12725 Outlook Ave					
21. City Forth Worth				22. State TX	
23. Province/Region QUEBEC PROVINCE		24. Country USA		25. Zip/Postal Code 76244	
26. Company Mailing Address 12725 Outlook Ave					
27. City Forth Worth				28. State TX	
29. Province/Region		30. Country USA		31. Zip/Postal Code 76244	
32. Telephone Number (817) 938-6843		33. Email kjordan@sanuvox.com			
34. Check if your company name or address changed <input type="checkbox"/>			35. Check if your Authorized Agent changed <input type="checkbox"/>		
36. Name of Company Official or Authorized Agent Vincent Gariepy			37. Title of Company Official or Authorized Agent Authorized Agent		
38. Signature of Company Official or Authorized Agent (Original signature required) _aa2a03a4-4ae7-434b-8fa4-6343bb1b29cf				39. Date 02/27/2020	
40. Telephone Number of Company Official or Authorized Agent 5143825823		41. Email Of Company Official or Authorized Agent mec@sanuvox.com			
EPA Reviewer	EPA Office	Post Mark Date	Date Entered	Comments	Reporting Year 2019

EPA Form 3540-16 (Rev. 11/19) Previous versions are obsolete.

1

ASEPT.2X CE/EC Certificate of Conformity

	
CE EC Certificate of Conformity	
Date of issue: 2017-01-09	
Sanuvox Technologies Inc. hereby declares that testing has been completed and report has been generated for:	
Product Name	: ASEPT.2X Mobile Sterilization Unit
Model Number	: ASEPT.2X
Applicant	: Sanuvox Technologies Inc. 146 Barr St-Laurent, QC H4T 1Y4 Canada
And, in accordance to the following applicable Directive: 2014/35/EC, 2011/65/EU, 2014/30/EU	
That this product has been tested according to the following Standard: EN 61010-1:2010, EN55011:2009+A1:2010 Class B, Group 1, IEC/EN 60598-2-1 (ballast), IEC/EN 61347-2-3 (ballast)	
Therefore, Sanuvox Technologies Inc. hereby acknowledges that the manufacturer may issue a DECLARATION OF CONFORMITY and apply the CE mark in accordance to European Union Rules	
Approved by:	Benoit Despatis, P. Eng
	<u>Benoit Despatis, P.Eng</u> <u>2017-01-09</u>
Signature	Date
The statement is based on a single evaluation of one sample of the above product.	



STUDY REPORT

STUDY TITLE

Evaluation of Antimicrobial Effectiveness of a UVC Generating Device
on Hard Nonporous Surfaces

Test Organisms:

Methicillin Resistant *Staphylococcus aureus* - MRSA (ATCC 33592)
Clostridium difficile - spore form (ATCC 43598)
Vancomycin Resistant *Enterococcus faecalis* – VRE (ATCC 51575)

PRODUCT IDENTITY

Aseptix

AUTHOR

Joshua Luedtke, M.S.
Study Director

STUDY COMPLETION DATE

January 17, 2014

PERFORMING LABORATORY

ATS Labs
1285 Corporate Center Drive, Suite 110
Eagan, MN 55121

SPONSOR

Sanuvox Technologies Inc.
146, rue Barr
Saint-Laurent, QC H4T 1Y4
Canada

PROJECT NUMBER

A15984

Page 1 of 13

Project No. A15984
TRF Number: SXT01121613.CUST

Sanuvax Technologies Inc.
Page 2 of 13



STUDY REPORT

GENERAL STUDY INFORMATION

Study Title: Evaluation of Antimicrobial Effectiveness of a UVC Generating Device on Hard Nonporous Surfaces
Project Number: A15984
TRF Number: SXT01121613.CUST

TEST SUBSTANCE IDENTITY

Test Device Name: Aseptix

STUDY DATES

Date Sample Received: December 27, 2013
Study Initiation Date: December 27, 2013
Experimental Start Date: January 8, 2014
Experimental End Date: January 13, 2014
Study Completion Date: January 17, 2014

Test Organism	ATCC #	Culture Medium	Incubation Parameters
Methicillin Resistant <i>Staphylococcus aureus</i> - MRSA	33592	Synthetic Broth	35-37°C, aerobic
<i>Clostridium difficile</i> - spore form	43598	CDC Anaerobic Blood Agar	35-37°C, anaerobic
Vancomycin Resistant <i>Enterococcus faecalis</i> - VRE	51575	Fluid Thioglycollate Medium	35-37°C, aerobic

The test organisms to be used in this study were obtained from the American Type Culture Collection (ATCC), Manassas, Virginia.

Exposure Times: 5 minutes, 10 minutes and 15 minutes
Exposure Temperature: Room temperature (20.3°C)
Number of Carriers Tested: 2 per organism per location
Soil Load Description: No organic soil load required
Neutralizer: Lethen Broth
Agar Plate Medium: Tryptic Soy Agar with 5% Sheep Blood (BAP) [for MRSA and VRE]
BHI-HT Agar [for *Clostridium difficile* - spore form]

Project No. A15984
TRF Number: SXT01121613.CUST

Sanuvox Technologies Inc.
Page 3 of 13



EXPERIMENTAL DESIGN

Glass carriers (1" x 3") inoculated with a dry film of the test organism were placed into the testing room and exposed to the UV generating device(s) for the Sponsor specified exposure times. Duplicate carriers per organism per location were placed around the testing area as indicated by the Sponsor. Briefly, two carriers per organism were placed on the bedrail of a hospital bed approximately 3 feet off the floor and approximately 2 feet from a UV device. The second set of carriers was located on the opposite side of the hospital bed, on a hospital table approximately 4 feet off the ground and approximately 5 feet from a UV device. Each test carrier was oriented so that the inoculated area of the carrier was perpendicular to the ground, parallel to the device and as vertical as possible. After exposure, the carriers were transferred to vessels containing subculture media and assayed for survivors. Appropriate culture purity, media sterility, carrier sterility, carrier quantitation, HCl resistance (for *Clostridium difficile*) and neutralization confirmation controls were performed.

Per Sponsor's direction, the study was not required to be conducted under US EPA 40 CFR Part 160 or US FDA 21 CFR Part 58.

STUDY RESULTS

TABLE 1: CONTROL RESULTS

The following results from controls confirmed study validity:

Type of Control	Results		
	Methicillin Resistant <i>Staphylococcus aureus</i> - MRSA (ATCC 33592)	<i>Clostridium difficile</i> - spore form (ATCC 43598)	Vancomycin Resistant <i>Enterococcus faecalis</i> - VRE (ATCC 51575)
Purity Control	Pure	Pure	Pure
Neutralizer Sterility Control	No Growth		
Carrier Sterility Control	No Growth		

Project No. A15984

TRF Number: SXT01121613.CUST

Sanuvox Technologies Inc.

Page 4 of 13



TABLE 2: NEUTRALIZATION CONFIRMATION CONTROL RESULTS

Test Device	Test Organism	Neutralization Confirmation (CFU/plate)		±1.0 Log ₁₀ Pass/Fail
		Numbers Control	Results	
Aseptix	Methicillin Resistant <i>Staphylococcus aureus</i> - MRSA (ATCC 33592)	20,13	21,14	-0.03 (Pass)
	<i>Clostridium difficile</i> - spore form (ATCC 43598)	40,36	40,34	0.01 (Pass)
	Vancomycin Resistant <i>Enterococcus faecalis</i> - VRE (ATCC 51575)	17,11	15,15	-0.03 (Pass)

CFU = Colony Forming Unit

TABLE 3: CARRIER QUANTITATION CONTROL RESULTS

Test Organism	Carrier #	Result	Average Log ₁₀	Geometric Mean
		CFU/Carrier (Log ₁₀)		
Methicillin Resistant <i>Staphylococcus aureus</i> - MRSA (ATCC 33592)	1	1.5 x 10 ⁶ (6.18)	6.18	1.51 x 10 ⁶
	2	1.5 x 10 ⁶ (6.18)		
<i>Clostridium difficile</i> - spore form (ATCC 43598)	1	9.0 x 10 ⁶ (6.95)	6.95	8.91 x 10 ⁶
	2	9.0 x 10 ⁶ (6.95)		
Vancomycin Resistant <i>Enterococcus faecalis</i> - VRE (ATCC 51575)	1	1.7 x 10 ⁵ (5.23)	5.17	1.48 x 10 ⁵
	2	1.3 x 10 ⁵ (5.11)		

CFU = Colony Forming Unit

Project No. A15984

TRF Number: SXT01121613.CUST

Sanuvox Technologies Inc.

Page 5 of 13



TABLE 4: EVALUATION OF TEST CARRIER DATA – 5 Minute Exposure

Test Device: Aseptix							
Carrier Location: Bedrail (approximately 3 feet off ground and 2 feet from device)							
Test Organism	Carrier #	Number of Survivors (CFU)					
		Dilution					
		Filtered 10 ⁰	10 ⁰ (1.00 mL)	10 ⁰	10 ⁻¹	10 ⁻²	10 ⁻³
Methicillin Resistant <i>Staphylococcus aureus</i> - MRSA (ATCC 33592)	1	0	0,0	0,0	0,0	0,0	0,0
	2	0	0,0	0,0	0,0	0,0	0,0
<i>Clostridium difficile</i> – spore form (ATCC 43598)	1	3	0,0	0,0	0,0	0,0	0,0
	2	52	3,8	0,0	0,0	0,0	0,0
Vancomycin Resistant <i>Enterococcus faecalis</i> – VRE (ATCC 51575)	1	0	0,0	0,0	0,0	0,0	0,0
	2	0	0,0	0,0	0,0	0,0	0,0
Carrier Location: Table (approximately 4 feet off ground and 5 feet from device)							
Test Organism	Carrier #	Number of Survivors (CFU)					
		Dilution					
		Filtered 10 ⁰	10 ⁰ (1.00 mL)	10 ⁰	10 ⁻¹	10 ⁻²	10 ⁻³
Methicillin Resistant <i>Staphylococcus aureus</i> - MRSA (ATCC 33592)	1	0	0,0	0,0	0,0	0,0	0,0
	2	0	0,0	0,0	0,0	0,0	0,0
<i>Clostridium difficile</i> – spore form (ATCC 43598)	1	TNTC	30,20	3,2	0,0	0,0	0,0
	2	TNTC	46,40	3,9	0,0	0,0	0,0
Vancomycin Resistant <i>Enterococcus faecalis</i> – VRE (ATCC 51575)	1	14	1,1	1,0	0,0	0,0	0,0
	2	2	1,0	0,0	0,0	0,0	0,0

CFU = Colony Forming Unit

TNTC = Too Numerous To Count

A value of <1 was used in place of zero for calculation purposes only.

Project No. A15984
TRF Number: SXT01121613.CUST

Sanuvax Technologies Inc. **ATS LABS**
Page 6 of 13

TABLE 5: EVALUATION OF TEST CARRIER DATA – 10 Minute Exposure

Test Device: Aseptix							
Carrier Location: Bedrail (approximately 3 feet off ground and 2 feet from device)							
Test Organism	Carrier #	Number of Survivors (CFU)					
		Dilution					
		Filtered 10⁰	10⁰ (1.00 mL)	10⁰	10⁻¹	10⁻²	10⁻³
Methicillin Resistant <i>Staphylococcus aureus</i> - MRSA (ATCC 33592)	1	0	0,0	0,0	0,0	0,0	0,0
	2	0	0,0	0,0	0,0	0,0	0,0
<i>Clostridium difficile</i> – spore form (ATCC 43598)	1	1	0,0	0,0	0,0	0,0	0,0
	2	15	0,0	0,0	0,0	0,0	0,0
Vancomycin Resistant <i>Enterococcus faecalis</i> – VRE (ATCC 51575)	1	0	0,0	0,0	0,0	0,0	0,0
	2	0	0,0	0,0	0,0	0,0	0,0
Carrier Location: Table (approximately 4 feet off ground and 5 feet from device)							
Test Organism	Carrier #	Number of Survivors (CFU)					
		Dilution					
		Filtered 10⁰	10⁰ (1.00 mL)	10⁰	10⁻¹	10⁻²	10⁻³
Methicillin Resistant <i>Staphylococcus aureus</i> - MRSA (ATCC 33592)	1	0	0,0	0,0	0,0	0,0	0,0
	2	0	0,0	0,0	0,0	0,0	0,0
<i>Clostridium difficile</i> – spore form (ATCC 43598)	1	80	1,8	0,0	0,0	0,0	0,0
	2	52	2,5	0,0	0,0	0,0	0,0
Vancomycin Resistant <i>Enterococcus faecalis</i> – VRE (ATCC 51575)	1	39	3,2	0,0	0,0	0,0	0,0
	2	TNTC	36,46	1,0	0,0	0,0	0,0

CFU = Colony Forming Unit

TNTC = Too Numerous To Count

A value of <1 was used in place of zero for calculation purposes only.

Project No. A15984
TRF Number: SXT01121613.CUST

Sanuvox Technologies Inc.
Page 7 of 13



TABLE 6: EVALUATION OF TEST CARRIER DATA – 15 Minute Exposure

Test Device: Aseptix							
Carrier Location: Bedrail (approximately 3 feet off ground and 2 feet from device)							
Test Organism	Carrier #	Number of Survivors (CFU)					
		Dilution					
		Filtered 10 ⁰	10 ⁰ (1.00 mL)	10 ⁰	10 ⁻¹	10 ⁻²	10 ⁻³
Methicillin Resistant <i>Staphylococcus aureus</i> - MRSA (ATCC 33592)	1	0	0,0	0,0	0,0	0,0	0,0
	2	0	0,0	0,0	0,0	0,0	0,0
<i>Clostridium difficile</i> – spore form (ATCC 43598)	1	0	0,0	0,0	0,0	0,0	0,0
	2	0	0,0	0,0	0,0	0,0	0,0
Vancomycin Resistant <i>Enterococcus faecalis</i> – VRE (ATCC 51575)	1	0	0,0	0,0	0,0	0,0	0,0
	2	0	0,0	0,0	0,0	0,0	0,0
Carrier Location: Table (approximately 4 feet off ground and 5 feet from device)							
Test Organism	Carrier #	Number of Survivors (CFU)					
		Dilution					
		Filtered 10 ⁰	10 ⁰ (1.00 mL)	10 ⁰	10 ⁻¹	10 ⁻²	10 ⁻³
Methicillin Resistant <i>Staphylococcus aureus</i> - MRSA (ATCC 33592)	1	0	0,0	0,0	0,0	0,0	0,0
	2	0	0,0	0,0	0,0	0,0	0,0
<i>Clostridium difficile</i> – spore form (ATCC 43598)	1	6	0,0	0,0	0,0	0,0	0,0
	2	18	0,0	0,0	0,0	0,0	0,0
Vancomycin Resistant <i>Enterococcus faecalis</i> – VRE (ATCC 51575)	1	0	0,0	0,0	0,0	0,0	0,0
	2	4	0,0	0,0	0,0	0,0	0,0

CFU = Colony Forming Unit

A value of <1 was used in place of zero for calculation purposes only.

Project No. A15984
TRF Number: SXT01121613.CUST

Sanuvox Technologies Inc.
Page 8 of 13



TABLE 7: CALCULATED VALUES – 5 Minute Exposure

Test Device: Aseptix					
Carrier Location: Bedrail (approximately 3 feet off ground and 2 feet from device)					
Test Organism	Carrier #	# Survivors/ Carrier (Log ₁₀)	Average Log ₁₀	Geometric Mean	Percent Reduction (Log ₁₀)
Methicillin Resistant <i>Staphylococcus aureus</i> - MRSA (ATCC 33592)	1	<1 (<0.00)	<0.00	<1	>99.9999% (>6.18)
	2	<1 (<0.00)			
<i>Clostridium difficile</i> - spore form (ATCC 43598)	1	4 (0.60)	1.20	1.58 x 10 ¹	>99.999% (5.75)
	2	6.2 x 10 ¹ (1.79)			
Vancomycin Resistant <i>Enterococcus faecalis</i> – VRE (ATCC 51575)	1	<1 (<0.00)	<0.00	<1	>99.999% (>5.17)
	2	<1 (<0.00)			
Carrier Location: Table (approximately 4 feet off ground and 5 feet from device)					
Test Organism	Carrier #	# Survivors/ Carrier (Log ₁₀)	Average Log ₁₀	Geometric Mean	Percent Reduction (Log ₁₀)
Methicillin Resistant <i>Staphylococcus aureus</i> - MRSA (ATCC 33592)	1	<1 (<0.00)	<0.00	<1	>99.9999% (>6.18)
	2	<1 (<0.00)			
<i>Clostridium difficile</i> - spore form (ATCC 43598)	1	5.0 x 10 ² (2.70)	2.82	6.60 x 10 ²	99.99% (4.13)
	2	8.6 x 10 ² (2.93)			
Vancomycin Resistant <i>Enterococcus faecalis</i> – VRE (ATCC 51575)	1	1.7 x 10 ¹ (1.23)	0.77	5.89 x 10 ⁰	>99.99% (4.40)
	2	2 (0.30)			

Project No. A15984

TRF Number: SXT01121613.CUST

Sanuvox Technologies Inc.

Page 9 of 13



TABLE 8: CALCULATED VALUES – 10 Minute Exposure

Test Device: Aseptix					
Carrier Location: Bedrail (approximately 3 feet off ground and 2 feet from device)					
Test Organism	Carrier #	# Survivors/ Carrier (Log ₁₀)	Average Log ₁₀	Geometric Mean	Percent Reduction (Log ₁₀)
Methicillin Resistant <i>Staphylococcus aureus</i> - MRSA (ATCC 33592)	1	<1 (<0.00)	<0.00	<1	>99.9999% (>6.18)
	2	<1 (<0.00)			
<i>Clostridium difficile</i> - spore form (ATCC 43598)	1	1 (0.00)	0.63	4.27 x 10 ⁰	>99.9999% (6.32)
	2	1.8 x 10 ¹ (1.26)			
Vancomycin Resistant <i>Enterococcus faecalis</i> – VRE (ATCC 51575)	1	<1 (<0.00)	<0.00	<1	>99.999% (>5.17)
	2	<1 (<0.00)			
Carrier Location: Table (approximately 4 feet off ground and 5 feet from device)					
Test Organism	Carrier #	# Survivors/ Carrier (Log ₁₀)	Average Log ₁₀	Geometric Mean	Percent Reduction (Log ₁₀)
Methicillin Resistant <i>Staphylococcus aureus</i> - MRSA (ATCC 33592)	1	<1 (<0.00)	<0.00	<1	>99.9999% (>6.18)
	2	<1 (<0.00)			
<i>Clostridium difficile</i> - spore form (ATCC 43598)	1	9.5 x 10 ¹ (1.98)	1.89	7.76 x 10 ¹	99.999% (5.06)
	2	6.2 x 10 ¹ (1.79)			
Vancomycin Resistant <i>Enterococcus faecalis</i> – VRE (ATCC 51575)	1	4.6 x 10 ¹ (1.66)	2.29	1.95 x 10 ²	>99.8% (2.88)
	2	8.2 x 10 ² (2.91)			

Project No. A15984

TRF Number: SXT01121613.CUST

Sanuvax Technologies Inc.

Page 10 of 13



TABLE 9: CALCULATED VALUES – 15 Minute Exposure

Test Device: Aseptix					
Carrier Location: Bedrail (approximately 3 feet off ground and 2 feet from device)					
Test Organism	Carrier #	# Survivors/ Carrier (Log ₁₀)	Average Log ₁₀	Geometric Mean	Percent Reduction (Log ₁₀)
Methicillin Resistant <i>Staphylococcus aureus</i> - MRSA (ATCC 33592)	1	<1 (<0.00)	<0.00	<1	>99.9999% (>6.18)
	2	<1 (<0.00)			
<i>Clostridium difficile</i> - spore form (ATCC 43598)	1	<1 (<0.00)	<0.00	<1	>99.9999% (>6.95)
	2	<1 (<0.00)			
Vancomycin Resistant <i>Enterococcus faecalis</i> – VRE (ATCC 51575)	1	<1 (<0.00)	<0.00	<1	>99.999% (>5.17)
	2	<1 (<0.00)			
Carrier Location: Table (approximately 4 feet off ground and 5 feet from device)					
Test Organism	Carrier #	# Survivors/ Carrier (Log ₁₀)	Average Log ₁₀	Geometric Mean	Percent Reduction (Log ₁₀)
Methicillin Resistant <i>Staphylococcus aureus</i> - MRSA (ATCC 33592)	1	<1 (<0.00)	<0.00	<1	>99.9999% (>6.18)
	2	<1 (<0.00)			
<i>Clostridium difficile</i> - spore form (ATCC 43598)	1	7 (0.85)	1.09	1.23 x 10 ¹	>99.999% (5.86)
	2	2.1 x 10 ¹ (1.32)			
Vancomycin Resistant <i>Enterococcus faecalis</i> – VRE (ATCC 51575)	1	<1 (<0.00)	<0.35	<2.24 x 10 ⁰	>99.99% (>4.82)
	2	5 (0.70)			

Project No. A15984

TRF Number: SXT01121613.CUST

Sanuvox Technologies Inc.

Page 11 of 13



TABLE 10: VERIFICATION OF ANTIBIOTIC RESISTANCE – *Staphylococcus aureus* - MRSA

Organism (ATCC)	Zone of Inhibition (mm)	CLSI* Resistant Range (mm)
Methicillin Resistant <i>Staphylococcus aureus</i> – MRSA (ATCC 33592)	6	≤ 10
Quality Control Organism (ATCC)	Zone of Inhibition (mm)	CLSI* Acceptable Range (mm)
<i>Staphylococcus aureus</i> (ATCC 25923)	20	18 - 24

*CLSI = Clinical and Laboratory Standards Institute

TABLE 11: VERIFICATION OF ANTIBIOTIC RESISTANCE – *Enterococcus faecalis* - VRE

Quality Control Organism	Zone of Inhibition (mm)	CLSI* Acceptable Range (mm)
Vancomycin Resistant <i>Enterococcus faecalis</i> – VRE (ATCC 51575)	10	≤14
Test Organism	Zone of Inhibition (mm)	CLSI* Resistant Range (mm)
<i>Staphylococcus aureus</i> (ATCC 25923)	17	17-21

*CLSI = Clinical and Laboratory Standards Institute

Interpretation of result and acceptable range are from the Clinical and Laboratory Standards Institute, Performance Standards for Antimicrobial Susceptibility Testing; Twenty-Second Information Supplement January 2012, Volume 31 Number 1, Approved Standard M02-A11 and M07-A9, Wayne, Pennsylvania.

Project No. A15584

TRF Number: SXT01121613.CUST

Sanuvax Technologies Inc.

Page 12 of 13



TABLE 12: HCL RESISTANCE VERIFICATION

Test Organism: <i>Clostridium difficile</i> – spore form (ATCC 43598)						
Exposure Time	10 ⁻²	10 ⁻³	10 ⁻⁴	CFU/mL (Log ₁₀)	Log ₁₀ Reduction from Control	Pass/Fail (≤ 2 log ₁₀ difference)
5 minutes (test)	184,204	27,27	6,1	1.94 x 10 ⁵ (5.29)	0.94	Not Applicable
10 minutes (test)	56,58	12,17	4,0	5.7 x 10 ⁴ (4.76)	1.47	Pass
20 minutes (test)	4,7	1,3	0,0	6 x 10 ³ (3.78)	2.45	Not Applicable
20 minutes (control)	T,T	142,200	29,37	1.71 x 10 ⁶ (6.23)	Not Applicable	Not Applicable

T = Too Numerous To Count (≥300 colonies)

CONTROL RESULTS

The results of controls run for purity, carrier sterility, neutralizer sterility, neutralization confirmation, HCl resistance control, antibiotic resistance and carrier quantitation were all acceptable.

Project No. A15984
TRF Number: SXT01121613.CUST

Sanuvox Technologies Inc.
Page 13 of 13



ANALYSIS

The UV light generating device, Aseptix, demonstrated a >99.999% (>5.75 log₁₀) reduction, >99.9999% (6.32 log₁₀) and >99.9999% (>6.95 log₁₀) reduction of *Clostridium difficile* – spore form on the test carriers located on the hospital bedrail following a 5, 10 and 15 minute exposure time, respectively, when tested at room temperature (20.3°C).

The UV light generating device, Aseptix, demonstrated a 99.99% (4.13 log₁₀) reduction, 99.999% (5.06 log₁₀) and >99.999% (5.86 log₁₀) reduction of *Clostridium difficile* – spore form on the test carriers located on the hospital table following a 5, 10 and 15 minute exposure time, respectively, when tested at room temperature (20.3°C).

The UV light generating device, Aseptix, demonstrated a >99.9999% (>6.18 log₁₀) reduction of Methicillin Resistant *Staphylococcus aureus* - MRSA on all test carriers following 5, 10 and 15 minute exposure times when tested at room temperature (20.3°C).

The UV light generating device, Aseptix, demonstrated a >99.999% (>5.17 log₁₀) reduction of Vancomycin Resistant *Enterococcus faecalis* – VRE (ATCC 51575) on the test carriers located on the hospital bedrail following 5, 10 and 15 minute exposure times when tested at room temperature (20.3°C).

The UV light generating device, Aseptix, demonstrated a >99.99% (4.40 log₁₀) reduction, >99.8% (2.88 log₁₀) reduction and a >99.99% (>4.82 log₁₀) reduction of Vancomycin Resistant *Enterococcus faecalis* – VRE (ATCC 51575) on the test carriers located on the hospital table following a 5, 10 and 15 minute exposure time, respectively, when tested at room temperature (20.3°C).

PREPARED BY:


Joshua Luedtke, M.S.
Microbiologist

1-17-14
Date

The use of the ATS Labs name, logo or any other representation of ATS Labs without the written approval of ATS Labs is prohibited. In addition, ATS Labs may not be referred to in any form of promotional materials, press releases, advertising or similar materials (whether by print, broadcast, communication or electronic means) without the express written permission of ATS Labs.

ASEPT.2X disinfection efficiency against SARS-CoV-2 virus



Sanuvox Technologies Inc.
146 Barr,
St-Laurent, Qc.,
H4T 1Y4

p. 1.888.726.8869
t. 1.888.582.6475
e. info@sanuvox.com

June 30th 2020,

Object: ASEPT.2X disinfection efficiency against SARS-CoV-2 virus

This letter is to confirm that based on specific test data performed in a P3 level lab, the ASEPT.2X disinfection unit provides at least 99.99% disinfection of SARS-CoV-2 in less than 2 minutes of exposure time within a radius of 3 meters around the unit. This allows a single ASEPT.2X to disinfect a room of 6m x 6m in 2 minutes.

Validations with UVC dosimeter indicators sensitive to 254 nm wavelength has also shown that for a disinfection cycle time of 5 minutes, the reach for the same disinfection level of 4 log extends up to a radius of 5.5 meters. Consequently, all exposed surfaces of a 11 m x 11m room can be effectively disinfected with a 5 minutes operating cycle. All the air inside the room also gets the same disinfection level.

Germicidal UV disinfection uses no chemicals and leaves no harmful residuals behind so that the disinfected room can be accessed immediately after.

For further information, please do not hesitate to email Sanuvox Technologies at info@sanuvox.com

Best regards,

Dr. Normand Brais, P. Eng, M.A.Sc., Ph.D.
VP Engineering & Founder